


FT-1R5090



LumiFlash Infinity Chemiluminescent Substrate, HRP

To detect femtogram of protein by Western Blot

Product Description

Catalog #: 1R5090, 2 x 50 ml 1R5091, 2 x 250 ml
Name : **LumiFlash Infinity Chemiluminescent Substrate, HRP**
Components : Solution A (Substrate Solution) 250mL
 Solution B (Peroxide Solution) 250mL
 User's manual 1 booklet
Storage : Substrate should be stored at 4°C and shielded from light.
 Please use up the product in in 1 year.
 For Research Use Only

Introduction

LumiFlash Infinity Chemiluminescent Substrate, HRP is an extremely sensitive ECL product by supplementing a proprietary luminol for chemiluminescent detection of immobilized proteins (Western blotting), conjugated with Horseradish Peroxidase (HRP) directly or indirectly. In Western blotting application, LumiFlash™ Infinity Chemiluminescent Substrate, HRP provide high signal and low background, which allows detection of protein targets at femtogram levels. This feature benefits the researchers with excellent low background results and without signal burn effects at the same time.

Technical and Scientific Information

- Duration >6 hours
- Detection Method X-ray film or imaging acquisition system
- Typical Antibody dilution Primary: 1:1,000–1:20,000
Secondary: 1:40,000–1:400,000
- Shelf life ≥1 year at time of receipt
- Recommended Initial Exposure Time 30 seconds (to X ray film)

Safety Information

Please wear gloves, lab coat and goggles while operating. Prevent contact product directly. In case of contacting, wash with large amount of water.

Materials needed but not provided

1. PVDF or nitrocellulose membrane.
2. Wash buffer: Phosphate-buffered saline (PBS) or Tris-buffered saline (TBS) containing 0.05–0.1% Tween®-20

P.1

FT-1R5090

PBS: 10 mM sodium phosphate, 150 mM NaCl, pH 7.2

TBS: 10 mM Tris, 150 mM NaCl, pH 7.4

3. Blocking buffer: 1–5% (w/v) blocking agent (e.g., casein, BSA, or gelatin) in wash buffer
4. Specific primary antibody for interested protein, diluted in blocking buffer
5. HRP-conjugated secondary antibody, specific for primary antibody, diluted in blocking buffer
6. X-ray film or chemiluminescence image acquisition systems

Instruction

A. Protein transfer

1. A Perform 1D or 2D electrophoresis for protein separation.
2. Move the electrophoretic gel into appropriate transfer buffer and equilibrate for 10 minutes.
3. Wet the PVDF or nitrocellulose membrane in transfer buffer. (For PVDF membrane, it is necessary to pre-wet it in methanol before moving into transfer buffer).
4. Assemble the transferring sandwich as the order of two filter papers, gel, membrane and two filter papers.
5. Transfer proteins according to blotting apparatus manufacturer's instruction.
6. Wash the transferred membrane of three times by Wash buffer for 5 minutes.

B. Antibody incubation

1. Add BSA, Skim milk based blocking buffer or BlockPRO™ Blocking buffer (BP01-1L) and incubate at room temperature for 30 minutes.
2. Prepare the primary antibody by diluting it with Blocking buffer according to the manufacturer's instruction or previous experience. (Due to the good sensitivity of chemiluminiscent detection, previous antibody dilution factor can be increased 2-5 folds for optimal signal to noise ratio.)
3. Add primary antibody and incubate at room temperature for at least 1 hour with gentle agitation. For more specific interaction between primary and antigen proteins, it is recommended to perform additional incubation at 4°C 8-12hours.
4. Decant the primary antibody solution thoroughly. Wash the membrane of at least three times with ample amount of fresh Wash buffer for 10 minutes.
5. Prepare the secondary antibody by diluting it with Blocking buffer according to the manufacturer's instruction. (Due to the good sensitivity of chemiluminiscent detection, secondary antibody dilution can start from 1:20000).
6. Add secondary antibody and incubate at room temperature for 1 hour with gentle agitation.
7. Decant the secondary antibody solution thoroughly. Wash the membrane of at least four times with ample amount of fresh Wash buffer for 10 minutes.

C. Chemiluminiscent detection

1. To prepare working HRP substrate, mix equal volume of Solution A and Solution B in a clean tube freshly. 0.1 mL of working HRP substrate is sufficient per 1 cm² membrane area.
2. Allow the HRP working substrate to stand at room temperature for 3 minutes with proper light shielding.
3. In the dark room or box, place the membrane side up in a clean box or plastic wrap. Add HRP working substrate onto the membrane.
4. Incubate the membrane at room temperature for 10 seconds.
5. Overlay plastic wrap or a transparency sheet on the wet membrane.

NOTE: Do not use filter papers to overdrain the HRP substrate. It will decrease the signal significantly. Keep the membrane wet while exposing!! (see next step)

6. Expose the membrane to appropriate X-ray film or by chemiluminiscent image acquisition system. It is recommended to use 30 seconds as the initial exposure time.

D. Stripping of PVDF membrane

The immunoblot of PVDF membrane can be stripped of antibodies, and then reprobed.

1. Incubate membrane in stripping buffer (62.5mM Tris-HCl pH 6.8, 100mM β-mercaptoethanol and 2% (w/v) SDS) for 30 minutes at 50-70°C.
2. Wash the membrane twice in Wash buffer for 10 minutes each.
3. To ensure complete removal of antibodies, incubate the membrane with LumiFlash Advance HRP working substrate and expose against X-ray film for 5 minutes. No signal should be observed for complete stripping.

FT-1R5090

Trouble Shooting

Problem	Possible cause	Remedy
No signal or weak signal	Poor transfer efficiency	Optimize the membrane transferring procedure.
	Insufficient antigen	Increase the amount of loaded antigen. Make sure the blot have been store correctly to avoid the degradation of target protein.
	The concentration of primary and secondary antibody is too low	Increase the concentration of the primary and/or the secondary antibody.
	Inappropriate storage/preparation of the ECL detection reagents	Use HRP or HRP conjugates to check the applicability of ECL reagents.
	Too short exposure time	Extend exposure time
Excessive signal	Antigen or antibody excess	Reduce the amount of loaded antigen. Dilute the primary antibody and/or the secondary antibody.
High Background	Antigen or antibody excess	Optimize the condition by reducing the amount of antigen, or the concentration of the primary antibody and/or secondary antibody. Initially, reduce the secondary antibody to 20% of the original usage.
	Inappropriate blocking	Try different blocking substrate such as gelatin, casein, skim milk or casein.
	Inadequate washing	Increase the concentration of Tween-20 in washing solution. Increase the washing steps between the hybridization procedures. Extend washing time.
	Overexposure to film	Shorten the exposure time.

References

- BioMed Research International, Volume 2014, Article ID 692061 (2014)
- Journal of Proteomics, Volume 120, Pages 204–214 (2015)
- Nature Communications 7, Article number: 10347 (2016)

Related products

- LuminolPen, HRP System, 1R7260, 1 pen
- LuminolPen Lite, HRP System, 1R7270, 1 pen

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.
Please inquire for higher quantities (availability, shipment conditions).

Please contact InterBioTech – Interchim for any other information
Hotline : +33(0)4 70 03 73 06 – Interbiotech@interchim.com

Disclaimer : Materials are sold for research use only, and are not intended for food, drug, household, or cosmetic use. Interchim® is not liable for any damage resulting from handling or contact with this product.