



FT-41391A




TMRE, TMRM

Product Information

Name :	TMRE Tetramethylrhodamine Ethyl Ester, perchlorate
Catalog Number :	FP-41391A , 25 mg
Structure :	C ₂₆ H ₂₇ ClN ₂ O ₇
Molecular Weight :	MW= 514.96
Solubility:	Soluble in DMSO, EtOH
Absorption / Emission :	$\lambda_{exc}/\lambda_{em}$ (MeOH) = 549/574 nm
EC (M⁻¹ cm⁻¹) :	109 000



Name :	TMRM Tetramethylrhodamine Methyl Ester, perchlorate
Catalog Number :	FP-21089A , 25 mg
Structure :	C ₂₅ H ₂₅ ClN ₂ O ₇
Molecular Weight :	MW= 500.93
Solubility:	Soluble in DMSO, MetOH
Absorption / Emission :	$\lambda_{exc}/\lambda_{em}$ (MeOH) = 548/574 nm
EC (M⁻¹ cm⁻¹) :	115 000



Storage: -20°C (M) Protect from light and moisture


Introduction

TMRE and TMRM, the methyl and ethyl esters of tetramethylrhodamine, are preferred as rapid membrane potential sensors for quantitative measurements using Nernst equation. They pass through the plasmatic membrane better than related rhodamine 123, and accumulate in mitochondria where fluorescence self quenches. The transmembrane distribution is directly related to the membrane potential, because they do not form aggregates in cell membranes or interact with membrane proteins. They are used to obtain unbiased images of potential-dependent dye distribution. Highly selective, potential-dependent staining of mitochondria is obtained by setting the extracellular K⁺ concentration close to intracellular values (~137 mM), thereby depolarizing the plasmatic membrane.

Their high fluorescence allows to work with a low concentration (avoid aggregation in cell membranes).

FT-41391A

Absorbance and emission spectra of TMR in pH7 buffer



Directions for use

Handling and Storage

TMRM, **TMRE**, were dissolved in methanol or ethanol and used directly. The final concentrations of methanol or ethanol were kept to <0.5% (v/v) in all incubations of mitochondria.

Some datas about use of TMRE/TMRA:

Cells	Source	Concentration	Temperature	Time	Notes
Rat hippocampi	Bindokas 1998	1µM TMRE/TMRM	+23°C	Maintained in dye sol. 100nM	Microfluorometer
Rat cardiomyocytes	Duchen 1998	3µM TMRE	Room temperature	20 min	Epifluorescence
NIH 3T3 fibroblast	Fink 1998	100 nM TMRE	+37°C	10 min	Confocal microscopy

Other protocol may be found in litterature.

Related products

- JC-1, [FP-52314A](#)
- JC-10, superior alternative to JC-1, [FP-CL0440](#)
- Dihydrorhodamine 123, [FP-AM352A](#)
- Valinomycin, [FP-09246B](#)
- CCCP, [091640](#)

References

- **Bindokas V.P.**, *et al.*, « Changes in Mitochondrial Function Resulting from Synaptic Activity in the Rat Hippocampal Slice », *The Journal of Neuroscience*, **18**, 4570(1998) [Article](#)
- **Duchen M.R.**, *et al.*, « Transient Mitochondrial Depolarizations Reflect Focal Sarcoplasmic Reticular Calcium Release in Single Rat Cardiomyocytes », *J. Cell Biol.*, **142**, 975(1998) [Article](#)
- **Fink C.**, *et al.*, « Intracellular Fluorescent Probe Concentrations by Confocal Microscopy », *Biophys J*, **75**, 1648 (1998) [Article](#)
- **Helen Everett**, *et al.*, « M11L: A Novel Mitochondria-localized Protein of Myxoma Virus That Blocks Apoptosis of Infected Leukocytes », *The Journal of Experimental Medicine*, **191**, 1487(2000) [Article](#)
- **Hüser J.**, *et al.*, « Imaging the Permeability Pore Transition in Single Mitochondria », *Biophys J*, **74**, 2129 (1998) [Article](#)
- **Jake Jacobson**, *et al.*, « Mitochondrial oxidative stress and cell death in astrocytes — requirement for stored Ca²⁺ and sustained opening of the permeability transition pore », *Journal of Cell Science* **115**, 1175 (2002)
- **Julian D.** *et al.*, Mitochondrial depolarization following hydrogen sulfide exposure in erythrocytes from a sulfide-tolerant marine invertebrate, *Journal of Experimental Biology* 208, 4109-4122 (2005) [Article](#)
- **O'Reilly C.M.**, *et al.*, « Quantitative Analysis of Spontaneous Mitochondrial Depolarizations », *Biophysical Journal* **85**, 3350 (2003) [Article](#)
- **Düßmann H.**, *et al.*, « Outer mitochondrial membrane permeabilization during apoptosis triggers caspase-independent mitochondrial and caspase-dependent plasma membrane potential depolarization: a single-cell analysis », *Journal of Cell Science*, **116**, 525 (2003) [Article](#)
- **Parra V.** *et al.*, Changes in mitochondrial dynamics during ceramide-induced cardiomyocyte early apoptosis, *Cardiovasc Res*, 77: 387 - 397 (2008) [Article](#)
- **Reichert S.A.**, *et al.*, « The Mitochondrial Permeability Transition Pore and Nitric Oxide Synthase Mediate Early Mitochondrial Depolarization in Astrocytes during Oxygen-Glucose Deprivation », *The Journal of Neuroscience*, **21**, 6608(2001) [Article](#)

FT-41391A

- **Scaduto R.C., et al.**, « Measurement of Mitochondrial Membrane Potential Using Fluorescent Rhodamine Derivatives », *Biophys J*, **76**, 469 (1999) [Article](#)
- **Schinder A.F., et al.**, « Mitochondrial Dysfunction Is a Primary Event in Glutamate Neurotoxicity », *The Journal of Neuroscience*, **16**, 6125(1996) [Article](#)
- **Slomovic S. & Schuster G.**, Stable PNPase RNAi silencing: Its effect on the processing and adenylation of human mitochondrial RNA, *RNA*, **14**:310-323 (2008) [Abstract](#)
- **Voronina S.G., et al.**, « Effects of Secretagogues and Bile Acids on Mitochondrial Membrane Potential of Pancreatic Acinar Cells COMPARISON OF DIFFERENT MODES OF EVALUATING », *J. Biol. Chem.*, **279**, 27327, (2004) [Article](#)
- **Ward M.W., et al.**, « Mitochondrial Membrane Potential and Glutamate Excitotoxicity in Cultured Cerebellar Granule Cells », *The Journal of Neuroscience*, **20**, 7208 (2000) [Article](#)
- **Xi Q. et al.**, Mitochondria-Derived Reactive Oxygen Species Dilate Cerebral Arteries by Activating Ca²⁺ Sparks, *Circulation Research*. **97**:354 (2005) [Article](#)

Ordering information

Catalog size quantities and prices may be found at <http://www.fluoprobes.com>
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.

Rev.F02VB