

# β-Glucuronidase/Arylsulfatase

obtained from *Helix pomatia*

β-D-Glucuronoside glucuronosohydrolase, EC 3.2.1.31

Aryl-sulfate sulfohydrolase, EC 3.1.6.1

Cat. No. 10 127 060 001 2 ml

Cat. No. 10 127 698 001 10 ml

**Version 08**  
Content version: March 2012

Store at +2 to +8°C


## Product overview

<b>Composition</b>	Enzyme mix of β-glucuronidase/arylsulfatase obtained from <i>Helix pomatia</i> in saline, stabilized with 0.02% sodium azide.
<b>pH optimum</b>	<ul style="list-style-type: none"> <li>The pH optimum value for β-glucuronidase activity is 4.5–5.0.</li> <li>The optimum for arylsulfatase activity is generally 6.2, but may be greater for some substrates in comparatively high concentrations; for instance, for 16.5 mM solutions of 2-hydroxy-5-nitrophenylhydrogen sulfate it is 7.2.</li> </ul>
<b>Inhibitors</b>	<ul style="list-style-type: none"> <li>β-Glucuronidase activity is inhibited by: <ul style="list-style-type: none"> <li>D-glucuronic acid</li> <li>D-galacturonic acid</li> <li>D-glucaro-1,4-lactone (saccharolactone found in urine)</li> </ul> </li> <li>Arylsulfatase activity is inhibited by phosphate.</li> </ul>
<b>Application</b>	<p>The enzyme preparation obtained from the Roman snail, <i>Helix pomatia</i>, which exhibits very strong β-glucuronidase and arylsulfatase activity, is widely used for the simultaneous hydrolysis of β-glucuronides (β-glucosiduronic acids) and sulfate esters in urine and other biological fluids (1, 2, 3, 4).</p> <ul style="list-style-type: none"> <li>Enzymatic hydrolysis of steroid β-glucuronides and sulfates</li> <li>Cell biology (removal of cell walls from yeasts in the preparation of protoplasts) (5)</li> <li>Enzyme immobilization studies (6)</li> <li>Determination of drugs in urine (7)</li> </ul>
<b>Storage/ Stability</b>	<p>Undiluted, the aqueous solution of β-glucuronidase/arylsulfatase is stable at +2 to +8°C until the expiration date printed on the label.</p> <p><b>Note:</b> Aliquots portions of the diluted preparation may be stored at -15 to -25°C; they should not be thawed and refrozen more than once or twice; and storage at lower temperatures does not lengthen their expiration date beyond that of the product kept at +2 to +8°C.</p>
<b>Working concentration</b>	<p>In many applications the product can be diluted with water immediately before use or used undiluted.</p> <p><b>Note:</b> This β-glucuronidase/arylsulfatase preparation is very concentrated and must be diluted for some applications. In particular for the preparation of protoplasts, the precise concentration to use for a given strain of yeast must be determined empirically.</p>

## Product description

<b>Specificity of β-Glucuronidase</b>	<p>The glycosides that β-D-glucuronic acid forms with a variety of compounds containing hydroxyl groups, hydrolyse readily in the presence of β-glucuronidase. Such compounds include:</p> <ul style="list-style-type: none"> <li>steroids, such as estriol (<math>K_m = 0.42</math> mM, pH 4.5), androsterone, pregnanediol, tetrahydrocortisone</li> <li>phenols, such as phenolphthalein (<math>K_m = 0.39</math> mM), 4-nitro-phenol, 4-methylumbelliferone</li> <li>drugs, such as chloramphenicol and tetrahydrocannabinols</li> <li>metabolites, such as thyroxine and bilirubin</li> </ul> <p>Polysaccharides that contain β-glucuronic acid residues, such as hyaluronic acid, are also hydrolyzed.</p> <p>β-Glucuronidase is highly specific for the carbohydrate moiety: neither α-glucosides nor β-glucosiduronic acids are hydrolyzed.</p>
<b>Specificity of Arylsulfatase</b>	<p>Sulfate esters of many phenols are hydrolyzed in the presence of arylsulfatase. Examples are steroid sulfates, such as estronesulfate, 4-nitrophenyl hydrogen sulfate (<math>K_m = 1.8</math> mM, pH 7.3), 4-nitro-pyrocatechol 2-sulfate (<math>K_m = 1.25</math> mM, pH 7.5), and phenolphthalein disulfate.</p>

## Principle of the standard test



At a wavelength of 405 nm, the molar absorption coefficient of 4-nitrophenol is  $1.85 \text{ mM}^{-1} \times \text{l} \times \text{cm}^{-1}$  at +25°C.

<b>Steroids in urine</b>	<p>The various steroids found in urine may be present in one or more of three forms:</p> <p>(a) the free compound (in minor or trace quantities and amounts)</p> <p>(b) the sulfate (predominant in some cases)</p> <p>(c) the β-glucuronide (the predominant form in most cases)</p>
--------------------------	---

**Tab 1:** Relative steroid proportions before hydrolysis are given in the following table:

Compound or category	free [%]	sulfate [%]	glucuronide [%]
7-Hydroxycorticosteroids	1	10–15	85–90
Pregnanediol	0	trace	~ 100
Pregnanetriol	trace	trace	~ 100
Estrone (O <sub>1</sub> )	1–3	10–15	85–89
Estradiols (O <sub>2</sub> )	1–3	5–10	90–95
Estriol (O <sub>3</sub> )	0–2	5–10	90–95
Androsterone	trace	20	80
Etiocholanolone	trace	10	90
Dehydroepiandrosterone (DHEA)	trace	~ 100	trace
Epiandrosterone	trace	~ 100	trace
11 β-Androsterone	trace	10	90
11 β-Etiocholanolone	trace	10	90
11 Ketoandrosterone	trace	trace	~ 100
11 Ketoetiocholanolone	trace	trace	~ 100

**Methods for hydrolysis**

Several methods of hydrolyzing steroid esters and glycosides are commonly used.

- For the sulfates of DHEA and androsterone, solvolysis is suitable; this involves treatment with excess organic solvent (*e.g.*, ethyl acetate, dioxan, or tetrahydrofuran) at a temperature of +38°C for 18–24 h.
- Acid hydrolysis at elevated temperatures is a more general method, but has two disadvantages: it can alter the structure and function of the steroids, and the resinified pigments formed need to be removed, because they are present in the extract.
- The third method, enzymatic hydrolysis (with β-glucuronidase and sulfatase), does not involve these drawbacks.

**Procedure for the hydrolysis of glucuronides and sulfates in urine**

**Protocol**

Please refer to the following table.

Step	Action
1	Adjust the pH of a portion of the sample (10 ml) to 5.5 by adding dilute acetic acid.
2	Add 1 ml acetate buffer (1 M pH 5.5) and 0.2 ml β-glucuronidase/ arylsulfatase solution.
3	Incubate at a temperature of +37°C for 16 h.
4	Cool and extract with an appropriate solvent ( <i>e.g.</i> , chloroform or dichloromethane) to isolate the hydrolysis products.

**Specific activity of Glucuronidase**

**Activity**

Generally, the β-glucuronidase activity of the preparation is not as high with respect to steroid β-glucuronides, as values obtained from the hydrolysis of synthetic phenyl β-glucuronides indicate. However, under certain conditions, results obtained with phenolphthalein β-glucuronide may be comparable with those given by steroid glycosides, such as estradiol β-glucuronide. For instance, at +37°C and pH 4.5, a β-glucuronidase/arylsulfatase preparation that promotes the hydrolysis of 300 μmol of phenolphthalein β-glucuronide in 1 h, also promotes the hydrolysis of 441 μmol of estradiol β-glucuronide in 1 h.

**Standard unit**

The standard unit of β-glucuronidase activity is the enzyme activity that increases the rate of release of 4-nitrophenol from 4-nitrophenyl β-D-glucosiduronic acid at a temperature of +25°C and pH 4.5 by 1 μM.

**Phenolphthalein unit**

The phenolphthalein unit of β-glucuronidase activity is the enzyme activity that increases the rate of release of phenolphthalein from phenolphthalein β-D-glucosiduronic acid at a temperature of +38°C by 1 μM.

Approximately 4.5 standard units are equivalent to 5.5 phenolphthalein units.

**Fishman unit**

The Fishman unit (8) of β-glucuronidase activity is the enzyme activity that increases the rate of release of phenolphthalein from phenolphthalein β-D-glucosiduronic acid at a temperature of +38°C by 1 μg.

Approximately 1 standard unit is equivalent to 22,000 Fishman units (1 phenolphthalein unit is equivalent to 19,000 Fishman units).

**Specific β-glucuronidase activity**

At +25°C and pH 4.5, the β-glucuronidase activity of 1 ml of the preparation is 4.5 standard units, equivalent to 5.5 phenolphthalein units or 100,000 Fishman units at +38°C.

**Specific activity of Arylsulfatase**

**Standard unit**

The standard unit of arylsulfatase activity is the enzyme activity that increases the rate of release of 4-nitrophenol from 4-nitrophenyl sulfate at a temperature of +25°C and pH 6.2 by 1 μM.

**Phenolphthalein unit**

The phenolphthalein unit of arylsulfatase activity is the enzyme activity that increases the rate of release of phenolphthalein from phenolphthalein disulfate at a temperature of +38°C and pH 6.2 by 1 μM.

Approximately 5.4 standard units are equivalent to 1 phenolphthalein unit.

**Roy unit**

The Roy unit of arylsulfatase activity is the enzyme activity that increases the rate of release of 4-nitropyrrocatechol from 2-hydroxy-5-nitrophenyl hydrogen sulfate (4-nitropyrrocatechol 2-sulfate) at a temperature of +38°C and pH 6.2 by 1 μg (9).

Approximately 1 standard unit is equivalent to 57,000 Roy units (1 phenolphthalein unit is equivalent to 308,000 Roy units).

**Specific aryl-sulfatase activity**

At +25°C and pH 6.2, the arylsulfatase activity of 1 ml of the preparation is 14 standard units, equivalent to 2.6 phenolphthalein units or 800,000 Roy units at +38°C.

**Changes to previous version**

Editorial changes

**Ordering Information**

Roche Applied Science provides a large selection of reagents and systems for life science research. For a complete overview of related products and manuals, please visit and bookmark our home page, [www.roche-applied-science.com](http://www.roche-applied-science.com)

**References**

- Bergmeyer, H.U., Grassel, M. & Walter, H.-E. (1983) in "Methods of Enzymatic Analysis" (Bergmeyer, H.U., eds.), 3. ed., Vol. 2, pp 206–207, VCH, Weinheim.
- Stahl, P.D. & Fishman, W.H. (1984) in "Methods of Enzymatic Analysis", (Bergmeyer, H.U., eds.), 3. eds., Vol. 4, pp 246–256, VCH, Weinheim.
- Tuite, M. et al. (1980) J. Biol. Chem. 255 8761.
- Jarrige, P. Yon, J. & Hayle, M.F. (1963) Bull. Soc. Chim. Biol. 45 783.
- von Hedenstroem, M. & Hoefer, M. (1974) Arch. Microbiol. 98 51.
- Rapatz, E. et al. (1988) "Studies on the Immobilization of Glucuronidase (Part I and II)", in "Applied Biochemistry and Biotechnology", Vol. 19, pp 223–242.
- Machata, G. & Kryspin-Eyner, K. (1970) Wiener klin. Wochenz. 47 849.
- Wakabayashi, M. & Fishman, W.H. (1961) J. Biol. Chem. 236 996.
- Roy, A.B. (1970) in "Chemical and Biological Aspects of Steroid Conjugation" (Bernstein, S. & Solomon, S., eds.), pp 106, Springer, New York.

**Trademarks**

All other product names and trademarks are the property of their respective owners

**Regulatory Disclaimer**

For life science research only. Not for use in diagnostic procedures.

**Contact and Support**

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support Site** at:

**[www.roche-applied-science.com/support](http://www.roche-applied-science.com/support)**

To call, write, fax, or email us, visit the Roche Applied Science home page, [www.roche-applied-science.com](http://www.roche-applied-science.com), and select your home country. Country-specific contact information will be displayed. Use the Product Search function to find Pack Inserts and Material Safety Data Sheets.



Roche Diagnostics GmbH  
Roche Applied Science  
68298 Mannheim  
Germany