



Biotin-X-dUTP

Biotin-X-dUTP can be enzymatically incorporated into DNA via nick translation, PCR random priming or 3'-end terminal labeling.

Description

Product name cat.number	MW (g·mol ⁻¹)	Comments
Biotin-11-dUTP, 1mM in pH 7.5 Tris-HCl buffer, 1 mM EDTA FP-AM539B, 100 µl	886.5	- Biotin-11-2'-deoxyuridine-5'-triphosphate, tetralithium salt - The length of linker "11" is optimal for the majority of applications.
Biotin-11-dUTP, lyophilized powder FP-AY489A, 50 µg	886.5	- Biotin-11-2'-deoxyuridine-5'-triphosphate, tetralithium salt - Lyophilized solid that is more suitable for long term storage. The product contains lyophilized TE buffer, so you only need to add an appropriate amount of deionized H ₂ O to reconstitute the solution.
Biotin-16-dUTP, 1mM in pH 7.5 Tris-HCl buffer FP-AM541A, 50 µl	971.5	- Biotin-16-2'-deoxyuridine-5'-triphosphate, tetralithium salt - The terminal deoxynucleotidyl transferase (TDT)-mediated biotin-dUTP nick end-labeling (TUNEL) method has been commonly used for apoptosis studies.
Biotin-16-dUTP, lyophilized powder FP-AM542A, 50 µg	971.5	- Biotin-16-2'-deoxyuridine-5'-triphosphate, tetralithium salt - Lyophilized solid that is more suitable for long term storage. The product contains lyophilized TE buffer, so you only need to add an appropriate amount of deionized H ₂ O to reconstitute the solution.
Biotin-20-dUTP, 1mM in pH 7.5 Tris-HCl buffer FP-AM543B, 50 µl	1020.5 4	- Biotin-20-2'-deoxyuridine-5'-triphosphate, tetralithium salt - Longer and water-soluble spacer group, which should facilitate interaction between the biotin group and avidin or streptavidin.
Biotin-20-dUTP, lyophilized powder FP-AM543A, 50 µg	1020.5 4	- Biotin-20-2'-deoxyuridine-5'-triphosphate, tetralithium salt - Lyophilized solid that is more suitable for long term storage. The product contains lyophilized TE buffer, so you only need to add an appropriate amount of deionized H ₂ O to reconstitute the solution.

Storage: -20°C >1 year Protect from light

Introduction

The number 'X' is the number of atoms in the linker between biotin and dUTP. The length of the linker affects the incorporation efficiency of the biotin-dUTP probe into DNA using DNA polymerases, and it also affects biotin/avidin or biotin/streptavidin. In general, the shorter the linker, the more efficiently the biotin-dUTP is

FT-595911

incorporated into DNA by DNA polymerases. On the other hand, the longer the linker, the better biotin can interact with avidin or streptavidin.

Standard protocol

The most popular approach for DNA PCR-labeling with labeled dUTP is based on the usage of dNTPs mixture which contains labeled dUTP and all 4 other dNTP in regular concentrations. The molar ratio dUTP/labeled dUTP (or dTTP/labeled dUTP) can vary from 3 :1 to 1 :1. The incorporation efficiency depends mainly on the usage of dTTP or dUTP (the incorporation efficiency of dTTP is slightly better than those for dUTP) and on the enzyme used for PCR. Regular Taq DNA polymerase incorporates dUTP (and especially labeled dUTP) less efficient than Taq DNA polymerase with modified active center.

In some special applications one may completely substitute dTTP by labeled dUTP and to get DNA with all « T » substituted to labeled dUTP. Meanwhile, this 100% labeled DNA will be quite different from regular DNA in terms of electrophoresis mobility, hydrophobic properties, denaturation behavior etc. If all these points can be neglected, one can completely substitute dTTP by labeled dUTP.

The regular protocol for DNA labeling by PCR use a ratio of 1:2 (labeled dUTP : non-labeled dTTP) :

Reagent	Final concentration	Quantity for 50 µl of reaction mixture
Sterile deionized water	-	Variable
10X PCR buffer	1X	5 µl
10 mM dNTP Mix	0.2 mM of each	1 µl
labeled dUTP, 1 mM	0.1 mM	5 µl
Primer I	0.1 – 1 µM	Variable
Primer II	0.1 – 1 µM	Variable
Taq DNA Polymerase	1.25 U – 2.5 U/50 µl	1.25 U – 2.5 U
100 mM MgCl ₂	1 – 4 mM	Variable
Template DNA	10 pg – 1 µg	Variable

PCR should be performed as optimized on the regular dNTPs – with the same MgCl₂ concentration, with the temperatures and cycles optimized for the particular template and primers.

Protocol with Biotin-dUTP, PCR grade

20 µl PCR labeling assay			
Component	Stock conc.	Amount	Final conc.
High yield buffer without MgCl ₂ (Cat.-No. PCR-201)	10x	2 µl	1x
MgCl ₂ stock solution	25 mM	1.6 µl	2 mM
dATP	1 mM	2 µl	100 µM
dCTP	1 mM	2 µl	100 µM
dGTP	1 mM	2 µl	100 µM
dTTP	1 mM	1 µl	50 µM
Biotin-dUTP-PCR	1 mM	1 µl ¹⁾	50 µM ¹⁾
forward Primer	10 µM	1 µl	500 nM
reverse Primer	10 µM	1 µl	500 nM
Template DNA		0.1-10 ng	5-500 pg/µl
Taq Pol (Cat.-No. PCR-201)	5 units/µl	0.2 µl (1 unit)	0.05 units/µl
PCR grade H ₂ O		Fill up to 20 µl	

FT-595911

Recommended cycling conditions

Initial denaturation	94 °C	2 min	1x
Denaturation	94 °C	30 sec	25-30x
Annealing ¹⁾	50-60 °C	30 sec	
Elongation ²⁾	72 °C	1 min	
Final elongation	72 °C	5 min	1x

- 1) The annealing temperature depends on the melting temperature of the primers used.
- 2) The elongation time depends on the length of the fragments to be amplified. A time of 2 min/kbp is recommended.

For optimal amplification results and high incorporation rates an individual optimization of the recommended PCR assay and cycling conditions may be necessary for each new primer-template pair.

Related products

Description	Cat. No.
UptiTherm DNA Polymerase (1000 Units) 5U/μl with Mg free Buffer + 50 mM MgCl ₂ buffer	UPS53921
Terminal deoxynucleotidyl transferase (1000 Units)	HP9020
2x Hot Start PCR Master Mix (500 Units)	CJ5361
Fast EvaGreen qPCR Master Mix (500 tests)	DV7221
PCR set 1 (dATP, dGTP, dCTP, dTTP 100mM each)	UP968640
PCR Mix 3 (200μl) (10 mM of each dA, dC, dG and dT)	UP984440
GelRed nucleic acid gel prestaining kit	FJ5570

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.