

Kinases inhibitors

Inhibitors based on specificity for different kinases :

Kinase target	Kinase inhibitors available (inquire)
PKA inhibitors :	H-89, HA-1004, H-7, H-8, HA-100, PKI, staurosporine
PKG inhibitors :	H-7, H-8, HA-1004, H-89, staurosporine
PKC inhibitors :	bisindolylmaleimide IX, chelerythrine, edelfosine, edelfosina, ET18OCH ₃ , H-7, HA-100, H89, HA-1004, Ro 31-8220, rottlerin, staurosporine, quercetin
ABL inhibitors :	PP1
Adenosine kinase inhibitors :	5-Iodotubercidin
Akt/PKB inhibitors :	Triciribine (TCN, Tricyclic Nucleoside, NSC-154020, API-2)
AMPK inhibitors :	indirubin-3'-oxime, 5-Iodotubercidin, H-89
BCR-ABL :	adaphostin, Tyrphostin AG 957
BTK inhibitors :	LFM-A13
CaMKII inhibitors :	KN-62, W-7, HA-1004, HA-1077, staurosporine
CDK1, CDK2, CDK5 inhibitors :	aloesine, aloesine A, alsterpaullone, indirubin-3'-oxime, kenpaullone, olomoucine, purvalanol A, roscovitine
CHK1 inhibitors :	H-89
CK1 inhibitors :	DRB
CK2 inhibitors :	TBB, DRB
EGF-R kinase inhibitors :	erbstatin analog, gefitinib, PD 153035, Tyrphostin AG 490, Tyrphostin AG 825, PP1, PP2
ERK 1/2 inhibitors :	ERK inhibitor
ERK kinase inhibitors :	PD 98059, SL327, olomoucine, 5-Iodotubercidin
FYN inhibitors :	PP1, PP2, SU6656
GSK3 inhibitors :	aloesine, aloesine A, indirubin-3'-oxime, kenpaullone
HCK inhibitors :	PP2
HER2 inhibitors :	Tyrphostin AG 825
IGF1-R PTK inhibitor :	picropodophyllin, Tyrphostin AG 1024
JAK inhibitors :	Tyrphostin AG 490, LFM-A13
JNK inhibitors :	aloesine A, SP600125
KIT inhibitors :	PP1
LCK inhibitors :	aminogestine, PP1, PP2, indirubin-3'-oxime, kenpaullone
LYN inhibitors :	SU6656, staurosporine
MAPKAP-K1b inhibitors :	HA-1077
MAPK inhibitors :	SB202190, SB203580
MAPKK inhibitors :	arctigenin, PD 98059, SL327, U0126
MEK inhibitors :	PD 98059, SL327, U0126
MKK inhibitors :	arctigenin, PD 98059, SL327, U0126
MLCK inhibitors :	ML-9, W-7, HA-100, HA-1004, H-7, H-8, HA-1077, staurosporine
MSK1 inhibitors :	H-89, ML-9, HA-1077
P38 MAPK inhibitors :	SB202190, SB203580
PDGF-R kinase inhibitors :	DMPQ
PDK1 inhibitors :	SB203580, NU6102
PI3K inhibitors :	LY294002, wortmannin, quercetin
PRK2 inhibitors :	HA-1077
RHO kinase inhibitors :	HA-1077, hydroxyfasudil
ROCK inhibitors :	Y-27632, H-89, ML-9, HA-1077
S6K1 inhibitors :	rapamycin, sirolimus, H-89, ML-9, HA-1077, staurosporine
SAPK inhibitors :	PD 169316, SB202190, SB203580
SGK inhibitors :	indirubin-3'-oxime, H-89
SRC family PTK inhibitors :	PP1, PP2, SU6656, staurosporine
SYK inhibitors :	piceatannol, staurosporine

FT-972971

TOR inhibitors :	rapamycin, sirolimus
PTK inhibitors :	aminogestrol, DMPQ, erastatin analog, gefitinib, LFM-A13, PD 153035, piceatannol, PP1, PP2, quercetin, SU 4312, SU6656, Tyrphostin AG 490, Tyrphostin AG 825, Tyrphostin AG 957
VEGF-R kinase inhibitors :	SU 4312
YES inhibitors :	SU6656
ZAP70 inhibitors :	piceatannol

Products Description


Product name cat.number	MW g·mol ⁻¹ (salt)	CAS		
Staurosporine (M) 74146E, 1mg 74146F, 5 mg	466.54	62996-74-1		Inhibitor of Protein Kinase (broad specificity: PKC, PKA, PKB...)
Chelerythrine Chloride (M) 972971, 1g	383.83			
Bisindolylmaleimide IX (M) S31322, 1mg S31323, 5mg	553.66			
DRB (M) 53944A	319.2 (553.66)	53-85-0		specific inhibitor of casein kinase (CH2 and CK1)
Edelfosine, ET180CH3 (L) CJ6031, 300mg	523.7			
H-7 (M) KV3842, 25mg	291.4 (+73/HCl)	108930-17-2		inhibitor of PKC, PKA and PKG and MLCK
H-8 (M) KV3851, 10mg	265.3 (+73/HCl)	84478-11-5		inhibitor of PKG and PKA and MLCK
H-89 (M) CQ7251, mg	446.4 (+73/HCl)	130964-39-5		potent and inhibitor of PKA, PKG ...
HA-100 (M) KV3871, 10mg	d (+73/HCl)			inhibitor of MLCK, PKA and PKC
HA-1004 (M) , mg	293.3 (+73/HCl)	92564-34-6		inhibitor of PKG and PKA and at higher ...
HA-1077 (M) KV3891, 10mg	291.4 (+73/HCl)	203911-27-7		Inhibitor of ROCKII, Rho kinase, PRK2 and MSK1
HydroxFasudil (M) KV3900, 5mg	343.8	155558-32-0		Rho kinase inhibitor.
5-IodoTubercidin (M) KV3920, 5µmol	392.2	24386-93-4		specific inhibitor of ERK2 and adenosine kinase

Storage: +4°C (or -20°C for long term) (L)
Protect from light and moisture

Technical and scientific information

Staurosporine #74146

Synonyms: antibiotic AM-2282, Staurosporine from *Streptomyces staurosporeus*
See the [technical sheet](#)




Bisindolylmaleimide IX #S3132

Product name: Bisindolylmaleimide IX, Methanesulfonate Salt

Synonyms: Ro 31-8220, Ro31-8220, Ro 318220, Ro318220, BIM IX, BIM-9, Bis IX

Appearance: Orange to red microcrystalline solid. Solubility: Soluble in DMSO, ethanol and water.

Bisindolylmaleimide IX (Ro 31-8220) is a selective and ATP-competitive PKC inhibitor with an IC₅₀ of 10 nM, whereas the IC₅₀ for the inhibition of PKA is 900 nM [Davis et al. 1992]. The IC₅₀ values for the inhibition of the group A PKCs range from 5-27 nM and that for PKCε is 24 nM [Wilkinson et al. 1993]. Ro 31-8220 inhibits growth factor-stimulated expression of MKP-1 and c-Fos but strongly stimulated c-Jun expression, even in the absence of growth factors. Ro 31-8220 strongly stimulates the stress-activated protein kinase JNK1 in a PKC-independent manner [Beltman et al. 1996].



References:

- Davis PD, Elliott LH, Harris W, Hill CH, Hurst SA, Keech E, Kumar MK, Lawton G, Nixon JS, Wilkinson SE (1992) "Inhibitors of protein kinase C. 2. Substituted bisindolylmaleimides with improved potency and selectivity" *J. Med. Chem.* 35(6):994-1001
- Wilkinson SE, Parker PJ, Nixon JS (1993) "Isoenzyme specificity of bisindolylmaleimides, selective inhibitors of protein kinase C" *Biochem. J.* 294 (Pt 2):335-7
- Bit RA, Davis PD, Elliott LH, Harris W, Hill CH, Keech E, Kumar H, Lawton G, Maw A, Nixon JS, et al. (1993) "Inhibitors of protein kinase C 3. Potent and highly selective bisindolylmaleimides by conformational restriction" *J. Med. Chem.* 36(1):21-9
- Dieter P, Fitzke E (1991) "RO 31-8220 and RO 31-7549 show improved selectivity for protein kinase C over staurosporine in macrophages" *Biochem. Biophys. Res. Commun.* 181 (1): 396-401
- Keller HU, Niggli V (1993) "The PKC-inhibitor Ro 31-8220 selectively suppresses PMA- and diacylglycerol-induced fluid pinocytosis and actin polymerization in PMNs" *Biochem. Biophys. Res. Commun.* 194 (3): 1111-6
- Beltman J, McCormick F, Cook SJ (1996) "The selective protein kinase C inhibitor, Ro-31-8220, inhibits mitogen-activated protein kinase phosphatase-1 (MKP-1) expression, induces c-Jun expression, and activates Jun N-terminal kinase" *J. Biol. Chem.* 271(43):27018-24
- Budworth J, Davies R, Malkhandi J, Gant TW, Ferry DR, Gescher A (1996) "Comparison of staurosporine and four analogues: their effects on growth, rhodamine 123 retention and binding to P-glycoprotein in multidrug-resistant MCF-7/Adr cells" *Br. J. Cancer.* 3 (9):1063-8
- Harris TE, Persaud SJ, Jones PM (1996) "Atypical isoforms of PKC and insulin secretion from pancreatic beta-cells: evidence using Go 6976 and Ro 31-8220 as PKC inhibitors" *Biochem. Biophys. Res. Commun.* 227 (3): 672-6
- Yeo EJ, Provost JJ, Exton JH (1997) "Dissociation of tyrosine phosphorylation and activation of phosphoinositide phospholipase C induced by the protein kinase C inhibitor Ro-31-8220 in Swiss 3T3 cells treated with platelet-derived growth factor" *Biochim. Biophys. Acta* 1356 (3): 308-20
- Begemann M, Kashimawo SA, Lunn RM, Delohery T, Choi YJ, Kim S, Heitjan DF, Santella RM, Schiff PB, Bruce JN, Weinstein IB (1998) "Growth inhibition induced by Ro 31-8220 and calphostin C in human glioblastoma cell lines is associated with apoptosis and inhibition of CDC2 kinase" *Anticancer Res.* 18(5A):3139-52
- Hers I, Tavaré JM, Denton RM (1999) "The protein kinase C inhibitors bisindolylmaleimide I (GF 109203x) and IX (Ro 31-8220) are potent inhibitors of glycogen synthase kinase-3 activity" *FEBS Lett.* 460(3):433-6.
- Standaert ML, Bandyopadhyay G, Antwi EK, Farese RV (1999) "RO 31-8220 activates c-Jun N-terminal kinase and glycogen synthase in rat adipocytes and L6 myotubes. Comparison to actions of insulin" *Endocrinology* 140(5):2145-51
- Lingameneni R, Vysotskaya TN, Duch DS, Hemmings HC Jr (2000) "Inhibition of voltage-dependent sodium channels by Ro 31-8220, a 'specific' protein kinase C inhibitor" *FEBS Lett.* 473(2):265-8

Chelerythrine #97297

Product name: Chelerythrine chloride

Chemical name: 1,2-Dimethoxy-N-methyl(1,3)benzodioxolo(5,6-c) phenanthridinium chloride, CAS: 3895-92-9

Molecular weight 383.8 g/mol. Purity: >98%. Appearance: Yellow to orange solid.

Solubility: soluble in DMSO and water. Rather insoluble in aqueous media that contain substantial saline concentrations.

Chelerythrine Chloride is a potent and selective protein kinase C inhibitor (PKC, IC₅₀ = 0.66 μM).

References:

Info@fluoprobes.com
Technical-support@fluoprobes.com
Order-online@fluoprobes.com

Contact your local distributor

FluoProbes®, powered by



213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

FT-972971

Herbert JM, Augereau JM, Gleye J, Maffrand JP (1990) "Chelerythrine is a potent and specific inhibitor of protein kinase C" *Biochem. Biophys. Res. Commun.* 172(3):993-9

Liu Y, Cohen MV, Downey JM (1994) "Chelerythrine, a highly selective protein kinase C inhibitor, blocks the anti-infarct effect of ischemic preconditioning in rabbit hearts" *Cardiovasc. Drugs Ther.* 8(6):881-2

Serrano PA, Rodriguez WA, Pope B, Bennett EL, Rosenzweig MR (1995) "Protein kinase C inhibitor chelerythrine disrupts memory formation in chicks" *Behav. Neurosci.* 109(2):278-84

Eckly-Michel AE, Le Bec A, Lugnier C (1997) "Chelerythrine, a protein kinase C inhibitor, interacts with cyclic nucleotide phosphodiesterases" *Eur. J. Pharmacol.* 324(1):85-8


Sacchetti B, Bielavska E (1998) "Chelerythrine, a specific PKC inhibitor, blocks acquisition but not consolidation and retrieval of conditioned taste aversion in rat" *Brain Res.* 799(1):84-90

Lundmark JL, Ramasamy R, Vulliet PR, Schaefer S (1999) "Chelerythrine increases Na-K-ATPase activity and limits ischemic injury in isolated rat hearts" *Am. J. Physiol.* 277(3 Pt 2):H999-H1006

Yu R, Mandelkar S, Tan TH, Kong AN (2000) "Activation of p38 and c-Jun N-terminal kinase pathways and induction of apoptosis by chelerythrine do not require inhibition of protein kinase C" *J. Biol. Chem.* 275(13):9612-9

DRB #53944

Product name: 5,6-Dichlorobenzimidazole riboside (DRB)
Synonyms: 5,6-Dichloro-1-β-D-ribofuranosyl benzimidazole; CAS [53-85-0]
Molecular weight 319.2 g/mol. **Chemical formula** C₁₂H₁₂Cl₂N₂O₄.
Purity: >98% (HPLC). **Appearance:** White solid.
Solubility: soluble in DMSO or ethanol (may require warming for complete dissolution).



DRB is a specific inhibitor of casein kinase 2 (CK2, IC₅₀ = 6 μM) and casein kinase 1 (CK1, IC₅₀ = 14 μM). It is also a specific inhibitor for RNA polymerase II transcription in vivo and in vitro whereat the effect on RNA polymerase II-specific transcription seems to be mediated by its inhibition of nuclear casein kinase 2. DRB is a competitive inhibitor in respect to its phosphate donor substrates ATP and GTP.

References:

Zandomeni R, Weinmann R. (1984) "Inhibitory effect of 5,6-dichloro 1-β-D ribofuranosyl benzimidazole on a protein kinase." *J. Biol. Chem.* 259(23):14804-11.

Zandomeni R, Zandomeni MC, Shugar D, Weinmann R (1986) "Casein kinase type II is involved in the inhibition by 5,6-dichloro-1-β-D-ribofuranosylbenzimidazole of specific RNA polymerase II transcription." *J. Biol. Chem.* 261(7):3414-9.

Chodosh LA, Fire A, Samuels M, Sharp PA. (1989) "5,6-Dichloro 1-β-D ribofuranosyl benzimidazole inhibits transcription elongation by RNA polymerase II in vitro." *J. Biol. Chem.* 264(4):2250-7

Zandomeni RO (1989) "Kinetics of inhibition by 5,6-dichloro 1-β-D ribofuranosyl benzimidazole on calf thymus casein kinase II." *Biochem. J.* 262(2):469-73.

Meggio F, Shugar D, Pinna LA (1990) "Ribofuranosyl-benzimidazole derivatives as inhibitors of casein kinase-2 and casein kinase-1." *Eur. J. Biochem.* 187(1):89-94

Blaydes JP, Hupp TR (1998) "DNA damage triggers DRB-resistant phosphorylation of human p53 at the CK2 site" *Oncogene* 17(8):1045-52


Egyhazi E, Ossoinak A, Filhol-Cochet O, Cochet C, Pigon A (1999) "The binding of the alpha subunit of protein kinase CK2 and RAP74 subunit of TFIIF to protein-coding genes in living cells is DRB sensitive" *Mol. Cell. Biochem.* 191(1-2):149-59

Johnston IM, Allison SJ, Morton JP, Schramm L, Scott PH, White RJ (2002) "CK2 forms a stable complex with TFIIB and activates RNA polymerase III transcription in human cells" *Mol. Cell. Biol.* 22(11):3757-68

Uhle S, Medalia O, Dubiel W et al. (2003) "Protein kinase CK2 and protein kinase D are associated with the COP9 signalosome" *EMBO J.* 22(6):1302-12

Edelfosine, ET18OCH3 #CJ603

Product name: edelfosine
Synonyms: edelfosin, edelfosinum, edelfosina, ET18OCH3
Chemical name: rac-2-methyl-1-octadecyl-glycero-(3)-phosphocholine
Chemical formula C₂₇H₅₈NO₆P. **Molecular weight** 523.7 g/mol.
Purity: >99.1% (GMP quality). **Appearance:** white, odourless powder.
Solubility: soluble in chloroform, ethanol and water.
 The substance may be carcinogenic or teratogenic.



Edelfosine or ET18OCH3 is a synthetic analog of lysophosphatidylcholine. It is a strong inhibitor of the PI3K pathway (anti-tumor lipid) and selectively inhibits protein kinase C (PKC) and phosphatidylinositol specific phospholipase C (PI PLC). Edelfosine induces expression of c-fos, c-jun and transcription factor AP-1, induces apoptosis in human promyelocytic and promonocytic leukemic cells and decreases multi drug resistance (MDR) when combined with other MDR modifiers and pH-lowering agents such as verapamil, tamoxifen or cyclosporin A.

References:

Zheng, B., et al. (1990). Inhibition of protein kinase C, (sodium plus potassium)-activated adenosine triphosphate, and sodium pump by synthetic phospholipid analogues. *Cancer Res.* 50: 3025-3031.

Info@fluoprobes.com

Technical-support@fluoprobes.com

Order-online@fluoprobes.com

FluoProbes®, powered by



213 Avenue J.F. Kennedy - BP 1140
 03103 Montluçon Cedex - France
 Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

P.4

Contact your local distributor

FT-972971

Powis, G., et al. (1992). Selective inhibition of phosphatidylinositol phospholipase C by cytotoxic ether lipid analogues., *Cancer Res.* 52: 2835-2840.

Mollinedo, F. et al. (1993). Early and selective induction of apoptosis in human leukemic cells by the alkyl-lysophospholipid ET-18-OCH₃. *Biochem Biophys Res Commun.*, 192(2), 603-609.

Mollinedo, F., et al. (1994). The ether lipid 1-octadecyl-2-methyl rac glycerol-3-phosphocholine induces expression of fos and jun proto-oncogenes and activates AP-1 transcription factor in human leukaemic cells. *Biochem J.* 302:325-329.

Mollinedo F., et al. (1997) Selective induction of apoptosis in cancer cells by the ether ET-18-OCH₃ (edelfosine): molecular structure requirements, cellular uptake, and protection by bcl-2 and bclxL. *Cancer Res.* 1997; 57: 1320-1328.

Harguindeguy, S. et al. (2000). Edelfosine, apoptosis, MDR and Na⁺/H⁺ exchanger: Induction mechanisms and treatment implications. *Apoptosis.* 5: 87-89.

H-7, #KV384


Product names: H-7. dihydrochloride

Syn.: H7, H 7, i-(5-Isoquinolinesulfonyl)-2-methylpiperazine . 2HCl

CAS number: 108930-17-2; Chemical formula: C₁₄H₁₇N₃O₂S . 2 HCl;

Molecular weight: 291.4 + 73.0 g/mol

Purity: > 97%; Appearance: powder; Solubility: soluble in water



H-7 is an inhibitor of protein kinase C (PKC), cAMP- and cGMP-dependent protein kinase (PKA and PKG) and also inhibits myosin light chain kinase (MLCK). It induces apoptotic DNA fragmentation and cell death and inhibits telomerase activity in treated NPC-076 cells. Molecular weight 291.4 + 73.0 g/mol. Solubility: soluble in water. Formula C₁₄H₁₇N₃O₂S . 2 HCl.

References:

Hidaka H, Inagaki M, Kawamoto S, Sasaki Y. (1984) "Isoquinolinesulfonamides, novel and potent inhibitors of cyclic nucleotide dependent protein kinase and protein kinase C." *Biochemistry.* 23(21):5036-41.

Kawamoto S, Hidaka H (1984) "1-(5-Isoquinolinesulfonyl)-2-methylpiperazine (H-7) is a selective inhibitor of protein kinase C in rabbit platelets" *Biochem. Biophys. Res. Commun.* 125(1):258-64

Schachtele C, Seifert R, Osswald H (1988) "Stimulus-dependent inhibition of platelet aggregation by the protein kinase C inhibitors polymyxin B, H-7 and staurosporine" *Biochem. Biophys. Res. Commun.* 151(1):542-7

Engl RA, Girod A, Kinzel V, Huber R, Bossemeyer D. (1996) "Crystal structures of catalytic subunit of cAMP-dependent protein kinase in complex with isoquinolinesulfonyl protein kinase inhibitors H7, H8, and H89. Structural implications for selectivity." *J. Biol. Chem.* 271(42):26157-64.

Ku WC, Cheng AJ, Wang TC (1997) "Inhibition of telomerase activity by PKC inhibitors in human nasopharyngeal cancer cells in culture" *Biochem. Biophys. Res. Commun.* 241(3):730-6

H-8, #KV385


Product name: H-8 . dihydrochloride, H8, H 8

N-[2-(Methylamino)ethyl]-5-isoquinolinesulfonamide . 2 HCl; CAS: 84478-11-5

Molecular weight 265.3 + 73.0 g/mol. Formula C₁₂H₁₅N₃O₂S . 2HCl.

Purity: > 98% Appearance: powder. Solubility: soluble in water and DMSO.

Potent inhibitor of cGMP-dependent protein kinase (PKG, IC₅₀ ~500nM) and cAMP-dependent protein kinase (PKA, IC₅₀ ~1200nM). Also inhibits myosin light chain kinase (MLCK).



References:

Hidaka H, Inagaki M, Kawamoto S, Sasaki Y. (1984) "Isoquinolinesulfonamides, novel and potent inhibitors of cyclic nucleotide dependent protein kinase and protein kinase C." *Biochemistry.* 23(21):5036-41.

Hagiwara M, Inagaki M, Hidaka H. (1987) "Specific binding of a novel compound, N-[2-(methylamino)ethyl]-5-isoquinolinesulfonamide (H-8) to the active site of cAMP-dependent protein kinase." *Mol Pharmacol.* 31(5):523-8

Engl RA, Girod A, Kinzel V, Huber R, Bossemeyer D. (1996) "Crystal structures of catalytic subunit of cAMP-dependent protein kinase in complex with isoquinolinesulfonyl protein kinase inhibitors H7, H8, and H89. Structural implications for selectivity." *J. Biol. Chem.* 271(42):26157-64.


H-89, #CQ7251

Product name: H-89 . dihydrochloride, H89, H 89

N-[2-(p-Bromocinnamylamino)ethyl]-5-isoquinolinesulfonamide . 2HCl; CAS: 130964-39-5

Molecular weight 446.4 + 73.0 g/mol. Formula C₂₀H₂₀BrN₃O₂S . 2HCl

Solubility: soluble in 100% ethanol, DMSO and water.



Info@fluoprobes.com

Technical-support@fluoprobes.com

Order-online@fluoprobes.com

Contact your local distributor

FluoProbes®, powered by



213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

P.5


FT-972971

Cell permeable potent and selective inhibitor of cAMP-dependent protein kinase (PKA, IC₅₀ ~50 nM), cGMP-dependent protein kinase (PKG) and protein kinase C_μ (PKC_μ). In contrast, most other protein kinase C (PKC) isotypes are much more weakly inhibited. Also inhibits ROCKII, Ca²⁺/calmodulin dependent protein kinase II (CaMKII), casein kinase I (CKI) and myosin light chain kinase (MLCK). Induces apoptosis. 1

References:

Chijiwa T, Mishima A, Hagiwara M, Sano M, Hayashi K, Inoue T, Naito K, Toshioka T, Hidaka H. (1990) "Inhibition of forskolin-induced neurite outgrowth and protein phosphorylation by a newly synthesized selective inhibitor of cyclic AMP dependent protein kinase, N-[2-(p-bromocinnamylamino) ethyl]- 5-isoquinoline sulfonamide (H-89), of PC12D pheochromocytoma cells." J. Biol. Chem. 265(9):5267-72

Geilen CC, Wiprecht M, Wieder T, Reutter W. "A selective inhibitor of cyclic AMP-dependent protein kinase, N-[2-bromocinnamyl (amino)ethyl] 5-isoquinoline sulfonamide (H-89), inhibits phosphatidylcholine biosynthesis in HeLa cells" (1992) FEBS Lett. 309(3):381-4.



HA-100, #KV3871


Product name: HA-100 . dihydrochloride, HA100, HA 100

Syn.: 1-(5-Isoquinolinesulfonyl)-piperazine dihydrochloride; CAS number: 105628-72-6
Molecular weight 277.3 + 73.0 = 350.27 g/mol. Formula C₁₃H₁₅N₃O₂S . 2HCl. Form solid.
Purity: > 97%. Form: solid. Solubility: soluble in water.

Inhibitor of myosin light chain kinase (MLCK), cAMP-dependent protein kinase (PKA) and protein kinase C (PKC)

References:

Ref.: Hagiwara M, Inagaki M, Watanabe M, Ito M, Onoda K, Tanaka T, Hidaka H. (1987)
"Selective modulation of calcium-dependent myosin phosphorylation by novel protein kinase inhibitors, isoquinolinesulfonamide derivatives" Mol Pharmacol. (1):7-12



HA-1004, #KV3881

Product name: HA-1004 . dihydrochloride, HA1004, HA 1004

N-(2'-Guanidinoethyl)-5-isoquinolinesulfonamide dihydrochloride; CAS number: 92564-34-6
Molecular weight 293.3 + 73.0 g/mol. Formula C₁₂H₁₅N₅O₂S . 2HCl. Form solid.
Purity: > 99%; Form: solid; Solubility: soluble in water.

Inhibitor of cAMP- and cGMP-dependent kinases (PKA/PKG, IC₅₀ ~1-2 μM). Also inhibits CaMKII, PKC and MLCK at higher concentrations. Intracellular Ca²⁺ antagonist with no effect on cardiac function. Causes selective pulmonary vasodilation during pulmonary hypertension.

References:

Crowley MR, Fineman JR, Soifer SJ. (1994) "HA1004, an intracellular calcium antagonist, selectively attenuates pulmonary hypertension in newborn lambs." J. Cardiovasc. Pharmacol. (5):806-13


HA-1077, fasudil, #KV389

Product name: HA-1077 . dihydrochloride, fasudil, HA1077, HA 1077, AT877
1-(5-Isoquinolinesulfonyl)-1H-hexahydro-1,4-diazepine dihydrochloride; CAS number: 203911-27-7
Molecular weight 291.4 + 73.0 g/mol. Purity: >98%. Formula C₁₄H₁₇N₃O₂S . 2HCl. Form solid.
Purity: >98%; Form: solid; Solubility: soluble in water.

Inhibitor of ROCKII, Rho kinase, PRK2 and MSK1, also inhibits PKA, PKG, S6K1 and MAPKAP-K1b, myosin light chain kinase (MLCK) and Ca²⁺/calmodulin dependent protein kinase II (CaMKII). Cell permeable Ca²⁺ antagonist with antivasospastic properties.

References:

Asano T, Suzuki T, Tsuchiya M, Satoh S, Ikegaki I, Shibuya M, Suzuki Y, Hidaka H. (1989) "Vasodilator actions of HA1077 in vitro and in vivo putatively mediated by the inhibition of protein kinase." Br J Pharmacol. 98(4):1091-100.
Seto M, Sasaki Y, Hidaka H, Sasaki Y. (1991) "Effects of HA1077, a protein kinase inhibitor, on myosin phosphorylation and tension in smooth muscle." Eur J Pharmacol. 195(2):267-72.



Info@fluoprobes.com

Technical-support@fluoprobes.com

Order-online@fluoprobes.com

Contact your local distributor

FluoProbes®, powered by




213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

P.6

FT-972971

- Takizawa S, Hori M, Ozaki H, Karaki H. (1993) "Effects of isoquinoline derivatives, HA1077 and H-7, on cytosolic Ca²⁺ level and contraction in vascular smooth muscle." *Eur J Pharmacol.* 250(3):431-7.
- Seto M, Shindo K, Ito K, Sasaki Y (1995) "Selective inhibition of myosin phosphorylation and tension of hyperplastic arteries by the kinase inhibitor HA1077" *Eur. J. Pharmacol.* 276(1-2):27-33
- Satoh S, Ikegaki I, Suzuki Y, Asano T, Shibuya M, Hidaka H (1996) "Neuroprotective properties of a protein kinase inhibitor against ischaemia-induced neuronal damage in rats and gerbils" *Br. J. Pharmacol.* 118(7):1592-6
- Negoro N, Hoshiga M, Seto M, Kohbayashi E, Ii M, Fukui R, Shibata N, Nakakoji T, Nishiguchi F, Sasaki Y, Ishihara T, Ohsawa N (1999) "The kinase inhibitor fasudil (HA-1077) reduces intimal hyperplasia through inhibiting migration and enhancing cell loss of vascular smooth muscle cells" *Biochem. Biophys. Res. Commun.* 262(1):211-5
- Sward K, Dreja K, Susnjak M, Hellstrand P, Hartshorne DJ, Walsh MP (2000) "Inhibition of Rho-associated kinase blocks agonist-induced Ca²⁺ sensitization of myosin phosphorylation and force in guinea-pig ileum" *J. Physiol.* 522 Pt 1:33-49
- Chrissobolis S, Sobey CG (2001) "Evidence that Rho-kinase activity contributes to cerebral vascular tone in vivo and is enhanced during chronic hypertension: comparison with protein kinase C" *Circ. Res.* 88(8):774-9



Hydroxyfasudil, #KV3900

Product name: Hydroxyfasudil

1-(1-Hydroxy-5-isoquinolinesulfonyl)homopiperazine hydrochloride; CAS number: 155558-32-0


Molecular weight 343.8 g/mol. Formula C₁₄H₁₇N₃O₃ · HCl.

Purity: >98%. Form solid. Solubility: soluble in water. Rho kinase inhibitor.

Active metabolite of fasudil (HA-1077). Exists in two tautomers. Hydroxyfasudil relaxes the rabbit basilar artery mainly by disinhibiting myosin light chain phosphatase through the inhibition of rho-associated kinase. The hydroxyfasudil-sensitive Rho-kinase-mediated pathway appears to mediate enhanced myosin light chain phosphorylations and plays a central role in the pathogenesis of coronary artery spasm.

References:

- Shimokawa H, Seto M, Takeshita A et al. (1999) "Rho-kinase-mediated pathway induces enhanced myosin light chain phosphorylations in a swine model of coronary artery spasm" *Cardiovasc Res.* 43(4):1029-39
- Utsunomiya T, Satoh S, Shimokawa H et al. (2001) "Antianginal effects of hydroxyfasudil, a Rho-kinase inhibitor, in a canine model of effort angina" *Br J Pharmacol.* 134(8):1724-30
- Nakamura K, Nishimura J, Kanaide H et al. (2001) "Hydroxyfasudil, an active metabolite of fasudil hydrochloride, relaxes the rabbit basilar artery by disinhibition of myosin light chain phosphatase" *J Cereb Blood Flow Metab.* 21(7):876-85
- Kandabashi T, Shimokawa H, Takeshita A et al. (2002) "Involvement of rho-kinase in agonists-induced contractions of arteriosclerotic human arteries" *Arterioscler Thromb Vasc Biol.* 22(2):243-8
- Takemoto M, Sun J, Hiroki J, Shimokawa H, Liao JK (2002) "Rho-kinase mediates hypoxia-induced downregulation of endothelial nitric oxide synthase" *Circulation* 106(1):57-62
- Yada T, Shimokawa H, Kajiji F et al. (2005) "Beneficial effect of hydroxyfasudil, a specific Rho-kinase inhibitor, on ischemia/reperfusion injury in canine coronary microcirculation in vivo" *J Am Coll Cardiol.* 45(4):599-607



Iodotubercidin, #KV392

Product name: Iodotubercidin, 5-Iodotubercidin

Chemical name: 4-Amino-5-iodo-7-(beta-D-ribofuranosyl) pyrrolo[2,3-d]-pyrimidine; CAS number: 24386-93-4

Alternate names: 5-I-dTu, 7-Deaza-7-iodoadenosine

Molecular weight 392.2 g/mol. Formula C₁₁H₁₃N₄O₄.

Appearance: light brown solid. Purity: >98% (HPLC). Solubility: soluble in DMSO. Handling: Protect from light.

5-Iodotubercidin is a potent and specific inhibitor of ERK2 (K_i ~0.5 μM) and adenosine kinase (K_i ~30 nM). It has also been reported to abolish the accumulation of AICAR-5'-monophosphate (ZMP) and thus inhibits the activation of AMPK.

References:

- Spychala J, Van den Berghe G (1987) "Adenine nucleotide metabolism in isolated chicken hepatocytes" *Biochem. J.* 242(2):551-8
- Cottam HB, Wasson DB, Shih HC, Raychaudhuri A, Di Pasquale G, Carson DA (1993) "New adenosine kinase inhibitors with oral antiinflammatory activity: synthesis and biological evaluation" *J. Med. Chem.* 36(22):3424-30
- Kroll K, Decking UK, Dreikorn K, Schrader J (1993) "Rapid turnover of the AMP-adenosine metabolic cycle in the guinea pig heart" *Circ Res.* 73(5):846-56
- Flückiger-Isler RE, Walter P (1993) "Stimulation of rat liver glycogen synthesis by the adenosine kinase inhibitor 5-iodotubercidin" *Biochem. J.* 292 (Pt 1):85-91
- Massillon D, Stalmans W, van de Werve G, Bollen M (1994) "Identification of the glycosidic compound 5-iodotubercidin as a general protein kinase inhibitor" *Biochem. J.* 299 (Pt 1):123-8
- Pak MA, Haas HL, Decking UK, Schrader J (1994) "Inhibition of adenosine kinase increases endogenous adenosine and depresses neuronal activity in hippocampal slices" *Neuropharmacology* 33(9):1049-53

Info@fluoprobes.com

Technical-support@fluoprobes.com

Order-online@fluoprobes.com

Contact your local distributor

FluoProbes®, powered by



213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

P.7

FT-972971

Parkinson FE, Geiger JD (1996) "Effects of iodotubercidin on adenosine kinase activity and nucleoside transport in DDT1 MF-2 smooth muscle cells" J. Pharmacol. Exp. Ther. 277(3):1397-401

Gu JG, Nath A, Geiger JD (1996) "Characterization of inhibitor-sensitive and -resistant adenosine transporters in cultured human fetal astrocytes" J. Neurochem. 67(3):972-7












Fox T, Coll JT, Xie X, Ford PJ, Germann UA, Porter MD, Pazhanisamy S, Fleming MA, Galullo V, Su MS, Wilson KP (1998) "A single amino acid substitution makes ERK2 susceptible to pyridinyl imidazole inhibitors of p38 MAP kinase" Protein Sci. 7(11):2249-55

Kowaluk EA, Kohlhaas KL, Bannon A, Gunther K, Lynch JJ 3rd, Jarvis MF (1999) "Characterization of the effects of adenosine kinase inhibitors on acute thermal nociception in mice" Pharmacol. Biochem. Behav. 63(1):83-9




Abbud W, Habinowski S, Zhang JZ, Kendrew J, Elkairi FS, Kemp BE, Witters LA, Ismail-Beigi F (2000) "Stimulation of AMP-activated protein kinase (AMPK) is associated with enhancement of Glut1-mediated glucose transport Arch. Biochem. Biophys. 380(2):347-52

Fleming I, Schulz C, Fichtlscherer B, Kemp BE, Fisslthaler B, Busse R (2003) "AMP-activated protein kinase (AMPK) regulates the insulin-induced activation of the nitric oxide synthase in human platelets" Thromb Haemost. 90(5):863-71


Other inhibitors:

Other Kinase Inhibitors - AOC (inquire)		
	#CJ5841, BMS-354825 (Dasatinib) CAS: 302962-49-8 	#CJ5851, CHIR-265 CAS: 405169-16-6 
#CJ5861, BIBX-1382 di HCl (Falinidamol) CAS: 196612-93-8 	#CJ5871, ZD-6474 (Zactima) CAS: 443913-73-3 	#CJ5881, TAK-715 (P38 MAPK Kinase Inhibitor) CAS: 303162-79-0 
#CJ5891, GW-2016 (Lapatinib; Tykerb) CAS: 231277-92-2 	#CJ5901, Iressa (Geftinib HCl) CAS: 184475-35-2 	#CJ5911, PTK-787 (Vatanalib Succinate) CAS: 212142-18-2 
#CJ5921, Sorafenib Tosylate (BAY-43-9006) CAS: 284461-73-0 	#CJ5931, BIBW-2992 (EGFR Inhibitor) CAS: 284461-73-0 	#CJ5941, Sunitinib Maleate (SU-11248; Sutent) CAS: 557795-19-4 


FT-972971




#CJ5961, Erlotinib HCl
(Tarceva)
CAS: 183319-69-9




#CJ5951, Doramapimod
(BIRB-0796)
CAS: 639089-54-6




#CJ5971, VX-745
(P38 MAPK Kinase Inhibitor)
CAS: 209410-46-8




#CJ5981, PKC-412
(Midostaurin)
CAS: 120685-11-2



#CJ5991, PI-103
(PI3-K Inhibitor)
CAS: 371935-74-9



#CJ5601, JNJ-7706621
(Aurora Inhibitor)
CAS: 443797-96-4



Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use.
FluoProbes® is not liable for any damage resulting from handling or contact with this product.

Rev.H06E

Info@fluoprobes.com
Technical-support@fluoprobes.com
Order-online@fluoprobes.com

Contact your local distributor

FluoProbes®, powered by



213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60