



CybLuc

mg

Product Description

Name :	CybLuc	
Catalog Number :	AXBRCA, 5 mg	AXBRCB, 10
Structure :	C15H15N3O2S2	
Molecular Weight :	333,42	
Properties:	Solubility: 10mM in DPBS	

erchim



Storage: -20°C

Introduction

In vivo imaging is one of the most promising areas of development of Luciferin-Luciferase bioluminescence (1–3). The group of Prof. Minyong Li recently developed CybLuc, a novel Luciferin derivative with longer signal duration and significant bathochromic (red) shift (λ max = 603 nm).

Application

Firefly luciferase bioluminescence is a powerful investigative tool for the understanding of fundamental biological and physiological questions (i.e. for studying tumor growth and inhibition).



211bis Avenue J.F. Kennedy - BP 1140 - 03103 Montluçon Cedex - France S.A. a conseil de surveillance 1.548.000 euro - RCS Montluçon 917 050 171 -APE 4669B

FT-AXBRCA

Úptima

CybLuc is particularly suitable for in vivo imaging, requiring low substrate loading: only 0.01 % of the typical dosage of native Luciferin (dLuc) is sufficient.

The half-life of CybLuc is superior to another red-shifted luciferin derivative, amino luciferin (aLuc), allowing longer circulation time in vivo. In addition, the increased lipophilicity of the compound allows better passage of the blood-brain barrier, which is particularly important for brain tissue imaging (Fig.1).

Advantages of CybLuc at a glance:

o Red-shift of light emission improves signal penetration depth in vivo

- o 20-fold higher bioluminescence signal than D-Luciferin in vivo
- o Stable luminescence signal for 1 hour (detectable as long as 12 hours) in vivo

Technical information

Light emission properties of CybLuc, CycLuc (a commercially available modified amino-luciferin), aLuc and dLuc:



Bathochromic effect of dLuc and aLuc derivatives. Comparison of bathochromic effect of the four Luciferins when exposed to Photinus pyralis luciferase: all three Luc derivatives have significant red-shift compared to native D-Luciferin (dLuc), in particular CybLuc shows an even slightly higher maximum wavelength shift.

Directions for use

The intensity of light emission was also tested in the following conditions: Tris HCl 50 mM with 10 mM MgCl 2 and 0.1 mM ZnCl 2, 1 mM ATP, pH 7.4, 10 μ g/mL recombinant P. pyralis luciferase, 1 μ M luciferines, showing that the emission intensity of CybLuc is in the same order of magnitude as CycLuc and significantly stronger than aLuc.

A. CybLuc for in vitro bioluminescent assays

Prepare aqueous solutions in ultrapure water (e.g. Milli-Q water).

(1.) Prepare a fresh 200X CybLuc stock solution (4 mM): Dissolve 1.3 mg/mL CybLuc in ultrapure water. Mix gently until CybLuc is completely dissolved.

 SCIENCES
 211bis Avenue J.F. Kennedy - BP 1140 - 03103 Montluçon Cedex - France
 S.A. a conseil de surveillance 1.548.000 euro - RCS Montluçon 917 050 171 - APE 4669B



FT-AXBRCA

(2.) Prepare freshly the CybLuc working solution: dilute 200X CybLuc solution 1:200 in pre-warmed cell culture medium or normal saline (NS) to give 20 μ M final concentration. e.g. to prepare 5 mL of substrate solution, dilute 25 μ L of 200X CybLuc stock solution in 5 mL of medium or NS.

(3.) Aspirate media from cultured cells.

(4.) Add an appropriate volume of CybLuc working solution to cells or tissue culture. e.g. on 96-well plates use 50 μ L per well

(5.) Measure the bioluminescence immediately after the addition of the substrate with an acquisition time of 1 or 20 s. It is recommended to use a recording system that is equipped with a cooled CCD camera for bioluminescent imaging at 37° C.

B. CybLuc for in vivo bioluminescent assays

Prepare all aqueous solutions in ultrapure water (e.g. Milli-Q water).

(1.) To prepare a fresh 10 mM CybLuc solution dissolve 3.3 mg CybLuc (C-8790) in 1mL DPBS (Dulbecco's phosphate-buffered saline, without Magnesium and Calcium) or in NS (normal saline).

Mix gently by inversion until CybLuc is completely dissolved. Filter sterilize through a 0.2 µm filter.

(2.) Inject 10 μ L of the 10 mM CybLuc Solution per gram of body weight. For example, inject intra-peritoneal (i.p.) 100 μ L for a 10 g mouse.

(3.) Bioluminescence images can be acquired 10 minutes after i.p. injection* with 1 second acquisition time.

* For signal optimization, CybLuc kinetic studies should be performed for each animal model to determine the peak signal time.

References

 Prescher J.A.; Steinhardt R.C.; Paley M.A.; McCutcheon D.C. J. Am. Chem. Soc. 2012, 134 (18), 7604-7607.
 Miller S.C.; Prescher J.A.; Aronin N.; Paley M.A.; Reddy G.R.; Adams Jr. S.T.; Chaurette J.P.; Evans M. S. Nat. Methods 2014, 11, 393-395.

3. Miller C.; Adams S.T. Jr.; Curr. Opin. Chem. Biol. 2014, 0, 112-120.

4. Wu W.; Su J.; TangC.; Xaixiu B.; Ma Z.; Zhang T.; Yuan Z.; Li Z.; Zou W.; Zhang H.; Liu Z.; Wang Y.; Zhou.; Du L.; Gu l.; Li M. Anal. Chem., 2017, 89,4808–4816.

Ordering information

Catalog size quantities and prices may be found at <u>http://www.interchim.com</u>. Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06 Order on-line or Contact your local distributor

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product.

