

Innovative and remarkable chemistries, conjugation methods, labeling and functionalisation (PH)

| <u>Hydrazone chemistry Introduction</u> | More info | <u>Conjugation kit</u> | <u>Peptide/Olig</u> | <u>Tagging</u> | <u>Immob.(beads)</u>| <u>Crosslinkers</u>

 Peptide/Oligos
 Antibody labeling
 WB RapidDirect

 Crosslinkers
 Organic Synthesis

■ Hydrazone Chemistry (Hydralink, Controlled Amine) - Introduction

The new technology to conjugate simply and efficiently your biomolecules

Advantages : efficient and flexible!

- **easy** operating (modify each molecule, mix),
- activate the biomolecules in advance (stable for months),
- **better control** of the coupling ratio,
- excellent yield of conjugation,
- no reduction or de-protection step,
- orientated heteroconjugation highly selective (no homoconjuguates),
- conjugation keeping the bioactivity of the components, and very stable,
- **flexibility** for many applications / samples.



ates),

This technique replaces advantageously the standard methods based on NHS/Maleimide, epoxides, glutaraldehyde, ... and is more flexible and satisfying for various applications, i.e. for conjugation in solution to solid phase synthesis. Its suits perfectly proteins and any biomolecule or support containing amines or derivatized by conventional biochemistry (proteins, peptides, oligonucleotides, cDNAs, carbohydrates, fluorophores, beads, glass, silica ...).

This technology is offered as <u>stand-alone reagents</u> for organic synthesis, peptides, antibodies, nucleic acids conjugation and labeling, but also in convenient conjugation or labeling kits:

All purpose controlled amine conjugation kit | Peptide-oligos conjugation kits | Antibody labeling kits (HRP, AP, RPE, APC)

Immobilisation/supports (agarose beads, magnetic beads)

label/conjugate proteins
 immobilisation on supports

■ Hydralink[™] ControlledAmine Protein conjugation kits

General use linking kits for proteins, peptides, aminated oligos, beads or surfaces. Everything included, no experience needed.

The S-HyNic Conjugation Kit can be used to covalently crosslink any two aminecontaining biomolecules. HydraLink[™] conjugation chemistry is based on the reaction of a HyNic linker with a 4FB linker to form a stable hydrazone bond. The bond created is both stable and UV-traceable. This unique covalent bond is created when the HyNic linker, incorporated into one type of biomolecule reacts with the 4FB linker, incorporated into the second biomolecule. S-HyNic and S-4FB react only with each other and not with other protein functional groups (selective, oriented and bioorthogonal reaction).

 All Purpose SANH Conjugation Kit (S-HyNic)
 BL150A, 1 kit

 Contains SANH and S-4FB activators, and all needed reagents for conjugating to amine-containing molecules. Technical sheet
 All Purpose SHTH Conjugation Kit
 BL152A, 1 kit

 Contains SHNH and S-4FB activators, and all needed reagents for conjugating to amine-containing molecules. Technical sheet
 BL152A, 1 kit

■ Hydralink[™] Protein-Oligos conjugation kits

Innovative conjugation chemistry to prepare proteins-oligonucleotide conjugates efficiently and easily.

- Easy protocol: activate, mix, clean-up
- >95% conversion
- UV_{354 nm}-traceable stable ligation
- The Antibody-Oligonucleotide

All-in-One Conjugation Kit includes buffers, spin columns, and a calculator to determine MSR.

No column chromatography is required.

Advantages

Simple and Easy to Use – Requires only µpipette, µcentrifuge, and UV spectrophotometer
Automated Calculations – Calculator with fully integrated input/output provided
High Yield – 30 - 95% yield based on starting protein (depends on kit/protein)
Faster kinetics for greater efficiency and yields thank to catalyzed conjugation
High Purity – >95% purity without chromatographic purification (kit FV9100 only)
High Stability – Conjugates are 10 times more stable than any other conjugation linker
Specificity – Two-linker method avoids homoconjugate formation
Quantifiable – Using a UV signature wavelength and simple UV scan

It generates high-purity conjugates virtually free of residual antibody or oligonucleotide (>95% conjugate).

Ask for White Paper: Antibody-Oligonucleotide Conjugate Preparation

How it works?

1. Amine-modified, 20 to 60-mer oligonucleotide is modified using an excess of the Sulfo-S-4FB linker2. Polyclonal or monoclonal antibody (100 μ g) is modified using the S-HyNic linker.

3. The 4FB and HyNic modified biomolecules are mixed together in the presence of the TurboLink[™] catalyst, leading to rapid and conjugation through formation of stable bis-arylhydrazone bonds, followed by magnetic-affinity, solid phase purification.

4. The antibody-oligonucleotide conjugate is ready for use in the most demanding and sensitive applications.

Application: Oligonucleotide-Protein conjugation

PAGE gel demonstrating a 5-aldehyde modified oligonucleotide (1 equivalent) was reacted with a 15mer peptide that was modified by C6-HNA at N-terminus). Simple addition of the hydrazine-modified peptide to the aldehyde-modified oligonucleotide (lane 1) directly yielded the peptide/oligonucleotide conjugate without the requirement of reducing reagents.



• The *Protein-oligos conjugation kit* is similar, and contains the material to synthesize two conjugation reactions for any immunoPCR or hybridization-type assay in just under 4 hours each, yielding between 40–60% conjugate after purification. 2 features facilitate easy purification and identification of the conjugate from the excess oligo using size exclusion purification methods such as FPLC or diafiltration: 1/the ligation consist in a **stable bond** (bis-arylhydrazone) with measurable **absorbance at 354 nm**. 2/it yields >95% conversion of protein to conjugate when more than 4 molar equivalents of oligo are added.

Antibody-Oligonucleotide All-in-One Conjugation Kit *:the kit conjugates 100ug of antibody; <u>Technical sheet</u> **Protein-Oligos All-in-One conjugation kit** *:the kit provide 2 conjugations reactions, <u>Technical sheet</u> FV9100, 1 kit*

Prices and technical sheets on-line

IV2490, 1 kit



Antibody labeling kits

■ HydraLink BioConjugation & Labeling kits (HRP, Fluorescein, PE,...)

The S-HyNic Conjugation kits are designed to conjugate proteins with pre-activated high-activity markers such as HRP, PE, Biotin. Any suitably pure and sized molecule containing amines (antibody, aminoallyl-oligo can also be conjugated and purified in ~4 hours (30 minutes hands-on time).

All-in-One Conjugation kits contain the activated reagent for HyNic activation of 2 x 100 µg of any user-supplied antibody, the activated marker, and dual purification tool (affinity beads).

One-Shot Conjugation kits contains same reagents but purification tool is Gelfiltration column, while conjugation level can be assessed.

These kits use the unique HydraLink[™] conjugation chemistry: the protein to label is activated as a HyNic group using SANH or SHTH reagents, while the marker is activated as an Aldehyde using SFB reagents. This allows more flexibility and control than conventional methods (glutaraldehyde, Mal/NHS, periodate...). Then HyNic and 4FB:

- react only with each other (no interference of chemical groups found in biomolecules) in mild conditions (gentle chemistry)

- create a unique covalent bond that is both stable and UV-traceable (@354nm).





HRP All-in-One Antibody Labeling Kit	FK1640, Kit for 2x100µg Ab
" (Large scale kit)	FK1641, Kit for 5 mg Ab
Contains all reagents to conjugate antibody with 4FB pre-activa	ated Horseradish Peroxidase*. Technical sheet
AP All-in-One Antibody Labeling Kit	FK8880, 1 kit (for 2x100µg Ab)
Contains all reagents to conjugate 2x100µg of antibody with pr	e-activated Alkaline Phosphatase. Technical sheet
R-PE All-in-One Antibody Labeling Kit	FK8900, 1 kit (for 2x100µg Ab)
Contains all reagents to conjugate 2x100µg of antibody with pr *: with affinity magnetic beads + spin filter (for AP Lab.Kit FK88	e-activated R-Phycoerythrin*. <u>Technical sheet</u> 880: 30K spin filter alone)r

R-PE Antibody Labeling Kit JO2380, 1 kit (for 2x500µg Ab) Contains all reagents to conjugate 2x~0.5mg of antibody with pre-activated R-Phycoerythrin and 0.5ml gelfiltration columns. Technical sheet APC Antibody Labeling Kit **FK8890**, **1** kit (for 2x500µg Ab) Contains all reagents to conjugate 2x~0.5mg of antibody with pre-activated Allophycocyanin and 0.5ml gelfiltration columns. Technical sheet Fluorescein One-Shot Antibody Labeling Kit RJ3201, 1 kit (for 2x100µg Ab) Contains all reagents to conjugate 2x~100µg of antibody with NHS pre-activated Fluorescein and 0.5ml gelfiltration columns. Technical sheet

Associated reagents:

RapidDirect[™] Primary Antibody polyHRP **1I6151** (for 2x100µg Ab) Contains all reagents to conjugate 2x100µg of antibody with pre-activated HRP See technical sheet Modify your primary antibody directly with polyHRP for superior sensitivity. No purification. No secondary antibodies needed, hence minimal background,, eliminate II Ab antibody incubation steps, thereby reduce the time required to complete the protocol, and eliminate any cross-species contamination. Ideal for IP/WB applications, and for those who use the same antibody for multiple western blots with the additional benefit of time savings

TurboLink Buffer	HT1820, 1.5ml
Catalyze the HyNic/4FB reaction to yield stable hyd	Irazonie link.
HRP - 4FB	Inquire
4FB-modified Horse Radish Peroxidase - Conjugat	es retain their high enzymatic activity.
R-PE - 4FB	Inquire
4FB-modified R-PhycoErythrin in 100mM Sodium F	Phosphate 150mM NaCL, pH6.0
Fluorescein-PEG ₃ -4FB	Inquire
Conjugates retain their high quantum yield (~0.98)	and are easy to purify using gel filtration methods





100% purified IgG-HRP conjugate

AP - 4FB Inquire 4FB-modified Alkaline Phosphatase - Conjugates retain their high enzymatic activity APC - 4FB Inquire

4FB-modified AlloPhycoCyanine in 100mM Sodium Phosphate 150mM NaCL, pH6.0.

211 bis Av. J.F. - BP 1140 03103 Montluçon Cedex - Tel. 33 (0) 4 70 03 88 55 - Fax 33 (0) 4 70 03 82 60 e-mail interchim@interchim.com - web www.interchim.com **Minterchim**

■ RapidDirectTM Western Blot Kits

The NEW RapidDirect[™] Western Blot Products, using patented linking technology, attaches multiple HRP proteins to any user-supplied primary antibody, offering several advantages over the use of secondary antibodies for western blot detection.

Publication of	quality gels	Solulink Direct Western Blot vs.	Classical 2° Western Blot
the firs	d lime!	Conjugate HRP to Antibody enough for 15 – 20 blots 4 h (30 min) time for up to 20 blots	
1 2	3 4	Run Gel Transfer to Men	ubrane. Blot Gel
		(step duration) 3 h	5 min) (hands on)
		1° Ab-polyHRP In	cubation
		1 h 🔶 (5 Wash	min)
			1 h 🔰 (5 min)
	and the second s	ECL Development	2° Ab-polyHRP Incubation
	Contraction of the second	Wot	30 min 🔰 (5 min)
	mining married	5 hours per bank)	Wash
and the second second second		Save 1. Hes hands on	15 min 🔰 (2 min)
	and the second se	(12 minute	ECL Development
		Total 4 h 15 min (22 min hands-on) 1 incubation	Total 5 h 45 min (34 min hands-on) 2 incubations
The second s	Contraction Contraction (5.0 pt 2.5 pt 1.25 pt 0.0625 pt	5.0 pl 2.5 pl 1.25 pl 0.0625 pl
The street suggest star			
	Contraction of the second s		
	Statements Statements	Sector Sector Sector Sector	And States and States
And the second second	And the second second second		
and the second second second	States and states	Heavy chain	-
and the second second second	Strength 20080950	■ ■ ■ ■ ■ ■ ■ ■ ■ ■ ■ ■ ■ ■ ■ ■ ■ ■ ■	
	<i>A A A</i>		
Solulink	Conventional	Light chain	
Western Blots	Western Blots		

- Publication quality gels the first time •
- Far less background does not detect any antibody contaminants on blot, such as light and heavy chains
- Sensitivity better or equal to the classical secondary antibody methods in ALL immunoassay formats •
- Significant time savings of 1.5 hours per blot or IHC tissue staining •

Select from Two Versions of Western Blot Kits

■ RapidDirectTM Primary Antibody polyHRP Western Blot Kit

For standard western blots and other immunoassay applications. Modify your primary antibody directly with HRP. No purification. No secondary antibodies needed for downstream applications. Ideal for those who use the same antibody for multiple western blot or immunoassay applications with the additional benefit of time savings.

RapidDirect Primary Antibody polyHRP<u>IP</u>/Western Blot Kits

IP/western techniques are typically used to detect low-copy antigens. The primary antibody-HRP/antigen complex is immobilized by α species IgG immobilized on NanoLink[™] magnetic beads (pull-dwon assay), then eluted wit hantigen and analyzed by WB.

Ask for more information about these kits

RapidDirect[™] Primary Antibody polyHRP

Contains all reagents to conjugate 2x100µg of antibody with pre-activated HRP See technical sheet Modify your primary antibody directly with polyHRP for superior sensitivity. No purification. No secondary antibodies needed, hence minimal background,, eliminate II Ab antibody incubation steps, thereby reduce the time required to complete the protocol, and eliminate any cross-species contamination

Ideal for IP/WB applications, and for those who use the same antibody for multiple western blots with the additional benefit of time savings

RapidDirect Primary Antibody polyHRP IP/Western Blot Kits / pull-down assay			
with goat anti-mouse	JO1990	(inquire)	
with goat anti-rabbit	JO2000	(inquire)	
with rabbit anti-goat	JO2010	(inquire)	
*Containe all reasons for 10.40 Immunopresinitations and Western blats analysis			

*Contains all reagents for 10-40 Immunoprecipitations and Western-blots analysi

Prices and technical sheets on-line

211 bis Av. J.F. - BP 1140 03103 Montluçon Cedex - Tel. 33 (0) 4 70 03 88 55 - Fax 33 (0) 4 70 03 82 60 e-mail interchim@interchim.com - web www.interchim.com **Minterchim**

116151 (for 2x100µg Ab)

Hydralink tagging

Hydralink technology* allows for tagging chemically any purified biomolecule, extract or supports to yield easy and complete labeling. Affinity tags such as c-Myc or TAT binding moieties are available activated by the hydralink technology. Bioactive peptides are also available on inquire at interbiotech@interchim.com.

* See the Hydrazone chemistry principle (HyNic and SFB reactions).

Digoxigenin ChromaLink Labeling

DIG label your Antibody with greater confidence, consistency and reproducibility

UV-Traceable Digoxigenin Linker allows simple UV measurement for fool-proof incorporation of Digoxigenin into any 1' or 2' antibody from any source (Rabbit, mouse, human, IgY)



ChromaLink Digoxigenin OneShot Ab Labeling Kit Complete kit to label 100%g of antibody in just 2 hours. Techn

ChromaLink Digoxigenin reagent

Inquire



Application: Multi-Color Immunofluorescence Technique using Primary Antibodies raised in the Same Host Species

Triple-labeling immunofluorescence in IHC was achieved using 3 different primary antibodies derived from a single host source:

-one primary antibody with biotin (ChromaLink Biotin #BT3614)

-a second primary antibody with DIG (ChromaLink Digoxigenin #DW133A).

Hydralink tagging

Tags such as c-Myc or TAT binding moieties are available activated by the hydralink technology*, to ensure easy and complete chemical labeling when the genetic expression tagging (cell systems) is not possible or suitable. They can be used for labeling peptides or other biomolecules to be detected using anti tag antibodies, for controls on cells expressing the tag, for preparing tag-affinity supports...

* See the Hydrazone chemistry principle (HyNic and SFB reactions).

•]	Penetrating	peptides
-----	-------------	----------

HyNic-TAT FK5491,	0.5 mg FK5492, 1	mg
HyNic-(Arg)ଃ FK5481,	0.5 mg FK5482, 1	mg

Protein expression reporter peptides:

HyNic-(His)₀ Tag HyNic-S-Tag	FK5471, 0.5 mg	FK5472, 1 mg FK5502, 1 mg		
HyNic c-Myc Tag	FK5511, 0.5 mg	FK5512, 1 mg HyNic-Flag Tag	FK5531, 0.5 mg	FK5532, 1 mg

Bioactive peptides HyNic or 4FB conjugated : (Ghrelin, Obestatin,...) **Custom PepLink Peptides** Inquire at interbiotech@interchim.com

Lipids 4FB-DSPF FK5671, 5x20 mg 4FB/disteroyl-phosphatidylethanolamide

Polated products: other tags			
Related products, other tags			
Biotin tag: see section 'biotin labeling'	, i.e.:	Chelate / poly-	His tag : see section 'chelate labeling', i.e.:
Biotin-PEO ₄ -NHS	R2027A, 50 mg	Maleimido-C ₃ -NTA	T3212A, 10 mg
Reactive with amino groups		Reactive with thiol groups.	
AB-NTA free acid	BE8210, 100mg		
Reactive with aldehydes and NHS groups			

Other tags: see the section 'protein expression' for transfection methods

Anti tags: see the section 'Secondary antibodies' for labeled anti tag antibodies.

tag detection reagents: see in proteins analysis section 'Fusion tags detection'.

Hydralink supports: resins and beads

Hydralink fuunctionnalized magnetic beads

NanoLink[™] Magnetic Beads are polymer encapsulated, 800 nanometer-sized, super-paramagnetic particles. MagnaLink[™] Magnetic Beads are polymer encapsulated, uniform, 2.8µm-sized, super-paramagnetic particles.

NanoLink and MagnaLink beads have the most consistent, ultra-high binding of any particle on the market (>12 nmol/mg), thanks to the proprietary coupling technology. The high surface area and lower non-specific binding of both Beads are ideal for immobilizations and Co-IP applications. These microspheres are particularly suited for high throughput robotic applications where high biotin loads must be removed or immobilized without the presence of an iron leachage.

Applications include capture and immobilization of biotinylated biomolecules, such as DNA, RNA, peptides, proteins and antibodies.



IO6430, 1mL NanoLink 4FB Magnetic Beads 0.8µm, at 10mg/mL NanoLink Amino Magnetic Beads 0.8µm, at 10mg/mL IO6440, 1mL CP5530, 1mL NanoLink Streptavidin Magnetic Beads 0.8µm, at 10mg/mL CP5531, 5ml CP5532, 10ml MagnaLink 4FB Magnetic Beads 2.8µm, at 10mg/mL IO3320, 1mL MagnaLink Amino Magnetic Beads 2.8µm, at 10mg/mL 106450, 1mL AWJ460, 1mL MagnaLink Streptavidin Magnetic Beads 2.8µm, at 10mg/mL AWJ461; 5ml AWJ462, 10ml Magnetic Stand (holds 1.5 mL tube for separation of magnetic particles) L77191, 1u DZ1381, 10 Slides

HyLink[™] Glass Slides, HyNic-activated

Hydralink immobilization on magnetic beads NanoLink BeadLink Kit MagnaLink BeadLink Kit

Streptavidin Agarose beads:



MagnaLinkTMand NanoLinkTM streptavidin beads:

106470, 1mL

IO6470, 1mL

Free-Biotin Binding Capacity



- •Highest biotin-binding capacity on the market
- •Binding capacity can be "dialed in"
- •Lower background noise
- •Available versions:
- Streptavidin - Amino
- 4FB

Hydralink immobilization on Agarose beads

Agarose Gel activated by the hydralink technology are available for easy immobilisation of HyNic derivatised ligands, for further capture, depletion or affinity purifications.

Agarose – 4 FB	IV2510, 1ml	
High binding of HyNic derivatised molect	cules	
Agarose – Streptavidin	IO6460, 2mL IO6461, 5mL	IO6462, 10mL
High biotin-binding streptavidin agarose	(>330nmol biotin/ml gel)	

More: see the section 'Purification by affinity' see the section 'Immobilization of biomolecules'

211 bis Av. J.F. - BP 1140 03103 Montluçon Cedex - Tel. 33 (0) 4 70 03 88 55 - Fax 33 (0) 4 70 03 82 60 e-mail interchim@interchim.com - web www.interchim.com **Pinterchim**

Hydrazone Crosslinkers – for biochemistry

(see <u>principle and benefits</u> of this chemistry).

_

-

_

All these reagents are packaged under dry argon, and to store at -20°C for long term. See the Technical sheet

SANH Succinimidyl 4-Hydrazinonicotinate Acetone Hydrazone; MW: 2	BL9270 10 mg 290.2	BL9271, 25 mg	Su0 ~~~ N ~~
Used to convert primary amines to hydr of the hydrazine is required. The protect hydrazone conjugate	azinopyridine moietio ting group leaves dur	es where protection ing formation of the	о́н с _м щ _м м ц
C6-SANH C6-Succinimidyl 4-Hydrazinonicotinate Acetone Hydrazone. M	BL9330 10 mg W: 403.4	BL9331, 25 mg	
Used to convert primary amines to hydr six carbon linker where protection of the group leaves during formation of the hydr	azinopyridine moietio e hydrazine is require drazoneconjugate.	es with an extended ed. The protecting	
SHNH Succinimidyl Hydraziniumnicotinate Hydrochloride: MW: 286.7	BL9360 10 mg	BL9361, 25 mg	
Used to convert primary amines to hydr chelates99mTc.	azinopyridine moieti	es. Also	
SHTH	BL9370 10 mg	BL9371, 25 mg	0
Succinimidyl 4-Hydrazidoterephthalate.Hydrochloride; MW: 313 Used to convert primary amines to arom	^{3.7} natic hydrazide moiet	ies.	
SATH	BL9390 25 mg		0
Used to incorporate 4-hydrazidoterephth amine-containing moieties.	halamide moieties on	proteins or other	
SFB	M11771 100mg		0
Used to convert primary amines to benz	aldehyde moieties.		SuO H
C6-SFB 360.4	BL9410, 25mg C6-Succir	nimidyl 4 –formylbenzoate; MW:	
Used to convert primary amines to benz carbon linker.	aldeyhde moieties wi	ith an extended six	
SUIfoSuccinimidyl 4-formylbenzoate; MW: = 327.3	BI1311 25mg		
The sulfo derivative of SFB, for improv for conversion of amino surfaces such a	ed water solubility. E s beads and plates.	Especially useful for	
PEO4-SFB (PEG4-SFB) SulfoSuccinimidyl 4-formylbenzoate; MW: = 563.5	CP3981 10mg		
This derivative of SFB contains an exter alleviates steric hindrance and increase	nded and hydrophilic the conjugate hydrop	PEO spacer, taht hilicity and stability.	
SS-SFB Succinimidyl-4-formylbezoate; MW 410.5	CP3991 10mg		
This derivative of SFB contains a S-S link in	n the spacer allowing th	iol mediated cleavage.	

MHPH BL9400 10 mg 5-Maleimido-2-hydraziniumpyridine hydrochloride; MW: 204.2

Used to convert thiol moieties to hydrazinopyridine moieties. Its advantages include thiol-reaction specificity, UV-traceability. 2-hydrazinopyridine-modified proteins then react to form stable conjugates in the presence of other (aromatic) aldehyde-modified proteins or biomolecules.

BZ0770 10mg

BL9420, 25 mg

Used to convert thiol moieties to 4FB (4-formylbenzamide) moieties. Possess a PEG3 linker for increased solubility of modified biomolecule.

Used to incorporate hydrazinopyridine moieties on silica or glass surfaces

Also available: Hydrazine Silane Coated Glass Slides #DZ1381, 10u			
2SBA (2-Sulfobenzaldehyde)	A42050, 100mg		
Used to quench or cap hydrazone	conjugation reactions		

TurboLink Catalist Buffer HT1820 Speeds up de reaction between 4FB and HyNic.

MTFB

MW: 503.5

MW: 396.4

Hydrazine-silane

Hydrazone linkers for organic synthesis

The following reagents are designed for modifying biomolecules and oriented conjugation according hydrazone chemistry. (see principle and benefits of this chemistry). See the Technical sheet

6-FMOC-HNA BL9740, 100mg BL9741, 500mg

6-FMOC-hydrazinonicotinic acid; MW: 375.2 Used to incorporate hydrazine moieties during solid or solution phase peptide synthesis.

6-FMOC-HNA-Osu BL9760, 100mg BL9761, 500mg

Succinimidyl 6-FMOC-hydrazinonicotinate; MW: 472.2 Used to incorporate hydrazine moieties during solid or solution phase peptide synthesis.

6-BOC-HNA BL9750, 100mg BL9751, 500mg

6-BOC-hydrazinonicotinic acid: MW: 253.1 Used to incorporate hydrazine moieties during solid or solution phase peptide synthesis.

6-BOC-HNA-OSu BL9770, 100mg BL9771, 500mg

Succinimidyl 6-BOC-hydrazinonicotinate; MW: = 350.3 Used to incorporate hydrazine moieties during solid or solution phase peptide synthesis.

C6-HNAA BL9780, 100mg BL9781, 500mg

6-hydrazinonicotinic acid acetone hydrazone; MW: 306.4 Used to incorporate protected hydrazine moieties with extended six carbon linker during peptide synthesis

- HNA BL9790, 100mg 6-Hydrazinonicotinic Acid; MW: 153.1 Precursor molecule.
- **BOC-HTA** BL9810, 100mg BL9811, 500mg (4-BOC-hydrazido) terephthalic acid; MW: 280.3
- Precursor molecule. **BOC-HTA-OSu** BL9820, 100mg BL9821, 500mg
- Succinimidyl 4-BOC-hydrazido) terephthalate; MW: 377.4 Precursor molecule.
- 4FB-Amidite

FV7031, 250mg Building block for nucleic acids synthesis / modification.







n



BL9401 25 mg

BZ0774 5x1mg

More about Hydrazone chemistry

Principle:

The method includes one (or 2 separate) activation step(s) of amine, thiols, or silanol to aromatic functionalities, aldehyde (4SFB) and hydrazine (HyNic) able to form an hydrazone bond. The level of activation is fully controllable using chemical groups quantitation reagents (4NBA #BL9650 for hydrazines, 2HP #019022 for aldehydes). Then the modified molecules are simply mixed to yield a stable conjugate. The hydrazone bound formed is fully stable, in contrast to the hydrazone formed by more commonly accessible hydrazides (unstable acyl hydrazones). It is the only known example of a stable Schiff base, which requires no additional steps (reduction) to stabilize the bond.



Applications:

Hydrazone chemistry provides additional benefits for surface modifications, compared to conventional methods such a glutaraldehyde (that is not oriented chemistry) or other aldehyde (4hours procedure), SMCC (Mal-NHS) (that has poor yield), CNBr activation (that yields a charged bond): it is specific and flexible, controllable activation steps, UV-traceability, and quick (2 hours).

HydraLink conjugation system is a privileged new method to conjugate and immobilize a variety of biomolecules, including peptides, peptides, carbohydrates and nucleic acids. It is covered by US patent 5.206.370. 5.420.285. 5.753.520 and 5.769.778. and EU Patent 0.384.769). The involved reaction is highly selective and mild, derivatized molecules are stable and non-susceptible to non-specific binding. Hydrazone link is kinetically and metabolically stable analog of a cysteine bridge. It is a superior alternative to step-wise methods, in which difficulties come from the need to separate the modified or conjugates molecules from excess modifiers, unavailable reactive groups requiring tedious activation steps or labor-intensive site specific engineering methods, from intramolecular undesired crosslinking, and from heterogeneous conjugation at the molecular level.

Comparison of conjugation methods:					
	Hydrazine/carbonyl	Avidin/Biotin	Maleimide/Thiol		
Stability of Activated biomolecules	+++	++	-		
High selectivity	+++	+++	+		
No reticulation	+++	+	-		
No undesirable covalent modification	+++	+++	-		
No non-specific binding of conjugate	++	+	++		
No need of reductant	+++	+++	-		
Covalent (stable) linkage	++	(non cov)	++		
Fast reaction kinetics	+++	+++	++		
Suitable to a variety of biomolecules&support	+++	+	+		
Amenable to solid phase synthesis	+++	++	+		
Reproducible/adjustable coupled ratio	+++	++	+		
Scalable	+++	-	++		
pH range (optimum)	3-(4.7)-7, 5-11	4-7.5			

The technology is ideal for:

-protein to protein conjugation (see <u>Conjugation Kit #BL150A</u>)

-peptide to oligonucleotides conjugation (see Conjugation Kit #IV2490)

-immobilization of peptides and nucleic acids onto microarrays, microplates. (see here)

-organic synthesis of peptides and nucleic acids (see here)

-other biomolecule conjugation (see application below)

-labeling biomolecules by enzymes, fluorophores (see All-In-One kits) or tags (see here).

Related products lines

Interbiotec - BioSciences innovation – proposes a complete range of products for protein biochemistry.

Innovative and remarkable chemistries, conjugation methods, labeling and functionalisation
 Standard Click Chemistry reagents
 (PH)

- •
- Copper-free Click Chemistry reagents (DBCO reagents) (PH) Staudinger reaction (effective conjugations/chemical modification) (PH) •
- **PEGylation reagents** (conjugation reagents, linkers and building blocks) ^(PH) **SDA reagents** (effective photo reactions) ^(PH) .
- •
- STELLA labeling (azocycloaddition reactions) (PH) •
- SAM reagents (Self-Assembled Monolayers for surface modification) (PH) •
- Gold nano-particules and materials
- Carbone nanotubes
- ITO slides (PH)
- FluoProbes labeling agents

Desalting tools – CelluSep tubings, SpectraPor tubings, GebaFlex, FloatALyser, SlideALyser,...

Products HighLights Overview

Information inquire

Reply by Fax : +33 (0) 4 70 03 82 60 or email at interbiotech@interchim.com

I wish to receive the complete documentation about:

Name:		2 nd name:		_ Position:	
Company/Institute:			Service, Lab:		
Adress:					
	Zip code:	Town:			
	Tel	Fax	Email:		