

Innovative biochemistry conjugation method from Interchim

Click Chemistry (Azide / alkyne reaction)

A versatile and reliable conjugation chemistry for linking covalently in very mild conditions

Click chemistry can be used for the synthesis of a variety of conjugates. Virtually any biomolecules can be involved, and labeling with small molecules, such as fluorescent dyes, biotin, and other groups can be readily achieved.

<u>Chemistry</u>

The Click chemistry involves the reaction between an azide and an alkyne (i.e. acetylene), forming a covalent chemical bond.



This Cu(I) catalyzed [3+2] azide-alkyne cycloaddition is also known as CuAAC reaction. It yelds the efficient formation of a non-toxic triazole from biological building blocks that have been modified with non-perturbing azides and unactivated alkynes.

This process has unprecedented tolerance and reliability. It is pH-independent, and it can be carried out in water at ambient temperature. It applies to conjugation and solid phase immobilization. *Note:* both azido and acetylenic groups are nearly never encountered in natural biomolecules. Hence, the reaction is **highly bioorthogonal and specific**.

The only needed things are azido- and alkyne-labeled biomolecules, dyes, and surfaces. Interchim BioScience provides a variety of reagents, notably in the Uptima range.

• <u>Standard Click Chemistry reagents (Alkyne reagents)</u> - CuAAC⁽⁻⁾

uses alkyne ligators and azide ligators, to activate i.e. amines (ZL5530 & ZL5540), or nucleic acids (with amidite-Alkynes) Alkyne and Azide react in presence of Copper(II)-TBTA complex <u>FY2780</u>

• Copper-free Click Chemistry reagents (DBCO & BCN reagents) - SPAAC (PH)

use cyclooctynes ligators (DBCO, BCN) and azide ligators, to activate i.e. amines or other biomolecules DBCO/BCN and azide partners react directly (SPAAC reaction).

Uptima and FluoProbes NHS reagents are useful to incorporate an alkyne, an azide or directly a label on a aminated molecule

Reagents (for CuAAC)

- Alkyne reagents: Fluorescent Fluor dyes⁰ Nucleotides⁰.
 Alkyne-activated ligands are used for CuAAC Click chemistry.
- Azide reagents: Fluorescent FluoProbes dyes⁰, Cy dyes⁰ Non Fluorescent labels (biotin) and ligands Nucleotides⁰.
 Azides are useful reagents for ligation of ligands via Click chemistries, <u>Hydrazone ligation</u>, or via <u>Staudinger ligation</u>.
- **Complementary** reagents for click chemistry: Azide ou Alkyne NHS⁰ to activate amines(proteins)
- Auxillary reagents for click chemistry: CLICK labeled molecules Capture Kit Catalyzers and buffers ⁰
- Please contact interbiotech@interchim;com for other molecules to bz custom-modified by acetylene or DBCO

■ Alkynes Click Chemistry reagents

• Fluorescent Dyes - Alkynes: see FT-DQP790

Alkyne-PEO ₄ -CR110 Fluor 488-Acetylene (5/6-CarboxyFluoresc	cein); Abs/Em = 501/525 nm	FP-DQP790, 1mg / 5mg	
Alkyne-PEO ₄ -CR6G	,, ,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	FP-DQP800, 1mg / 5mg	
Fluor525-Acetylene (5/6-Carboxyrhodamii Alkyne-PFO ₄ -TAMRA	ne 6G); Abs/Em = 522/544 nm	FP-DQP810, 1mg / 5mg	
Fluor545-Acetylene (TMRA-PEO4-Alkyne); MW:644.73; Abs/Em = 546/565 nr	n	
Alkyne-PEO4-SRB		FP-DQP820, 1mg / 5mg	
Alkvne-PEO4-SR101		FP-DQP30. 1mg / 5mg	
Fluor585 –Acetylene; Abs/Em = 584/603 r	nm		
* (with a sulfo-prop	ovl substituent):		
Alkyne-TriSulfo-CY _{anine} 3	,	FP-1C8830, 1mg / 5mg / 25mg / 10)0mg
Tri-SulfoCy3 –Acetylene; MW:761.92; Abs	s/Em = 550/570nm	ED 109940 1mg / Emg / 25mg / 10	00mg
Tri-SulfoCy5 –Acetylene; MW:787.96; Abs	s/Em = 647/663nm	FF-108840, 111g/ 511g/ 2511g/ 10	long
Alkyne-TriSulfo-CY _{anine} 7		FP-11220, 1mg / 5mg / 25mg / 100)mg
Tri-Sulfo-CY7-Acetylene; MW:1010.22; At	os/Em = 753/775nm		
Alkyne-DiSulfo-CY _{anine} 3	intuent).	FP-I QV030, 1mg	
DiSulfo-Eth-CY7-Acetylene, CF ₃ CO ²⁻ salt;	MW:781.86; Abs/Em.=555/565nm;	Soluble in DMSO	
Alkyne-DiSulfo-CYanine5	MM/807 00: Abc/Em =640/66pm; S	FP-LQV090, 1mg	
Alkvne-TetraSulfo-CYanine5.5	WW.007.90, ADS/EIII049/001111, 3	FP-LQV320. 1mg	
TetraSulfo-Eth-CY5.5-Acetylen, CF3CO2-	salt; MW:1068.14; Abs/Em.=578/70	1nm; Soluble in DMSO.	
Alkyne-DiSulfo-CYanine7 DiSulfo-Eth-CY7-Acetylene_CEaCO ²⁻ salt:	MW/833 93. Abs/Em =749/776m. S	FP-LQV250, 1mg	
* (with a methyl sul	bstituent):		
Alkyne-SulfoCYanine3	,	FP-1C4620, 1mg / 5mg / 25mg / 10	l0mg
MonoCy3 – Acetylene; MW:573.75; Abs/El	m = 550/567nm; EC:96800; QY:0.1	5 ED 144630, 1mg / 5mg / 25mg / 10	10mg
DiSulfoCy3 – Acetylene Na salt; MW:675.7	79; Abs/Em = 548/567nm; EC:16200	10; QY:0.15	long
Alkyne-DiSulfo-CY _{anine} 5		FP-005590, 1mg / 5mg / 25mg / 10	00mg
VISulfoCy5 – Acetylene Na salt; MW: 701.8	33; Abs/em.=646/664nm; EC:2/1000 hstituent):	J; QY:0.28	
Alkvne-DiSulfo-CY _{anine} 3	Ustituent).	FP-0B8390. 1mg / 5mg / 25mg / 10	10ma
DiSulfo-Cy3 –Acetylene K salt; MW:691.9	(+653.2); Abs/Em = 548/563nm; EC	:162000; QY:0.1; CF260:0.03; CF280:0.06 - S	olubility is good in Water, DMF, DMSO
Alkyne-DiSulfo-CY _{anine} 3.5		Inquire	
Alkyne-DiSulfo-CY _{anine} 5		FP-SJI060, 1mg / 5mg / 25mg / 100	Dmg
DiSulfo-Cy5 –Acetylene Ksalt; MW:717.94	; Abs/em.=649/662nm; EC:271000;	QY:0.28; CF260:0.04; CF280:0.04 – Solubility	r is very good in Water, DMF, DMSO
AIKyne-IetraSulto-CY anine 3.3 TetraSulfo-Cv5.5 – Acetvlene tri-K salt: MV	V:1054.36: Abs/Em = 673/691nm: E	C:195000: QY:-: CF260:0.09: CF280:0.11 – So	umg plubility is good in Water. DMF. DMSO
Alkyne-DiSulfo-CY _{anine} 7		FP-0B8280, 1mg / 5mg / 25mg / 10	10mg
DiSulfo-Cy7 – Acetylene K salt; MW:745.3	; Abs/Em = 750/773nm; EC:240600	QY:-; CF260:0.04; CF280:0.04 – Solubility is	good in DMF, DMSO, significant in water
TetraSulfo-CY7.5-Acetylene tri-K salt; MW	/:1120.46(+1005.2); Abs/Em = 778/	797nm; EC:222000; QY:-; CF260:0.09; CF280:	0.09 – Solubility is good in Water, DMF, DMSO
* (no sulfo):			
Alkyne-CYanine3	FF/F70 FO:4F0000, OV:0.04, OF	FP-1A6320, 1mg / 5mg / 25mg / 10	10mg
Acetylene; MW:530.14; Abs/Em = 5 Alkvne-CYanine 3.5	55/570nm; EC:150000; QY:0.31; CF	280:0.04; CF280:0.09	
Cy3.5 –Acetylene			
Alkyne-CYanine5	6/662pm; EC:250000; OV:0 2; CE2	FP-OO5590, 1mg / 5mg / 25mg / 10	00mg
Alkyne-CY _{anine} 5.5	0/0021111, E0.200000, Q1.0.2, Of 20	FP-SJH910, 1mg / 5mg / 25mg / 10	00mg
Cy5.5 –Acetylene; MW:656.30; Abs/Em =	684/710nm; EC:209000; QY:0.2; C	F260:0.02; CF280:0.03	00
AIKyITe-C Tanine / Cy7 –Acetylene; MW:622.38: Abs/Em = 7:	50/773nm; EC:199000: QY:0.23: CF	rr-wzeizu, img / omg / 2omg / 1 260:0.022; CF280:0.029	uung
Alkyne-CY _{anine} 7.5	, , ,	FP-WZE141, 1mg / 5mg / 25mg / 1	00mg
CY7.5-Acetylene; MW:722.40(+686.4); Ab	os/Em = 788n808nm; EC:223000		
See technical sheets:	FT-1A6320 (FluoProb	es CYanine-Alkvnes)	FT-DOP790 (FluoProbes CYanine-DBCO)
	FT-HO7250 (FluoProt	bes CYanine-Azide)	2200)
See also	PH- <mark>BB060k</mark> Clickable	Fluorescent dyes	



Biotin Alkynes:
Acetylene-PEO(-Biotir

Acetylene-PEO₄ -Biotin MW:457.58; (M) DQP650, 25mg / 100mg / 1g

• PEO spacer(PEGs) - Alkynes:

Acetylene-PEO₄ -Amine, - Acid, -NHS, -Maleimide: inquire +DQP61/3/4

• Nucleotides - Alkynes: +:

Propargyl-CEP-**oligonucleotides** are used for attachment of azide-containing reporter groups such as biotin or fluorescent dyes by click chemistry. This allows the synthesis of highly-modified DNA strands carrying multiple labels in a density that is not achieved by classic labeling techniques.

propargyl is an alkyl functional group of 2-propynyl (structure: HC=C-CH₂-).

CEPs (Cyanoethyl-N,N-diisopropyl Phosphoramidites) functional group is used to synthesize oligonucleotides by automated standard protocols.

5-Propargyloxy-dU CEP	More description here.
5-Octadiynyl-dU CEP	More description here
Alkynyl-modifier-C6-dT CEP	More description here .
5-(Propargyloxy)-2'-deoxyuridine	More description here
5-(1,7-Octadiyn-1-yl)-2'-deoxyuridine	More description here.
5'-O-(Dimethoxytrityl)-5-(propargyloxy)-2'-deoxyuridine	More description here
5-Octadiynyl-TMS-dU CEP	More description here.
5-OctadiynyI-TMS-dC CEP	More description here
5-Octadiynyl-dC CEP	More description here
5-Octadiynyl-TIPS-dU CEP	More description here
More reagents: inquire +BA0174; +DOI62,DOI57,Ml	M982,ZC68/9 DOP20/1

More reagents: inquire +BA0174; +DQI62,DQI57,MM982,ZC68/9DQP20/1/ZC667/8|X028|DQ071/2/3/5 | DQP21 2'-O- and 3'-O-Propargyl CEPs and other alkyne-containing CEPs (C8- TIPS- and TMS- versions)



• Azides Click Chemistry reagents

• Fluorescent dyes - Azides:

Standard Dyes such as Cy3, Tamra or Te xas Re d have been thoroughly selected to cover the whole UV-Vis spectrum. Novel Alternative Dyes are available as well.

• FluoProbes - Azides are superior fluorescent dyes (see caracteristics, protocol) +

FluoProbe488 – Azide (494/519nm)	B38820	see All FluoProbes – Azides
FluoProbe565A – Azide (563/592nm)	YE4990	
FluoProbe647N – Azide (653/674nm)	YE5020	
FluoProbe594 –Azide (591/617nm)	FJ0800	
FluoProbe782 –Azide (738/800nm)	JV1940	
FluoProbe831 –Azide (844/nm)	JV7790	

• **CYanine – Azides :** See technical sheet:

FT-HO7250 (FluoProbes CYanine-Azide)

• Classic fluorescent dyes+: See PH-<u>BB060k</u>: incl. (Picolyl-) Azide – Fluorescent dyes. See <u>FT-DQP790</u>: FluoProbes CR110, CR6G, TAMRA, SRB, SR101, Cy3/5

6-Carboxyfluorescein-PEO azide 6-Carboxyfluorescein-dipivalate PEO azide BBQ-650[®]-PEO azide 6-TET-PEO-azide DANSYL-TEG-N₃

More description here . More description here . More description here . More description here . FJ5740

• Biotin Azides Biotin-azide	ZC6710.	
⊬ Desthiobiotin-PEO azide Biotin-PEO azide	More description here More description here	

• PEO spacer (PEGs) - Azides:

PEO₂ to PEO₈-Azides linkers and crosslinkers: inquire +*DQP22/3*, *ZC684 to ZC689* +*DQP66* rem: PEO is also called PEG and TEG

For example: Aminooxy-PEO-azide Amino-PEO azide

More description here . More description here .

• Nucleotides - Azides: inquire +

Azide-CEP-**oligonucleotides** are used for attachment of azide-containing reporter groups such as biotin or fluorescent dyes by click chemistry. This allows the synthesis of highly-modified DNA strands carrying multiple labels in a density that is not achieved by classic labeling techniques.

Azide is a functional group for click reactino with alhynes CEPs (Cyanoethyl-N,N-diisopropyl Phosphoramidites) functional group is used to synthesize oligonucleotides by automated standard protocols. Inquire for a products list

Other Azides:
Folate-PEO3 azide
Tocopherol-PEO azide
Water soluble dansyl-PEO azide
PQQ-PEO azide
More description here
Cholesteryl-PEO azide
More description here



Activators for Click Chemistry

• NHS activated Alkyne and Azide - to activate amines(proteins) Alkyne Amine-Activator #ZL5532, Azide Amine-Activator #ZL5542

Minterchim

Auxillary reagents for Click Chemistry: catalyzers; buffers:

• Copper sources, reduction and stabilization

The efficiency of a copper (Cu(I))-catalyzed Azide-Alkyne click chemistry reaction (CuAAC) strongly depends on the presence of copper ions in the +1 oxidation state (Cu(I)).

Different copper catalyst sources, reduction reagents and Cu(I) stabilizing ligands are available however, for most bioconjugation applications the combination of CuSO4 as copper catalyst source, sodium ascorbate as a reduction reagent and a water-soluble Cu(I) stabilizing ligand such as THPTA[1,2] or BTTAA[3,4] is recommended (Tab. 1).

An optimal balance between reaction speed and Cu(I) concentration can be achieved using THPTA or BTTAA in combination with Picolyl-Azide detection reagents (Picolyl-Azides of Biotin or Picolyl-Azides of fluorescent dyes) that contain an additional internal copper chelating moiety[4].

Presolski et. al.[1] (Download pdf) and Hong et. al.[2] provide a general protocol for CuAAC reactions that may be used as a starting point for the set up and optimization of individual assays. Table 1: Overview of available reaction reagents for Cu(I) catalyzed Azide-Alkyne click chemistry (CuAAC).

Available reaction reagents for Cu(I) catalyzed Azide-Alkyne click chemistry (CuAAC):

Copper (Cu) catalyst	Reduction reagent	Cu(I) stabilizing Ligand		
CuSO4	Sodium Ascorbate	ТНРТА		
		BTTAA		
CuBr		TBTA		

Protein Reaction Buffer Kit

Kit Contents: 2.5 ml Reaction Buffer (250 mM), 0.6 mg Additive 1, 11 mg Additive 2, 0.25 ml Copper (II) Sulfate (20 mM), Reducing Agent

Click Solvent – Click chemistry grade	ZC6950, 1 ml	ZC6951, 10x 1 ml	
CuSO4 - click chemistry grade – Click chemistry grade	1H3690, 10 mg	1H3691, 5x 10 mg	
Cupper(II) Guphate) CuBr - click chemistry grade	ZC6930, 5mg	ZC6931, 10x 5 mg	
Na-Ascorbate – Click chemistry grade	10487F, 10 mg	10487G, 10mg	
THPTA (Tris(3-hydroxypropyltriazolylmethyl)amine)	MRU020, 100mg	MRU021, 500mg	MRU022, 1g
(Tris[(1-benzyl-1H-1,2,3-triazol-4-yl)methyl]amine)	ZC6940, 5mg	ZC6941, 10mg	
BIIAA			

MRU690

2-(4-((bis((1-(tert-butyl)-1H-1,2,3-triazol-4-yl)methyl)amino)methyl)-1H-1,2,3-triazol-1-yl)acetic acid

Copper(II)-TBTA complex #FY2780

See <u>FT-*FY*2780</u>

AzidoAniline Hydro. ZC6900

Selected references:

[1] Presolski et al. (2011) Copper-Catalyzed Azide-Alkyne Click Chemistry for Bioconjugation. Current Protocols in Chemical Biology 3:153.

[2] Hong et al. (2011) Analysis and Optimization of Copper-Catalyzed Azide-Alkyne Cycloaddition for Bioconjugation. Angew. Chem. Int. Ed. 48:9879.
[3] Besanceney-Webler et al. (2011) Increasing the Efficiacy of Bioorthogonal Click Reactions for Bioconjugation: A Comparative Study. Angew. Chem. Int.

Ed. 50:8051

[4] Uttamapinant et al. (2012) Fast, Cell-Compatible Click Chemistry with Copper-Chelating Azides for Biomolecular Labeling. Angew. Chem. Int. Ed. 51:5852.

•Fluorimetric assessment of Click Chemistry reaction ⁰.

The efficiency of Cu(I)-catalyzed Azide-terminal Alkyne Click Chemistry reactions (CuAAC) strongly depends on the presence of the copper catalyst in the +1 oxidation state (Cu(I)). While different copper sources, reduction reagents and Cu(I) stabilizing ligands are available, the Cu(II) salt $CuSO_4$ – in combination with sodium ascorbate (reduction reagent) and THPTA (Cu(I) stabilizing ligand) – has become the catalyst of choice for most biomolecule labeling applications^[1,2].

Biomolecule labeling however, often **requires optimization of reaction conditions** (e.g. type of copper source, ligand or copper source/ligand ratio...). **Assay performance (i.e. rate of product formation) can be conveniently monitored in real time using the pro-fluorogenic dye 3-Azido-7-hydroxycoumarin.** The non-fluorescent Azide form of **3-Azido-7-hydroxycoumarin strongly increases its fluorescence upon Cu(I)-catalyzed conjugation to a terminal Alkyne** that is triggered by formation of a triazole moiety^[3] (Fig. 1). Removal of unreacted dye before measurements is therefore not required.

When setting up a new CuAAC labeling assay, the general protocol provided by Presolski *et. al.*^[1] and Hong *et. al.*^[2] may be used as a starting point for subsequent optimization.



Figure 1: The pro-fluorogenic dye **3-Azido-7-hydroxycoumarin** becomes highly fluorescent upon Cu(I)-catalyzed conjugation to a terminal Alkyne moiety. Progress of product formation can be monitored by measuring the fluorescence increase at 477 nm.

3-Azido-7-hydroxycoumarin

IOK790, 1mg IOK791, 5mg

3-Azido-7-hydroxy-chromen-2-one, MW: 203.15 g/mol Abs/Em = 404/477 nm, Soluble in DMSO, DMF, MeOH, MeCN (rem: solution in DMSO and DMF are not suitable for longterm storage); (M)

Selected References:

Presolski *et al.* (2011) Copper-Catalyzed Azide-Alkyne Click Chemistry for Bioconjugation. *Current Protocols in Chemical Biology* 3:153.
Stanislav *et al.* (2009) Analysis and Optimization of Copper-Catalyzed Azide-Alkyne Cycloaddition for Bioconjugation. *Angew. Chem. Int. Ed.* 48:9879.

[3] Sivakumar et al. (2004) A fluorogenic 1,3-dipolar cycloaddition reaction of 3-azidocoumarins and acetylenes. Org Lett. 6 (24):4603.

CLICK labeled molecules Capture Kit

Click Chemistry Capture Kit

WXS271, 1kit

The Click Chemistry Capture Kit provides all necessary reagents to covalently capture specific sub-classes of proteins by a Cu(I)-catalyzed azide-alkyne cycloaddition reaction (CuAAC). The proteins of interest need to be metabolically, enzymatically or chemically azido- or alkyne-tagged. Subsequently, the resin containing the covalently attached proteins can be washed with high stringency, virtually eliminating any non-specifically bound proteins. Upon protease digestion, this yields a highly pure peptide pool that is ideal for mass spectroscopy (e.g. LC MS/MS) based analysis.

Kit Contents:

- 7 ml Lysis buffer store at 4 °C
- 4.8 g Urea store at ambient temperature
- 1.5 ml Additive 1 store at 4 °C

0.5 ml Copper (II) Sulfate (100 mM) - store at ambient temperature

- 400 mg Additive 2 store at ambient temperature
- 7.7 g Agarose wash buffer - store at ambient temperature
- 10 Empty spin columns store at ambient temperature



Detailled technical information

• Nucleotides – Alkynes⁰.

Alkyne-bearing nucleoside phosphoramidites enable click chemistry. An attractive strategy for nucleic acid conjugation involves the click reaction of alkyne-bearing oligonucleotides with azide-bearing species to join them via a triazole linkage^[r]. Interchim Biosciences offers options: Ethynyl-dU CEP,8 5-Octadiynyl-dU CEP,7-10 Alkynyl-Modifier-C6-dT CEP, and other new compounds.

Use of **5-Propargyloxy-dU CEP**: Employ acetonitrile diluent at the concentration recommended by the synthesizer manufacturer. Use standard coupling protocols; extended coupling times are not required. Cleavage from the solid support may be carried out by standard procedures. Nucleobase deprotection should be done at 55 °C for 8-12 hours. Extended heating is not recommended due to by-product formation.

CEPs Nucleic acids	Product #	
5-Octadiynyl-dU CEP IV9770 CAS938186-76-6; MW: 879.98 – Store -20°C 3'-O-[(Diisopropylamino)(2-cyanoethoxy)phosphino]-5'- O-(4.4'-dimethoxytrityl)-5-[N-(hex-5-ynyl)-3-(E)-acrylamido]-2'- deoxyuridine A much superior reagent for the installation of an alkyne-bearing nucleoside into an oligonucleotide when click reactions are contemplated. The terminal alkyne is more flexible and extends farther away from the pyrimidine ring, thereby allowing more efficient ligation via click chemistry, even when multiple ligations to a single oligo are desired.	.BA0308	
5-Propargyloxy-dU CEP IV9780 CAS:; MW: 784.83 – Store -20°C, dry A modified nucleotides for introducing an alkyne into an oligonucleotide internally or at the 5'-terminus.	.BA0174	
Alkynyl-modifier-C6-dT CEP CAS:; MW: 879.98 – Store -20°C, dry 3'-O-[(Diisopropylamino)(2-cyanoethoxy)phosphino]-5'-O-(4,4'- dimethoxytrityl)-5-[N-(hex-5-ynyl)-3-(E)-acrylamido]-2'-deoxyuridine This phosphoramidite features a well- established linker that places the alkyne at a sufficient distance from the oligonucleotide to allow efficient click conjugation.	.BA0361	
5-(Propargyloxy)-2'-deoxyuridine CAS[65367-85-3]; MW: 282.25 – Store at -20°C	.PY7712	
5-(1,7-Octadiyn-1-yl)-2'-deoxyuridine CAS NO: 909398-18-1; MW: 332.35 – Store at -20°C	.PY7713	



•Azides for click chemistry ⁰.

Azide reagents		
Desthiobiotin-PEO azide Desthiobiotin-TEG azide CAS[1306615-47-3]; MW: 414.50 – Store at +4°C for ligation of desthiobiotin via Saudinger ligation or click chemistry.	.BT1075	
Biotin-PEO azide Biotin-TEG azide, N-(11-Azido-3,6,9-trioxa-undecylamine)-D-(+)- biotinamide CAS[875770-34-6]; MW: 444.55 – Store at +4°C for ligation of biotin via Staudinger ligation or click chemistry.	.BT1085	
See also : Biotin-Hydrazide <u>FJ6741</u> , Biotin- PEG-Azides <u>78631A</u>		



Amino-PEO ₃ -Azide	.LK4310	$\sim 10^{-100}$ NH ₂
Amino-TEG-Aide ; 1-Amino-11-azido-3,6,9-trioxaundecane		
A bi functional linkor with potential		
applications in oligonucleotide ligations.		
Aminooxy-PEO ₃ -Azide	FZ8700	$H_{\nu}N^{\prime}O^{\prime}O^{\prime}O^{\prime}O^{\prime}O^{\prime}N_{3}$
AMINOOXY-TEG-AZIDE CASL1: MW: 234 25 - Store at +4°C		Conjugation by azide can be used to introduce a PEO-
A bi-functional linker Aminooxy-PEO-azide		aminooyu functionality. Aculation of the experime and
A bi-functional miker, Ammooxy-i EO-azide.		animoxy functionality. Acylation of the oxyanime end
		yields a hydroxamic acid that mimics carboxylic acid.
Other AminoOxy-PEOs: <u>F1-JV2290</u>		Alternatively, the oxyamine end can condense with an
		aldehyde, affording an oxime that bears a PEO-linked
		azide group.
Folate- PEO ₃ -Azide	070600	0 CO2H
Folate-TEG-Azide		
CAS[1313026-32-2]; MW: 641.64 – Store at +4°C		
For the introduction by Click or Staudinger		
ligations of the folate tag for recognition by		H ₂ N ^ N ^ N
overexpressed folate receptors in cancer cells		N3 🗸
Tocopherol- PEO ₃ -Azide	.FC8160	
Tocopherol-TEG-Azide		
CAS[-]; MW: 631.93 – Store at -20°C		
Lipophilic carrier tag for Click and		j l
Staudinger ligations collection.		
		O Na
Water caluble demand DEO Arida	ED4200	N(CH)
dansyl-TEG-Azide : 4-Dansyl-14-aza-7 10 13-trioxa-15-	.FD1300	10(01)2
azido)pentadecylsulfonate		
C23H34N5NaO8S2; CAS[-]; MW: 595.67 - Store at +4°C, in the dark		
		ΥΎ
		0=\$=0
6-FAM-PFO2-Azide	F.10011	H0. A .0. A .0H
6-Carboxyfluorescein-TEG azide	100011	
CAS[412319-45-0]; MW: 576.55 – Store at +4°C, in the dark		
fluorescein azide for fluorescein		
incorporation via Staudinger ligation or click		
chemistry. See also FP488-Azide FP-YE4970		
6-FAM-dinivalate- PFO2-Azide	FF6120	
6-Carboxyfluorescein-dipivalate TEG azide		$\times 0.$
MW: 744.79 – Store at -20°C, dark		
		$\gamma \geq \rho$
		\sim° \sim°
BBQ-650®-PEO ₃ -Azide	.BL3030	0
BlackBerry® Quencher 650 - TEG - azide		
A useful long weyelength guencher for Click		H III
A useful long wavelength quencher for Click		
Chemistry and Staudinger ligation		H ₃ CO
Distance Operation to the state of a second sector day UC Determined		
#7.879.986. Berry & Associates. Inc.		
"BlackBerry is a trademark of Berry & Associates, Inc. Products		
derived from BlackBerry [®] Quencher reagents are sold		0 ₂ N
exclusively for research and development use by the purchaser. They may not be used for clinical or diagnostic purposes without		
prior agreement and consent of Berry & Associates, Inc."		
6-TET-PEO ₃ -Azide	.FF6130	HO
6-Carboxy-2',4,7,7'-tetrachlorofluorescein-TEG azide		
CAS[-]; MW: 714.33 – Store at -20°C, dark		
for tetrachloro-fluorescein incorporation via		
Staudinger ligation or click chemistry		
		NH A
		GI C S Na
		$\sim \sim$





Related products lines

Interbiotec - BioSciences innovation - proposes a complete range of products for protein biochemistry.

Innovative and remarkable chemistries, conjugation methods, labeling and functionalisation (PH)

- Standard Click Chemistry reagents (PH)
- Copper-free Click Chemistry reagents (DBCO reagents) (PH)
- Staudinger reaction (effective conjugations/chemical modification) (PH)
- **PEGylation reagents** (conjugation reagents, linkers and building blocks) (PH)
- SAM reagents (Self-Assembled Monolayers for surface modification) (PH)
- **SDA reagents** (effective photo reactiions) ^(PH)
- STELLA labeling (azocycloaddition reactions) (PH)
- Gold nano-particules and materials
- Carbone nanotubes
- ITO slides ^(PH)

•

• <u>FluoProbes labeling agents</u>

Desalting tools - CelluSep tubings, SpectraPor tubings, GebaFlex, FloatALyser, SlideALyser,...

Products HighLights Overview

Information inquire

Reply by Fax : +33 (0) 4 70 03 82 60 or email at <u>interbiotech@interchim.com</u>

☐ I wish to receive the complete documentation about: _

Name:		2 nd name:		Position:
Company/Institute:			Service, Lab:	
Adress:				
	Zip code:	Town:		