


Affinity Chromatography of sugar derivatives

Catch and Release of Diols with Boronic Acid Agarose

Boronic acids bind diol compounds with high selectivity.^[1] Since many biomolecules are **diol-decorated sugar derivatives**, **immobilized boronic acids** permit the single step enrichment of

- Carbohydrates^[2]
- Nucleotides^[3]
- Cofactors^[4]
- Entire Glycoproteins^[5] and RNAs^[6,7]

To this end, **agarose resins** stand out for their large exclusion limit (1×10^4 - 4×10^5 Da) and are therefore particularly well suited for the complexation of **large biomolecules**.^[8,9] Capitalizing on the **reversibility** of boronate ester formation, captured compounds can be **easily retrieved under mild conditions** for further applications or analysis.



Scheme 1: Depending on the *pH* of the solution, **boronic acids coexist in equilibrium with their corresponding boronates** ($pK_a \sim 9$). While boronic esters are prone to fast hydrolysis, their tetragonal boronate counterparts are stable and allow for the efficient immobilization of **cis-diol** functionalities. Trapped diols are then either released by **a) lowering the *pH*** or by **b) replacement with sorbitol**.

Ordering information:

30 µmol m-Aminophenylboronic acid/ml gel;

Beads/Particle size: 45 - 165 µm (for 1 x 10⁴ - 4 x 10⁵); Recommended linear flow rate 30 cm/h; Maximum pressure: 0.25 bar (3.6 psi)

Stable at pH4-9 (short term) – 7.5 for long term); and to all solutions commonly used in gel filtration including 8 M urea and 6 M guanidine hydrochloride Not stable in organic solvents! ; Not autoclavable!; Store at +4°C (H)

Selected references:

- [1] Siegel (2012) Applications of reversible covalent chemistry in analytical sample preparation. *Analyst*. 137:5457.
- [2] Weith et al. (1970) Synthesis of cellulose derivatives containing the dihydroxyboryl group and a study of their capacity to form specific complexes with sugars and nucleic acid components. *Biochemistry*. 9:4396.
- [3] Moore et al. (1974) Separation of ribonucleotides and deoxyribonucleotides on columns of borate covalently linked to cellulose. Application to the assay of ribonucleoside diphosphate reductase. *Biochemistry*. 13:2904.
- [4] Gellekink et al. (2005) Stable-isotope dilution liquid chromatography-electrospray injection tandem mass spectrometry method for fast, selective measurement of S-adenosylmethionine and S-adenosylhomocysteine in plasma. *Clin. Chem*. 51:1487.
- [5] Mallia et al. (1981) Preparation and use of a boronic acid affinity support for the separation and quantitation of glycosylated hemoglobins. *Anal. Lett.* 14:649.
- [6] Schott et al. (1973) Dihydroxyboryl-substituted methacrylic polymer for the column chromatographic separation of mononucleotides, oligonucleotides, and transfer ribonucleic acid. *Biochemistry*. 12:932
- [7] Nübel et al. (2017) Boronate affinity electrophoresis for the purification and analysis of cofactor-modified RNAs. *Methods*. 117:14.
- [8] Chen et al. (2017) Novel boronate material affords efficient enrichment of glycopeptides by synergized hydrophilic and affinity interactions. *Anal. Bioanal. Chem*. 409:519.
- [9] Grundy et al. (2016) PARP3 is a sensor of nicked nucleosomes and monoribosylates histone H2B(Glu2). *Nat. Commun*. 7:12404.

Applications using boronate affinity media.

Application	Application Buffer	Molecules in Vo	Retained Molecules	Elution Buffer	Support Used	References
Adenylate cyclase assay	HEPES pH 7.5 + MgCl ₂	cAMP	ATP, AMP, and adenosine	0.05 M NaOAc	Bio-Gel P-150 boronate gel	Hageman and Kuehn 1977
Isolation of catecholamines from urine	0.1 M phosphate pH 7.0 + EDTA	Other urine components and DOPA	Norepinephrine, epinephrine	0.025 N HCl	Boric acid gel ¹	Higa et al. 1977
Modified nucleosides in urine	0.25 M NH ₄ OAc, pH 8.8	Thymidine, adenine	Pseudouridine	0.1 M HCOOH	Bio-Gel P-2 gel, 200–400 mesh, boronate	Davis et al. 1977; Uziel et al. 1976
Separation of mono- and oligonucleotides	Triethyl ammonium	Deoxyribo-nucleotides	Ribonucleotides	H ₂ O and others	Dihydroxyboryl methacrylate	Schott et al. 1973
Sugars	0.05 M N-methyl-morpholinium-Cl, pH 7.5 + 1 M NaOAc	Erythritol, adonitol, sucrose, and D-glucose	L-arabinitol, xylitol, D-mannitol, dulcitol, sorbitol, and D-fructose*	Isocratic run, elution buffer same as application buffer	Dihydroxyboryl cellulose	Weith et al., 1970
Assay of ribonucleotide reductase	Tris, MgCl ₂	dUMP, dCMP	CDP	Citrate	Dihydroxyboryl cellulose Dowex 1(AG) resin	Moore et al. 1974

Regeneration and Storage

Regeneration of the Boronic gel is facilitated by washing with starting buffer if a low-pH elution buffer has been used. If boric acid buffer or diol solutions have been used for elution, the gel must first be washed with 0.1 M acetic acid before reequilibrating with starting buffer

Typical Procedure

Stepwise elution is typical for separation of a ribonucleotide (such as GMP) from a cyclic nucleotide (cGMP):

1. Prepare a 1 x 4 cm column of Affi-Gel boronate in 0.1 M HEPES, pH 8.5.
2. Load a 200 µl solution containing 1 µmol each cGMP and GMP in the starting buffer onto the column.
3. Wash the column with 5 column volumes of starting buffer at a flow rate of 1 ml/min and collect the cyclic nucleotide.
4. Elute the ribonucleotide with 0.1 M NaPO₄, pH 6.0.

Related products/documents

• Other affinity resins for the purification of nucleotide-binding, phosphorylated or click-tagged proteins.

Boronic acid on magnetic beads, 1µm, 2.5% (w/v) #7A2610, 2ml

#

• Other Boronic-based conjugates [\[PH-BB032n\]](#), e.g. labeled probes (fluorescent, biotinylated, tagged).

Search at www.interchim.com or ask at interbiotech@interchim.com

See [Products HighLights Overview](#)

Information inquire

Reply by Fax : +33 (0) 4 70 03 82 60 or email at interbiotech@interchim.com

☐ I would like to receive further information on: _____

Title : _____ First name: _____ Surname: _____ Position: _____

Company/Institute: _____ Service, Lab: _____

Adress: _____

Postcode: _____ Town: _____

Tel _____ Fax _____ Email: _____