

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### Label-nucleotide spacers

The number "11" is the number of atoms in the linker between biotin and dUTP. The linker length affects the incorporation efficiency of the biotin-dUTP probe into DNA using DNA polymerases, and it also affects biotin/avidin or biotin/streptavidin interaction. In general, the shorter the linker, the more efficiently the biotin-dUTP is incorporated into DNA by DNA polymerases. On the other hand, the longer the linker, the better biotin can interact with avidin or streptavidin.

## Biotinylated nucleotides

Biotin is a versatile label, allowing color choice and high sensitivity through the use of various signal amplification system. Biotinylated nucleotides can be incorporated into nucleic materials by Polymerase Chain Reaction (PCR). dUTP serves for DNA amplifications. UTP serves for RNA amplifications. The resulting biotin-containing DNA or RNA can be subsequently labeled with a fluorescent labeled (strept)avidins or detected with anti-biotin antibodies.

Biotin-dUTP can be enzymatically incorporated into DNA via nick translation, random priming, or 3'-end terminal labeling.

### Biotin-11-dUTP

(Biotin-11-2'-deoxyuridine-5'-triphosphate, tetralithium salt)

1 mM in pH 7.5 Tris-HCl buffer

MW : 886.4

| Description                       | Cat.#     | Qty   |
|-----------------------------------|-----------|-------|
| Biotin-11-dUTP                    | FP-AM539A | 50 µl |
| Biotin-11-dUTP Lyophilized powder | FP-AY489A | 50 µg |

### Biotin-16-dUTP

(Biotin-16-2'-deoxyuridine-5'-triphosphate, tetralithium salt)

1 mM in pH 7.5 Tris-HCl buffer

MW : 971.5

| Description                        | Cat.#     | Qty   |
|------------------------------------|-----------|-------|
| Biotin-16-dUTP                     | FP-AM541A | 50 µl |
| Biotin-16-dUTP, Lyophilized powder | FP-AM542A | 50 µg |

### Biotin-20-dUTP

(Biotin-20-2'-deoxyuridine-5'-triphosphate, tetralithium salt)

MW : 1020.54

| Description    | Cat.#     | Qty   |
|----------------|-----------|-------|
| Biotin-20-dUTP | FP-AM543A | 50 µg |

### Biotin-11-UTP

(Biotin-11-uridine-5'-triphosphate, tetralithium salt)

1 mM in pH 7.5 Tris-HCl buffer

MW : 902.5

| Description   | Cat.#     | Qty   |
|---------------|-----------|-------|
| Biotin-11-UTP | FP-AM550A | 50 µl |

### Biotin-16-UTP

(Biotin-16-uridine-5'-triphosphate, tetralithium salt)

1 mM in pH 7.5 Tris-HCl buffer

MW : 987.52

| Description   | Cat.#     | Qty   |
|---------------|-----------|-------|
| Biotin-16-UTP | FP-AM551A | 50 µl |

### Biotin-EDA-ATP

2'/3'-O-(2-(Biotinyl-butyl-carbamoyl)-ethyl-carbamoyl)-Adenosine-5'-triphosphate, Na salt

MW : 815.53 (Anion)

| Description    | Cat.#  | Qty  |
|----------------|--------|------|
| Biotin-EDA-ATP | FQ5630 | 20 u |
|                | FQ5331 | 10 u |

### Biotin-EDA-AppNHp (Biotin-EDA-AMPPNP)

2'/3'-O-(2-(Biotinyl-butyl-carbamoyl)-ethyl-carbamoyl)-Adenosine-5'-[(β,γ)-imido]triphosphate, Na salt

MW : 814.55 (Anion)

| Description                           | Cat.#  | Qty  |
|---------------------------------------|--------|------|
| Biotin-EDA-AppNHp (Biotin-EDA-AMPPNP) | FQ2831 | 10 u |
|                                       | FQ2830 | 50 u |

Related Products :  
Aminated nucleotides page D.137

### Amine labeled nucleotides

Amine labeled nucleotides are used to link nucleotides to amine-reactive probes or to stationary phases for affinity chromatography.

#### N<sup>6</sup>-(4-Amino)butyl-ATP

N<sup>6</sup>-(4-Amino)butyl-adenosine-5'-triphosphate, Sodium salt  
MW : 575.28 (Anion)

| Description                        | Cat.#  | Qty   |
|------------------------------------|--------|-------|
| N <sup>6</sup> -(4-Amino)butyl-ATP | FQ2960 | 50 u  |
|                                    | FQ2961 | 250 u |

#### N<sup>6</sup>-(6-Amino)hexyl-ATP

N<sup>6</sup>-(6-Amino)hexyl-adenosine-5'-triphosphate, Sodium salt  
MW : 603.33 (Anion)

| Description                        | Cat.#  | Qty   |
|------------------------------------|--------|-------|
| N <sup>6</sup> -(6-Amino)hexyl-ATP | FQ2970 | 50 u  |
|                                    | FQ2971 | 250 u |

Ito, *et al.*, "Modified nucleic acid for systematic evolution of RNA ligands by exponential enrichment.", *J. Bioact. Compat. Pol.*, 13 (2), 114 (1998)

Trayer, "Affinity chromatography of some adenosine phosphate-requiring systems", *Biochem. Soc. T.*, 2 (6), 1302 (1974)

#### 8-[(4-Amino)butyl]-amino-ATP

8-[(4-Amino)butyl]-amino-adenosine-5'-triphosphate, Sodium salt  
MW : 590.32 (Anion)

| Description                  | Cat.#  | Qty   |
|------------------------------|--------|-------|
| 8-[(4-Amino)butyl]-amino-ATP | FQ2630 | 50 u  |
|                              | FQ2631 | 250 u |

#### 8-[(6-Amino)hexyl]-amino-ATP

8-[(6-Amino)hexyl]-amino-adenosine-5'-triphosphate, Sodium salt  
MW : 618.35 (Anion)

| Description                  | Cat.#  | Qty   |
|------------------------------|--------|-------|
| 8-[(6-Amino)hexyl]-amino-ATP | FQ2640 | 50 u  |
|                              | FQ2641 | 250 u |

#### EDA-ADP

2'/3'-O-(2-Aminoethyl-carbamoyl)-Adenosine-5'-diphosphate, Sodium salt  
MW: 511.28 (Anion)

| Description | Cat.#  | Qty   |
|-------------|--------|-------|
| EDA-ADP     | FQ4800 | 30 u  |
|             | FQ4801 | 150 u |

Oiwa, *et al.*, "Comparative Single-Molecule and Ensemble Myosin Enzymology: Sulfoindocyanine ATP and ADP Derivatives.", *Biophys. J.* 78:3048 (2000)

Jameson, *et al.*, "Fluorescent analogs: Synthesis and Applications.", *Methods in Enzymology*, 278, 363 (1997)

#### EDA-ATP

2'/3'-O-(2-Aminoethyl-carbamoyl)-Adenosine-5'-triphosphate, Sodium salt  
MW : 590.25 (Anion)

| Description | Cat.#  | Qty   |
|-------------|--------|-------|
| EDA-ATP     | FQ4850 | 30 u  |
|             | FQ4851 | 150 u |

Oiwa, *et al.*, "Comparative Single-Molecule and Ensemble Myosin Enzymology: Sulfoindocyanine ATP and ADP Derivatives.", *Biophys. J.* 78:3048 (2000)

Jameson, *et al.*, "Fluorescent analogs: Synthesis and Applications.", *Methods in Enzymology*, 278, 363 (1997)

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### $\gamma$ -Aminophenyl-ATP

Adenosine-5'-[  $\gamma$ -(4-aminophenyl)]triphosphate, Sodium salt

MW : 595.27 (Anion)

| Description               | Cat.#         | Qty   |
|---------------------------|---------------|-------|
| $\gamma$ -Aminophenyl-ATP | <b>FQ5120</b> | 3 mg  |
|                           | <b>FQ5121</b> | 15 mg |

Haystead, *et al.*, "Gamma-phosphate-linked ATP-Sepharose for the affinity purification of protein-kinases -rapid purification to homogeneity of skeletal-muscle mitogen-activated protein-kinase kinase.", *Eur. J. Biochem.*, 214 (2), 459 (1993)

Trayer, *et al.*, "Preparation of adenosine nucleotide derivatives suitable for affinity chromatography.", *Biochem. J.* 139 (3):609 (1974)

### 5-Propargylamino-Dctp

5-Propargylamino-2'-deoxy-cytidine-5'-triphosphate, Sodium salt

MW : 517.19 (Anion)

| Description           | Cat.#         | Qty   |
|-----------------------|---------------|-------|
| 5-Propargylamino-Dctp | <b>FQ4580</b> | 30 u  |
|                       | <b>FQ4581</b> | 150 u |

Hobbs, "Palladium-catalyzed synthesis of alkynylamino nucleosides - a universal linker for nucleic-acids.", *J. Org.Chem.*, 54 (14), 3420 (1989)

Related products :

| Description                | Cat.#         | Qty  |
|----------------------------|---------------|------|
| EDA-AppNHp<br>(EDA-AMPPNP) | <b>FQ5690</b> | 20 u |

## Aminoallyl nucleotides

Aminoallyl-nucleotides can be enzymatically incorporated into nucleic materials, by Polymerase Chain Reaction (PCR). dUTP serves for DNA amplifications. The resulting amine-containing DNA or RNA can be subsequently labeled with a fluorescent dye, biotin or other haptens via conventional peptide coupling methods. This two-step method for labeling nucleic acids is considerably more economical than the one-step method using a prelabeled nucleotide.

### 5-Aminoallyl-DUTP

Sodium salt

(5-(3-aminoallyl)-2'-deoxyuridine-5'- triphosphate, trisodium salt, AA-dUTP)

MW : 589.17

| Description       | Cat.#            | Qty                              |
|-------------------|------------------|----------------------------------|
| 5-Aminoallyl-dUTP | <b>FP-AK218A</b> | 100 $\mu$ l (10 mM in TE buffer) |
|                   | <b>FP-AY490A</b> | 1 mg (Powder)                    |

### 5-Aminoallyl-UTP

Sodium salt

(5-(3-aminoallyl)uridine-5'-triphosphate, trisodium salt, AA-UTP)

MW : 605.17

| Description      | Cat.#            | Qty                              |
|------------------|------------------|----------------------------------|
| 5-Aminoallyl-UTP | <b>FP-AM555A</b> | 100 $\mu$ l (10 mM in TE buffer) |
|                  | <b>FP-AM558A</b> | 1 mg (Powder)                    |

**D.138** Proc Natl Acad Sci U S A 90, 4206(1993)  
J. Clin Microbiol 29, 583(1991)

### Mant-labeled oligonucleotides

Mant fluorophore (N-Methyl-anthraniloyl) is compact so, when attached on ribose ring it causes minimal perturbation of nucleotide-protein interactions. It is however sensitive to environment, notably to conformational changes of nucleotide-binding proteins. Mant-labeled nucleotides offer valuable tools to investigate the structure, protein-ligand interactions, and enzymatic activity of nucleotide-binding proteins (as myosin with ATP, P21 with GTP).

Literature for MAnt-ATP (also available for other Mant nucleotides).

-Cheng JQ, Jiang W, Hackney DD 1998 ; Interaction of mant-adenosine nucleotides and magnesium with kinesin ; *Biochemistry* 37, 5288-5295

-Churchich 1997 ; Conformational changes at the nucleotide binding of GroEL induced by binding of protein substrates – Luminescence studies ; *J.Biol. Chem.*


272(32) : 19645

-Malnasi-Csizmadia A, Woolley RJ, Bagshaw CR. 2000; Resolution of conformational states of dictyostelium myosin II motor domain using tryptophan (W501)

mutant s: implications for the open-closed transition identified by crystallography. ; *Biochemistry* 39, 16135-16146

-Shaffer J, Sun G, Adams JA. 2001 ; Nucleotide release and associated conformational changes regulate function in the COOH-terminal Src kinase, Csk; *Biochemistry*

40, 11149-11155



NPE-caged-ATP releases ATP generally on a millisecond time scale upon flash photolysis with near-UV light.

### N<sup>6</sup>-[4-(Mant-Amino)]butyl-ATP

N<sup>6</sup>-[4-((N-methyl-anthraniloyl)-amino)]butyl-adenosine-5'-triphosphate, Sodium salt

MW : 708.43 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 335 / 440 nm

$\lambda_{max}$  : 280/335 nm ;  $\epsilon$  : 18 000/2 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description                               | Cat.#     | Qty   |
|-------------------------------------------|-----------|-------|
| N <sup>6</sup> -[4-(Mant-Amino)]butyl-ATP | FP-FQ6611 | 20 u' |
|                                           | FP-FQ6610 | 100 u |

### N<sup>6</sup>-[6-(Mant-Amino)]hexyl-ATP

N<sup>6</sup>-[6-((N-methyl-anthraniloyl)-amino)]hexyl-adenosine-5'-triphosphate, Sodium salt

MW : 736.48 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 335 / 440 nm

$\lambda_{max}$  : 280/335 nm ;  $\epsilon$  : 18 000/2 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description                               | Cat.#     | Qty   |
|-------------------------------------------|-----------|-------|
| N <sup>6</sup> -[6-(Mant-Amino)]hexyl-ATP | FP-FQ6621 | 20 u  |
|                                           | FP-FQ6620 | 100 u |

### 8-[4-(Mant-Amino)]butyl-ATP

8-[(4-(N-methyl-anthraniloyl)-amino)butyl]-amino-adenosine-5'-triphosphate,

Sodium salt

MW : 723.44 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 335 / 440 nm

$\lambda_{max}$  : 280/335 nm ;  $\epsilon$  : 18 000/2 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description                 | Cat.#     | Qty   |
|-----------------------------|-----------|-------|
| 8-[4-(Mant-Amino)]butyl-ATP | FP-FQ6461 | 20 u  |
| (MABA-ATP)                  | FP-FQ6460 | 100 u |

Groemping, et al., "Regulation of ATPase and chaperone cycle of DnaK from *Thermus thermophilus* by the nucleotide exchange factor GrpE.", *J. Mol. Biol.*, 305 (5):1173(2001)

Weikl, et al., "C-terminal regions of Hsp90 are important for trapping the nucleotide during the ATPase cycle.", *J. Mol. Biol.* 303 (4):583(2000)

### 8-[6-(Mant-Amino)]hexyl-ATP

8-[(6-(N-methyl-anthraniloyl)-amino)hexyl]-amino-adenosine-5'-triphosphate,

Sodium salt

MW : 751.50 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 335 / 440 nm

$\lambda_{max}$  : 280/335 nm ;  $\epsilon$  : 18 000/2 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description                 | Cat.#     | Qty   |
|-----------------------------|-----------|-------|
| 8-[6-(Mant-Amino)]hexyl-ATP | FP-FQ6471 | 20 u  |
| (MAHA-ATP)                  | FP-FQ6470 | 100 u |

1 u = 1 µl at 10 mM

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### Mant-EDA-ATP

2'/3'-[(2-(N-methyl-anthraniloyl)-amino)ethyl-carbamoyl]-adenosine-5'-triphosphate, Sodium salt

MW : 723.40 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 255/355 nm ;  $\epsilon$  : 23 300/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description  | Cat.#     | Qty   |
|--------------|-----------|-------|
| Mant-EDA-ATP | FP-FQ5811 | 20 u  |
|              | FP-FQ5810 | 100 u |

### Mant-ADP

2'/3'-O-(N-Methyl-anthraniloyl)-adenosine-5'-diphosphate, Triethylammonium salt

MW : 558.32 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 255/355 nm ;  $\epsilon$  : 23 300/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description | Cat.#     | Qty   |
|-------------|-----------|-------|
| Mant-ADP    | FP-FQ2061 | 150 u |
|             | FP-FQ2060 | 750 u |

Bujalowski et al, "Kinetic mechanism of nucleotide cofactor binding to Escherichia coli replicative helicase DnaB protein. stopped-flow kinetic studies using fluorescent, ribose-, and base-modified nucleotide analogues.", Biochemistry 39:2106(2000)  
Cheng et al, "Interaction of mant-adenosine nucleotides and magnesium with kinesin.", Biochemistry 37:5288 (1998)

### Mant-ATP

2'/3'-(N-Methyl-anthraniloyl)-adenosine-5'-triphosphate, Triethylammonium salt

MW : 637.30 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 255/355 nm ;  $\epsilon$  : 23 300/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description | Cat.#     | Qty   |
|-------------|-----------|-------|
| Mant-ATP    | FP-FQ2071 | 150 u |

Tung-Chung Mou, et al, "Structural Basis for the Inhibition of Mammalian Membrane Adenylyl Cyclase by 2'(3')-O-(NMethylanthraniloyl) guanosine 5'-Triphosphate.", J. Biol. Chem. 280 (8):7253(2005)  
Booth, et al, "Analysis of the properties of the N-terminal nucleotide-binding domain of human P-glycoprotein.", Biochemistry 39:5518(2000)  
Thoenges, et al., "Tight binding of bulky fluorescent derivatives of adenosine to the low affinity E2ATP site leads to inhibition of Na<sup>+</sup>/K<sup>+</sup>-ATPase. Analysis of structural requirements of fluorescent ATP derivatives with a Koshland-Nemethy-Filmer model of two interacting ATP sites.", J. Biol. Chem. 274:1971(1999)

### Mant-dATP

3'-(N-Methyl-anthraniloyl)-2'-deoxy-adenosine-5'-triphosphate, Triethylammonium salt

MW : 621.31 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 255/355 nm ;  $\epsilon$  : 23 300/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description | Cat.#     | Qty   |
|-------------|-----------|-------|
| Mant-dATP   | FP-FQ2901 | 50 u  |
|             | FP-FQ2900 | 250 u |

Tung-Chung Mou, et al, "Structural Basis for the Inhibition of Mammalian Membrane Adenylyl Cyclase by 2'(3')-O-(NMethylanthraniloyl)- guanosine 5'-Triphosphate.", J. Biol. Chem. 280 (8):7253(2005)  
Churchich, "Binding of a fluorescent nucleotide analog to hsc70 - the effect of peptide protein interactions on the luminescence properties of the probe.", Eur. J. Biochem. 231 (3):736(1995)

### Mant-GDP

2'/3'-O-(N-Methyl-anthraniloyl)-guanosine-5'-diphosphate, Triethylammonium salt

MW : 574.32 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 252/355 nm ;  $\epsilon$  : 22 600/5 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description | Cat.#     | Qty   |
|-------------|-----------|-------|
| Mant-GDP    | FP-FQ2081 | 150 u |
|             | FP-FQ2080 | 750 u |

Thanbichler, et al. "Kinetics of the interaction of translation factor SelB from Escherichia coli with guanosine nucleotides and selenocysteine insertion sequence RNA.", J. Biol. Chem. 275:20458(2000)  
Murthy, et al, "Nucleotide binding by the erythrocyte transglutaminase/Gh protein, probed with fluorescent analogs of GTP and GDP.", Proc. Natl. Acad. Sci. USA 97:7744 (2000)

### Mant-dGDP

3'-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-diphosphate, Triethylammonium salt

MW : 558.33 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 252/355 nm ;  $\epsilon$  : 22 600/5 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description | Cat.#            | Qty   |
|-------------|------------------|-------|
| Mant-dGDP   | <b>FP-FQ4901</b> | 30 u  |
|             | <b>FP-FQ4900</b> | 150 u |

Nomanbhoy, et al., "Investigation of the GTP-binding GTPase cycle of Cdc42Hs using extrinsic reporter group fluorescence.", Biochemistry 35 (14):4602 (1996).

Leonard, et al., "Investigation of the GTP-binding GTPase cycle of cdc42hs using fluorescence spectroscopy.", Biochemistry 33 (40):12323. (1994)

### Mant-GTP

2'/3'-O-(N-methyl-anthraniloyl)-guanosine-5'-triphosphate, Triethylammonium salt

MW : 653.30 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 255/355 nm ;  $\epsilon$  : 23 300/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description | Cat.#            | Qty   |
|-------------|------------------|-------|
| Mant-GTP    | <b>FP-FQ2091</b> | 150 u |
|             | <b>FP-FQ2090</b> | 750 u |

Tung-Chung Mou, et al, "Structural Basis for the Inhibition of Mammalian Membrane Adenylyl Cyclase by 2'(3')-O-(NMethylanthraniloyl)-guanosine 5'-Triphosphate.", J. Biol. Chem. 280 (8):7253 (2005) : the mechanism of discrimination between guanosine and adenosine nucleotides.", Biochemistry 34:593(1995)

### Mant-dGTP

3'-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-triphosphate, Triethylammonium salt

MW : 637.30

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 252/355 nm ;  $\epsilon$  : 22 600/5 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description | Cat.#            | Qty  |
|-------------|------------------|------|
| Mant-dGTP   | <b>FP-FQ2911</b> | 10 u |
|             | <b>FP-FQ2910</b> | 50 u |

### Mant-GppNHp

2'/3'-O-(N-Methyl-anthraniloyl)-guanosine-5'-[( $\beta$ ,  $\gamma$ )-imido]triphosphate, Triethylammonium salt

MW : 651.31 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 252/355 nm ;  $\epsilon$  : 22 600/5 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description               | Cat.#            | Qty  |
|---------------------------|------------------|------|
| Mant-GppNHp (MAnt-GMPPNP) | <b>FP-FQ2551</b> | 10 u |
|                           | <b>FP-FQ2550</b> | 50 u |

Diebold, et al. "Molecular basis for Rac2 regulation of phagocyte NADPH oxidase", Nature Immunol. 2:211(2001)

Graham, et al., "The conserved arginine in rho-GTPase-activating protein is essential for efficient catalysis but not for complex formation with Rho.GDP and aluminum fluoride.", Biochemistry 38:985(1999)

### Mant-dGppNHp

3'-O-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-[( $\gamma$ ,  $\beta$ )-imido]triphosphate, Triethylammonium salt

MW : 636.32 (Anion)

$\lambda_{exc} / \lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  252/355 nm ;  $\epsilon$  22 600/5 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description                 | Cat.#            | Qty  |
|-----------------------------|------------------|------|
| Mant-dGppNHp (MAnt-dGMPPNP) | <b>FP-FQ2541</b> | 10 u |
|                             | <b>FP-FQ2540</b> | 50 u |

Scheidig, et al., "X-ray crystal-structure analysis of the catalytic domain of the oncogene product p21(H-Ras) complexed with caged GTP and mant dGppNHp.", J. Mol. Biol. 253 (1):132.(1995)

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### Mant-GTPgammaS

2'/3'-O-(N-Methyl-anthraniloyl)-guanosine-5'-( $\gamma$ -thio)-triphosphate, Triethylammonium salt

MW : 669.36 (Anion)

$\lambda_{exc}$  /  $\lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  252/355 nm ;  $\epsilon$  22 600/5 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description    | Cat.#     | Qty  |
|----------------|-----------|------|
| Mant-GTPgammaS | FP-FQ2921 | 10 u |
|                | FP-FQ2920 | 50 u |

Remmers, "Detection and quantitation of heterotrimeric G proteins by fluorescence resonance energy transfer.", Anal. Biochem. 257 (1):89(1998)

Remmers, et al. "Interdomain interactions regulate GDP release from heterotrimeric G proteins.", Biochemistry 38(42):13795(1999)

### Mant-XDP

2'/3'-O-(N-Methyl-anthraniloyl)-xanthosine-5'-diphosphate, Triethylammonium salt

MW : 575.32 (Anion)

$\lambda_{exc}$  /  $\lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  254/355 nm ;  $\epsilon$  19 000/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description | Cat.#     | Qty  |
|-------------|-----------|------|
| Mant-XDP    | FP-FQ2941 | 10 u |
|             | FP-FQ2940 | 50 u |

### Mant-XppNHp

2'/3'-O-(N-Methyl-anthraniloyl)-xanthosine-5'-[( $\gamma$ ,  $\beta$ )-imido]triphosphate, Triethylammonium salt

MW : 653.30 (Anion)

$\lambda_{exc}$  /  $\lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  254/355 nm ;  $\epsilon$  19 000/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description               | Cat.#     | Qty  |
|---------------------------|-----------|------|
| Mant-XppNHp (mant-XMPPNP) | FP-FQ5181 | 5 u  |
|                           | FP-FQ5180 | 25 u |

Rodnina, et al., "Codon-dependent conformational change of elongation-factor tu preceding GTP hydrolysis on the Ribosome", EMBO J. 14 (11):2613(1995)

Hazlett, et al., "Solution dynamics of p21(Ras) proteins bound with fluorescent nucleotides - a time-resolved fluorescence study.", Biochemistry-US 32 (49):13575.(1993)

### Mant-XTP

2'/3'-O-(N-Methyl-anthraniloyl)-xanthosine-5'-triphosphate, Triethylammonium salt

MW : 654.28 (Anion)

$\lambda_{exc}$  /  $\lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  254/355 nm ;  $\epsilon$  19 000/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description | Cat.#     | Qty  |
|-------------|-----------|------|
| Mant-XTP    | FP-FQ2951 | 10 u |
|             | FP-FQ2950 | 50 u |

Tung-Chung Mou, et al., "Structural Basis for the Inhibition of Mammalian Membrane Adenylyl Cyclase by 2'(3')-O-(NMethylanthraniloyl)- guanosine 5'-Triphosphate.", J. Biol. Chem. 280 (8):7253(2005)

### Mant-AppNHp

2'/3'-O-(N-Methyl-anthraniloyl)-Adenosine-5'-[( $\beta$ , $\gamma$ )-imido]triphosphate, Triethylammonium salt

MW : 636.32 (Anion)

$\lambda_{exc}$  /  $\lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  255/355 nm ;  $\epsilon$  23 300/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description               | Cat.#     | Qty  |
|---------------------------|-----------|------|
| Mant-AppNHp (Mant-AMPPNP) | FP-FQ2531 | 10 u |
|                           | FP-FQ2530 | 50 u |

Bujalowski, et al., "Kinetic mechanism of nucleotide cofactor binding to Escherichia coli replicative helicase DnaB protein. Stopped-flow kinetic studies using fluorescent, ribose-, and base-modified nucleotide analogues.", Biochemistry 39:2106(2000)

Moore, et al. "Kinetic mechanism of adenine nucleotide binding to and hydrolysis by the Escherichia coli Rep monomer. 1. Use of fluorescent nucleotide analogues.", Biochemistry 33:14550(1994)

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### Mant-ITP $\gamma$ S

2'/3'-O-(N-Methyl-anthraniloyl)-inosine-5'-( $\gamma$ -thio)-triphosphate, Sodium salt

MW : 654.35 (Anion)

| Description         | Cat.#     | Qty  |
|---------------------|-----------|------|
| Mant-ITP $\gamma$ S | FP-FQ2931 | 10 u |
|                     | FP-FQ2930 | 50 u |

### Mant-N6-Methyl-ATP

2'/3'-O-(N-Methyl-anthraniloyl)-N6-methyl-adenosine-5'-triphosphate, Sodium salt

MW : 650.32 (Anion)

$\lambda_{exc}$  /  $\lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 255/355 nm ;  $\epsilon$  : 23 300/5 800 mol<sup>-1</sup> cm<sup>-1</sup>

| Description        | Cat.#     | Qty   |
|--------------------|-----------|-------|
| Mant-N6-Methyl-ATP | FP-FQ4910 | 30 u  |
|                    | FP-FQ4911 | 150 u |

### NPE-caged-Mant-dGTP

3'-O-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-triphosphate, P3-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt

MW : 786.45 (Anion)

$\lambda_{exc}$  /  $\lambda_{em}$  : 355 / 448 nm

$\lambda_{max}$  : 252/355 nm ;  $\epsilon$  : 22 600/5 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description         | Cat.#     | Qty  |
|---------------------|-----------|------|
| NPE-caged-Mant-dGTP | FP-FQ3021 | 10 u |
|                     | FP-FQ3020 | 50 u |

# Nucleic acid labeling and modification

## Labeled and modified nucleotides


### Technical tip

"Caged compounds" serve as photo-induced precursors of biological effector compounds, i.e. nucleotides and  $\text{Ca}^{2+}$  chelators. Combined with reactive fluorescent probes they have been used to label biological molecules ; by FT-IR coupled with flash photolysis and by fast time-resolved UV-Visible spectroscopy coupled with laser flash photolysis. Applications :

- ◆ Kinetic studies of release and export exchanges of nucleotides (mitochondria<sup>3</sup>)
- ◆ Phosphorylation activation of enzymes ( $\text{Ca}^{2+}$  ATPase<sup>2</sup>)

1) Electrical currents associated with nucleotide transport by the reconstituted mitochondrial ADP/ATP carrier ; Nickolay B. et al ; Proc. Natl. Acad. Sci. USA, Vol. 93, pp. 664-668 (1996) Article.

2) ATP-Induced phosphorylation of the sarcoplasmic reticulum  $\text{Ca}^{2+}$  ATPase : Molecular interpretation of infrared difference spectra ; A. Barth, et al (1998) ; Biophys. J. 75 538-544 / Abstract.



Capacitive currents generated by the reconstituted AAC under UV flashes with loaded (A and C) and unloaded (B and D) vesicles in the presence of caged ATP (c-ATP) or caged ADP (c-ADP), and inhibitor of CAT and BKA.

### Caged nucleotides

Caged nucleotides contain a photolabile group that allows monitoring of fast kinetics especially in time-resolved X-ray crystallography : Flash photolysis of the caging group by ultraviolet light leads to a rapid and highly localized release of the biologically active nucleotide from the inactive caged compound.

Caged nucleotides are increasingly useful for :

- ◆ Identification of drug-targets and its binding site
- ◆ Determination of the affinity and selectivity of the drug-target interaction.
- ◆ Determination of the conformation of the binding pocket in structure-based drug design.

#### NPE-caged-ATP

Adenosine-5'-triphosphate, P<sup>3</sup>-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt  
MW : 653.30 (Anion)

| Description   | Cat.#            | Qty   |
|---------------|------------------|-------|
| NPE-caged-ATP | <b>FP-FQ2121</b> | 150 u |
|               | <b>FP-FQ2120</b> | 750 u |

Scheirlinckx, et al., "Monitoring of secondary and tertiary structure changes in the gastric H<sup>+</sup>/K<sup>+</sup>-ATPase by infrared spectroscopy", Eur. J. Biochem. 268 (13):3644(2001)

Barth, et al., "Substrate binding and enzyme function investigated by infrared spectroscopy.", FEBS Lett. 477:151(2000)

#### DMB-caged-ATP

Adenosine-5'-triphosphate, P<sup>3</sup>-(1-(3',5'-dimethoxyphenyl)-2-oxo-2-phenyl-ethyl)-ester, Triethylammonium salt  
MW : 758.44 (Anion)

| Description   | Cat.#            | Qty  |
|---------------|------------------|------|
| DMB-caged-ATP | <b>FP-FQ3001</b> | 10 u |
|               | <b>FP-FQ3000</b> | 50 u |

Sokolov, et al., «Fast transient currents in Na,K-ATPase induced by ATP concentration jumps from the P-3-[1-(3',5'-dimethoxyphenyl)-2-phenyl-2-oxoethyl ester of ATP.», Biophys. J. 74 (5):2285(1998)

Sokolov, et al., «Fast transient currents in the Na,K-ATPase induced by ATP concentration jump experiments from DMB caged ATP.», Biophys. J. 72 (2):242.

#### NPE-caged-GTP

Guanosine-5'-triphosphate, P<sup>3</sup>-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt  
MW : 669.30

| Description   | Cat.#            | Qty   |
|---------------|------------------|-------|
| NPE-caged-GTP | <b>FP-FQ2421</b> | 100 u |
|               | <b>FP-FQ2420</b> | 500 u |

Allin, et al., "Ras catalyzes CTP hydrolysis by shifting negative charges from gamma- to beta-phosphate as revealed by time-resolved FTIR difference spectroscopy.", Biochemistry-US 40 (10):3037(2001)

Scheidig, et al., "The pre-hydrolysis state of p21(ras) in complex with GTP : new insights into the role of water molecules in the GTP hydrolysis reaction of ras-like proteins.", Struct. Fold Des. 7 (11):1311(1999)

#### NPE-caged-mant-dGTP

3'-O-(N-Methyl-anthraniloyl)-2'-deoxy-guanosine-5'-triphosphate, P<sup>3</sup>-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt  
MW : 786.45 (Anion)

$\lambda_{\text{exc}} / \lambda_{\text{em}}$  : 355 / 448 nm

$\lambda_{\text{max}}$  : 252/355 nm ;  $\epsilon$  : 22 600/5 700 mol<sup>-1</sup> cm<sup>-1</sup>

| Description         | Cat.#            | Qty  |
|---------------------|------------------|------|
| NPE-caged-mant-dGTP | <b>FP-FQ3021</b> | 10 u |
|                     | <b>FP-FQ3020</b> | 50 u |

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### NPE-caged-AppNHp

Adenosine-5'-[( $\beta,\gamma$ )-imido]triphosphate, P<sup>3</sup>-(1-(2-nitrophenyl)-ethyl)-ester, triethylammonium salt  
MW : 652.32 (Anion)

| Description                         | Cat.#     | Qty   |
|-------------------------------------|-----------|-------|
| NPE-caged-AppNHp (NPE-caged-AMPPNP) | FP-FQ2411 | 100 u |
|                                     | FP-FQ2410 | 500 u |

### NPE-caged-GpCpp

Guanosine-5'-[( $\beta,\gamma$ )-methylene]triphosphate, P<sup>3</sup>-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt  
MW : 667.33 (Anion)

| Description                        | Cat.#     | Qty  |
|------------------------------------|-----------|------|
| NPE-caged-GpCpp (NPE-caged-GMPCPP) | FP-FQ3011 | 10 u |
|                                    | FP-FQ3010 | 50 u |

### NPE-caged-XDP

Xanthosine-5'-diphosphate, P<sup>3</sup>-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt  
MW : 591.32 (Anion)

| Description   | Cat.#     | Qty  |
|---------------|-----------|------|
| NPE-caged-XDP | FP-FQ3031 | 10 u |
|               | FP-FQ3030 | 50 u |

### NPE-caged-XppNHp

Xanthosine-5'-[( $\beta,\gamma$ )-imido]triphosphate, P<sup>3</sup>-(1-(2-nitrophenyl)-ethyl)-ester, Triethylammonium salt  
MW : 669.30 (Anion)

| Description                         | Cat.#     | Qty  |
|-------------------------------------|-----------|------|
| NPE-caged-XppNHp (NPE-caged-XMPPNP) | FP-FQ5201 | 5 u  |
|                                     | FP-FQ5200 | 25 u |

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### BrDU nucleotides

BrDU is an analog of thymidine, used alternatively to radioactive thymidine to track DNA synthesis by (sub)populations of cells. It is added to the culture at high concentration. Newly synthesized DNA thymidine is (partially) replaced by BrdU. After fixation and permeabilization of the cells the BrdU can be shown using a fluorochrome labeled antibody against the BrdU. Using DNase during the staining period the BrdU becomes more available to the antibody without influencing the fluorescence or the structure of the cell.

#### BrdU

5-Bromo-2'-deoxyuridine

MW : 307.1

BrdU can be incorporated into DNA during cell division and subsequently detected by a anti-BrdU antibody. The probe can be used to study cell-cycle kinetics.

| Description                    | Cat.#            | Qty    |
|--------------------------------|------------------|--------|
| BrdU (5-Bromo-2'-deoxyuridine) | <b>FP-18538A</b> | 100 mg |

Methods Cell Biol 41, 297(1994).

#### BrUTP

5-Bromouridine-5'-triphosphate

10 mM in TE buffer

MW : 628.99

BrUTP can be enzymatically incorporated into RNA.

| Description | Cat.#            | Qty   |
|-------------|------------------|-------|
| BrUTP       | <b>FP-T8200A</b> | 25 µl |

Literature : Nakayama C, et al., "Utilizations of various uridine 5'-triphosphate analogues by DNA-dependent RNA polymerases I and II purified from liver nuclei of the cherry salmon (*Oncorhynchus masou*)", J Biochem (Tokyo) 96, 1501(1984)

#### Br-dUTP

5-Bromo-2'-deoxyuridine-5'-triphosphate,

10 mM in TE buffer

5-Bromo-dUTP is widely used in TUNEL assay to detect apoptosis cells(1,2) and it is also a good substrate for reverse transcriptase.

| Description | Cat.#            | Qty   |
|-------------|------------------|-------|
| Br-dUTP     | <b>FP-T8199A</b> | 25 µl |

Li X., et al., "Detection of apoptosis and DNA replication by differential labeling of DNA strand breaks with fluorochromes of different color.", Exp Cell Res, 222, 28(1996). Li X., et al., "Labelling DNA strand breaks with BrdUTP. Detection of apoptosis and cell proliferation.", Cell Prolif, 28, 571 (1995). Rytting AS, et al., "Colorimetric capture assay for human-immunodeficiency-virus-I reverse transcriptase activity", Biotechnol Appl Biochem, 29, 241(1999)

Ekstrand DH, et al., "A sensitive assay for the quantification of reverse transcriptase activity based on the use of carrier-bound template and non-radioactive-product detection, with special reference to human-immunodeficiency-virus isolation", Biotechnol Appl Biochem, 23, 95(1996)

#### 5Br-dUDP

5-Bromo-2'-deoxy-uridine-5'-diphosphate, Sodium salt

MW : 465.04 (Anion)

10 mM solution

| Description | Cat.#         | Qty    |
|-------------|---------------|--------|
| 5Br-dUDP    | <b>FQ4570</b> | 30 µl  |
|             | <b>FQ4571</b> | 150 µl |

### Cyclic nucleotides

#### Guanosine-3',5'-cyclic monophosphate

Guanosine-3',5'-cyclic monophosphate, Sodium salt

MW : 343.19 (Anion)

| Description | Cat.#         | Qty    |
|-------------|---------------|--------|
| cGMP        | <b>FQ2840</b> | 50 mg  |
|             | <b>FQ2841</b> | 250 mg |

Nisoli, et al. «Mitochondrial biogenesis in mammals: the role of endogenous nitric oxide.», Science, 299 (5608), 896 (200). Rybalkin, et al, «PDE5 is converted to an activated state upon cGMP binding to the GAF A domain.», EMBO J., 22 (3), 469 (2003). Coates, et al., «Antagonistic pathways in neurons exposed to body fluid regulate social feeding in *Caenorhabditis elegans*.», Nature, 419(6910), 925 (2002)

#### Related products :

| Description              | Cat.#         |
|--------------------------|---------------|
| PCR Cleanup Column Kit   | <b>T66380</b> |
| Proteinase K             | <b>718960</b> |
| TrueBlue PCR cloning kit | <b>R58140</b> |



### Nucleoside bisphosphates

#### Adenosine-3',5'-bisphosphate

Adenosine-3',5'-bisphosphate, Triethylammonium salt

MW : 425.18 (Anion)

| Description                                                                                                                                                                                                                               | Cat.#         | Qty   |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|
| Adenosine-3',5'-bisphosphate (pAp)                                                                                                                                                                                                        | <b>FQ2281</b> | 100 u |
| Jacobson, et al., "Ribose modified nucleosides and nucleotides as ligands for purine receptors.", Nucleos. Nucleot. Nucl., 20 (4-7), 333 (2001)                                                                                           |               |       |
| Jung, et al., "Platelet collagen receptor integrin alpha(2)beta(1) activation involves differential participation of ADP-receptor subtypes P2Y1 and P2Y12 but not intracellular calcium change.", Eur. J. Biochem., 268 (12), 3513 (2001) |               |       |
| Sheng, et al., "Bacterial expression, purification, and characterization of rat hydroxysteroid sulfotransferase STa.", Protein Expres. Purif, 21 (1), 235 (2001)                                                                          |               |       |

#### Adenosine-2',5'-bisphosphate

Adenosine-2',5'-bisphosphate, Triethylammonium salt

MW : 425.18 (Anion)

| Description                                                                                                                                                                                                                           | Cat.#         | Qty   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|
| Adenosine-2',5'-bisphosphate                                                                                                                                                                                                          | <b>FQ2271</b> | 100 u |
|                                                                                                                                                                                                                                       | <b>FQ2270</b> | 500 u |
| Toth-Zsomboki, et al., "The P2Y(1) receptor antagonist adenosine-2',5'-diphosphate non-selectively antagonizes the platelet P2X(1) ion channel.", Thromb. Haemostasis, 86 (5), 1338 (2001)                                            |               |       |
| Leonidas, et al., "Mapping the ribonucleolytic active site of eosinophil-derived neurotoxin (EDN) - High resolution crystal structures of EDN complexes with adenylic nucleotide inhibitors.", J. Biol. Chem., 276 (18), 15009 (2001) |               |       |
| Gu, et al., "Crystal structures of the complexes of trichosanthin with four substrate analogs and catalytic mechanism of RNA N-glycosidase.", Proteins, 39 (1), 37 (2000)                                                             |               |       |

1 unit = 1 µl at 10 mM solution

#### 2'OMe-Adenosine-3',5'-bisphosphate

2'-O-Methyl-adenosine-3',5'-bisphosphate, Sodium salt

MW : 439.21 (Anion)

| Description                                    | Cat.#         | Qty   |
|------------------------------------------------|---------------|-------|
| 2'OMe-Adenosine-3',5'-bisphosphate (2'OMe-pAp) | <b>FQ6430</b> | 100 u |

#### Guanosine-3',5'-bisphosphate

Guanosine-3',5'-bisphosphate, Triethylammonium salt

MW : 441.18 (Anion)

| Description                                                                                                                                                                                                  | Cat.#         | Qty   |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|-------|
| Guanosine-3',5'-bisphosphate (pGp)                                                                                                                                                                           | <b>FQ2371</b> | 100 u |
|                                                                                                                                                                                                              | <b>FQ2370</b> | 500 u |
| Kvint, et al., "Emergency derepression: stringency allows RNA polymerase to override negative control by an active repressor.", Mol. Microbiol., 35 (2), 435 (2000)                                          |               |       |
| Acharya, et al., "The transmission of the electronic character of guanine-9-yl drives the sugar-phosphate backbone torsions in guanosine 3',5'-bisphosphate.", Angew. Chem. Int. Edit., 38 (24), 3645 (1999) |               |       |

#### Guanosine-2',5'-bisphosphate

Guanosine-2',5'-bisphosphate, Triethylammonium salt

MW: 441.18 (Anion)

| Description                  | Cat.#         | Qty   |
|------------------------------|---------------|-------|
| Guanosine-2',5'-bisphosphate | <b>FQ2361</b> | 100 u |
|                              | <b>FQ2360</b> | 500 u |

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### Xanthosine and Inosine nucleotides

Xanthosine nucleotides are used as G-nucleotide analogs in signal transduction research, molecular and cell biology.

Inosine is an analog of purine, and inosine nucleotides are widely used for mechanistic studies on ATP (cAMP) or GTP (cGMP) binding proteins.

#### XDP

Xanthosine-5'-diphosphate, Triethylammonium salt

MW : 442.17 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| XDP         | <b>FQ2401</b> | 100 u |
|             | <b>FQ2400</b> | 500 u |

Legate, et al., "Nucleotide-dependent binding of the GTPase domain of the signal recognition particle receptor betasubunit to the alpha-subunit.", *J. Biol. Chem.*, 275 (35), 27439 (2000). Yu, et al., "Interaction of the xanthine nucleotide binding Go alpha mutant with G protein-coupled receptors.", *J. Biol.Chem.* 273 (46), 30183 (1998)

#### XTP

Xanthosine-5'-triphosphate, Triethylammonium salt

MW : 521.14 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| XTP         | <b>FQ2111</b> | 150 u |
|             | <b>FQ2110</b> | 750 u |

Fulga, et al., "SR beta coordinates signal sequence release from SRP with ribosome binding to the translocon.", *EMBO J.*, 20 (9), 2338(2001). Legate, et al., "Nucleotide-dependent binding of the GTPase domain of the signal recognition particle receptor betasubunit to the alpha-subunit.", *J. Biol. Chem.* 275 (35), 27439 (2000)

#### IMP

Inosine-5'-monophosphate, Triethylammonium salt

MW : 347.20 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| IMP         | <b>FQ2041</b> | 150 u |
|             | <b>FQ2040</b> | 750 u |

Bhattacharya, et al., "Regiospecificity of nucleotide-amino acid mating vs. water dynamics: a key to protein-nucleic acid assemblies: structure of undecahydrated inosine-5'-monophosphate and L-glutamic acid.", *J. Chem. Crystallogr.* 30 (10), 655 (2000)  
Heroux, et al., "Crystal structures of the Toxoplasma gondii hypoxanthine-guanine phosphoribosyltransferase-GMP and -IMP complexes : Comparison of purine binding interactions with the XMP complex.", *Biochemistry-US* 38 (44), 14485 (1999)

#### IDP

Inosine-5'-diphosphate, Triethylammonium salt

MW : 426.17 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| IDP         | <b>FQ2381</b> | 100 u |
|             | <b>FQ2380</b> | 500 u |

AlAli, et al., "Effects of metal ions on the activity of cytosolic phosphoenolpyruvate carboxykinase from camel kidney.", *Arab. Gulf J. Sci. Res.* 14 (3), 535 (1996). Vial, et al., "Purification, partial kinetic characterization and reactive sulfhydryl-groups of the phosphoenolpyruvate carboxykinase from peromyscus-purpuratus adductor muscle.", *Comp. Biochem. Phys. B.* 112 (3):451 (1995)

#### ITP

Inosine-5'-triphosphate, Triethylammonium salt

MW : 505.14

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| ITP         | <b>FQ2051</b> | 150 u |
|             | <b>FQ2050</b> | 750 u |

Noji, et al., "Purine but not pyrimidine nucleotides support rotation of F(1)-ATPase.", *J. Biol. Chem.*, 276 (27), 25480 (2001)  
Bianchi, et al., "Intramolecular equilibria in metal ion complexes of guanosine 5'-triphosphate (GTP(4-)) and inosine 5'-triphosphate (ITP(4-)) in aqueous solution.", *J. Inorg. Biochem.*, 86 (1), 148 (2001). Chakrabarti, et al., "Nucleoside triphosphate specificity of tubulin.", *Biochemistry*, 39 (33), 10269 (2000)

#### cIMP

Inosine-3',5'-cyclic monophosphate, Triethylammonium salt

MW : 329.18 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| cIMP        | <b>FQ2850</b> | 50 u  |
|             | <b>FQ2851</b> | 250 u |

Wang, et al., "RNA polymerase-cNMP-ligated cAMP receptor protein (CRP) mutant interactions in the enhancement of transcription by CRP mutants.", *J. Biol. Chem.*, 275 (43), 33457 (2000). Shapiro, et al., "Structural basis for ligand selectivity of heteromeric olfactory cyclic nucleotide-gated channels.", *Biophys. J.*, 78 (5), 2307 (2000). Sunderman, et al., "Sequence of events underlying the allosteric transition of rod cyclic nucleotide-gated channels.", *J.Gen. Physiol.*, 113 (5), 621 (1999)

1 u = 1 µl at 10 mM

Related products :

| Description                         | Cat.#            | Qty  |
|-------------------------------------|------------------|------|
| XppNHp (XMPPNP)                     | <b>FQ2990</b>    | 50 u |
| XTPγS,                              | <b>FQ3640</b>    | 5 u  |
| Mant-XDP                            | <b>FP-FQ2941</b> | 10 u |
| Mant-XTP                            | <b>FP-FQ2951</b> | 10 u |
| Mant-XppNHp (mant-XMPPNP)           | <b>FP-FQ5181</b> | 5 u  |
| NPE-caged-XDP                       | <b>FP-FQ3031</b> | 10 u |
| NPE-caged-XTP                       | <b>FP-FQ3041</b> | 10 u |
| NPE-caged-XppNHp (NPE-caged-XMPPNP) | <b>FP-FQ5201</b> | 5 u  |
| IppNHp (IMPPNP)                     | <b>FQ2890</b>    | 50 u |
| Mant-EDA-ATP                        | <b>FP-FQ5811</b> | 20 u |
| Mant-ITPgS                          | <b>FP-FQ2931</b> | 10 u |

### Halogen-containing nucleotides

Halogen-containing nucleotides are used as rational phasing tools for protein crystallography, as specific substrates for studies on nucleotide binding proteins or as starting material for several nucleotide modifications.

#### 5I-dUTP

5-Iodo-2'-deoxy-uridine-5'-triphosphate, Sodium salt

MW : 591.01 (Anion)

| Description | Cat.#  | Qty   |
|-------------|--------|-------|
| 5I-dUTP     | FQ2600 | 50 u  |
|             | FQ2601 | 250 u |

#### 5I-dCTP

5-Iodo-2'-deoxy-cytidine-5'-triphosphate, Sodium salt

MW : 590.03 (Anion)

| Description | Cat.#  | Qty   |
|-------------|--------|-------|
| 5I-dCTP     | FQ2590 | 50 u  |
|             | FQ2591 | 250 u |

Related products :

| Description | Cat.#     | Qty   |
|-------------|-----------|-------|
| 5Br-dUDP    | FQ4570    | 30 u  |
| 5Br-dUTP    | FP-T8199A | 25 µl |

### Bisubstrate inhibitors

Bisubstrate inhibitors are used for mechanistic studies on nucleotide kinases.

#### AP<sub>3</sub>A

P<sup>1</sup>-(5'-Adenosyl) P<sup>3</sup>-(5'-adenosyl) triphosphate, Sodium salt

MW : 753.38 (Anion)

| Description       | Cat.#  | Qty   |
|-------------------|--------|-------|
| AP <sub>3</sub> A | FQ2730 | 50 u  |
|                   | FQ2731 | 250 u |

Guranowski, *et al.*, "Selective degradation of 2'-adenylated diadenosine tri- and tetraphosphates, Ap(3)A and Ap(4)A, by two specific human dinucleoside polyphosphate hydrolases.", *Arch. Biochem. Biophys.*, 373, 218 (2000)

Luo, *et al.*, "Identification and characterization of diadenosine 5',5''-P<sub>1</sub>,P<sub>2</sub>-diphosphate and diadenosine 5',5''-P<sub>1</sub>,P<sub>3</sub>-triphosphate in human myocardial tissue.", *FASEB J.*, 13, 695 (1999)

Luthje, *et al.*, "Catabolism of Ap<sub>4</sub>A and Ap<sub>3</sub>A in whole blood. The dinucleotides are long-lived signal molecules in the blood ending up as intracellular ATP in the erythrocytes.", *Eur. J. Biochem.*, 173, 241 (1988)

#### AP<sub>4</sub>A

P<sup>1</sup>-(5'-Adenosyl) P<sup>4</sup>-(5'-adenosyl) tetraphosphate, Sodium salt

MW : 832.36 (Anion)

| Description       | Cat.#  | Qty   |
|-------------------|--------|-------|
| AP <sub>4</sub> A | FQ2301 | 100 u |
|                   | FQ2300 | 500 u |

Vartanian, *et al.*, "Ap<sub>4</sub>A induces apoptosis in human cultured cells.", *FEBS Lett.*, 456, 175 (1999)

Campbell, *et al.*, "Characterization of P<sub>1</sub>,P<sub>4</sub>-diadenosine 5'-tetraphosphate binding on bovine aortic endothelial cells.", *Arch. Biochem. Biophys.*, 364, 280 (1999)

Guedon, *et al.*, "Effect of diadenosine tetraphosphate microinjection on heat shock protein synthesis in *Xenopus laevis* oocytes.", *EMBO J.*, 4, 3743 (1985)

#### AP<sub>5</sub>A

P<sup>1</sup>-(5'-Adenosyl) P<sup>5</sup>-(5'-adenosyl) pentaphosphate, Sodium salt

MW : 911.83 (Anion)

| Description       | Cat.#  | Qty   |
|-------------------|--------|-------|
| AP <sub>5</sub> A | FQ2760 | 50 u  |
|                   | FQ2761 | 250 u |

Wildman, *et al.*, "Selectivity of diadenosine polyphosphates for rat P<sub>2</sub>X receptor subunits.", *Eur. J. Pharmacol.*, 367, 119 (1999)

Pintor, *et al.*, "Presence of dinucleotide and ATP receptors in human cerebrocortical synaptic terminals.", *Eur. J. Pharmacol.* 366, 159 (1999)

1 u = 1 µl at 10 mM

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### AP<sub>6</sub>A

P<sup>1</sup>-(5'-Adenosyl) P<sup>6</sup>-(5'-adenosyl) hexaphosphate, Sodium salt

MW : 990.30 (Anion)

| Description       | Cat.#  | Qty   |
|-------------------|--------|-------|
| AP <sub>6</sub> A | FQ5600 | 20 u  |
|                   | FQ5601 | 100 u |

Morii, *et al.*, "Adenosine(5')hexaphospho(5')adenosine stimulation of a Ca(2+)-induced Ca(2+)-release channel from skeletal muscle sarcoplasmic reticulum.", *Eur. J. Biochem.*, 205, 979 (1992)

Luo, *et al.*, "Identification of diadenosine hexaphosphate in human erythrocytes.", *Hypertension*, 34, 872 (1999)

### AP<sub>5</sub>U

P<sup>1</sup>-(5'-Adenosyl) P<sup>5</sup>-(5'-uridyl) pentaphosphate, Sodium salt

MW : 888.29 (Anion)

| Description       | Cat.#  | Qty   |
|-------------------|--------|-------|
| AP <sub>5</sub> U | FQ2790 | 50 u  |
|                   | FQ2791 | 250 u |

Cheng, *et al.*, "Homogeneous uridine kinase from Ehrlich ascites tumor: substrate specificity and inhibition by bisubstrate analogs.", *Mol Pharmacol.*, 30 (2), 159 (1986)

### AP<sub>4</sub>dT

P<sup>1</sup>-(5'-adenosyl) P<sup>4</sup>-(5'-(2'-deoxy-thymidyl)) tetraphosphate, Triethylammonium salt

MW : 807.34 (Anion)

| Description        | Cat.#  | Qty   |
|--------------------|--------|-------|
| AP <sub>4</sub> dT | FQ2740 | 50 u  |
|                    | FQ2741 | 250 u |

### AP<sub>5</sub>dT

P<sup>1</sup>-(5'-adenosyl) P<sup>5</sup>-(5'-(2'-deoxy-thymidyl)) pentaphosphate, Triethylammonium salt

MW : 886.32 (Anion)

| Description        | Cat.#  | Qty   |
|--------------------|--------|-------|
| AP <sub>5</sub> dT | FQ2770 | 50 u  |
|                    | FQ2771 | 250 u |

Lavie, *et al.*, "Structural basis for efficient phosphorylation of 3'-azidothymidine monophosphate by Escherichia coli thymidylate kinase.", *Proc. Natl. Acad. Sci., USA* 95 (24), 14045 (1998)

Lavie, *et al.*, "Crystal structure of yeast thymidylate kinase complexed with the bisubstrate inhibitor P-1-(5'-adenosyl) P-5-(5'-thymidyl) pentaphosphate (TP(5)A) at 2.0 angstrom resolution: Implications for catalysis and AZT activation.", *Biochemistry*, 37, 3677 (1998)

### AP<sub>4</sub>G

P<sup>1</sup>-(5'-Adenosyl) P<sup>4</sup>-(5'-guanosyl) tetraphosphate, Triethylammonium salt

MW : 843.36 (Anion)

| Description       | Cat.#  | Qty   |
|-------------------|--------|-------|
| AP <sub>4</sub> G | FQ2750 | 50 u  |
|                   | FQ2751 | 250 u |

Ortiz, *et al.*, "specific synthesis of adenosine (5')tetraphospho(5')nucleoside and adenosine- (5')oligophospho(5')adenosine (n-greaterthan-4) catalyzed by firefly luciferase.", *Eur. J. Biochem.*, 212 (1), 263 (1993)

Palfi, *et al.*, "Alterations in the accumulation of adenylylated nucleotides in heavy-metal-ion-stressed and heat-stressed *Synechococcus* sp strain pcc-6301, a cyanobacterium, in light and dark.", *Biochem. J.* 276:487 (1991)

### AP<sub>5</sub>G

P<sup>1</sup>-(5'-Adenosyl) P<sup>5</sup>-(5'-guanosyl) pentaphosphate, Triethylammonium salt

MW : 927.33 (Anion)

| Description       | Cat.#  | Qty   |
|-------------------|--------|-------|
| AP <sub>5</sub> G | FQ2780 | 50 u  |
|                   | FQ2781 | 250 u |

Prinz, *et al.*, "Binding of nucleotides to guanylate kinase, p21(ras), and nucleoside-diphosphate kinase studied by nano-electrospray mass spectrometry.", *J. Biol. Chem.*, 274 (50), 35337 (1999)

Ortiz, *et al.*, "specific synthesis of adenosine (5')tetraphospho(5')nucleoside and adenosine-(5')oligophospho(5')adenosine (n-greaterthan-4) catalyzed by firefly luciferase.", *Eur. J. Biochem.*, 212 (1), 263 (1993)

1 u = 1 µl at 10 mM

### Antiviral nucleotides

#### AzTMP

3'-Azido-2',3'-dideoxy-thymidine-5'-monophosphate, Sodium salt

MW : 346.21 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| AzTMP       | <b>FQ5620</b> | 20 u  |
|             | <b>FQ5621</b> | 100 u |

Cruchaga, *et al.*, "Inhibition of Phosphorolysis Catalyzed by HIV-1 Reverse Transcriptase Is Responsible for the Synergy Found in Combinations of 3'-Azido-3'-deoxythymidine with Nonnucleoside Inhibitors", *Biochemistry*, 44 (9), 3535 (2005)

Chenal-Francisque, *et al.*, "The highly similar TMP kinases of *Yersinia pestis* and *Escherichia coli* differ markedly in their AZTMP phosphorylating activity.", *Eur. J. Biochem.* 265, 112 (1999)

#### AzTTP

3'-Azido-2',3'-dideoxy-thymidine-5'-triphosphate, Sodium salt

MW : 504.16 (Anion)

| Description | Cat.#         | Qty  |
|-------------|---------------|------|
| AzTTP       | <b>FQ2810</b> | 10 u |
|             | <b>FQ2811</b> | 50 u |

Cruchaga, *et al.*, « Inhibition of Phosphorolysis Catalyzed by HIV-1 Reverse Transcriptase Is Responsible for the Synergy Found in Combinations of 3'-Azido-3'-deoxythymidine with Nonnucleoside Inhibitors », *Biochemistry*, 44 (9), 3535 (2005)

Akeb, *et al.*, «The production and evaluation of antibodies for enzyme immunoassay of AZTTP.», *Nucleosides Nucleotides Nucleic Acids*. 20, 243 (2001) Faraj, *et al.*, « Effects of beta-L-3'-azido-3'-deoxythymidine 5'-triphosphate on host and viral DNA polymerases.», *Antiviral Res.* 47, 97 (2000)

#### d4TMP

2',3'-Didehydro-2',3'-dideoxy-thymidine-5'-monophosphate, Sodium salt

MW : 303.18 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| d4TMP       | <b>FQ5640</b> | 20 u  |
|             | <b>FQ5641</b> | 100 u |

De Clercq, "Highlights in the development of new antiviral agents.", *Mini Rev. Med. Chem.*, 2, 163 (2002)

Mas, *et al.*, "Multidrug-resistant HIV-1 reverse transcriptase: involvement of ribonucleotide-dependent phosphorolysis in cross-resistance to nucleoside analogue inhibitors.", *J. Mol. Biol.*, 267, 181 (2002)

Selmi, *et al.*, "The valine-to-threonine 75 substitution in human immunodeficiency virus type 1 reverse transcriptase and its relation with stavudine resistance.", *J. Biol. Chem.*, 276, 13965 (2001)

#### d4TTP

2',3'-Didehydro-2',3'-dideoxy-thymidine-5'-triphosphate, Sodium salt

MW : 461.13 (Anion)

| Description | Cat.#         | Qty  |
|-------------|---------------|------|
| d4TTP       | <b>FQ2861</b> | 10 u |
|             | <b>FQ2860</b> | 50 u |

Ray, *et al.*, "Insights into the molecular mechanism of inhibition and drug resistance for HIV-1 RT with carbovir triphosphate.", *Biochemistry*, 41, 5150 (2002)

Vaccaro, *et al.*, "Mechanism of inhibition of the human immunodeficiency virus type 1 reverse transcriptase by d4TTP: an equivalent incorporation efficiency relative to the natural substrate dTTP.", *Antimicrob. Agents. Chemother.*, 44, 217 (2000)

#### 3TCMP

2', 3'-Dideoxy-3'-thia-cytidine-5'-monophosphate, Sodium salt (L isomer)

MW : 308.22 (Anion)

| Description | Cat.#         | Qty  |
|-------------|---------------|------|
| 3TCMP       | <b>FQ2561</b> | 10 u |
|             | <b>FQ2560</b> | 50 u |

Mas, *et al.*, "Multidrug-resistant HIV-1 reverse transcriptase: involvement of ribonucleotide-dependent phosphorolysis in cross-resistance to nucleoside analogue inhibitors.", *J. Mol. Biol.*, 323, 181 (2002)

#### 3TCTP

2', 3'-Dideoxy-3'-thia-cytidine-5'-triphosphate, Sodium salt (L isomer)

MW : 466.16 (Anion)

| Description | Cat.#         | Qty  |
|-------------|---------------|------|
| 3TCTP       | <b>FQ3650</b> | 5 u  |
|             | <b>FQ3651</b> | 25 u |

Boyer, *et al.*, "YADD mutants of human immunodeficiency virus type 1 and Moloney murine leukemia virus reverse transcriptase are resistant to lamivudine triphosphate (3TCTP) in vitro.", *J. Virol.*, 75, 6321 (2001).

Kewn, *et al.*, "Development of enzymatic assays for quantification of intracellular lamivudine and carbovir triphosphate levels in peripheral blood mononuclear cells from human immunodeficiency virus-infected patients.", *Antimicrob. Agents Chemother.*, 46, 135 (2002)

1 u = 1 µl at 10 mM

# Nucleic acid labeling and modification

## Labeled and modified nucleotides

### Ribavirin-5'-triphosphate

1-β-D-Ribofuranosyl-1,2,4-triazole-3-carboxamide-5'-triphosphate, Sodium salt  
MW: 481.12 (Anion)

| Description               | Cat.#         | Qty   |
|---------------------------|---------------|-------|
| Ribavirin-5'-triphosphate | <b>FQ5940</b> | 20 u  |
|                           | <b>FQ5941</b> | 100 u |

Maag, *et al.*, "Hepatitis C virus RNA-dependent RNA polymerase (NS5B) as a mediator of the antiviral activity of ribavirin.", *J. Biol. Chem.*, 276 (49), 46094 (2001)

Lanford, *et al.*, "Ribavirin induces error-prone replication of GB virus B in primary tamarin hepatocytes.", *J. Virol.*, 75 (17), 8074 (2001)

Crotty, *et al.*, "The broad-spectrum antiviral ribonucleoside ribavirin is an RNA virus mutagen.", *Nat. Med.*, 6 (12), 1375 (2000)

## Other nucleotides analogs

### N6-Methyl-ATP

N6-Methyl-adenosine-5'-triphosphate, Potassium salt  
MW : 518.8 (Anion)

| Description   | Cat.#         | Qty   |
|---------------|---------------|-------|
| N6-Methyl-ATP | <b>FQ2101</b> | 150 u |
|               | <b>FQ2100</b> | 750 u |

### AP4

Adenosine-5'-tetraphosphate, Triethylammonium salt  
MW : 583.13 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| AP4         | <b>FQ2291</b> | 100 u |
|             | <b>FQ2290</b> | 500 u |

### ara-ATP

Adenine-arabinofuranoside-5'-triphosphate, Sodium salt  
MW : 504.16 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| ara-ATP     | <b>FQ2800</b> | 50 u  |
|             | <b>FQ2801</b> | 250 u |

Genini, *et al.*, "Nucleotide requirements for the in vitro activation of the apoptosis protein-activating factor-1-mediated caspase pathway.", *J. Biol. Chem.*, 275 (1), 29 (2000)

Ghoshal, *et al.*, "Ara-ATP impairs 3'-end processing of pre-messenger-RNAs by inhibiting both cleavage and polyadenylation.", *Nucleic Acids Res.* 19 (21):5871 (1991)

### ε-ATP

1,N<sup>6</sup>-Etheno-adenosine-5'-triphosphate, Sodium salt  
MW : 528.18 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| ε-ATP       | <b>FQ5780</b> | 200 u |
|             | <b>FQ5781</b> | 100 u |

Literature : Aguilar, *et al.*, "Ectoenzymatic breakdown of diadenosine polyphosphates by *Xenopus laevis* oocytes.", *Eur. J. Biochem.*, 268 (5), 1289 (2001)

Churchich, *et al.*, "A catalytic site of protein disulfide isomerase probed with adenosine-5'-triphosphate analogs.", *BBA-Protein Struct.M.*, 1479 (1-2), 293 (2000)

Gualix, *et al.*, "Studies of chromaffin granule functioning by flow cytometry : Transport of fluorescent epsilon-ATP and granular size increase induced by ATP.", *Receptor Channel*, 6 (6), 449 (1999)

### 8-Oxo-GTP

8-Oxo-guanosine-5'-triphosphate, Sodium salt  
MW : 536.15 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| 8-Oxo-GTP   | <b>FQ2710</b> | 50 u  |
|             | <b>FQ2711</b> | 250 u |

Kuryavii, *et al.*, "Kinetics of inhibition of escherichia-coli RNA polymerase-catalyzed synthesis of the dinucleotide pppApU by 8-oxo-GTP and 8-BrGTP on the A1 promoter of Bacteriophage-t7 delta-d111 DNA using a limited set of substrates.", *Mol. Biol.* + 23 (3), 648 (1989)

1 u = 1 μl at 10 mM



### 8-Oxo-dGTP

8-Oxo-2'-deoxy-guanosine-5'-triphosphate, Sodium salt

MW : 519.15 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| 8-Oxo-dGTP  | <b>FQ4790</b> | 30 u  |
|             | <b>FQ4791</b> | 150 u |

Kai, et al., "An oxidized nucleotide affects DNA replication through activation of protein kinases in *Xenopus* egg lysates.", *Nucl. Acids Res.*, 30 (2), 569 (2002). Canitrot, et al., "Nucleotide excision repair DNA synthesis by excess DNA polymerase beta: a potential source of genetic instability in cancer cells.", *FASEB J.* 14 (12), 1765 (2000)  
Nampalli, et al., "Efficient synthesis of 8-oxo-dGTP: A mutagenic nucleotide.", *Bioorg. Med. Chem. Lett.*, 10 (15), 1677 (2000)

### 6-Thio-GTP

6-Thio-guanosine-5'-triphosphate, Sodium salt

MW : 536.22 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| 6-Thio-GTP  | <b>FQ1991</b> | 150 u |
|             | <b>FQ1990</b> | 750 u |

Tiede, et al., "CD28-dependent Rac1 activation is the molecular target of azathioprine in primary human CD4 + T lymphocytes.", *J. Clin. Invest.*, 111, 1133 (2003). Poland, et al., "Entrapment of 6-thiophosphoryl-IMP in the active site of crystalline adenylosuccinate synthetase from *Escherichia coli*.", *J. Biol. Chem.*, 272 (24), 15200 (1997)

### 5F-UTP

5-Fluoro-uridine-5'-triphosphate, Sodium salt

MW : 499.11 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| 5F-UTP      | <b>FQ2221</b> | 100 u |
|             | <b>FQ2220</b> | 500 u |

Au, et al., "Reversed-phase ion-pair high-performance liquid-chromatographic assay of 5-fluorouracil, 5'-deoxy-5- fluorouridine, their nucleosides, monophosphate, diphosphate, and triphosphate nucleotides with a mixture of quaternary ammonium-ions.", *J. Chromatogr.*, 228, 245 (1982). Glazer, et al., "The effect of 5-fluorouridine 5'-triphosphate on RNA transcribed in isolated-nuclei in vitro.", *Mol.Pharmacol.* 17 (2), 279 (1980)

### 5-Aza-dCTP

5-Aza-2'-deoxy-cytidine-5'-triphosphate, Sodium salt

MW : 465.12 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| 5-Aza-dCTP  | <b>FQ5270</b> | 20 u  |
|             | <b>FQ5271</b> | 100 u |

### 6-Chlorpurine-riboside-5'-triphosphate

6-Chlorpurine-riboside-5'-triphosphate, Sodium salt

MW : 523.99 (Anion)

| Description                            | Cat.#         | Qty   |
|----------------------------------------|---------------|-------|
| 6-Chlorpurine-riboside-5'-triphosphate | <b>FQ2231</b> | 100 u |
|                                        | <b>FQ2230</b> | 500 u |

Patzelt, et al., "Study of (Na+K)-ATPase with 6-chlorpurine triphosphate, 6-mercaptapurine triphosphate and dinitrophenyl-6-thiopurine triphosphate.", *H-S Z. Physiol. Chem.*, 355 (10), 1237 (1974)

### 6-Mercaptapurine-riboside-5'-triphosphate

6-Mercaptapurine-riboside-5'-triphosphate, Sodium salt

MW : 521.20 (Anion)

| Description                               | Cat.#         | Qty   |
|-------------------------------------------|---------------|-------|
| 6-Mercaptapurine-riboside-5'-triphosphate | <b>FQ2241</b> | 100 u |
|                                           | <b>FQ2240</b> | 500 u |

Patzelt, et al., "Study of (Na+K)-ATPase with 6-chlorpurine triphosphate, 6-mercaptapurine triphosphate and dinitrophenyl-6-thiopurine triphosphate.", *H-S Z. Physiol. Chem.*, 355 (10), 1237 (1974)

### dPTP

6H,8H-3,4-Dihydro-pyrimido[4,5-c][1,2]oxazin-7-one-8-β-D-2'-deoxy-ribofuranosid-5'-triphosphate, Sodium salt

MW : 506.17 (Anion)

| Description | Cat.#         | Qty   |
|-------------|---------------|-------|
| dPTP        | <b>FQ4520</b> | 30 u  |
|             | <b>FQ4521</b> | 150 u |

Cummins, et al., "Synthesis and study of the fluorescein conjugate of the nucleotide dPTP.", *Nucleos. Nucleot.* 20 (4-7), 1049 (2001) ; Zaccolo, et al., "An approach to random mutagenesis of DNA using mixtures of Triphosphate derivatives of Nucleoside analogues.", *J. Mol. Biol.* 255 (4), 589 (1996)

1 u = 1 µl at 10 mM

Related products :

| Description               | Cat.#         |
|---------------------------|---------------|
| ε-AppNHp (ε-AMPPNP)       | <b>FQ3660</b> |
| Ribavirin-5'-triphosphate | <b>FQ5940</b> |

# Nucleic acid labeling and modification

## Fish labeling

### Label IT® Fluorescence In Situ Hybridization Kits (FISH)

- ◆ Enables convenient, high efficiency labeling of any probe.
- ◆ Obtain consistent and controlled labeling with the Label IT® non-enzymatic covalent attachment method.
- ◆ Optimized protocols for use in centromeric and PAINT analysis.
- ◆ Provides sensitive hybridization detection.

Fluorescence In Situ hybridization (FISH) uses fluorescent DNA probes to vividly paint genes or chromosomes, revealing the localization of specific nucleic acid sequences. This valuable technique has proven useful for gene mapping and the identification of chromosomal abnormalities. FISH entails the preparation of short sequences of fluorescently labeled DNA (probes) which are complementary to the target DNA sequences. Metaphase chromosome spreads are prepared, RNase treated, dehydrated and denatured according to the protocol provided. The slides are then hybridized with the fluorescently labeled DNA probe. Post hybridization, the slides are washed, counterstained with a DAPI antifade mixture, and analyzed using fluorescent microscopy.

The optimized Label IT® FISH Kits utilize Label IT® technology to prepare superior labeled probes, by design, and hybridize them to targeted interphase/metaphase chromosome spreads.

The Label IT® FISH Kits are supplied with labeling reagents, hybridization buffer, and an optimized protocol, which can be used to generate fluorescently labeled probes for a variety of FISH applications, including :


- ◆ Alpha-satellite centromeric probes - generated from repetitive sequences found at the centromeres of specific chromosomes. Researchers use this technique to enumerate specific chromosomes.
- ◆ Whole chromosome PAINT probes (Protocol Aberration Identification and Nomenclature Terminology) - collections of small probe sequences that hybridize to different regions along the length of the same chromosome. These probes are useful for examining chromosomal abnormalities.

Each FISH kit includes the appropriate Cy™3, TM-Rhodamine, Fluorescein or Biotin labeling reagent, labeling buffers and hybridization buffer. The Label IT® FISH Biotin Kit does not contain Biotin detection reagents.

| Description                              | Cat.#  | Qty     |
|------------------------------------------|--------|---------|
| Label IT® FISH Cy™3 Labeling Kit         | AO0901 | 10 µg * |
|                                          | AO0900 | 2 µg ** |
| Label IT® FISH TM-Rhodamine Labeling Kit | AO0891 | 10 µg   |
|                                          | AO0890 | 2 µg    |
| Label IT® FISH Fluorescein Labeling Kit  | AO0991 | 10 µg   |
|                                          | AO0990 | 2 µg    |
| Label IT® FISH Biotin Labeling Kit       | AO0921 | 10 µg   |
|                                          | AO0920 | 2 µg    |

\* Each kit contains sufficient reagents to label 10 µg of DNA and perform 40 PAINT hybridizations or 200 centromeric hybridizations.

\*\* Each Kit contains sufficient reagents to label 2 µg of DNA and perform 8 PAINT hybridizations or 40 centromeric hybridizations.



Multicolor chromosome PAINT analysis using Label IT® FISH Biotin labeled human chromosome 7 (Streptavidin-Cy™3 conjugate detection) and Label IT® FISH Fluorescein labeled human chromosome 6 on Jurkat metaphase spread.



Multicolor chromosome centromeric analysis using Label IT® FISH Cy™3 labeled human chromosome 7 and Label IT® FISH Fluorescein labeled human chromosome 6 on Jurkat metaphase spread.

#### Related products

| Description          | Cat.#     | Qty   |
|----------------------|-----------|-------|
| Acridine Orange      | FP-05920D | 50 g  |
| 7-Aminoactinomycin D | FP-132303 | 1 mg  |
| DAPI                 | FP-371867 | 10 mg |
| Hoechst 33342        | FP-BB1340 | 5 ml  |

D.154

