


FT-BE928A



Casein

Protease assay substrate

Product Information

Name : Casein-TAMRA
Catalog Number : FP-BE928A, 5 mg
Solubility: Dionized water, aqueous buffer or 1 M NaOH
Absorption / Emission : $\lambda_{exc} \lambda_{em} = 546/574$ nm

Name : Casein-Fluorescein
Catalog Number : FP-65802B, 5 mg
Solubility: Dionized water, aqueous buffer or 1 M NaOH
Absorption / Emission : $\lambda_{exc} \lambda_{em} = 494/521$ nm

Storage: $-20^{\circ}\text{C} > 1$ year

Introduction

Casein is a phosphoprotein found in milk. This product contains α -, β -, γ -, and κ -casein subunits. In the intact substrate, casein is heavily labeled with TAMRA, resulting in almost total quenching of the conjugate's fluorescence. Compared to the fluoresceinated casein, TAMRA-labeled casein demonstrates pH-insensitive fluorescence. It is readily used for continuous fluorometric measurement of protease activity. Protease-catalyzed hydrolysis relieves this quenching conjugate, yielding brightly green fluorescent dye-labeled peptides. The increase in fluorescence intensity is directly proportional to protease activity. The lightly TAMRA-labeled casein may be useful for a continuous assay if monitored by fluorescence polarization.

Directions for use

For a proteolytic assay¹, incubate the fluorescent conjugated casein with proteases in appropriate assay buffer. The protease activity is demonstrated by the increment of fluorescence. A detail protocol can be found in Universal Fluorimetric Protease Assay Kit (Cat# BK962A).

Sample Protocol For Trypsin

1. Make a 5-10mg/mL Casein, fluorescent conjugated stock solution in PBS buffer. Unused stock solution can be divided into single use aliquots and stored at -20°C , and avoid exposure to light.
2. Prepare 2X assay working solution by diluting the fluorescent conjugated stock solution into 50-100 mM Tris buffer (pH 7.4) at 100-400 $\mu\text{g/mL}$



FT-BE928A

Note1: The 2X Assay working solution is designed for detecting the activity of chymotrypsin, trypsin, thermolysin, proteinase K, protease XIV, and human leukocyte elastase. For other proteases, please refer to Appendix I for the appropriate assay buffer formula.

Note2: The optimum concentration of the assay working solution should be determined experimentally for individual proteases.

- Mix equal volume of the trypsin standards or samples with 2X Assay working solution.
- Monitor the fluorescence increase at Ex/Em = 490/525 nm for Fluorescein or 546/574 nm for TAMRA. For kinetic reading: Immediately start measuring fluorescence intensity continuously and record data every 5 minutes for 30 minutes. For end-point reading: Incubate the reaction at a desired temperature for 30 to 60 minutes, protected from light. Then measure the fluorescence intensity.

Appendix I

<i>Protease</i>	<i>1X Assay Buffer*</i>
Cathepsin D	20 mM Sodium Citrate, pH 3.0
Papain	20 mM sodium acetate, 20 mM cysteine, 2 mM EDTA, pH 6.5
PAE	20 mM sodium phosphate, pH 8.0
Pepsin	10 mM HCl, pH 2.0
Porcine pancreas elastase	10 mM Tris-HCl, pH 8.8
Subtilisin	20 mM potassium phosphate buffer, pH 7.6, 150 mM NaCl

References

- Twining SS** (1984). Fluorescein isothiocyanate-labeled casein assay for proteolytic enzymes. Anal Biochem 143, 30-4.
- Arthur JS and Mykles DL** (2000). Calpain zymography with casein or fluorescein isothiocyanate casein. Methods Mol Biol. 144, 109-16
- Farmer WH and Yuan ZY** (1991). A continuous fluorescent assay for measuring protease activity using natural protein substrate. Anal Biochem 197, 347-52

Related Product

- Universal Fluorimetric Protease Assay Kit (Green), [BK962A](#)
- Universal Fluorimetric Protease Assay Kit (Red), [BK963A](#)

Ordering information

Catalog size quantities and prices may be found at <http://www.fluoprobes.com>
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.