

Multi-Functional CrossLinkers

Products Information

Multifunctional crosslinkers, for the preparation of multimeric conjugate:

/ Homo-tri-&quadri- functional Amine-reactive homo multi / Sulfhydryl-reactive /

/ Hetero-tri-functional Amino and sulfhydryl reactive crosslinking reagents /

Applications:


Preparation of multivalent conjugates, e.g. antigenic carrier-peptides

Oriented coupling for probe design

- *Amine-reactive homo-multi-functional crosslinkers.*


Catalog #: L7962A, 250mg
Name: **NHS-3 (TSAT)**
 tris-Succinimidyl aminotriacetate
MW: 482.36 (K)
 Spacer Arm length: 4.2 Å.

- Amino reactive trifunctional crosslinking reagent
- Amino core for preparation of dendritic or molecular aggregates.
- Can be selectively aminated to generate mixed trimer




Catalog #: BU243A, 50mg
Name: **lc-NHS-3 (lc-TSAT)**
 tris-Succinimidyl (6-aminocaproyl)aminotriacetate
MW: 821.83 (K)

- TSAT analog with long chain spacer




Catalog #: BU248A, 50mg
Name: **NHS-4 (TKSPE)**
 tetrakis-(N-succinimidylcarboxypropyl)pentaerythritol
MW= 812.69 (K)
 Mp: 158-159 °C

- Reacts with 4 amines with its Succinimidyl groups
- Useful in forming multivalent ligands.



- *Sulfhydryl-reactive homo-multi-functional crosslinkers.*

Catalog #: 86685A, 50 mg
Name: **Mal-3 (TMEA)**
 tris-(2-Maleimidoethyl)amine
MW: 386.36 (K)
 Spacer length: 10.3 Å




FT-BU245A

Catalog #: BU247B, 25mg
Name: **MAL-4 (TKMPE)**

tetrakis-(3-Maleimidopropyl)pentaerythritol
 Mp: 90-91 °C

MW= 684.70 (K)

- Reacts with 4 sulfhydryl with its Maleimide groups
- Useful in forming multivalent ligands.



Store at +4°C (or -20°C for long term) (K)


- *Hetero-tri-functional Amino and sulfhydryl reactive crosslinking reagents.*

Catalog #: BU245A, 100mg
Name: **SDMB (Mal₂-NHS)**

Succinimidyl-3,5-dimaleimidophenyl benzoate
 Mp: 198-199 °C

- Reacts with one sulfhydryl with its Maleimide group, and 2 amino groups with its Succinimidyl groups.
- Useful for generating multivalent antitumor antibody fragments.¹

1. Schott, M.E., et. al. (1993) Bioconjugate Chem. 4, 153-165.




MW= 409.31

Store at +4°C (or -20°C for long term) (K)

Catalog #: [BU246A](#), 50mg
Name: **MDSI (Mal-NHS₂)**

Maleimido-3,5-disuccinimidyl isophthalate
 Mp: 247 °C (Decomp)

- Reacts with 2 sulfhydryl with its Maleimide groups, and one amino group with Succinimidyl group.
- Useful in forming mixed, trimeric aggregates of ligands or therapeutics.



MW= 455.34

Store at +4°C (or -20°C for long term) (K)

Introduction

Cross-linkers are chemical reagents used to conjugate molecules together by a covalent bound (forming a 'spacer arm') to associate the characteristics and biological activities of each component. In the extensive range of Uptima reagents for this purpose, this sheet presents **multi-functional cross-linkers**, that is to say they have distinct reactivities (toward amine, and sulfhydryls), allowing the conjugation of molecules in an oriented manner, and/or multiple same reactive groups. Multimeric conjugates can be prepared, that is taken to good account for the preparation of antigen in immunisation procedures, and tri- or quadri-valent probes.

Technical and Scientific Information

- The chemical group **N-hydroxysuccinimydyl (NHS)** reacts in aqueous phase on primary (-NH₂) and secondary amines (=NH) (in fact on its deprotonated form), optimally at neutral pH or higher : amines present in proteins (Lys aminoacid) and in a lower proportion on NH₂ located in terminal peptidic chains. The reaction competes with hydrolysis, that increases with pH, and with the high dilutions of the molecule that should be derivatized. In usual conditions, one should start with a ratio of 5-40 moles of maleimide per mole of protein.

FT-BU245A

- The **maleimide** group reacts very specifically with sulfhydryls –SH at neutral pH 6.5-7.5, forming a thioester, in comparison to iodoacetamides which react also with tyr, his, met. The reaction is rapid (a few minutes for cystein), but usually performed for 30-60min at RT or +4°C. In the absence of –SH, maleimide is well stable, but hydrolysis forming maleimic acid may compete significantly when pH goes up 8.0, and the reactivity specificity is lower. At pH higher than 8.0-8.5, a nucleophilic reaction with an adjacent amine can ring-open the crosslinker to yield crosslinked products with poor stability.
In usual conditions, one should start with a ratio of 1-20 moles of maleimide per mole of protein. With SH-peptides in optimal conditions, a molar 1:1 incubation ratio allows almost 1:1 coupling.

Related / associated products

- NHS reaction catalyzer: **EDAC** # [UP52005A](#) or DCC # [01202A](#)
- Quenchers for amine conjugations: **Hydroxylamine** # [13072](#), NHS-Acetate [UP68723A](#)
- 6-(N-trifluoroacetyl)caproic acid NHS (**TFCs**) # [L7727B](#) to protect amino group that can then be unmasked at pH 7.8-8.1
- 2-Iminothiolane (**Traut's reagent**) # [UP42425A](#) to convert amino group to un-protected sulfhydryl group
- SATA** # [84235A](#) to converts amino group to protected sulfhydryl group
- Reductants: DTT # [UP284250](#) or TCEP # [UP242210](#)
- DTNB (Elman's reagent)** ([01566H](#)) (SH dosage)
- PDA(Pyridine dithioethylamine hydrochloride) # [B11321](#) to converts sulfhydryl group in amino group.
- HomobiFunctional crosslinkers: i.e. DSS / BS3#[UP54940A](#)
- HeterobiFunctional crosslinkers: i.e. (sulfo)SMCC # [UP17412A](#)
- Hydrazone chemistry: [Conjugation kit \(BL1501\)](#) and [HynNic crosslinkers](#) (SANH #[BL9270](#), MHPH #[BL9401](#) SH-reactive)
- Desalting tools, i.e [CelluSep dialysis tubings](#), Desalting gel filtration columns # [UP84874](#)
- Buffers: PBS(Phosphate Buffer Saline) # [68723A](#), TBS(Tris Buffer Saline) # [74004A](#), MES Buffer [GS2960](#)
- Preservatives: AEBSF # [401070](#) and other protease inhibitors, SodiumAzide # [08112A](#)

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. Uptima is not liable for any damage resulting from handling or contact with this product.

Rev.T09E-M09E-J05E-H05E