

β-Galactosidase substrates regulating cell growth

Products Description

β-Galactosidase substrates regulating cell growth

Substrate name cat.number	MW (g·mol ⁻¹)	$\lambda_{exc}\lambda_{em}$ max. (nm) [a]	mol. abs. (M ⁻¹ cm ⁻¹) [a]	mol. abs. (M ⁻¹ cm ⁻¹)	
Inhibitors					
Fluorouridine-Gal (M) FP-BU7751 , 2mg	424.33	[218 / - nm]	.	H ₂ O, DMSO	5-Fluorouridine 5-O-β-D-Galactopyranoside
Tetracycline-Gal (M) FP-BU7761 , 5mg	604.6	[268 : - nm] [355 / - nm]	18 000 [13 000]	H ₂ O, DMSO	Tetracycline 10-O-β-D-Galactopyranoside
Chloramphenicol-Gal (M) FP-BU7721 , 5mg	485.27	[278 / - nm]	.	H ₂ O, DMSO	Chloramphenicol β-D-Galactoside
Stimulators					
	MW				Comment
BU7771, 5mg	506.64		1-Oleoyl-2-acetyl-3-β-D-Galactopyranosyl rac-glycerol		Synthetic diglyceride analog that activates protein kinase C. Soluble : H ₂ O, DMSO
BU7781, 5mg	560.71		1-Oleoyl-2-acetyl-3-β-D-Galactopyranosyl sn-glycerol		Synthetic diglyceride analog that activates protein kinase C. Soluble : H ₂ O, DMSO
BU7791, 10mg	331.33		Pyridoxine Galactoside		Vitamin B6 analog (pyridoxine). Soluble : DMSO, H ₂ O
BU7801, 2mg	315.28		Pyridoxal Galactoside		Vitamin B6 analog (pyridoxine), useful media component for lacZ positive cell selection; Soluble : DMSO, abs. EtOH
BU7811, 5mg	342.30		Myo-inositol Galactoside		Derivative of myo-inositol, a component of membrane phospholipids Soluble : H ₂ O, DMSO
BU7831, 2mg	538.50		Riboflavin Galactoside		Vitamin B2 analog (riboflavin), for transfected lacZ cell selection; Soluble : H ₂ O, DMSO, DMF
BU7841, 10mg	543.53		Pantothenic acid 2,4-di-O-β-D-Galactopyranoside		Vitamin B5 analog (pantothenic acid), for transfected lacZ cell selection. Soluble : DMSO, H ₂ O
BU7861, 2mg	786.50		Thiamine Galactoside		Vitamin B1 analog (thiamine), for transfected lacZ cell selection. Soluble: DMSO, H ₂ O
BU7921, 5mg	477.5		7-N-Benzoyl-Cephalosporanic acid L-glutamate ester		Cephalosporin C conjugate, released upon ampicillinase activity, Abs: 280nm

[a] in brackets: values upon β-galactosidase cleavage.

β-Glucuronidase substrates regulating cell growth

Others / Stimulators	MW				Comment
Bromoxynil-Glc (M) FP-BU7891 , 2mg	467.06			DMF, DMSO, EtOH, H ₂ O	Bromoxynil glucuronic acid methyl ester
Bromoxynil-Glc TAM (M) FP-BU7871 , 2mg	593.17			DMF, DMSO, EtOH	Bromoxynil glucuronic acid triacetate methyl ester
PET GLU (M) FP-BM8770 , 50mg			The hydrophobic thio-glycoside can be used as a inhibitor of β-Glucuronidase activity		Phenethyl 1-thio-β-D-Glucopyranosiduronic acid (

Storage: -20°C >1 year. (M)

Protect from light and moisture

Info@fluoprobes.com
Technical-support@fluoprobes.com
Order-online@fluoprobes.com

Contact your local distributor

FluoProbes®, powered by



213 Avenue J.F. Kennedy - BP 1140
 03103 Montluçon Cedex - France
 Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

Introduction

One of the most common reporter genes used in molecular biology is the E.coli lac Z gene that coded for an active subunit of β -galactosidase in vivo, due to several features: The E. coli β -galactosidase enzyme (EC 3.2.1.23) is generally absent in normal mammalian, yeast, some bacteria and even plant cells. Also, the enzyme has a wide substrate specificity (hydrolyzing D-galactose from various labeled β -galactosides), and it has a high turnover rate. As a result, monitoring lacZ expression has become routine to the point of detection of a few as 5 copies of β -galactosidase per cell by FACS analysis. This reporter system is commonly used for monitoring transfection efficiency in mammalian, yeast, and bacterial cells and identifying expression of recombinant fusion genes, as well as co-expressed genes or promoter efficiency.

This technical sheet presents our high quality β -galactosidase substrates that can be used to regulate cell growth, inhibitors, and stimulators.

[Fluorouridine-Gal](#)

[Chloramphenicol-Gal](#)

[Tetracycline-Gal](#)

[Bromoxynil-GLUc](#)


[Associated products and other gene reporter substrates](#)

#BU775 technical & scientific information

5-Fluorouridine 5-O- β -D-Galactopyranoside

This analog of 5-fluorouridine (5-FUR) is useful for lacZ specific release of 5-FUR in recombinant cells or tissues. 5-FUR is a potent anti-metabolite, with mode of action by inhibition of thymidylate synthetase or by incorporation into RNA in vivo.

NOTE: For research purposes only. Not for human or drug use!!



references – BU775/BU776

Watanabe KA, Matsuda A, Halat MJ, Hollenberg DH, Nisselbaum JS, Fox JJ. (1981) "Nucleosides. 114. 5-O-glucuronides of 5-fluorouridine and 5-fluorocytidine. Masked precursors of anticancer nucleosides." *J Med. Chem.* 24: 893-897.

Kanzawa F., Hoshi A., Kuretani K., (1980), "Differences between 5-fluoro-2'-deoxyuridine and 5-fluorouridine in their cytotoxic effect on growth of murine lymphoma L5178Y cells in in vivo and in vitro systems." *European Journal of Cancer* 16(8):1087-92.

Iapalucci-Espinoza S., Franze-Fernandez M.T., (1982), "Regulation of rRNA synthesis and processing in animal cells. Effect of nucleoside analogues." *The Biochemical Journal* 202(2):325-32.

Abraham, R., Aman, N., von Borstel, R., Darsley, M., Kamireddy, B., Kenten, J., Morris, G., Titmas, R., (1994), "Conjugates of COL-1 monoclonal antibody and β -d-galactosidase can specifically kill tumor cells by generation of 5-fluorouridine from the prodrug β -d-galactosyl-5-fluorouridine." *Cell Biophysics* 24/25:127-33.


FT-BU7721

#BU772 technical & scientific information

Chloramphenicol β -D-Galactoside

Upon enzymatic or chemical hydrolysis of the galactoside group, chloramphenicol, an antibiotic (bacteriostatic), is produced.

NOTE: For research purposes only. Not for human or drug use!!




#BU776 technical & scientific information

Tetracycline 10-O- β -D-Galactopyranoside

Tetracycline analog useful as an antibiotic for specific release in lacZ transfected cells and tissues.

NOTE: Extinction = 268nm(ϵ =18K), 355nm(13K) ; Absorption = 268nm(ϵ =18K), 355nm(13K)

NOTE: For research purposes only. Not for human or drug use!!



#BM787/BM789 technical & scientific information



Bromoxynil [1689-84-5] is a nitrile herbicide that is used for post-emergent control of annual broadleaf weeds. It is especially effective in the control of weeds in cereal, corn, sorghum, onions, flax and mint. This bromoxynil derivative can be used for selective release in GUS-positive plant systems, in a tissue specific manner. The compound works by inhibiting photosynthesis in the target plants.

Bromoxynil glucuronic acid triacetate methyl ester

The #BM787 beta-glucuronide derivative lacks activity until acted upon by intracellular beta-glucuronidase (cloned with the GUS marker gene). Additional ester groups help with cell permeation of this herbicide derivative.

The #BM789 beta-glucuronide derivative lacks activity until acted upon by intracellular beta-glucuronidase (cloned with the GUS marker gene). The methyl ester helps with cell permeation in vivo. Hence, the compound can selectively ablate GUS-positive cells and tissues in plant systems.

NOTE: For research purposes only. Not for human or drug use!!



references – BM787/BM789

Schafer, D.E., Chilcote, D. O., "Translocation and degradation of bromoxynil in a resistant and a susceptible species." Weed Science (1970) 18(6):729-32.

Davis, William H. "Production of genetically-controlled herbicide resistance in cotton plants in the absence of genetic engineering." U.S. Pat. Appl. Publ.(2003) US 20030024015.

FT-BU7721

Heim, Ute; Hillerbrand, Helke; Kunze, Irene; Herbers, Karin; Sonnewald, Uwe; Glickmann, Eric; Lein, Wolfgang; Hell, Ruediger; Jost, Ricarda., "Constitutive promoters of Arabidopsis and their use in expression of foreign genes in transgenic plants." PCT Int. Appl.(2003): WO 0306660.

Cessna, A. J., (1980) "Simultaneous extraction and detection of residues of (2,4-dichlorophenoxy)acetic acid and bromoxynil from wheat." J. Agricultural and Food Chemistry 28(6):1229-32.

Related / associated products and documents

See [BioSciences Innovations catalogue](#) and [e-search tool](#).

- IPTG Lac z inducer ([84853C](#))
- X-Gal ([40534M](#)), Blue lacZ β -Galactosidase Detection Kit (X-Gal based), [FP-BM8410](#))
- substrates for β -Glucuronidase (MUGlcU [FP-37744](#)) and β -Glucosidase (X-GLU [193325](#))
- Fluorescent Galactoside derivatives (MUGal/CUGal/FDGal/FMGal/TFMU-Gal/Res-Gal: [FP-BM8400](#))

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>

Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.

Rev.I03E