

CF™488A Dye



A highly photostable and bright green dye for 488 nm laser

Technical Summary

Abs/Em Maxima:	490/515 nm
Extinction coefficient:	70,000
Molecular weight:	~910
Flow cytometry laser line:	488 nm
Microscopy laser line:	488 nm
Direct replacement for:	FITC, FAM, Cy™2

Advantages

- Yields biologically more specific antibody conjugates and has less “spill-over” fluorescence in the red channel than other dye
- Extremely photostable
- Highly water-soluble and pH-insensitive




Figure 1. Absorption and emission spectra of CF™488A (blue), other dye (green) and FITC (red) conjugated to goat anti-mouse IgG in PBS.




Figure 2. Jurkat cells were stained with intracellular CD3 or isotype control (BD Biosciences) followed by goat anti-mouse IgG conjugates from the manufacturer's shown above. The inset graph depicts the median fluorescence of the population as analyzed on a Beckman Coulter FC-500. The histogram depicts an overlay of CD3-stained cells with CF488A (green), other dyes (as analyzed on a BD FACS Calibur in the FL1 channel).

CF™488A is a green fluorescent dye optimally excitable by the 488 nm argon laser line. Under common detection conditions, CF™488A is at least as bright as other 488 dye. However, a major advantage of CF™488A over other dye is that antibody conjugates prepared from the former are biologically more specific. Other dye carries multiple negative charges, which can significantly change the isoelectric point of the proteins the dye labels and consequently alter the specificity of the protein conjugates. CF™488A, on the other hand, is minimally charged. Thus, antibody conjugates prepared from the dye ensure biological detection with high signal-to-noise ratio. Another feature of CF™488A is that the emission peak wavelength is about 10 nm shorter than that of other dye and 15 nm shorter than that of the traditional green dye FITC (or FAM). The shorter wavelength of CF™488A offers the advantage of less fluorescence “spill-over” in the red channel in multi-color detection applications.

A list of CF™488A-based products are shown in Table 1. A full selection of secondary antibodies, antibody labeling kits, and other bio-conjugates including phalloidins, annexin V and α -bungarotoxin are also available for many CF™ dyes. Please visit the Biotium website at www.biotium.com for details.




Figure 3. Cryosections (6 μ m) of human control heart sections stained with anti-fibronectin followed by goat anti-rabbit IgG labeled with CF™488A, Cy™2 and other dye, respectively. Courtesy of Dr. Sawa Kostin at Max-Planck-Institute für Herz- und Lungenforschung (W.G. Kerckhoff-Institut) in Bad Nauheim (Hessen)-Germany.

If you are looking for an antibody conjugate not listed in our catalog, please let us know. We might be able to add it as a new product, or perform a custom conjugation for you.



CF™488A fluorescent reagents

Table 1. CF™488 Product List

Product Name	Size	Cat No.
CF™488A-Labeled Secondary Antibody Conjugates		
Donkey Anti-Goat IgG (H+L) whole antibody , 2 mg/mL (min X Chicken, Guinea Pig, Syrian Hamster, Horse, Human, Mouse, Rabbit, and Rat)	0.5 mL	20016
Donkey Anti-Mouse IgG (H+L) whole antibody , 2 mg/mL (min X Bovine, Chicken, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Rabbit, and Sheep)	0.5 mL	20014
Donkey Anti-Rabbit IgG (H+L) whole antibody , 2 mg/mL (min X Bovine, Chicken, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Mouse, Rat, and Sheep)	0.5 mL	20015
Donkey Anti-Rat IgG (H+L) whole antibody , 2 mg/mL (min X Bovine, Chicken, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Rabbit, and Sheep)	0.5 mL	20027
Donkey Anti-Sheep IgG (H+L) whole antibody , 2 mg/mL (min X Chicken, Guinea Pig, Syrian Hamster, Horse, Human, Mouse, Rabbit, and Rat)	0.5 mL	20024
Goat Anti-Chicken IgY (IgG) (H+L) whole antibody , 2 mg/mL (min X Bovine, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Mouse, Rabbit, Rat, and Sheep)	0.5 mL	20020
Goat Anti-Guinea Pig IgG (H+L) whole antibody , 2 mg/mL	0.5 mL	20017
Goat Anti-Human IgG (H+L) whole antibody , 2 mg/mL (min X Bovine, Horse, and Mouse)	0.5 mL	20022
Goat Anti-Mouse IgG (H+L) whole antibody , 2 mg/mL	0.5 mL	20010
Goat Anti-Mouse IgG (H+L) whole antibody , 2 mg/mL (min x Human, Bovine, Horse, Rabbit, and Swine)	0.5 mL	20018
Goat Anti-Mouse IgG (H+L), F(ab')₂ fragment , 2 mg/mL	0.25 mL	20011
Goat Anti-Rabbit IgG (H+L) whole antibody , 2 mg/mL	0.5 mL	20012
Goat Anti-Rabbit IgG (H+L) whole antibody , 2 mg/mL (min X Human, Mouse, and Rat)	0.5 mL	20019
Goat Anti-Rabbit IgG (H+L), F(ab')₂ fragment , 2 mg/mL	0.25 mL	20013
Goat Anti-Rat IgG (H+L) whole antibody , 2 mg/mL (min X Human, Bovine, Horse, and Rabbit)	0.5 mL	20023
Goat Anti-Swine IgG (H+L) whole antibody , 2 mg/mL	0.5 mL	20028
Rabbit Anti-Goat IgG (H+L) whole antibody , 2 mg/mL	0.5 mL	20021
Rabbit Anti-Mouse IgG (H+L) whole antibody , 2 mg/mL (min X Human)	0.5 mL	20026
Rabbit Anti-Rat IgG (H+L) whole antibody , 2 mg/mL (min X Human)	0.5 mL	20025
Other CF™488A-Labeled Products		
Annexin V , 50 µg/mL	0.5 mL	29005
α-Bungarotoxin	0.5 mg	00005
Phalloidin	300 U	00042
Streptavidin	1 mg	29034
CF™488A Reactive Dyes and Labeling Kits		
CF™488A, aminoxy	1 mg	92051
CF™488A hydrazide	1 mg	92152
CF™488A maleimide	1 µmole	92022
CF™488A succinimidyl ester	1 µmole	92120
CF™488A SE protein labeling kit	3 labelings (for 1 mg protein each)	92213
CF™488A SE microscale protein labeling kit	3 labelings (for 20-100 µg protein each)	92233

Listed products are for research use only. Not for use in diagnostic or therapeutic procedures.
CF is a trademark of Biotium; CF dye technologies are covered by pending US and international patents.



Interchim - 211 bis Avenue Kennedy - BP 1140 - 03103 Montluçon - France
Tél. (33) (0)4 70 03 73 01 - Fax (33) (0)4 70 03 82 60 - www.interchim.com