



Marker Gene Technologies, Inc

Instruction Manual

pNosdcGUS Plant GUS Expression Vector

Version 1.01

© Marker Gene Technologies, Inc. January, 2009

A Limited Label License covers this product (see Purchaser Notification). By use of this product, you accept the terms and conditions of the Limited Label License.

A plasmid expression vector for cloning and expression of GUS into plants.

Catalog no. M1390

Table of Contents

Table of Contents.....	2
Important Information.....	2
Shipping and Storage	2
Contents	2
Quality Control	2
Accessory Products	3
Methods.....	4
Using pNosdcGUS GUS Expression Vector.....	5
Transfection.....	6
Expression and Analysis.....	7
Creating Stable Cell Lines.....	7
Appendix.....	9
Recipes	9
Map and Features of pNosdcGUS GUS Expression Vecotor.....	10
Technical Service.....	19
Purchaser Notification	19
References.....	20

Important Information:

MSDS Sheets and product safety information are available by request from Marker Gene Technologies, Inc. and by accessing our web site at www.markergene.com.

Shipping and Storage

The pNosdcGUS plant GUS expression vector is shipped at room temperature. Store all samples of the vector at -20°C or below, once resuspended. Products are guaranteed for six months from date of shipment when stored properly.

Contents

Item	Concentration
pNosdcGUS plant GUS Expression Vector, lyophilized in TE buffer, pH 8.0	20 μg

Quality Control

The pNosdcGUS plant GUS expression vector has been qualified by restriction endonuclease digestion. pNosdcGUS is further qualified by transformation using an appropriate *Agrobacterium* strain in culture into *Arabidopsis thaliana* plant species and verified for activity.

Accessory Products

Additional products that may be used with the pNosdcGUS Plant GUS Expression Vector are now available from Marker Gene.

Ordering information is provided below.

Product	Unit Size	Catalog no.
----------------	------------------	--------------------

Expression of your recombinant fusion protein can be detected using:

MarkerGene™ β -Glucuronidase (GUS) Reporter Gene Activity Detection Kit	1 kit	M0877
4-Methylumbelliferyl β -D-Glucuronide (MUGlcU)	25 mg	M0240
Carboxyumbelliferyl β -D-Glucuronide (CUGlcU)	10 mg	M0256
Fluorescein di- β -D-Glucuronide, di-methyl ester	5 mg	M0969

Inhibition of cloned GUS activity can be afforded using:

Phenethyl-1-thio- β -D-Glucopyranosiduronic Acid (PETGU)	10 mg	M0261
---	-------	-------

Methods

Overview

Description

This high copy number plant vector, pNosdcGUS, expresses the GUS gene under the control of the nopaline synthase (pNos) promoter. This vector may be used to monitor virus production and transfection efficiency. GUS expression can be easily detected by addition of X-Gluc to the medium or by using our MarkerGene™ β -Glucuronidase (GUS) Reporter Gene Activity Detection Kit (M0877).

- pNosdcGUS expression vector also contains the kanamycin gene, which acts as a selection marker (@ 75 μ g/mL kanamycin resistance) in *Agrobacterium* host.
- The GUS gene can be excised using the 5' *Xba*I and 3' *Hind*III sites to allow the insertion of other genes to be expressed under the same regulatory elements in plants.
- For a map of pNosdcGUS, see page 10.

The pNosdcGUS Plant GUS Expression Vector System

The pNosdcGUS vector is a cloning vector that takes advantage of the site-specific recombination properties of bacteriophage lambda (Landy, 1989) to provide a rapid and highly efficient way to move your gene of interest into multiple vector systems. To express your gene of interest simply:

1. Clone your gene of interest into an entry vector to create an entry clone.
2. Perform an LR recombination reaction between the entry clone and a destination vector (e.g. pCMV β , pGPTV-KAN, etc.) to generate the expression vector.
3. Transfect your expression clone into the cell line of choice for stable expression of your gene of interest.

Using pNosdcGUS Plant GUS Expression Vector

The pNosdcGUS plant GUS expression vector is supplied as a supercoiled plasmid. Although Marker Gene has previously recommended using a linearized destination vector for more efficient recombination, further testing has found that linearization of this vector is NOT required to obtain optimal results for a downstream application.

Propagating pNosdcGUS :

If you wish to propagate and maintain pNosdcGUS, we recommend using *Agrobacterium tumefaciens* for transformation.

Entry Clone:

To recombine your gene into pNosdcGUS, you should have an entry clone containing your gene of interest.

Points to consider before recombining:

Your insert should contain a Kozak consensus sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak, 1991; Kozak, 1990). Other sequences are also possible, but the G or A at position -3 and the G at position +4 are the most critical for functional expression. If you wish to include the V5 epitope and 6xHis tag, your gene in the entry clone should not contain a stop codon. The gene should also be designed to be in frame with the C-terminal epitope tag after recombination. If you do NOT wish to include the V5 epitope and 6xHis tag, please be sure that your gene contains a stop codon in the entry clone.

Transfection

Introduction: This section provides general information for transfecting your expression clone into the plant of choice. We recommend that you include a positive control vector and a mock transfection (negative control) in your experiments to evaluate your results.

Methods of Transformation and Transfection:

Preparation of competent cells, transformation and recovery and floral dip with *Agrobacterium tumefaciens* protocols can be found in Weigel and Glazebrook's "Arabidopsis: A Laboratory Manual" (2002). For more information contact our Technical Assistance Staff (www.markergene.com or techservice@markergene.com).

Preparation of competent *Agrobacterium tumefaciens* cells:

1. Inoculate 250 mL of LB medium with 750 μ L *Agrobacterium tumefaciens* (for example GV3101::pMP90).
2. Incubate at 28°C (Dubenoff shaking water bath) with vigorous shaking for 11 hours until the OD is approximately 0.75.
3. Pellet the cells by centrifugation (5000 rpm).
4. Wash the pellet with 2 mL sterile TE.
5. Centrifuge at room temperature for 10 min at 5000 rpm.
6. Resuspend the pellet in 20 mL LB medium. Aliquot 250 μ L of this suspension into 1.5 mL Eppendorf microcentrifuge tubes.

Transformation and recovery of *Agrobacterium tumefaciens* cells:

1. Add 20 μ L of the **M1390** pNosdcGUS Plant GUS Expression Vector (5 μ g) to the competent GV3101::pMP90 *Agrobacterium tumefaciens* sample (250 μ L) in LB medium from above.
2. Incubate on ice (0°C) for 5 minutes.
3. Incubate in liquid nitrogen (-80°C) for 5 min.
4. Incubate at 37°C (water bath) for 5 min.
5. Add LB medium (1 mL) to each vial and incubate at room temperature (with rotation) for 4 hours.
6. Streak the bacteria onto LB Agar plates containing Kanamycin (50 μ g/mL) and incubate at 28°C for 3 days.
7. Pick one colony and re-streak onto an LB Agar plate containing Kanamycin (50 μ g/mL) and incubate at 28°C for 2 days.

Floral dip with *Agrobacterium tumefaciens*:

1. Prepare Infiltration Medium (250 mL) (see recipe below in Appendix)

2. Prepare MS Agar Gel for plant growth (see recipe below in Appendix)
3. Pipette 15 – 20 mL into 25 x 250 mm sterile test tubes for each plant to be grown after dip.
4. Grow up one colony of M1390 transformed bacteria (See step 7 from above) in LB media contain Kanamycin (50 ug/mL) (2 mL) overnight at 28°C (Dubenoff shaking water bath) with vigorous shaking.
5. Inoculate 225 mL of LB medium with 1.25 mL of preculture from step 4 above. Incubate overnight at 28°C (Dubenoff shaking water bath) with vigorous shaking.
6. Pellet the bacteria by centrifugation at 6000 rpm at 25°C for 10 min.
7. Resuspend the bacteria in Infiltration Medium (50 mL) and swirl to ensure complete mixture of the bacteria.
8. Transfer the bacterial suspension to a sterile dipping box (the lid of a pipette tip box works well). This process is best done in a laminar flow hood under sterile conditions.
9. Dip the roots of plants in the bacterial suspension for 30 seconds.
10. After dipping, place the plants on agar plates (MS medium) and seal the lids with parafilm. Allow to set overnight at room temperature.
11. After overnight storage, rinse the plant roots with sterile water and transfer the transformed plants to sterile tubes containing the MS Agar Gel for growth.
12. Grow the plants under long-day light conditions (16 hour daylight, 8 hour darkness) for several weeks until flowering and seed pod production occurs.
13. Remove seed pods and dry seeds.
14. Replant seeds on MS Agar Gel plates and collect seeds as above. Regrowth of these plants will provide stably transformed GUS plants and seeds for use.

Positive Control:

We recommend the use of a positive control vector for plant transfection and expression which may be used to optimize recombinant protein expression levels in your particular plant. A vector that allows expression of a C-terminally tagged β -glucuronidase fusion protein that may be detected by Western blot or functional assay provides the easiest way to measure protein expression levels. Consult our technical assistance for more information about C-terminal fusion protein expression systems.

To propagate and maintain the plasmid:

1. Resuspend the vector in 20 μ l sterile water to prepare a 1 μ g/ μ L stock solution and store at -20°C. Use the stock solution to transform *Agrobacterium tumefaciens*, or equivalent vehicle.
2. Select transformants on LB agar plates containing 50-100 μ g/ml kanamycin.
3. Prepare a glycerol stock of a transformant containing plasmid for long-term storage.

Expression and Analysis

Expression of your gene of interest from the expression clone can be performed in transiently transfected cells or stable cell lines.

GUS Assay using Cell Lysates

Cell lysate from above step can be assayed for GUS activity using the MarkerGene™ β -Glucuronidase (GUS) Reporter Gene Activity Detection Kit, M0877. For more detailed description of the MarkerGene™ Live Cell GUS Assay Kit please visit our website: <http://www.markergene.com/kits.php>

Detecting Recombinant Fusion Proteins:

To detect expression of your recombinant fusion protein by Western blot analysis, you may use the Anti-V5 antibodies or the Anti-His(C-term) antibodies available from Invitrogen or Amersham Biosciences or an antibody to your protein of interest.

Assay for β -glucuronidase:

If you use a positive control vector, you may assay for β -glucuronidase expression by Western blot analysis or activity assay using cell-free lysates (Miller, 1972). X-GlcU staining, or fluorescence detection are common methods of analysis. Marker Gene offers several reagents and kits for fast and easy detection of β -glucuronidase expression. See Accessory Products (page 3) for more information about these products.

Purification of Recombinant Fusion Proteins:

The presence of the C-terminal polyhistidine (6xHis) tag in your recombinant fusion protein allows for purification using a metal-chelating resin (available from Invitrogen). Note: Other purification methods may also be suitable including affinity chromatography.

Creating Stable Cell Lines:

The neomycin resistance gene can be cloned into the pNosdcGUS vector to allow selection of stable cell lines. If you wish to create stable cell lines, transfect your construct into the plant cell line of choice and select for foci using Geneticin®. General guidelines are provided below.

To obtain stable transfectants, we recommend that you linearize your pNosdcGUS construct before transfection. While linearizing the vector may not improve the efficiency of transfection, it increases the chances that the vector does not integrate in a way that disrupts elements necessary for expression in plant cells. To linearize your construct, cut at a unique site that is neither located within a critical element nor within your gene of interest.

Geneticin® (G418) blocks protein synthesis in plant cells by interfering with ribosomal function. It is an aminoglycoside, similar in structure to neomycin, gentamycin, and kanamycin. Expression in plant cells of the bacterial aminoglycoside phosphotransferase gene (APH), derived from *Tn5*, results in detoxification of Geneticin® (Southern and Berg, 1982).

Geneticin® Selection Guidelines:

Geneticin® is available from GIBCO (Catalog no. 11811-023). Use as follows:

1. Prepare Geneticin® in a buffered solution (e.g. 100 mM HEPES, pH 7.3).
2. Use 100 to 1000 μ g/ml of Geneticin® in complete medium.
3. Calculate concentration based on the amount of active drug.
4. Test varying concentrations of Geneticin® on your cell line to determine the concentration that kills your cells (kill curve). Cells differ in their susceptibility to Geneticin®. Cells will divide once or twice in the presence of lethal doses of Geneticin® Selective Antibiotic, so the effects of the drug

takes several days to become apparent. Complete selection for positive clones of cells can take up to 2 to 3 weeks of growth in selection medium.

Appendix: Recipes

LB (Luria-Bertani) Medium and Plates Composition:

1.0% Tryptone

0.5% Yeast Extract

1.0% NaCl

pH 7.0

1. For 1 liter, dissolve 10 g tryptone, 5 g yeast extract, and 10 g NaCl in 950 ml deionized water.
2. Adjust the pH of the solution to 7.0 with NaOH and bring the volume up to 1 liter.
3. Autoclave on liquid cycle for 20 minutes at 15 psi. Allow solution to cool to 55°C and add antibiotic if needed. Store at room temperature or at +4°C.

LB agar plates:

1. Prepare LB medium as above, but add 15 g/L agar before autoclaving.
2. Autoclave on liquid cycle for 20 minutes at 15 psi.
3. After autoclaving, cool to ~55°C, add antibiotic if needed, and pour into 10 cm plates.
4. Let harden, then invert and store at +4°C.

4X SDS-PAGE Sample Buffer:

1. Combine the following reagents:

0.5 M Tris-HCl, pH 6.8, 5 ml

Glycerol (100%), 4 ml

β -mercaptoethanol, 0.8 ml

Bromophenol Blue, 0.04 g

SDS, 0.8 g

2. Bring the volume to 10 ml with sterile water.
3. Aliquot and freeze at -20°C until needed.

Infiltration Medium:

$\frac{1}{2}$ MS salts (MS salts: 1650.00 mg/L ammonium nitrate, 332.02 mg/L calcium chloride anhydrous, 180.70 mg/L magnesium sulfate anhydrous, 1900.00 mg/L potassium nitrate, 170.00 mg/L potassium phosphate monobasic, 6.20 mg/L boric acid, 0.025 mg/L cobalt chloride·6H₂O, 0.025 mg/L cupric sulfate·5H₂O, 37.26 mg/L Na² EDTA, 16.90 mg/L manganese sulfate·H₂O, 0.250 mg/L molybdc acid sodium salt, 0.83 mg/L potassium iodide, 27.80 mg/L ferrous sulfate·7H₂O, 8.60 mg/L zinc sulfate·7H₂O) (Sigma, St. Louis, M0654)

1X Gamborg's B5 Vitamin (Sigma, St. Louis, G1019)

5% Sucrose (w/v)

0.044 μ M 6-benzylaminopurine (BAP) (stock solution 1mg/mL DMSO)

0.05% Silwet L77 (Lehle Seeds)


1. For 250mL Infiltration media, mixed 12.5mL 10X MS salts, 0.25mL Gamborg's B5 Vitamin 1000X, 12.5mg Sucrose, 2.5 μ L 1mg/mL stock solution of Benzylaminopurine (BAP), 12.5 μ L Silwet L77, 224.73mL dH₂O. Store at +4°C.

MS Agar Plant Growth Medium:

1. Autoclave 10.6g Murashige & Skoog (MS) Medium and 250 mL dH₂O on the liquid cycle for 20 minutes at 15 psi.
2. After autoclaving, cool to ~55°C, adjust to pH 5.7 with 1N NaOH.
3. Autoclave solution, again, on the liquid cycle for 20 minutes at 15 psi.
4. After autoclaving, cool to ~55°C, and pour into 10 cm plates.
5. Let harden, store at 25°C.

Map and Features of pNosdcGUS:

The map below shows the elements of pNosdcGUS. The GUS gene can be excised using the 5' *XhoI* and 3' *NotI* sites to allow the insertion of other genes to be expressed under the same regulatory elements in plant cells.



Restriction Digests of pNosdcGUS give the following sizes (kb):

EcoRI 13.6

BglII 8.2, 3.2, 2.7

HindIII 13.6

Technical Service

Please visit our website at www.markergene.com
Or write to techservice@markergene.com

Marker Gene is committed to providing our customers with premium-quality goods and services. We want to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about our products or service, please contact us at:

Marker Gene Technologies, Inc.
1850 Millrace Dr.
Eugene, OR 97403-1992
Toll Free: 1-888-218-4062
Local: 1-541-342-3760
Fax: 1-541-342-1960

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). No license is conveyed under the foregoing patents to use this product with any recombination sites other than those purchased from Marker Gene Technologies or its authorized distributor. The buyer cannot modify the recombination sequence(s) contained in this product for any purpose. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for the commercial purposes of the buyer, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for commercial purposes. Commercial purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. The user warrants to Marker Gene Technologies, Inc. that it will not use the purchase products or any part thereof for any purpose that is in contravention of any applicable law or regulation or in a manner that will infringe the copyright, trademarks, service marks or other intellectual property rights of third parties or violate the privacy, publicity or other personal rights of others or in any defamatory, obscene, threatening, abusive or hateful manner. The user shall indemnify Marker Gene Technologies, Inc., its employees and agents, harmless against all claims, liability, losses, damages and expenses including, without limitation, legal fees and costs arising out of or incurred as a result of any claims made, or litigation brought, against Marker Gene Technologies, Inc., its employees and agents, as a result of the use by the user of the purchased products or any part thereof for whatever purpose, to the fullest extent permitted by law.

References

- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). Current Protocols in Molecular Biology (New York: Greene Publishing Associates and Wiley-Interscience).
- Becker D, Kemper E, Schell J, Masterson R. (1992) "New plant binary vectors with selectable markers located proximal to the left T-DNA border." *Plant Mol. Biol.* 20: 1195-1197.
- Chen, C., and Okayama, H. (1987). "High-Efficiency Transformation of Plant Cells by Plasmid DNA." *Molec. Cell. Biol.* 7: 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). "Electroporation for the Efficient Transfection of Plant Cells with DNA." *Nucleic Acids Res.* 15: 1311-1326.
- de Wet JR, Wood KV, Helinski DR, DeLuca M (1985) "Cloning of firefly GUS cDNA and the expression of active GUS in Escherichia coli." *Proc Natl Acad Sci U S A.* 82(23):7870-7873.
- Felgner, P. L., Holm, M., and Chan, H. (1989). "Cationic Liposome Mediated Transfection." *Proc. West. Pharmacol. Soc.* 32: 115-121.
- Felgner, P. L. a., and Ringold, G. M. (1989). "Cationic Liposome-Mediated Transfection." *Nature* 337: 387-388.
- Kajiyama, N, Nakano E, (1991) "Isolation and characterization of mutants of firefly GUS which produce different colors of light." *Protein Eng.* 4: 691.
- Kajiyama N, Nakano E, (1993) "Thermostabilization of Firefly GUS by a Single Amino Acid Substitution at Position 217." *Biochemistry* 32: 13795-1 3799.
- Landy, A. (1989). 'Dynamic, Structural, and Regulatory Aspects of Lambda Site-specific Recombination.' *Annu. Rev. Biochem.* 58: 913-949.
- MacGregor, G.R. & Caskey, C.T., (1989). "Construction of Plasmids that Express E. Coli Beta-galactosidase in Plant Cells." *Nucleic Acids Res.* 17: 2365.
- Mamaev SV, Laikhter AL, Arslan T, Hecht SM, (1996) "Firefly GUS: Alteration of the Color of Emitted Light Resulting from Substitutions at Position 286."
- Masuda T, Tatsumi H, Nakano E (1989) "Cloning and sequence analysis of cDNA for GUS of a Japanese firefly, *Luciola cruciata*." *Gene.* 77(2):265-270.
- Miller, J. H. (1972). Experiments in Molecular Genetics (Cold Spring Harbor, New York: Cold Spring Harbor Laboratory).
- Southern, P. J., and Berg, P. (1982). "Transformation of Plant Cells to Antibiotic Resistance with a Bacterial Gene Under Control of the SV40 Early Region Promoter." *J. Molec. Appl. Gen.* 1: 327-339.
- Tatsumi H, Masuda T, Kajiyama N, Nakano E (1989) "GUS cDNA from Japanese firefly, *Luciola cruciata*: cloning, structure and expression in *Escherichia coli*." *J. Biolumin. Chemilumin.* 3(2):75-78.