

## Creatine Kinase-MB Procedure No. FT702

For the Quantitative Determination  
of Serum Creatine Kinase-MB (CK-MB)

### Summary and Principle

Creatine Kinase are dimeric molecules composed of M and B subunits and exist as the isoenzymes MM, MB, and BB. The subunits M and B are immunologically distinct. CK-MM and CK-MB are distributed primarily in the skeletal muscle and heart muscle, respectively, while CK-BB is present mainly in the brain and in tissues composed of smooth muscle.<sup>2</sup> Following acute myocardial infarction, CK-MB activity increases significantly and this elevation is highly specific for the laboratory diagnosis of myocardial infarction.<sup>3,4</sup> Although total CK activity usually increases following myocardial infarction, in some patients only the CK-MB activity increases, while the total CK remains in the normal range.<sup>5</sup> In this procedure CK activity is measured in the presence of an antibody to CK-M monomer. This antibody completely inhibits the activity of CKMM and half of the activity of CK-MB, while not affecting the B subunit activity of CK-MB and CK-BB. Due to negligible concentrations of CK-BB in the circulation, the remaining activity, multiplied by a factor of 2, represents the activity of the CK-MB isoenzyme.

### Reagents

**CK-MB Reagent (powder), Cat. No. FT702a:** 20vials for 6.5 ml each

ADP	2.0 mmol/L
G6PDH	2000 U/L
Creatine Phosphate	20 mmol/L
NAD	2.0 mmol/L
Hexokinase (yeast)	3000 U/L
D-Glucose	20 mmol/L

Anti-human CK-M Antibody (Goat) - Sufficient amount to inhibit up to 1500 U/L of CK-MM at 37° C. Also contains buffers, activators, surfactants and AK inhibitors.

**Precautions:** For in Vitro Diagnostic Use

Contains sodium azide. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Reagent Preparation: Reconstitute reagent with distilled water as indicated on vial label. Invert gently and swirl to dissolve completely.

Reagent Storage and Stability: CK-MB reagent is stable at refrigerated temperature (2-8° C) until the expiration date on the label. After reconstitution, the reagent is stable for 8 hours at room temperature or 30 days when stored at 2-8° C.

### Materials Required But Not Provided

Spectrophotometer capable of absorbance readings at 340 nm.  
Accurate pipetting devices.

Interval timer      Cuvets      Test Tubes  
Constant temperature bath, or block, 37° C or temperature controlled cuvet well.  
Mixer (Vortex type)

### Specimen Collection and Preparation

Clear unhemolyzed serum is the recommended specimen sample. No special additives or preservatives are required.

**Sample Stability:** Serum CK appears stable for 3 days at 2-8° C. It is recommended that specimens be assayed immediately after collection. Reconstituted control sera may show a decrease in CK activity in only a few hours.

**Interfering Substances:** Extremely hemolyzed samples are not suitable for testing since they may contain high levels of adenylate kinase, ATP and glucose-6-phosphate which interfere with the assay to yield false results. Drugs and other substances which may interfere with the determination of CK activity have been listed by Young et al.<sup>6</sup> This procedure may overestimate CK-MB values if CK-BB activity in the serum is high. CK-BB activity is usually absent in sera from normal individuals and patients with myocardial infarction.<sup>7</sup> The presence of a macro form of CK-BB in the specimen should be suspected if the CK-MB activity measured by this procedure represents more than 20% of the total CK activity.

### Automated Analyzer

**Parameters:**

Wavelength	340 nm
Reaction Type	Kinetic
Reaction Direction	Increasing
Reaction Temperature	37° C
Sample/Reagent Ratio	1:20
Equilibration Time	20 Seconds
Read Time	120 Seconds
Lag Time	300 Seconds
Blank Absorbance Limit	0.800A
High Absorbance Change/Min	0.375A/Min.
Factor	3376
Low Normal	0 IU/L
High Normal	24 IU/L
Linearity	1500 IU/L

Above parameters should be employed in programming automated analyzers for CK-MB. Consult your instrument manual for programming instructions. Specific programming applications for most automated analyzers are available from Interchim Customer Service Department.

### Manual Procedure

- For each sample add 1.0 mL CK-MB reagent into a cuvet/test tube and warm for approximately 5 minutes at 37° C.
- Add 50 uL of Sample to its respective tube, mix well and incubate for 5 minutes at 37° C.
- Set the wavelength of the instrument at 340 nm. Zero the instrument with distilled water.
- After the 5 minutes has elapsed, read and record the absorbance.
- Record the increase in absorbance at 60 second intervals for the next 2 minutes. (ΔA/Min.). The rate of change should be constant.

**NOTE:** If cuvet/tube is not temperature controlled, incubate samples at 37° C between readings.

**Quality Control:** Two levels of control material with known CK/CK-MB levels determined by this method should be analyzed each day of testing. Interchim recommends the use of its Cardiac Isoenzyme Control Set; No. 0985, (Level I & II) for this purpose.

### Results

A) **Total CK Activity:** Determine Total CK activity in serum according to the directions provided in the package insert for the CK reagent, Cat. No. FT7000.

B) **CK-B Activity:** Values are derived based on the "absorptivity micromolar extinction coefficient" of NADH at 340 nm (0.00622). A unit per liter (U/L) of CK-B activity is that amount of enzyme which oxidizes one umol/L of NADH per minute.

$$\text{CK-B activity (U/L)} = \frac{\Delta A/\text{Min.}}{\text{Absorptivity}} \times \frac{\text{Total Volume}}{\text{Sample Volume}}$$

$$\text{CK-B activity (U/L)} = \frac{\Delta A/\text{Min.}}{0.00622} \times \frac{1.050}{0.050}$$

$$\text{CK-B activity (U/L)} = \Delta A/\text{Min.} \times 3376$$

$$\text{C) CK-MB Activity (U/L)} = \text{CK-B activity (U/L)} \times 2$$

$$\% \text{ CK-MB Activity} = \frac{\text{CK-MB activity (U/L)} \times 100}{\text{Total CK activity (U/L)}}$$

**NOTE:** Samples having CK-MB values greater than 1500 U/L are diluted 1:2 (1 + 1) with distilled water, the assay repeated and results multiplied by the dilution factor of 2.

### Temperature Conversion Factor

To convert CK-MB activity at 37° C to 30° C value, multiply the result by 0.60.

### Expected Values<sup>8</sup>

0 - 24 U/L (37° C)
0 - 14 U/L (30° C)
% CK-MB < 5.5 %

This range should serve only as a guideline. It is recommended that each laboratory establish its own range of expected values, since differences exist between instruments, laboratories and local populations.

### Performance Characteristics

**Precision:** Two pools were subjected to analysis over five consecutive days to determine total precision. Another two pools were subjected to triplicate analysis daily over a five day period to determine within run data.

The data follows:

Pool	Total Precision:		Within Run:	
	1	2	1	2
X	32	1228	X	34 132
SD	3.1	9.2	SD	2.8 9.9
CV	9.8%	7.4%	CV	8.2% 7.5%

**Sensitivity:** Based on an instrument resolution of A = 0.001, this procedure has a sensitivity of 4 U/L.

**Correlation:** Serum CK-MB determinations by the procedure described (y) and by a similar method (Sigma) (x) on 40 sera, showed a correlation coefficient (r) of 0.991 and a regression equation of y = 0.98x - 0.823.

**Linearity:** This procedure is linear to 1500 U/L. Samples having values greater than 1500 U/L should be diluted and reassayed.

### References

- Dawson, DM, et al., Biochem Biophys. Res. Comm 21: 346 (1965)
- Neumeir D: Tissue Specific Distribution of Creatine Kinase Isoenzyme, Lang, Editor, Springer Verlag, New York, 1981, pp 85-109
- Wagner, et al, Circulation: 47 263 (1973)
- Bais R., Crit. Rev. Clin. Lab Sci. 18: 291 (1982)
- D'Souza JP et al, Clin. Biochem. 11: 204 (1978)
- Young, D.S., et al.: Clin Chem 21: 10 (1975)
- Lang, H. et al, Clin Chem 28: 1439 (1982)
- Wu AHB, Bowers CN Jr: Evaluation and comparison of immunoinhibition and immunoprecipitation methods for differentiating MB from BB and macro forms of creatine kinase isoenzymes in patients and healthy individuals. Clin Chem: 2017, 1982

For any information (technical service call, ordering,...), please contact Interchim 213 av. JF Kennedy, BP 1140 - 03 103 Montlucon (France)  
Phone: +33 4 70 03 73 06 (hot line) ; Fax +33 4 70 03 82 60  
e-mail: [interbiotech@interchim.com](mailto:interbiotech@interchim.com) ; web : <http://www.interchim.com>