

Amplite™ Luminometric Peroxidase Assay Kit

Ordering Information	Storage Conditions	Instrument Platform
Product Number: 11559 (500 assays)	Keep in freezer Avoid exposure to light	Luminescence microplate readers

Introduction

Enhanced chemiluminescence is a common technique for a variety of detection assays in biology. A horseradish peroxidase enzyme (HRP) is tethered to the molecule of interest (usually through labeling an immunoglobulin that specifically recognizes the molecule). This enzyme complex catalyzes the conversion of the enhanced chemiluminescent substrate into a sensitized reagent in the vicinity of the molecule of interest. The further oxidation of the substrate by hydrogen peroxide produces an excited molecule which emits light.

This kit uses our Amplite™ luminometric HRP substrate to quantify peroxidase in solutions. It provides an optimized “mix and read” assay protocol. Our Amplite™ Luminometric Peroxidase Assay Kit can detect as low as 100 $\mu\text{U}/\text{mL}$ HRP. The assay can be performed in a convenient 96-well or 384-well microtiter-plate format. Its signal can be easily read by a luminescence microplate reader. The kit can be used for ELISAs, characterizing kinetics of enzyme reaction and high throughput screenings.

Kit Key Features

Broad Application:	Can be used for quantifying HRP activities in solutions and solid surfaces (e.g. ELISA)
Sensitive:	Detect as low as 100 $\mu\text{U}/\text{mL}$ HRP in solution.
Continuous:	Easily adapted to automation without a separation step.
Convenient:	Formulated to have minimal hands-on time. No wash is required.
Non-Radioactive:	No special requirements for waste treatment.

Kit Components

Components	Amount
Component A: Assay Buffer	1 bottle (25 mL)
Component B: H_2O_2	1 vial (3% stabilized solution, 200 μL)
Component C: Horseradish Peroxidase	1 vial (20 units)

Assay Protocol for One 96-Well Plate

Brief Summary

**Prepare HRP reaction mixture (50 μL) → Add HRP standards and/or test samples (50 μL) →
Incubate at room temperature for 30 minutes to 2 hours → Monitor luminescent intensity**

Note: Thaw all the kit components at room temperature before starting the experiment.

1. Prepare HRP reaction mixture:

Add 30 μL of 3% stabilized H_2O_2 solution (Component B) into 5 mL of Assay Buffer (Component A), and keep from light.

Note: The HRP reaction mixture is stable at room temperature for at least 8 hours without loose activity if kept from light.

Data Analysis

The luminescence in blank wells with PBS and 0.1% BSA is used as a control, and subtracted from the values for those wells with HRP reactions. A HRP standard curve is shown in Figure 1.




Figure 1. HRP dose response was measured with the Amplitude™ Luminometric Peroxidase Assay Kit in a black 384-well plate using a NOVOstar plate reader (BMG Labtech). As low as 150 μ U/mL peroxidase can be detected with 30 minutes incubation (n=3).

References

1. Didenko VV, Baskin DS. (2006) Horseradish peroxidase-driven fluorescent labeling of nanotubes with quantum dots. *Biotechniques*, 40, 295.
2. Almeida LE, Imasato H, Tabak M. (2006) Enzymatic oxidation of dipyrromethole in homogeneous and micellar solutions in the horseradish peroxidase-hydrogen peroxide system. *Biochim Biophys Acta*, 1760, 216.
3. Krieg R, Halhuber KJ. (2003) Recent advances in catalytic peroxidase histochemistry. *Cell Mol Biol (Noisy-le-grand)*, 49, 547.
4. Matsui T, Nakayama H, Yoshida K, Shinmyo A. (2003) Vesicular transport route of horseradish C1a peroxidase is regulated by N- and C-terminal propeptides in tobacco cells. *Appl Microbiol Biotechnol*, 62, 517.
5. Wu TP, Zheng L, Ruan KC. (1998) Effect of Calcium Ion on Conformation of Horseradish Peroxidase Isoenzyme C. *Sheng Wu Hua Xue Yu Sheng Wu Wu Li Xue Bao (Shanghai)*, 30, 510.
6. Greco O, Folkes LK, Wardman P, Tozer GM, Dachs GU. (2000) Development of a novel enzyme/prodrug combination for gene therapy of cancer: horseradish peroxidase/indole-3-acetic acid. *Cancer Gene Ther*, 7, 1414.

Warning: This kit is only sold to end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest. Chemical analysis of the kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at info@aatbio.com if you have any questions.