

Cell Meter™ Phosphatidylserine Apoptosis Assay Kit

Green Fluorescence, Optimized for Microplate Readers

Ordering Information	Storage Conditions	Instrument Platform
Product Number: 22791 (100 assays)	Keep at 4 °C and protect from light	Fluorescence microscopes Fluorescence microplate readers

Introduction

Our Cell Meter™ assay kits are a set of tools for monitoring cellular functions. There are a variety of parameters that can be used for monitoring cell viability. This particular kit is designed to monitor cell apoptosis by measuring the translocation of phosphatidylserine (PS). In apoptosis, PS is transferred to the outer leaflet of the plasma membrane. The appearance of phosphatidylserine on the cell surface is a universal indicator of the initial/intermediate stages of cell apoptosis and can be detected before morphological changes can be observed.

Compared to Annexin V, the Apopxin™ PS sensor used in this kit is reengineered to have a much higher affinity to PS with $K_d < 10$ nM. The PS sensor used in this kit has green fluorescence upon binding to membrane PS. The fluorescent sensor has Ex/Em = 490/520 nm, making assay readily run with FITC filter set. Due to its highly enhanced affinity to PS, this kit is more robust than other commercial Annexin V-based apoptosis kits that are only used with either microscope or flow cytometry platform. This kit is optimized for a fluorescence microplate reader besides the microscope platform. It has been used for HTS applications.

Kit Key Features

Non-Radioactive:	No special requirements for waste treatment.
Continuous:	Easily adapted to automation with minimal hands on time.
Convenient:	All essential assay components are included.
Optimized Performance:	Optimized for detecting the translocation of phosphatidylserine (PS).
Enhanced value:	Less expensive than the sum of individual components.

Kit Components

Components	Amount
Component A: Apopxin™ Green (100X Stock Solution)	1 vial (100 µL/vial)
Component B: Assay Buffer	10 mL

Materials and Instruments Required (but not provided)

- 96 or 384-well microplates: Tissue culture microplates with black wall and clear bottom are recommended.
- A fluorescence microplate reader: Capable of detecting excitation and emission at 490/525 nm with bottom read mode.
- A fluorescence microscope

Assay Protocol for One 96-Well Plate

Brief Summary

Prepare cells with test compounds (100 µL/well/96-well plate or 25 µL/well/384-well plate) → Add equal volume of Apopxin™ Green assay solution → Incubate at room temperature for 1 hour
→ Monitor fluorescence intensity at Ex/Em = 490/525 nm (bottom read mode)

1. Prepare cells:

- 1.1 **For adherent cells:** Plate cells overnight in growth medium at 20,000 cells/well/90 μ L for a 96-well plate or 5,000 cells/well/20 μ L for a 384-well plate.
- 1.2 **For non-adherent cells:** Centrifuge the cells from the culture medium and then suspend the cell pellet in culture medium at 80,000 to 200,000 cells/well/90 μ L for a 96-well poly-D lysine plate or 20,000 to 50,000 cells/well/20 μ L for a 384-well poly-D lysine plate. Centrifuge the plate at 800 rpm for 2 minutes with brake off prior to the experiments.

Note: Each cell line should be evaluated on an individual basis to determine the optimal cell density for apoptosis induction.

2. Prepare Apopxin™ Green assay working solution:

- 2.1 Warm Kit Component B at room temperature before use.
- 2.2 Add 10 μ L of Apopxin™ Green (Component A) into 1 mL of Assay Buffer (Component B), and mix them well.

Note: 100 μ L of Apopxin™ Green assay working solution is enough for one well. Prepare fresh before use.

3. Run apoptosis assay:

- 3.1 Treat cells with test compounds by adding 10 μ L/well (96-well plate) or 2.5 μ L/well (384-well plate) of 10X test compound stock solution into PBS or the desired buffer. For blank wells (medium without the cells), add the same amount of compound buffer.
- 3.2 Incubate the cell plate in a 5% CO₂, 37 °C incubator for a desired period of time (4-6 hours for Jurkat cells treated with camptothecin) to induce apoptosis.
- 3.3 Add 100 μ L/well (96-well plate) or 25 μ L/well (384-well plate) of 2X Apopxin™ Green assay working solution (from Step 2.2) into each well.
- 3.4 Incubate the cell plate at room temperature for at least 1 hour, protected from light.
- 3.5 Centrifuge cell plate (especially for the non-adherent cells) at 800 rpm for 2 minutes (brake off).
- 3.6 Monitor the fluorescence intensity at Ex/Em = 490/525 nm (cut off at 510 nm) by using a fluorescent microplate reader (bottom read mode) or using a fluorescent microscope (FITC channel).

Data Analysis

The fluorescence in blank wells with the growth medium is subtracted from the values for those wells with the cells. The background fluorescence of the blank wells may vary depending on the sources of the growth media or the microtiter plates.




Figure 1. Detection of Apopxin™ Green-PS binding activity in Jurkat cells. Jurkat cells were seeded on the same day at 200,000 cells/90 μ L/well in a Costar black wall/clear bottom 96-well plate. The cells were treated with different doses of camptothecin for 5 hours as indicated. The Apopxin™ Green assay solution (100 μ L/well) was added and incubated at room temperature for 1 hour. The fluorescence intensity was measured at Ex/Em = 490/525 nm with NOVOstar instrument (from BMG Labtech) using bottom read mode.




Figure 2. Images of Jurkat cells stained with the Cell Meter™ Phosphatidylserine Apoptosis Assay Kit in a Costar black wall/clear bottom 96-well plate. A: Untreated control cells. B: Cells treated with 20 μ M camptothecin for 5 hours.

References

1. Hanshaw RG, Lakshmi C, Lambert TN, Johnson JR, Smith BD. (2005) Fluorescent detection of apoptotic cells by using zinc coordination complexes with a selective affinity for membrane surfaces enriched with phosphatidylserine. *Chembiochem*, 6, 2214.
2. Koulov AV, Stucker KA, Lakshmi C, Robinson JP, Smith BD. (2003) Detection of apoptotic cells using a synthetic fluorescent sensor for membrane surfaces that contain phosphatidylserine. *Cell Death Differ*, 10, 1357.
3. Hall MP, Burson KK, Huestis WH. (1998) Interactions of a vesicular stomatitis virus G protein fragment with phosphatidylserine: NMR and fluorescence studies. *Biochim Biophys Acta*, 1415, 101.
4. Saurel O, Cezanne L, Milon A, Tocanne JF, Demange P. (1998) Influence of annexin V on the structure and dynamics of phosphatidylcholine/phosphatidylserine bilayers: a fluorescence and NMR study. *Biochemistry*, 37, 1403.
5. Hanada K, Pagano RE. (1995) A Chinese hamster ovary cell mutant defective in the nonendocytic uptake of fluorescent analogs of phosphatidylserine: isolation using a cytosol acidification protocol. *J Cell Biol*, 128, 793.
6. Mancheno JM, Gasset M, Lacadena J, Ramon F, Martinez del Pozo A, Onaderra M, Gavilanes JG. (1994) Kinetic study of the aggregation and lipid mixing produced by alpha-sarcin on phosphatidylglycerol and phosphatidylserine vesicles: stopped-flow light scattering and fluorescence energy transfer measurements. *Biophys J*, 67, 1117.
7. Arnhold J, Wiegel D, Hussler O, Arnold K. (1994) Quenching and dequenching of octadecyl Rhodamine B chloride fluorescence in Ca^{2+} -induced fusion of phosphatidylserine vesicles: effects of poly(ethylene glycol). *Biochim Biophys Acta*, 1191, 375.
8. Walter A, Siegel DP. (1993) Divalent cation-induced lipid mixing between phosphatidylserine liposomes studied by stopped-flow fluorescence measurements: effects of temperature, comparison of barium and calcium, and perturbation by DPX. *Biochemistry*, 32, 3271.
9. Kobayashi T, Arakawa Y. (1991) Transport of exogenous fluorescent phosphatidylserine analogue to the Golgi apparatus in cultured fibroblasts. *J Cell Biol*, 113, 235.
10. Wrzosek A, Famulski KS, Lehotsky J, Pikula S. (1989) Conformational changes of $(\text{Ca}^{2+}\text{-Mg}^{2+})\text{-ATPase}$ of erythrocyte plasma membrane caused by calmodulin and phosphatidylserine as revealed by circular dichroism and fluorescence studies. *Biochim Biophys Acta*, 986, 263.

Warning: This kit is only sold to end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest. Chemical analysis of the kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at info@aatbio.com if you have any questions.