

Soil Total RNA Purification Kit
Product # 27700

Product Insert

Norgen’s Soil Total RNA Purification Kit provides a convenient and rapid method to purify total RNA from small amounts of soil samples. All types of soil samples can be processed with this kit, including common soil samples and difficult soil samples with high humic acid content such as compost and manure. The kit removes all traces of humic acid using the provided Bead Tubes and a combination of chemical and physical homogenization and lysis. A simple and rapid spin column procedure is then used to further purify the RNA. The kit purifies all sizes of RNA, from large mRNA and ribosomal RNA down to microRNA and small interfering RNA. The protocol does not rely on the use of phenol or chloroform, thereby providing a user friendly procedure and allowing high-throughput analysis on the lab bench. The purified RNA is of the highest integrity, and can be used in a number of downstream applications including real time PCR and reverse transcription PCR for gene expression analysis.

Norgen’s Purification Technology

Purification is based on spin column chromatography using Norgen’s proprietary resin as the separation matrix. The process involves first adding the soil sample and Lysis Solution to a provided Bead Tube, and the tube is vortexed for 5 minutes in order to efficiently and rapidly homogenize the sample and extract the RNA. The sample is then centrifuged, and the supernatant is transferred to a RNase-free microcentrifuge tube. Binding Solution I and II are added sequentially and mixed by inversion, and the lysate is incubated for 5 minutes on ice to remove all humic acids. The lysate is then spun to pellet any protein or soil particles, the supernatant is collected, an equal volume of ethanol is added to the lysate and the solution is loaded onto a spin-column. Norgen’s resin binds RNA in a manner that depends on ionic concentrations. Thus only the RNA will bind to the column, while the contaminating proteins will be removed in the flowthrough or retained on the top of the resin. An optional on-column DNase treatment can be performed at this point also to remove any residual DNA. The bound RNA is then washed with the provided Wash Solution in order to remove any remaining impurities, and the purified total RNA is eluted with the Elution Solution. The purified RNA is of the highest integrity, and can be used in a number of downstream applications.

Specifications

Kit Specifications	
Suggested Soil Input (Clay, loam, sand, feces, compost)	500 mg
Type of Soil Processed	All types, including common soil, compost and manure
Maximum Column Binding Capacity	50 µg
Maximum Column Loading Volume	650 µL
Time to Complete 10 Purifications	30 minutes

Advantages

- No phenol or chloroform extractions
- Fast and easy processing (less than 30 minutes for purification)
- Isolate high quality total RNA from a variety of soil samples
- Process all types of soil, including common soil, compost and manure

Kit Components

Component	Product # 26500 (25 preps)
Lysis Solution	20 mL
Binding Solution I	2 mL
Binding Solution II	2 mL
Wash Solution	9 mL
Elution Solution	2.5 mL
Bead Tubes	25
Mini Spin Columns	25
Collection Tubes	25
Elution tubes (1.7 mL)	25
Product Insert	1

Storage Conditions and Product Stability

All solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for at least 1 year in their unopened containers.

Precautions and Disclaimers

This kit is designed for research purposes only. It is not intended for human or diagnostic use.

Ensure that a suitable lab coat, disposable gloves and protective goggles are worn when working with chemicals. For more information, please consult the appropriate Material Safety Data Sheets (MSDSs). These are available as convenient PDF files online at www.norgenbiotek.com.

Customer-Supplied Reagents and Equipment

You must have the following in order to use Norgen's Soil Total RNA Purification Kit:

- Benchtop microcentrifuge
- RNase-free microcentrifuge tubes
- Flat bed vortex or bead beater equipment
- 95-100% ethanol
- 70% ethanol

Working with RNA

RNases are very stable and robust enzymes that degrade RNA. Autoclaving solutions and glassware is not always sufficient to actively remove these enzymes. The first step when preparing to work with RNA is to create an RNase-free environment. The following precautions are recommended as your best defense against these enzymes.

- The RNA area should be located away from microbiological work stations
- Clean, disposable gloves should be worn at all times when handling reagents, samples, pipettes, disposable tubes, etc. It is recommended that gloves are changed frequently to avoid contamination
- There should be designated solutions, tips, tubes, lab coats, pipettes, etc. for RNA only
- All RNA solutions should be prepared using at least 0.05% DEPC-treated autoclaved water or molecular biology grade nuclease-free water
- Clean all surfaces with commercially available RNase decontamination solutions
- When working with purified RNA samples, ensure that they remain on ice during downstream applications

Flow Chart

Procedure for Purifying Soil Total RNA using Norgen's Soil Total RNA Purification Kit

Add soil sample and Lysis Solution to Bead Tube



Vortex for 5 minutes.
Centrifuge. Transfer lysate.



Add Binding Solution I and
Binding Solution II
Incubate for 5 minutes on ice.


SPIN



Transfer lysate



Add Ethanol




Bind to column

SPIN




Wash two times with
Wash Solution



SPIN



Elute RNA with
Elution Solution



SPIN



Purified Total RNA

Procedure

All centrifugation steps are carried out in a benchtop microcentrifuge. Various speeds are required for different steps, so please check your microcentrifuge specifications to ensure that it is capable of the proper speeds. All centrifugation steps are performed at room temperature. The correct rpm can be calculated using the formula:

$$\text{RPM} = \sqrt{\frac{\text{RCF}}{(1.118 \times 10^{-5}) (r)}}$$

where *RCF* = required gravitational acceleration (relative centrifugal force in units of g); *r* = radius of the rotor in cm; and *RPM* = the number of revolutions per minute required to achieve the necessary *g*-force.

Notes Prior to Use

- All centrifugation steps are carried out in a benchtop microcentrifuge at **14,000 x g (~14,000 RPM)** except where noted. All centrifugation steps are performed at room temperature.
- A variable speed centrifuge should be used for maximum kit performance. If a variable speed centrifuge is not available a fixed speed centrifuge can be used, however reduced yields may be observed.
- Ensure that all solutions are at room temperature prior to use.
- Prepare a working concentration of **Wash Solution** by adding 21 mL of 95 - 100 % ethanol (provided by the user) to the supplied bottle containing the concentrated **Wash Solution**. This will give a final volume of 30 mL. The label on the bottle has a box that may be checked to indicate that the ethanol has been added.
- The maximum suggested soil input is 500 mg.

1. Lysate Preparation

- a. Add 500 mg of soil sample (maximum input varies – please see Notes Prior to Use) to a Bead Tube and add 700 µL of **Lysis Solution**. Vortex briefly to mix soil and Lysis Solution.

Note: In case of a wet soil sample, transfer the sample to a clean 1.7 mL microcentrifuge tube and centrifugation for 30 seconds at **14000 x g (~14,000 RPM)**. Remove the water carefully using a pipette, and resuspend the soil pellet in 700 µL of **Lysis Solution**. Transfer the soil to a Bead Tube using a pipette. **Proceed to Step 1b.**

- b. Secure tube horizontally on a flat-bed vortex pad with tape, or secure the tube in any commercially available bead beater equipment (e.g. Scientific Industries' Disruptor Genie™). Vortex for 5 minutes at maximum speed or optimize the conditions depending on the machine.
- c. Centrifuge the tube for 1 minute at **14000 x g (~14,000 RPM)**.
- d. Transfer up to 500 µL of supernatant to an RNAase-free microcentrifuge tube (not provided).
- e. Add 50 µL of **Binding Solution I** and add 50 µL **Binding Solution II** sequentially. Mix by inverting the tube a few times.
- f. Incubate for 5 minutes on ice.
- g. Spin the lysate for 1 minute at **14000 x g (~14,000 RPM)** to pellet any protein and soil particles.
- h. Using a pipette, transfer up to 400 µL of supernatant into a RNAase-free microcentrifuge tube (not provided).

Note: Avoid any contact with the pellet when collecting the supernatant. Also, depending on the soil type, some residue may be present on top of the supernatant. It is important to avoid collection of this residue while collecting the supernatant.

- i. Add an equal volume of 70% ethanol (provided by the user) to the lysate collected above (100 μ L of ethanol is added to every 100 μ L of lysate). Vortex to mix. **Proceed to Step 2.**

2. Binding to Column

- a. Assemble a spin column with one of the provided collection tubes.
- b. Apply up to 600 μ L of the clarified lysate with ethanol onto the column and centrifuge for 1 minute at **14000 \times g (~14,000 RPM)**. Discard the flowthrough and reassemble the spin column with the collection tube.

Note: Ensure the entire lysate volume has passed through into the collection tube by inspecting the column. If the entire lysate volume has not passed, spin for an additional minute.

- c. Depending on your lysate volume, repeat step **2b** if necessary.

Optional Step:

Norgen's Soil Total RNA Purification Kit isolates total RNA with minimal amounts of genomic DNA contamination. However, an optional **On-Column DNA Removal Protocol** is provided in Appendix A for maximum removal of residual DNA that may affect sensitive downstream applications. It is recommended that Norgen's RNase-Free DNase I Kit (Product # 25710) be used for this step. This step should be performed at this point in the protocol

3. Column Wash

- a. Apply 500 μ L of **Wash Solution** to the column and centrifuge for 1 minute.

Note: Ensure the entire wash solution has passed through into the collection tube by inspecting the column. If the entire wash volume has not passed, spin for an additional minute.

- b. Discard the flowthrough and reassemble the spin column with its collection tube.
- c. Repeat **3a** and **3b**.
- d. Spin the column for 2 minutes in order to thoroughly dry the resin. Discard the collection tube.

4. RNA Elution

- a. Place the column into a fresh 1.7 mL Elution tube provided with the kit.
- b. Add 50 μ L of **Elution Buffer** to the column.
- c. Centrifuge for 2 minutes at **200 \times g (~2,000 RPM)**, followed by a 1 minute spin at **14,000 \times g (~14,000 RPM)**. Note the volume eluted from the column. If the entire volume has not been eluted, spin the column at **14,000 \times g (~14,000 RPM)** for 1 additional minute.

5. Storage of RNA

The purified RNA may be stored at $-20\text{ }^{\circ}\text{C}$ for a few days. It is recommended that samples be placed at $-70\text{ }^{\circ}\text{C}$ for long term storage.

Appendix A

Protocol for Optional On-Column DNA Removal

Norgen's Soil Total RNA Purification Kit isolates total RNA with minimal amounts of DNA contamination. However, an optional protocol is provided below for maximum removal of residual DNA that may affect sensitive downstream applications. It is recommended that Norgen's RNase-Free DNase I Kit (Product # 25710) be used for this step.

1. For every on-column reaction to be performed, prepare a mix of $15\text{ }\mu\text{L}$ of **DNase I** and $100\text{ }\mu\text{L}$ of **Enzyme Incubation Buffer** using Norgen's RNase-Free DNase I Kit (Product # 25710). Mix gently by inverting the tube a few times. **DO NOT VORTEX**.

Note: If using an alternative DNase I, prepare a working stock of $0.25\text{ Kunitz unit}/\mu\text{L}$ RNase-free DNase I solution according to the manufacturer's instructions. A $100\text{ }\mu\text{L}$ aliquot is required for each column to be treated.

2. Perform the appropriate Total RNA Isolation Procedure for your starting material up to and including "**Binding to Column**"
3. Apply $400\text{ }\mu\text{L}$ of **Wash Solution** to the column and centrifuge for 2 minute. Discard the flowthrough. Reassemble the spin column with its collection tube.
4. Apply $100\text{ }\mu\text{L}$ of the RNase-free DNase I solution prepared in Step 1 to the column and centrifuge at $14,000\times g$ ($\sim 14,000\text{ RPM}$) for 1 minute.

Note: Ensure that the entire DNase I solution passes through the column. If needed, spin at $14,000\times g$ ($\sim 14,000\text{ RPM}$) for an additional minute.

5. After the centrifugation in Step 4, pipette the flowthrough that is present in the collection tube back onto the top of the column.

Note: Ensure Step 5 is performed in order to ensure maximum DNase activity and to obtain maximum yields of RNA, in particular for small RNA species.

6. Incubate the column assembly at $25 - 30^{\circ}\text{C}$ for 15 minutes.
7. Without any further centrifugation, proceed directly to the "**Column Wash**" procedure, (Step 3).

Troubleshooting Guide

Problem	Possible Cause	Solution and Explanation
Poor RNA Recovery	Homogenization was incomplete	Depending on the type of soil, further vortexing with the flat bed vortex or bead beater equipment may be required. However, it is not recommended to increase the vortex time to longer than 10 minutes at maximum speed.
	Lysis Additive was not added to the lysate	Ensure that the provided Lysis Additive is added to separate humic acid and increase DNA yield.
	Column has become clogged	Do not exceed the recommended input amount of 500mg mg soil. The amount of starting material may need to be decreased if the column shows clogging below the recommended levels. See also “Clogged Column” below.
	An alternative elution solution was used	It is recommended that the Elution Solution supplied with this kit be used for maximum RNA recovery.
	Ethanol was not added to the lysate	Ensure that the appropriate amount of ethanol is added to the lysate before binding to the column.
	Ethanol was not added to the Wash Solution	Ensure that 21 mL of 95% ethanol is added to the supplied Wash Solution prior to use.
RNA is Degraded	RNase contamination	RNases may be introduced during the use of the kit. Ensure proper procedures are followed when working with RNA. Please refer to “ <i>Working with RNA</i> ” at the beginning of this user guide.
	Improper storage of the purified RNA	For short term storage RNA samples may be stored at -20°C for a few days. It is recommended that samples be stored at -70°C for longer term storage.
	DNase used may not be RNase-free	Ensure that the DNase being used for the optional On-Column DNA Removal step RNase-free, in order to prevent possible problems with RNA degradation.

Problem	Possible Cause	Solution and Explanation
RNA does not perform well in downstream applications	RNA was not washed 2 times with the provided Wash Solution	Traces of salt from the binding step may remain in the sample if the column is not washed 2 times with Wash Solution. Salt may interfere with downstream applications, and thus must be washed from the column.
	Ethanol carryover	Ensure that the dry spin under the Column Wash procedure is performed, in order to remove traces of ethanol prior to elution. Ethanol is known to interfere with many downstream applications.

Related Products	Product #
RNase-Free DNase I Kit	25710
Soil DNA Isolation kit	26500
Plant/Fungi DNA Isolation kit	26200
Plant/Fungi RNA Isolation kit	25800
Plant RNA/DNA Purification Kit	26400
Water RNA/DNA Purification Kit	24400
Total RNA Purification Kit	17200
1kb RNA Ladder	15003
HighRanger 1kb DNA Ladder	11900
UltraRanger 1kb DNA Ladder	12100

Technical Support

Contact our Technical Support Team between the hours of 8:30 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362.

Technical support can also be obtained from our website (www.norgenbiotek.com) or through email at techsupport@norgenbiotek.com.



211 bis Avenue Kennedy - BP 1140
03103 Montluçon - France
33 (0) 4 70 03 88 55
Fax 33 (0) 4 70 03 82 60
e-mail interchim@interchim.com

Agence Paris - Normandie
33 (0) 1 41 32 34 40
Fax 33 (0) 1 47 91 23 90
e-mail interchim.paris@interchim.com