

TMB ELISA Stop solutions

Products Description

High quality Stop solutions to stop peroxidase/TMB reaction in microwell immunoassays.

The following Stop Reagents are propriety formulations intended for Microwell immunoassays , to stop Peroxydase/TMB substrate reactions in the blue or yellow spectral range before OpticalDensity reading:

•**650 nm Stop Reagents** stop TMB microwell substrate reactions without causing a color or absorbance change. The blue color is read in at 650 nm. It is available as a dry Blend, and a ready to use solution.

•**450 nm Stop Reagents** stop TMB microwell substrate reactions causing a color change from blue to yellow. The yellow color is read in at 450 nm. It is available as a ready to use solution, and a unique formulation (NovaStop):

Nova-Stop Solution for TMB MicrowellSubstrates is a patented acidic formulation used to stop TMB substrate reactions and change the chromogen from blue to yellow in color. Compared to sulfuric or hydrochloric acids, 450 nmLiquid Nova-Stop Solution produces a more stable signal. It is also non-corrosive to skin and eyes. U.S. Patent No.8,927,226

Product name:

Purchase number ,size

TMB Stop Solution, 650nm reading, Dry Blend Powder

S5168F BSTP-0100-01, 100ml *S5168G* BSTP-1000-01, 1L

TMB Stop Solution, 650nm reading, Liquid

S51681 LBSP-0100-01, pack for 100ml *S51682* LBSP-1000-01, pack for 1L

TMB Stop Solution, 450nm reading, Dry Blend Powder

S2959C STPR-0100-01, 100ml *S2959D* STPR -1000-01, 1L

TMB Stop Solution, 450nm reading, Liquid

S2959F LSTP-0100-01, 100ml *S2959G* LSTP-1000-01, 1L

TMB Stop Solution, 450nm reading, Liquid Nova-Stop

0C3440 NSTP-0100-01, 100ml *0C3441* NSTP-1000-01, 1L

Products stability: 3 years from the date of manufacture.

Storage: Storage recommendations are 2-25°C. Avoid exposure to moisture.

Specification: Stop HRP/TMB reaction and coloration:

Blue color is read at 620-650 nm.

Yellow color is read at 450 nm.

Directions for use

* **Stop TMB/peroxidase reactions:**

Add equal volumes of substrate and ready-to-use stop solution

(i.g. 100µl + 100µL stop solution per ELISA MiroWells)

* **For dry Blends, prepare ready-to-use solution (LBSP, STPR):**

1) **Add distilled or deionized water** to dissolve: 100 mL or 1L of water depending on pack size

Add 100 mL of distilled or deionized water to dissolve BSTP - 0100 - 01 .

Add 1 L of distilled or deionized water to dissolve BSTP - 1000 - 01.

2) Mix until dissolved.

This may take 5 minutes (BSTP) or 15-20min (STPR) to 15-20 minutes based on mixing rate.

3) The Stop solution is ready to stop the TMB reaction

Note: Excess solution may be stored on or few days in the fridge.

However it is recommended to prepare only require quantity.

FT-S51681

Related products & Categories

See the [BioSciences catalogue](#)

TMB Substrates for ELISA: see pages A147-148

Other immunoassay reagents: see pages A54-A193 including other HR Peroxidase substrates (colorimetric, fluorimetric, luminescent: pages A136-A139), coating solutions, saturating agents, buffers, diluants,...)

Ordering information

Please contact InterBioTech – Interchim for any other information

Hotline : +33(0)4 70 03 73 06 – Interbiotech@interchim.com

Catalog size quantities and prices may be found at <http://www.interchim.com>.

Please inquire for higher quantities (availability, shipment conditions).

Disclaimer : Materials are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. Interchim® is not liable for any damage resulting from handling or contact with this product.

Rev.R09E

