

FT-U47711



Fluorescent Cholesterol

Products Information

Name :

25-NBD Cholesterol

25-[N-[(7-nitro-2-1,3-benzoxadiazol-4-yl)methyl]amino]-27-norcholesterol

CAS [105539-27-3]

Catalog Number :


FP-U4771A, 1mg

Molecular Weight :

564.75,

Solubility:

Ethanol >99%pur



Name :

3-NBD-C6 Cholesterol

5-cholesten-3β-ol 6-[(7-nitro-2-1,3-benzoxadiazol-4-yl)amino]caproate

CAS [201731-19-3]

Catalog Number :


FP-1H2141, 1mg

Molecular Weight :

MW=662.90

Solubility:

Ethanol >99%pur



Name :

3-NBD-C12 Cholesterol

3-dodecanoyl-NBD Cholesterol; 3-C12-NBD Cholesterol conjugate

5-cholesten-3β-ol 12-[(7-nitro-2-1,3-benzoxadiazol-4-yl)amino]dodecanoate ; 12-[(7-nitro-2-1,3-benzoxadiazol-4-yl)amino]-(3S,10R,13R)-3-methoxy-10,13-dimethyl-17-((R)-6-methylheptan-2-yl)-tetradecahydro-1H-cyclopenta[α]phenanthrene

CAS [1246303-05-8]

Catalog Number :

FP-AYJ651, 1mg


Molecular Weight :

MW=747.06

Solubility:

Ethanol >99%pur

Exc/Em.= 464/512 to 497/551nm



Name :

Cholesterol – BrDIPY FL

5,5-difluoro-7-(3-((2-((3S,10R,13S,17R)-3-hydroxy-10,13-dimethyl-2,3,4,7,8,9,10,11,12,13,14,15,16,17-tetradecahydro-1H-cyclopenta[α]phenanthren-17-yl)propyl)amino)-3-oxopropyl)-1,3-dimethyl-5H-dipyrrolo[1,2-c:2',1'-f][1,3,2]diazaborinin-4-ium-5-uide

Catalog Number :

FP-05619A, 1mg


Molecular Weight :

MW=605.6

Solubility:

DMSO, DMF, Acetonitrile and Chloroform

Exc/Em.= 500/ 506, EC: 88 000



Storage:

-20°C

Protect from light and moisture

LIFE SCIENCES


interchim

211 bis, avenue JF Kennedy
BP 1140 - 03100 Montluçon
Fax +33 4 70 03 82 60


Hotline +33 4 70 03 73 06 • interbiotech@interchim.com

Technical information

25-NBD Cholesterol closely resembles the structure of native cholesterol. It incorporates into membranes and more effectively probes cholesterol containing domains, an obvious advantage over short-chain analogues.

NBD-C6-3-Cholesterol and **NBD-C12-3-Cholesterol** are fluorescently-tagged cholesterol with the hydrophilic NBD fluorophore attached to the hydrophilic end of cholesterol, separated by a 6- and 12-carbon spacer. This design allows the cholesterol to properly orient in membrane bilayers while the fluorescent tag is presented outside of the bilayer. This should model the behavior of cholesterol in membranes better than the previously-used 25-NBD cholesterol, which positions NBD directly on the 25th carbon of cholesterol at the hydrophobic terminus. Fluorescently-tagged lipids have been used to study their interactions with proteins, their utilization by cells and liposomes, and for the development of assays for lipid metabolism.

NBD fluorescent spectra



BDP fluorescence

The **BDP dye** is based on the 4,4-Difluoro-4-bora-3a,4a-diaza-s-indacene (or **DIPYrometheneBORon difluoride**) core. This dye fluoresce with max Exc/Em. = 500/ 506, is stable in physiological conditions, is relatively insensitive to variations in pH, and generally have high quantum yields and sharp emission peaks.

3-NBD Cholesterol technical info

3-dodecanoyl-NBD Cholesterol is supplied as a solution in ethanol. To change the solvent, simply evaporate the ethanol under a gentle stream of nitrogen and immediately add the solvent of choice. Solvents such as DMSO and dimethyl formamide (DMF) purged with an inert gas can be used.

The solubility of 3-dodecanoyl-NBD cholesterol in DMSO is approximately 0.25 mg/ml and approximately 1 mg/ml in DMF. If aqueous stock solutions are required for biological experiments, they can best be prepared by diluting the organic solvent into aqueous buffers or isotonic saline. We do not recommend storing the aqueous solution for more than one day.

Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations.

References

1. Kumagai, K., Yasuda, S., Okemoto, K., et al . CERT mediates intermembrane transfer of various molecular species of ceramides. *J. Biol. Chem.* 280(8) , 6488-6495 (2005).
2. Luo, M., Jones, S.M., Peters-Golden, M., et al . Nuclear localization of 5-lipoxygenase as a determinant of leukotriene B 4 synthetic capacity. *Proc. Natl. Acad. Sci. USA* 100(21) , 12165-12170 (2003).
3. Moreno, M.J., Estronca, L.M.B.B., and Vaz, W.L.C. Translocation of phospholipids and dithionite permeability in liquid-ordered and liquid-disordered membranes. *Biophys. J.* 91 , 873-881 (2006).
4. Loidl, A., Claus, R., Daigner, H.P., et al . High-precision fluorescence assay for sphingomyelinase activity of isolated enzymes and cell lysates. *J. Lipid Res.* 43 , 815-823 (2002).
5. Tani, M., Okino, N., Mitsutake, S., et al . Specific and sensitive assay for alkaline and neutral ceramidases involving C12-NBD-ceramide. *J. Biochem.* 125 , 746-749 (1999).

FT-U47711

Related products

ROS Probes, i.e.
Dihydrorhodamine 123, FP-83775A),

H2DCFDA (FP-467312)
Fluorescent Phalloidin (cytoskeleton probes) [AZ0330](#)

Ordering information

Catalog size quantities and prices may be found at <http://www.fluoprobes.com>
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold for research use only, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.

Rev.P03E

