



## Fluorescent dyes conjugated Tetrazine conjugates – for Click chemistry

### Products Description



CY<sub>anine</sub>3/5/5.5/7 dyes are classic fluorescent dyes used in a variety of detection techniques. They are provided activated by **Tetrazine** functional group to perform easy conjugation of biomolecules (proteins, peptides, amino-modified DNAs and oligos).

Tetrazines demonstrate exceptionally fast cycloaddition kinetics (up to 30 000 M<sup>-1</sup> s<sup>-1</sup>) with *trans*-cyclooctenes (TCO) as the dienophile, the fastest kinetics ever reported for any bioorthogonal reaction. In addition, inverse-Electron-Demand Diels-Alder reaction of tetrazines *trans*-cyclooctene forms a stable covalent bond and does not require Cu-catalyst or elevated temperatures. In applications such as in vivo cancer imaging or pre-targeted cell labeling studies where rapid reaction kinetics are desired, a faster hydrogen substituted CY<sub>anine</sub> Tetrazine probe would be a probe of choice.







Actually, tetrazines react not only with *trans*-cyclooctenes, but also with (methyl)cyclopropenes and strained non-benzo-fused cycloalkynes.

[Directions for use](#)
[Scientific and technical information](#)
[Related products](#)



- **Tetrazine** for conjugation to *trans*-cyclooctenes (TCO) modified biomolecules via Click Chemistry :

Fluorescent dye conjugate	Cat.number	MW	$\lambda_{Abs./Em.}$ (nm)	EC	QY	Store
<b>•TriSulfo-CY<sub>anine</sub>3-Tetrazine</b> 	FP-AW2NI0, 1mg	894.05	550/570 nm	150 000		(M)
Spectrally similar to AF555, A555, CF555, DL555 Soluble in Water, DMSO, DMF						
<b>•DiSulfo-CY<sub>anine</sub>3-Tetrazine</b> 	FP-WXS720, 0.5mg	799.96	550/570 nm	150 000		(J)
3-(p-Benzylamino)-1,2,4,5-tetrazine - Cy3						

FT-WXS720

Fluorescent dye conjugate	Cat.number	MW	$\lambda_{Abs./Em.}$ (nm)	EC	QY	Store
<b>•DiSulfoCY<sub>aniline</sub>3-Tetrazine, K salt</b> 	FP-0B8300, 0.5mg	838.05				(M)
	MW:838.05 () A/E: 548/563 nm; EC:162 000; QY:0.1; CF <sub>260</sub> : 0.03 ; CF <sub>280</sub> : 0.06					
<b>•CY<sub>aniline</sub>3-Tetrazine</b> 	FP-AW2ND0, 0.5mg	727.64				(M)
	MW:727.64 (611.4 M+ increment) A/E: 550/570 nm; EC:150 000; QY:0.31; CF <sub>260</sub> : 0.04 ; CF <sub>280</sub> : 0.09					
<b>•TriSulfo-CY<sub>aniline</sub>5-Tetrazine</b> 	FP-AS4K70, 1mg	919.27	649/670nm	250 000		(M)
	Spectrally similar to AF647, A647, CF647, DL649 Soluble in Water, DMSO, DMF					
<b>•DiSulfo-CY<sub>aniline</sub>5-Tetrazine</b> 	FP-0B1770, 0.5mg	826.00	649/670nm	250 000		(J)
	3-(p-Benzylamino)-1,2,4,5-tetrazine – Cy5					
<b>•DiSulfoCY<sub>aniline</sub>5-Tetrazine, K salt</b> 	FP-AW2NH0, 0.5mg	838.05				(M)
	MW:838.05 () A/E: 646/662 nm; EC:271 000; QY:0.28; CF <sub>260</sub> : 0.04 ; CF <sub>280</sub> : 0.04					
<b>•CY<sub>aniline</sub>5-Tetrazine</b> 	FP-AW2NE0, 51mg	811.84				(M)
	MW:811.84 (637.4 M+ increment) A/E: 646/662 nm; EC:250 000; QY:0.2; CF <sub>260</sub> : 0.03 ; CF <sub>280</sub> : 0.04					

FT-WXS720

Fluorescent dye conjugate	Cat.number	MW	$\lambda_{Abs./Em.}$ (nm)	EC	QY	Store
•Tri-Sulfo-CY <sub>anine</sub> 5.5-Tetrazine	Inquire					(M)
•DiSulfo-CY <sub>anine</sub> 5.5-Tetrazine	Inquire					(J)
•CY <sub>anine</sub> 5.5-Tetrazine	FP-AW2NF0, 5mg	911.96	684/710	209 000		(M)
 <p>MW:911.96 (734.4 M+ increment) A/E: 684/710 nm; EC:209 000; QY:0.2-; CF<sub>260</sub>: 0.07; CF<sub>280</sub>: 0.03</p>						
•Tri-Sulfo-CY <sub>anine</sub> 7-Tetrazine	Inquire					
•DiSulfo-CY <sub>anine</sub> 7-Tetrazine	Inquire					
•CY <sub>anine</sub> 7-Tetrazine	FP-AW2NG0, 5mg					
 <p>MW:877.94 (703.4 M+ increment) A/E: 750/773 nm; EC:199 000; QY:-; CF<sub>260</sub>: 0.022 ; CF<sub>280</sub>: 0.029 C47H54N7F6OP</p>						

Inquire for other size (1mg size(0), 25mg(2); bulk)

## Technical information

Tetrazine – CYanin products can not be stored in solution.

Dissolve the product in methanol, split into aliquots of sufficient amounts for one experiment and dry the product e.g. with a speedvac. Dried aliquots can be stored at -20°C and freshly dissolved immediately before use e.g. with DMSO.

### Information on fluorescent labels:

• **CY<sub>anine</sub>3** (Cy3) can replace Fluoresceins, A546 and DyLight 549&550 dyes. Alternative superior dye is [FluoProbes547H](#).

Cy3 is one of the most broadly used fluorophore which can be detected by various fluorometers, imagers, and microscopes with a wavelength range of 530-555 nm. Due to inherently high extinction coefficient, this dye is also easily detected by naked eye on gels, and in solution.

• **CY<sub>anine</sub>5** (Cy5) can replace TMR, A647, DyLight649&650, CF647, C645A, PF647. Alternative superior dye is [FluoProbes647H](#).

Cy5 is excited maximally at 650 nm to about 98% of maximum with a krypton/argon laser (647 nm line) or to about 63% of maximum with a helium/neon laser (633 nm line). It fluoresce maximally at 670 nm, that has a lower autofluorescence of biological specimens than shorter wavelengths. Cy5 can be used with a variety of other fluorophores for multiple labeling due to a wide separation of its emission from that of shorter wavelength-emitting fluorophores. However, emission cannot be seen well by eye, and Cy5 cannot be excited optimally with a mercury lamp. Therefore, this dye is not recommended for use with conventional epifluorescence microscopes. It is most

FT-WXS720

commonly visualized with a confocal microscope equipped with an appropriate laser for excitation and a far-red detector. Cy5 and FP547H is less expensive and equally bright alternative to Allophycocyanin conjugates for flow cytometry.

## Protocols

See the literature

## Related / associated products and documents

\*Other fluorescent reagents:

- Superior **FluoProbes** dyes, activated by [-Azide \(protocol\)](#), i.e. FP488-Azide [#YE4970](#)
- **Cyanine Dyes** (Cy3, Cy5, Cy5.5, Cy7...), activated by [NHS](#)(BB7493) , Maleimide(JO6660), Azide([EV0910](#)), Alkyne (A6320), Amine (CY3AM), Carboxyl(CY3CA), dCTP(ZD0210),...
- **Conventional dyes** activated by Azide, i.e. FAM-Azide [FP-EV0920](#), TAMRA-Azide,...

\*Other Click Chemistry reagents:

- DBCO reagents [#DQP580](#)
- Pentynoic acid activator [ZL5530](#) to modify amine containing biomolecules with an alkyne residue
- Copper(II)-TBTA complex [FY2780](#) to catalyse the conjugation reaction by click chemistry.
- Alkyne Amidite, 5'-terminal [ZL5500](#)
- Alkyne Amidite, hydroxyprolinol [ZL5510](#)

\*See [BioSciences Innovations catalogue](#) and [e-search tool](#).

## Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.

Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

**Disclaimer :** Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product.

Rev.S09E