

# Biochromatography - Affinity

## Introduction /Modified Agarose [Affarose™]


Interchim offer a range of affinity based products that combine a number of biological sample binding methods with a variety of housing formats. These formats allow the affinity supports to utilise a range of techniques to assist with the biological purification process i.e.

Affarose™ : Affinity agarose gels  
 UptiSpin™ : Affinity /Centrifugation  
 UptiSpin™ : Affinity /Vacuum  
 Ami R Gels : Amine specific affinity agarose  
 UptiBeads™ : Affinity /Magnetic

### Immobilized Metal Affinity Chromatography (IMAC)

Affarose™ : Affinity IMAC gels  
 UptiSpin™ : Affinity IMAC /Centrifugation  
 Cellufine™ : Cellulitic Affinity Purification/ Immobilizaton

### Protein A, G, L - Affarose™



Antibody purification was popularized by Protein A affinity chromatography. Interchim offers high-quality supports with Ig-binding proteins i.e. rProtein A, rProtein G, rProtein L supporting a range of requirements from benchtop R&D to production (ask for bulk sizes).

Support selection is dependant upon a number of parameters such as ligand specificity, stability, leaching and flow rate.

Protein A is particularly suited to the purification of monoclonal antibodies from FCS supplemented cell cultures and is often the affinity of choice for many other Ig species. Protein G is generally considered as an alternative when cross-reactivity to several species is desired, or when Protein A fails. Protein L is dedicated to specific applications.

Our recombinant grade Protein A, G & L provide unsurpassed binding capacity.

Affarose™ format : Ligands are covalently bound to the modified agarose matrix called Affarose™. Affarose™ ensures a minimal leaching of ligand and displays excellent stability, allowing for support re-use.

UptiSpin™ format : Designed for simple, rapid, multiple antibody purification from serum, ascites and tissue culture supernatant. Re-usable.

	Protein A (Staphylococcus)	Protein G (streptococcus)	Protein L (Peptostreptococcus)
(native) Origin	E. coli	E. coli	E. coli
Nb binding sites	(4) 4	(4) 2	4
Ig binding site	Fc	Fc	K light chain
Albumin binding	(-) -	(+) -	+
Optimal pH binding	8 - 9	5 - 6	7.5

	Protein A	Protein G	Protein L
Hu IgGs	+++	+++	+++
Hu IgG1	+++	+++	+++
Hu IgG2	+++	+++	+++
Hu IgG3	+	+++	+++
Hu IgG4	+++	+++	+++
Hu IgM	+	-	+++
Hu IgE	++	-	+++
Hu IgD	-	-	+++
Hu IgA	+	-	+++
Ms IgGs	+++	+++	+++
Ms IgG1	+	++	+++
Ms IgG2a	+++	+++	+++
Ms IgG2b	+++	+++	+++
Ms IgG2c	+++	+++	+++
Ms IgG3	+++	+++	+++
Ms IgM	-	-	++
Rt IgGs	+	++	+++
Rt IgG1	+	++	++
Rt IgG2a	+++	+++	+++
Rt IgG2b	+++	+++	+++
Rt IgG3	+++	+++	+++
Rb IgGs	+++	+++	+
Gp IgGs	+++	+	+/-
Hs IgGs	++		+/-
Dog IgGs	+++	+	
Cat IgGs	+++	+	
Mk IgGs	+++	+++	
Pig IgGs	+++	+	+++
Gt IgGs	+ / ++	+++	-
Sh IgGs	+ / ++	+++	+
Bv IgGs	+ / ++	+++	+ / -
Dk IgGs	++	+++	
Hs IgGs	+	+++	
Ck IgGs	-	-	

Bv : Bovine, Ck : Chicken, Ct : Cat, Dg : Dog, Gp : Guinea-pig, Gt : Goat, Hm : Hamster Hs : Horse  
 Hu : Human, Mk : Monkey, Ms : Mouse, Rt : Rat, Rb : Rabbit, Sh : Sheep

\* binding capacities for proteins L are given for K subclass only / +++ : strong binding / - : no binding

# Biochromatography - Affinity

## Modified Agarose [Affarose™]

### Affarose™ rProtein A, XTrem

This is our highest capacity affinity Protein A support. The protein A is immobilised onto a 6% cross-linked agarose. Supplied as a 50% suspension.

- Capacity : 30-40 mg of human IgG per ml of wet gel (static)
- Excellent stability and minimal leaching of ligand
- High flow rate, up to 400 cm /hour


#### Applications

- Purification of monoclonal antibodies, notably from FCS supplemented culture media (no Bovine IgG binding)
- Purification of polyclonal antibodies
- Immunoprecipitation
- IgG immobilization (affinity column preparation)
- Pre-adsorption, IgG depletion
- Isolation of immune complexes

Cycle	Capacity (% initial)	Protein A in eluate
1	100%	2.5 ng/ml*
2	100%	2.4 ng/ml*
16	100%	1.5 ng/ml*
32	>95%	6.7 ng/ml*
100	>85%	3.3 ng/ml*

\*representative data. - Eluate : polyclonal human IgG 7.4 - 8.5mg /ml

Product	P/N	Qty
Affarose™ rProtein A, Xtrem	UP904670	5 ml
Affarose™ rProtein A, Xtrem	UP904672	20 ml



# Biochromatography - Affinity

## Modified Agarose [Affarose™]

### Affarose™ rProtein A

Interchim's classical Protein A affinity support suitable for most applications. Protein A is immobilised on a 6% modified agarose and supplied as a 50% suspension.


- Minimal leaching (<3ng Protein /mL)
- Capacity : 18-22 mg of human IgG per ml of wet gel
- High flow rate - up to 300cm /hour

### Applications

- Purification of monoclonal antibodies, notably from FCS supplemented culture media (no Bovine IgGs binding)
- Purification of polyclonal antibodies
- Immunoprecipitation
- IgG immobilisation (affinity column preparation)
- Preadsorption, IgG depletion
- Isolation of immune complexes

Product	P/N	Qty
Affarose™ rProtein A	<b>UP49981G</b>	2 ml
Affarose™ rProtein A	<b>UP49981H</b>	5 ml

Uptima Protein A is produced by recombinant technology in E.coli, with culture media exempt of any animal proteins. Highly controlled procedures result in a material free from toxic bacterial contaminant (Staphylococcus endotoxins and hemolysin), found normally in native Protein A.



# Biochromatography - Affinity

## Modified Agarose [Affarose™]

### Affarose™ rProtein L


Interchim Protein L is immobilised on 4% modified agarose and supplied as a 50% suspension. It is particularly referenced as an excellent support for the purification of Ig fragments

- Capacity : 3-5 mg of human IgG per ml of wet gel
- Ig binding does not interfere with antigen binding site
- Excellent stability

#### Applications

- Ideal for the purification of F(ab')<sub>2</sub>, Fab and scFV fragments containing kappa light chains
- Purification of kappa monoclonal antibodies directly from media supplemented with FCS
- Purification of IgG, IgM, IgA, IgE, IgD containing kappa light chains from serum, ascites and hybridoma supernatant
- Antibodies sensitive to alkaline pH

Product	P/N	Qty
Affarose rProtein L	UP52746A	2 ml
	UP52746C	0.5 ml



3D structure of one sub-unit of rProtein L Drs T. Wan, B.J. Sutton et m. G. Gore, University of Southampton, GB.

# Biochromatography - Affinity

## Centrifuge Columns

### UptiSpin™ Protein A and UptiSpin™ Protein G

UptiSpin™ centrifugal columns are one of the simplest means to purify antibodies from serum, ascites and tissue culture supernatant. UptiSpin columns achieve high recovery of pure antibodies in minutes.

- Fast, convenient and easy to use
- Cost effective : less than 1 Euro to purify 1 mg
- Efficient : Proteins are sufficiently pure for further applications.

UptiSpin™ Protein A & G kits are designed for simple, rapid and multiple antibody purification from serum, ascites & tissue culture supernatant derived from static cultures & bioreactors. Working volume range for mini : 650 µl and midi : 20 ml

Product	P/N	Kit
Protein A UptiSpin™		
Mini		
Protein A Mini Sample Kit : 2 A Plugs (includes buffers)	UPBB9720	1
Protein A Mini Kit – 16 A Plugs (includes buffers)	UPBB9730	1
Protein A Mini Bulk Pack – 48 A plugs (does not include buffers)	UPBB9740	1
Midi		
Protein A Midi Kit – 4 Midi A Plugs (includes buffers)	UPBB9750	1
Protein A Midi Bulk Pack – 12 A Plugs (does not include buffers)	UPBB9760	1
Protein G UptiSpin™		
Mini		
Protein G Mini Sample Kit : 2 G Plug (includes buffers)	UPBB9770	1
Protein G Mini Kit – 16 G Plugs (includes buffers)	UPBA7770	1
Protein G Mini Bulk Pack – 48 G plugs (does not include buffers)	UPBB9780	1
Midi		
Protein G Midi Kit – 4 Midi G Plugs (includes buffers)	UPBB9790	1
Protein G Midi Bulk Pack – 12 G Plugs (does not include buffers)	UPBB9800	1
Starter pack A and G (Mini Only)		
Starter Pack: 2 A and 2 G Mini Plugs (does not include buffers)	UPBB9820	1

### Features

#### Key Features

Centrifugal spin column

Uses centrifugal g force driving force

Many samples can be centrifuged in parallel

Economical

Working range  
(1 -to- 100's of mg)

#### Benefits

Simple to use. Only a centrifuge required. No HPLC or format FPLC equipment has to be used.

Extremely fast separation times  
No hold up volume or spluttering.  
High recovery during fast elution step.

Multiple samples can be purified at the same time.

Low cost purification.  
(less than 1 Euro per mg IgG)\*  
(\*5 times re-use of Midi).

Variable volume and amounts can be loaded and recovered with ease.  
Typical amount of rabbit IgG purified :  
Midi kit : 400 mg (4 plugs)  
> 20 mg rabbit IgG /plug in one run ;  
> 10 mg rat IgG /plug /run).

Loading the plug into the spin column



Place the plug into the spin column with the recessed end uppermost



Push the plug FULLY into the tapered end of the spin column using the plug insertion tool.



It is now ready for pre-equilibration with binding buffer followed by centrifugation

# Biochromatography - Affinity

## Vacuum based affinity

### UptiSpin™ Protein A and G Vacuum kits

UptiSpin™ Protein A & G Vacuum kits are designed to concentrate and purify antibodies, at high speed, from large volumes of dilute cell culture supernatants (up to 2 L) using vacuum as the driving force.


#### Features

- Uses universal 50 ml collection tubes and standard glass laboratory bottles up to 2 l
- Works under low or high vacuum (capture efficiency is enhanced at lower vacuum)
- Simple to set up and easy to use
- Ultra-filtration concentrators to desalt and concentrate the purified monoclonal antibody
- Can be performed at room temperature or in cold room

#### Applications

- High speed, hands-free concentration and purification of large volumes of dilute cell culture supernatant without the need of a centrifuge.

Description	P/N	No. of purifications
UptiSpin™ Protein A Midi Vacuum Kit with Vacuum Head (4 columns, UF concentrators, Buffers)	UPBM8140	20
UptiSpin™ Protein G Midi Vacuum Kit with Vacuum Head (4 columns, UF concentrators, Buffers)	UPBM8150	20
UptiSpin™ Protein A Midi Vacuum Kit without Vacuum Head (4 columns, UF concentrators, Buffers)	UPBM8220	20
UptiSpin™ Protein G Midi Vacuum Kit without Vacuum Head (4 columns, UF concentrators, Buffers)	UPBM8240	20
UptiSpin™ Protein A Bulk Pack (12 columns only)	UPBM8260	60
UptiSpin™ Protein G Bulk Pack (12 columns only)	UPBM8270	60
UptiSpin™ Vacuum Accessory pack containing two complete sets of O-rings (6 x O-rings) plus 6 mm and 8 mm hose barb connectors	UPBM8280	1 unit
Protein A Buffer Pack (Buffers only)	UPBM8360	1 unit
Protein G Buffer Pack (Buffers only)	UPBM8370	1 unit
UptiSpin™ Vacuum Head and a sample loading sealing cap	UPBM8190	1 unit



# Biochromatography - Affinity

## Amine specific affinity

### Ami. R. Gel & Ami. Fast Flow Gel

Ami R and Ami fast flow affinity gels provide an excellent means to immobilize antibodies, proteins and peptides through available amines for affinity purifications.

These supports preserve protein activity compared to conventional coupling chemistries such as Glutaraldehyde, cyanogen bromide, epoxy, succinimide, etc... and provide an efficient and ease to use support for bench top R & D applications.

Efficient high binding capacities combined with ease of use, favour diagnostic applications. Ami.R. Gel and Ami.R. Gel Fast Flow are excellent alternatives to conventional activated gels in pharmaceutical processing thanks to their extreme stability and regeneration features.

#### Features

- High binding capacity (up to 60 mg/ml BSA)
- Higher activity of immobilized antibodies than with "oriented" coupling
- Amine specificity (couple on NH<sub>2</sub>-bearing molecules)
- High coupling efficiency (85-99%), irrespective of MW
- Coupling compatible at pH 3-10, with temperatures 0-40°C
- Fast coupling (from 20 min to 2 hr)
- No end-capping necessary
- Undetectable leaching of ligand (<0.1 ppm)
- High flow rate (up to 250 cm/h) (200cm/hr for Fast Flow gel)
- Sanitizable with NaOH, autoclavable

Product	P/N	Qty
Ami. R. Gel	UP564089	2 ml
	UP56408A	5 ml
	UP56408B	10 ml
Ami. R. Fast Flow Gel	UPR2289A	5 ml
	UPR2289B	10 ml

# Biochromatography - Affinity

## Magnetic based affinity

### UptiBeads™ Magnetic Beads

UptiBeads™ are magnetized iron core polystyrene microspheres for affinity /magnetic purification applications. The excellent rate of capture of our para-magnetic beads is due to the high diameter /ligand load ratio, and very low dispersal rates. UptiBeads™ are stable between pH 2 and 12, and are compatible with sonication applications. They are supplied ready-to-use in a 1% PBS suspension in 0.3 and 0.86 µm diameters.

The grafted ligands are high affinity anti IgG (all isotype), anti IgM, Streptavidin and anti Biotin. Other sizes and ligands are available upon request.

- High binding capacity
- Monodisperse
- Excellent colloidal stability
- Stable
- Excellent specific surface /magnetic capacity

Diameter (µm)	% ferrite	Density (g/cm <sup>3</sup> )	Surface area (cm <sup>2</sup> )	Number of beads (nb/ml)
0.3 µm	40%	1.54	3.0 x 10 <sup>-9</sup>	4.2 x 10 <sup>10</sup>
0.86 µm	55%	1.54	3.1 x 10 <sup>-8</sup>	1.6 x 10 <sup>10</sup>

Specificity	Diameter (µm)	Binding capacity (mg IgG/g particules)	P/N 2 ml	P/N 10 ml
Anti- Human IgG	0.3	0.3-0.6 mg	Please contact us for P/N information	
Anti-Mouse IgG	0.3	0.4-0.6 mg		
Anti-Rabbit IgG	0.3	0.3-0.6 mg		
Anti-Rat IgG	0.3	0.3-0.6 mg		
Anti-Human IgM	0.3	0.75-1 mg		
Anti-Mouse IgM	0.3	2-3 mg		
Anti-Rabbit IgM	0.3	nd		
Streptavidine	0.3	3-4 mg		
Anti Biotin	0.3	3-4 mg		
Anti- Human IgG	0.86	0.1-0.3		
Anti-Mouse IgG	0.86	0.1-0.3		
Anti-Rabbit IgG	0.86	0.5-1		
Anti-Rat IgG	0.86	0.2-0.4 mg		
Anti-Human IgM	0.86	0.1-0.3 mg		
Anti-Mouse IgM	0.86	0.75-1.25 mg		
Anti-Rabbit IgM	0.86	nd		
Streptavidine	0.86	2-3 mg		
Anti Biotin	0.86	2-3 mg		



# Biochromatography - Affinity

## Immobilized Metal Affinity Chromatography - IMAC

**IMAC (Immobilized Metal Affinity Chromatography)** is based on the interaction between certain superficial protein residues (histidines, cysteines and to a lesser extent tryptophans), and transition metal cations forming chelates ( $\text{Ni}^{2+}$ ,  $\text{Co}^{2+}$ ,  $\text{Zn}^{2+}$ ,  $\text{Cu}^{2+}$ ). This technique has proven particularly useful for the purification of recombinant proteins.

Affarose™ format :

Ligands are covalently bound to a modified agarose matrix called Affarose, this ensures minimal leaching of ligand, excellent stability and can be re-used.

UptiSpin™ centrifugal format :

Ready -to-use kits designed for simple, rapid and multiple recombinant protein purification utilizing centrifugation.

### Affarose™ / IMAC

Our Affarose™ IMAC affinity columns are activated with different metal ions in a range of densities for His-Tag recombinant protein purification. Ions show varying affinity & loading capacities with poly-His residues. The management of these dual features offers flexibility for purification applications.

Range :

- Four different ions
- Two different ion densities

Product	P/N	Qty	Product	P/N	Qty
<b>Chelate Affarose Beads</b>			<b>Chelate Test Kits</b>		
High Density IDA-Affarose Nickel charged	BG7020	5 ml	Nickel Chelate kit	BG7080	1 kit
	BG7021	10 ml	Includes 2ml high Density IDA-Affarose Ni Charged		
	BG7022	25 ml	2ml Low Density IDA-Affarose Ni Charged		
Low Density IDA-Affarose Nickel charged	BG7030	5 ml	Zinc Chelate kit	BG7090	1 kit
	BG7031	10 ml	Includes 2ml high Density IDA-Affarose Zn Charged		
	BG7032	25 ml	2ml Low Density IDA-Affarose Zn Charged		
High Density IDA-Affarose Zinc Charged	BG7040	5 ml	High Density Chelate kit	BG7100	1 kit
	BG7041	10 ml	Includes 2ml high Density IDA-Affarose Metal Free		
	BG7042	25 ml	2ml High Density IDA-Affarose Ni Charged		
Low Density IDA-Affarose Zinc Charged	BG7050	5 ml	2ml High Density IDA-Affarose Zn Charged		
	BG7051	10 ml	2ml High Density IDA-Affarose Co Charged		
	BG7052	25 ml			
High Density IDA-Affarose Cobalt charged	BG7060	5 ml	Low Density Chelate kit	BG7110	1 kit
	BG7061	10 ml	Includes 2ml Low Density IDA-Affarose Metal Free		
	BG7062	25 ml	2ml Low Density IDA-Affarose Ni Charged		
Low Density IDA-Affarose Cooper charged	BG7070	5 ml	2ml Low Density IDA-Affarose Zn Charged		
	BG7071	10 ml	2ml Low Density IDA-Affarose Cu Charged		
	BG7072	25 ml			

### Selectivity / loading

Cobalt displays the highest selectivity towards a poly-His sequence but the lowest loading capacity : it is the ion of choice when your recombinant protein is valuable and in limited quantities. Conversely, Copper shows the highest loading capacity but the lowest selectivity.

Nickel and Zinc are often more useful as they show intermediate selectivity and loading capacity.


High density :

The higher the density of ions within the support, the greater the purification of recombinant proteins. However, unwanted proteins within the sample will also be bounded.

Low density :

This will enhance the qualitative purification of recombinant protein but with lower target recoveries.

Interchim offer chelate test kits as screening tools. These enable you to identify the right ion and density for your application.



# Biochromatography - Affinity

## Immobilized Metal Affinity Chromatography - IMAC

### UptiSpin™ / IMAC

UptiSpin™ Metal Chelate Mini and Midi kits contain all necessary components for effective 6x His-tagged recombinant protein purification from bacterial cells, insect vectors, mammalian cells and yeast under native or denaturing conditions. They are ideal for low -to- medium throughput expression scouting, and semi-preparative /preparative purification.

UptiSpin™ spin columns can purify and concentrate protein of interest in a single step and with high reproducibility (low % CV). Proteins are sufficiently pure for structural and activity-based studies or for raising antibodies.

- Speed : Isolate /purify /concentrate sample in less than 45 minutes for the Midi units and under 20 min for Mini units. Simultaneous purifications can be performed for rapid screening of recombinant protein expression.
- Performance : Bind 0.1 mg - 20 mg of 6x His-tagged protein, elute with > 90% recovery & high purity and yield. Devices are re-usable.
- Convenience : Only a centrifuge is required.
- Price : Spin columns are substantially cheaper than research grade FPLC column.

Product	P/N
Mini	
Metal Chelate Mini Sample Kit - 8 Mini MC Plugs (includes buffers)	UPBB9580
Metal Chelate Mini Kit – 24 MC Plugs (includes buffers)	UPBB9600
Metal Chelate Mini Pack – 24 MC plugs (includes buffers, does not include UF concentrators)	UPBB9620
Metal Chelate Mini Bulk Pack (does not include buffers or UF concentrators)	UPBB9630
Midi	
Metal Chelate Midi Kit - 8 MC plugs (includes buffers)	UPBB9640
Metal Chelate Midi Pack - 8 MC plugs (includes buffers, does not include UF concentrators)	UPBB9650
Metal Chelate Midi Bulk Pack – 24 MC Plugs (does not include buffers or UF concentrators)	UPBB9660

### Specifications

Feature	Mini	Midi
Max. sample volume per load :	0.65 ml, fixed angle roto	20 ml, swing bucket rotor
Charged metal ion	Ni <sup>2+</sup>	Ni <sup>2+</sup>
Typical binding capacity per use :	1 mg 6x His-tagged protein	10 mg 6x His-tagged protein
Min. number of uses per plug :	2 (assuming no further Ni <sup>2+</sup> charging)	2 (assuming no further Ni <sup>2+</sup> charging)
Supporting Proteus matrix :	Covalently coupled to agarose resin	Covalently coupled to agarose resin
Resin bed volume :	0.23 ml	1.6 ml
Bead size range :	45-165 µm	45-165 µm
Storage temp. for resin plugs :	2 -to- 8 °C	2 -to- 8 °C
Recommended working pH :	pH 2.0-12.0	pH 2.0-12.0
Chemical stability	High	High
Plastic construction :	Polypropylene	Polypropylene
Purifications time :	15 min	48 min
Color coded end-caps :	Black	Black



# Biochromatography - Affinity

## Immobilized Metal Affinity Chromatography - IMAC

### UptiSpin™ Micro / IMAC

Ready to use spin micro columns filled in with activated agarose beads. Ni<sup>2+</sup>, Zn<sup>2+</sup>, Cu<sup>2+</sup> or Co<sup>2+</sup> metal ions are bound on the resin at different densities.

- For fast, small scale purification of histidine-tagged proteins using standard microcentrifuge.
- Saves time
- Single use spin columns
- Native and denaturing conditions

Product	P/N	Qty
UptiSpin™ Micro IMAC Ni high density	CE1130	25
	CE1131	50
UptiSpin™ Micro IMAC Ni Low density	CE1140	25
	CE1141	50
UptiSpin™ Micro IMAC Zn high density	CE1150	25
	CE1151	50
UptiSpin™ Micro IMAC Zn Low density	CE1160	25
	CE1161	50
UptiSpin™ Micro IMAC Co high density	CE1170	25
	CE1171	50
UptiSpin™ Micro IMAC Cu Low density	CE1180	25
	CE1181	50



### Specifications

Column materials : Polypropylene columns, PE frits  
Bed volume : 100 µl  
Sample Volume : 200 µl  
Bead structure : 6% agarose  
Storage : 20% Ethanol, 2-8 °C

# Biochromatography - Affinity

## Detoxication

### Cellufine™ Sulfate


For concentration, purification and depyrogenation of virus, viral/microbial antigens and heparin-binding proteins

#### Features

- Affinity for wide range of live, dead or disrupted viruses, viral or microbial antigens & heparin-binding proteins
- Closed column operation assures safety and product sterility
- Endotoxins do not bind, allowing a rapid and contaminant free depyrogenation
- Animal extract free
- Rigid, high-strength beads
- Autoclavable

#### Benefits

- More effective than ultracentrifugation at removing contaminants from culture media and host cells
- Avoids excessive product handling and safety concerns, particularly with viral preparations
- Simultaneous concentration /purification improves yields and reduces processing steps.
- Gentle binding and elution conditions provide high capacity and product yield
- Resists compression, providing rapid flow for high-speed processing, even in large columns, making it easily scalable
- Resistant to chemical depyrogenation with base. Chemically sterilizable with formalin



Partial Structure of Cellufine Sulfate

### Cellufine™ ET clean S and L : Endotoxin removal


The Cellufine™ ET clean is poly(ε-lysine) immobilized Cellufine™. The beads bind and remove endotoxin from solution.

#### Features

- Resistant against the cleanup solutions, including 0.2 M sodium hydroxide and 2 M sodium chloride
- Affinity for a wide spectrum of endotoxin

Cellufine ET clean S (small beads) : diameter : 45-105 μm

Cellufine ET clean L (large beads) : diameter : 53-125 μm



Product	P/N	Qty
Cellufine Sulfate	676943324	10 ml
	19845	50 ml
Endotoxin removal		
ET Clean L 53-125 μm	681984324	10 ml
	681984326	50 ml
ET Clean S 44-105 μm	682985324	10 ml
	682985326	50 ml

# Biochromatography - Affinity

## Protein Immobilization


### Cellufine™ Amino and Formyl

Activated supports for Immobilization of Antibodies, Antigens, Affinity Ligands and Enzymes


#### Features

- High flow rates in laboratory and process columns for high throughput
- Low ligand leakage due to exceptionally stable coupling chemistry and support matrix
- Excellent mechanical, chemical and environmental resistance
- High ligand loading capacity
- Compatible with high molecular weight ligands and target proteins due to structural similarity to 4 % cross-linked agarose media
- Unreacted formyl groups easily converted during reduction to neutral hydroxyls for low non-specific adsorption
- Built in hydrophilic spacer arms for maximum ligand accessibility and low non-specific adsorption
- No media damage or fines generation with extended mixing to allow use of simple coupling apparatus
- Ligand coupling occurs under mild conditions in short reaction times
- Thermal stability of media allows high temperature reactions
- Long shelf-life of unreacted media

Product	P/N	Qty
Cellufine formyl	676944324	10 ml
	19853	50 ml
Cellufine Amino	676945324	10 ml
	19856	50 ml



Cellufine Formyl



Cellufine Amino