

# Product Information

## PMA™ Real-Time PCR Bacterial Viability Kit - *Salmonella enterica* (invA)

### Catalog Numbers:

31033 (kit containing PMA)

31033-X (kit containing PMAxx)

**Unit Size:** 1 kit (200 PCR reactions)

### Kit Contents

Component	31033	31033-X
40019: PMA™ dye, 20 mM in H <sub>2</sub> O	1 X 100 uL	
40069: PMAxx™ dye, 20 mM in H <sub>2</sub> O		1 X 100 uL
31038: PMA Enhancer for Gram Negative Bacteria, 5X Solution	1 X 16 mL	1 X 16 mL
31003A: Fast EvaGreen® qPCR Master Mix	2 X 1 mL	2 X 1 mL
99939: 10X ROX Reference Dye	1 X 1 mL	1 X 1 mL
31033A: invA primer mix, 5 uM each primer For: 5'-ATCTGGTACTAATGGTGATGATC-3' Rev: 5'-GCCAGGCTATCGCCAATAAC-3'	1 X 400 uL	1 X 400 uL

### Storage and Handling

Store kit at -20 °C. After first thaw, PMA Enhancer should be stored at 4°C. Store PMA, PMAxx and Fast EvaGreen Master Mix protected from light. Protect PMA and PMAxx from light during use. Components are stable for at least 6 months when stored as recommended. Kit components are stable for several freeze/thaw cycles.

### Spectral Properties

**PMA and PMAxx:**  $\lambda_{abs} = 464$  nm (before photolysis);

$\lambda_{abs}/\lambda_{em} = \sim 510/-610$  nm (following photolysis and reaction with DNA/RNA)

**EvaGreen:**  $\lambda_{abs} = 471$  nm (without DNA)

$\lambda_{abs}/\lambda_{em} = 500/530$  nm (with DNA)

### Product Description

Viability-PCR kits are designed for selective detection of viable bacteria by real-time PCR. Each kit contains a viability dye (PMA™ or PMAxx™), Fast EvaGreen® qPCR Master Mix, and PCR primers for detection of a specific strain of bacteria.

This kit contains primers for amplification within the *Salmonella enterica* invA gene, with reagents sufficient to treat 80 bacterial cultures and perform 200 PCR reactions. The number of samples that can be treated with PMA or PMAxx using the kit may vary depending on sample type.


PMA is a photoreactive DNA binding dye developed by Biotium. It is cell membrane-impermeable and so selectively binds to DNA from dead cells with compromised membrane integrity, while leaving DNA from viable cells intact. Upon photolysis, the dye forms a stable covalent bond, resulting in permanent DNA modification. PMA inhibits PCR amplification of modified DNA templates by a combination of removal of modified DNA during purification and inhibition of template amplification by DNA polymerases (1). Thus the dye is useful in the selective detection of viable pathogenic cells by real-time qPCR (Figure 1).

PMAxx was developed by Biotium as an improved version of our popular PMA dye. In experiments with laboratory bacterial strains, PMAxx increases the difference between live and dead a further 3-7 Ct compared to PMA. Therefore viability PCR with PMAxx is more effective at discriminating between live and dead bacteria. Because PMAxx works the same way as PMA, it can directly replace PMA in your current PMA-PCR protocol.

PMA Enhancer for Gram Negative Bacteria is designed for use with PMA or PMAxx. When Enhancer is added to gram negative bacteria before treatment with PMA, dead cell DNA levels are further decreased, and thus live-dead cell discrimination is improved. The amount of improvement varies depending on such factors as the bacterial strain and the way that the bacterial were killed. PMA Enhancer gives the most improvement when bacteria are dead but their membranes are not completely disrupted, as occurs after mild heat treatment.

Fast EvaGreen Master Mix is a ready-to-use hot-start mix for qPCR and DNA melt curve analysis of PCR amplicons. It is formulated for qPCR using a fast cycling protocol, but can also be used for qPCR using regular cycling protocols. Fast EvaGreen Master Mix contains Cheetah™ Taq, Biotium's fast-activating chemically-modified hot-start Taq polymerase, which is particularly suitable for fast PCR cycling protocols.

*Salmonella enterica* is a gram-negative bacteria that causes the food-borne illness salmonellosis. PMA-based viability PCR for *Salmonella enterica* has been reported using the primers provided in the kit (2), and these primers have been validated at Biotium for real-time qPCR using EvaGreen Master Mix (Figures 2-5).



**Figure 1.** The cell membrane-impermeable PMA dye selectively and covalently modifies DNA from dead bacteria with compromised membranes. Subsequent PCR amplification of PMA-modified DNA templates is inhibited, allowing selective quantitation of viable bacteria.

### Quick guide protocol

(Detailed protocol on following page)

1. Aliquot 400 uL cell culture or sample into tubes. If desired, prepare live and dead cell controls.
2. Add 100 uL 5X PMA Enhancer to a final concentration of 1X. See detailed protocol for more information.
3. Working in dim light, add 25-50 uM PMA or PMAxx viability dye to tubes. Include no-dye controls.
4. Incubate for 10 min, rocking, protected from light.
5. Expose samples to light to crosslink dye to DNA. We recommend 15 min in the PMA-Lite.
6. Isolate DNA using a commercial kit or other protocol.
7. Set up qPCR reactions, using 2 uL of each isolated DNA sample as templates. Do not normalize the DNA concentrations.
8. Compare the amount of total and live-cell-derived DNA in your sample by calculating the dCt ( $dCt = Ct$  with viability dye -  $Ct$  without viability dye). See detailed protocol for more information.

### References

1. Nocker A., et al. Comparison of propidium monoazide with ethidium monoazide for differentiation of live vs. dead bacteria by selective removal of DNA from dead cells. *J. Microbiol. Meth.* 67(2), 310-320 (2006).
2. Nocker A., et al. Molecular monitoring of disinfection efficacy using propidium monoazide in combination with quantitative PCR. *J. Microbiol. Meth.* 70(2), 252-260 (2007).

## Detailed protocol for treating gram-negative bacteria with PMA or PMAxx plus Enhancer for qPCR

The following is a protocol for treating cultured laboratory strains of gram-negative bacteria with PMA or PMAxx. Treatment of complex biological or environmental samples such as feces or soil may require optimization of sample dilution for dye and light treatment. PMA Enhancer generally improves the activity of PMA and PMAxx on gram-negative bacteria, but has a detrimental effect on gram-positive bacteria. However, you may want to test whether it is beneficial in your assay of interest. If both gram-negative and gram-positive bacteria are to be treated in one sample, Enhancer should not be used.

- Inoculate an appropriate media broth with bacteria (volume is dependent on size of experiment).
- Shake cultures at 200 RPM at 37°C overnight.
- Continuing culturing bacteria until the OD<sub>600</sub> of the culture is approximately 1.
- For dead cell control samples, heat inactivate bacteria at 58°C for 3 hours or 95°C for 5 min. To confirm killing of bacteria, plate 10 uL of heat inactivated bacteria on the appropriate media plate, and 10 uL of a 1:100 dilution of control bacteria on another plate. Place the plate at 37°C and check for colony growth after 24-48 hours.
- Pipette 400 uL aliquots of bacterial culture into clear microcentrifuge tubes.
- [Optional]: Add 100 uL of 5X Enhancer to each tube, for a 1X final Enhancer concentration.
- Working quickly and in low light, thaw the 20 mM PMA or PMAxx stock and prepare a working stock by diluting to 5 mM in water. If using Enhancer, add 2.5 uL of working stock to each tube for a final concentration of 25 uM. If Enhancer is not used, we recommend adding 4 uL of PMA working stock to 400 uL of sample for a final concentration of 50 uM. 25 uM PMAxx should be sufficient for most assays.
- Incubate tubes in the dark for 10 minutes at room temperature. Flick tubes occasionally to mix, or incubate on a rocker covered with aluminum foil.
- Expose samples to light to cross-link PMA or PMAxx to DNA.
  - For best results, we recommend that the photo-crosslinking be carried out on Biotium's PMA-Lite LED Photolysis Device see next page for more information). 15 min exposure should be sufficient for complete PMA or PMAxx activation.
  - Commercial halogen lamps (>600 W) for home use have been employed for photoactivating PMA in some publications, though results have not been consistent due to inevitable variation in the set-up configurations. If you decide to use a halogen lamp, we recommend that you lay tubes on a block of ice set 20 cm from the light source. The ice block should be in a clear tray with a piece of aluminum foil under the clear tray to reflect the light upward. Set the lamp so that the light source is pointing directly downward onto the samples. Expose samples to light for 5-15 min.
- Pellet cells by centrifuging at 5,000 x g for 10 minutes.
- Extract genomic DNA using a standard protocol or commercially available kit. Use an appropriate protocol or kit for DNA extraction from complex biological or environmental samples (e.g., feces or soil).

- Perform qPCR using the invA primers included with this kit to detect *Salmonella enterica*. See reaction setup and fast cycling parameters below.
- Data analysis: Compare the amount of total and live-cell-derived DNA in your sample by calculating the dCt (dCt = Ct with viability dye - Ct without viability dye). The dCt of a control sample of killed cells can be calculated to determine the maximum inhibition that can be achieved by PMA or PMAxx in your sample, and the dCt of control live cells can be calculated to control for false negatives that may arise from dye getting into live cells.

**Note 1:** Part of the proposed mechanism of action of PMA is the removal of PMA-bound DNA from samples via precipitation; therefore the amount of input DNA in each sample should not be normalized between samples. Instead, PCR should be performed using equal volumes of gDNA eluate from each sample. For a positive control, 1 ng of live cell gDNA per reaction should be sufficient for achieving good signal. For gDNA extracted from bacterial cultures using a commercial extraction kit, 2 uL of eluted DNA can be used as a starting point for optimization.

### PCR Reaction Setup

Add reaction components to each PCR tube or well according to the table below:

Reaction component	Amount per 20 uL reaction	Final concentration
2X Fast EvaGreen Master Mix	10 uL	1X
invA primer mix, 5 uM	2 uL	0.5 uM each
Template	x uL See Notes 1&2	See Note 1
ROX	Optional	See Table 1
dH <sub>2</sub> O	Add to 20 uL	

**Note 2:** Template volume should not exceed 10% of final reaction volume.


### Fast-cycling parameters for invA real-time PCR on *S. enterica* gDNA

Hold	
95 °C for 2-10 minutes (see Note 3)	
Cycling	
95 °C for 5 seconds	Cycle 40 times
60 °C for 30 seconds (acquire data)	
Melt	
57 °C to 99 °C	


**Note 3 -** Activation of Cheetah™ Taq DNA Polymerase requires only 2 minutes at 95 °C, but genomic DNA can take longer to fully denature. If you observe high background fluorescence during initial amplification cycles, try increasing the hold time.

**Table 1. Recommended ROX Concentration for PCR Instruments**


PCR Instrument	Recommended Rox Concentration	Amount of 10X ROX per 20 uL reaction
BioRad: iCycler, MyiQ, MiQ 2, iQ 5, CFX-96, CFX-384, MJ Opticon, Option2, Chromo4, MiniOpticon Qiagen: Rotor-Gene Q, Rotor-Gene3000, Rotor-Gene 6000 Eppendorf: Mastercycler realplex Illumina: Eco RealTime PCR System Cepheid: SmartCycler Roche: LightCycler 480, LightCycler 2.0	No ROX	None
ABI: 7500, 7500 Fast, Viia 7 Stratagene: MX4000P, MX3000P, MX3005P	Low ROX 0.05-0.1X final	Dilute 10X ROX 1:10 with dH <sub>2</sub> O to obtain 1X ROX; add 1 to 2 uL of 1X ROX per 20 uL reaction
ABI: 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast, StepOne, StepOne plus	High ROX 1X final	2 uL of 10X ROX per 20 uL reaction



**Figure 2.** Fast EvaGreen® qPCR Master Mix was used to amplify a fragment of *invA* from 10 ng, 1 ng, 100 pg, 10 pg, or 1 pg of *S. enterica* gDNA (ATCC). The real-time PCR was performed on a RotorGeneQ (Qiagen).




**Figure 3.** Melt curve analysis of the *invA* real-time PCR product generated in Figure 2, from 1 ng of *S. enterica* gDNA input.



**Figure 4.** Reaction product from real-time PCR amplification of *invA* (288 bp fragment) from 1 ng *S. enterica* gDNA input. Biotium's 100 bp DNA ladder was run in the first lane. The 1% agarose 1X TBE gel was post-stained with 3X GelRed in water and imaged on a UVP GelDoc-iT using UV illumination and an ethidium bromide filter (3 second exposure).

### Light sources for photoactivation

Biotium offers the PMA-Lite™ LED Photolysis Device for light-induced cross-linking of PMA to dsDNA. The PMA-Lite™ LED Photolysis Device is a thermally-stable blue LED light source that provides even illumination to all samples. It contains a cooling unit to prevent sample overheating as well as several timer settings to allow for precisely timed light treatment.



**Figure 5.** Live and heat-killed *S. enterica* were treated with 25  $\mu$ M PMAxx, followed by photoactivation. Real-time PCR was performed with the *invA* primers from this kit and Fast EvaGreen Master Mix, on a RotorGeneQ (Qiagen). PMAxx reduced the dead cell signal by ~8 Ct, but had no effect on the live cell signal.

### Related Products

Catalog number	Product
E90002	PMA-Lite™ LED Photolysis Device
40013	PMA™ dye, 1 mg
40019	PMA™ dye, 20 mM in dH <sub>2</sub> O, 100 $\mu$ L
40069	PMAxx™ dye, 20 mM in dH <sub>2</sub> O, 100 $\mu$ L
31038	PMA Enhancer for Gram Negative Bacteria, 5X Solution
31034	PMA Real-Time PCR Bacterial Viability Kit - Mycobacterium tuberculosis ( <i>groEL2</i> )
31035	PMA Real-Time PCR Bacterial Viability Kit - Staphylococcus aureus ( <i>nuc</i> )
31036	PMA Real-Time PCR Bacterial Viability Kit - Staphylococcus aureus ( <i>mecA</i> )
31037	PMA Real-Time PCR Bacterial Viability Kit - E. coli 0157:H7 ( <i>stx1</i> )
31050	PMA Real-Time PCR Bacterial Viability Kit - E. coli ( <i>uidA</i> )
31051	PMA Real-Time PCR Bacterial Viability Kit - Listeria monocytogenes ( <i>hly</i> )
31053	PMA Real-Time PCR Bacterial Viability Kit - Legionella pneumophila ( <i>mip</i> )
31003	Fast EvaGreen® qPCR Master Mix (200 rxn), 2 x 1 mL
31021	1 kb DNA Ladder (100ng/ $\mu$ L), 300 $\mu$ g/300 $\mu$ L
31022	Ready-to-Use 1 kb DNA Ladder, 150 applications (1.5 mL)
31031	100 bp DNA Ladder, 30 $\mu$ g/300 $\mu$ L
31032	Ready-to-Use 100 bp DNA Ladder, 150 applications (1.5 mL)
41003	GelRed™ Nucleic Acid Gel Stain, 10,000X in water, 0.5 mL
32000-1	Live Bacterial Gram Stain Kit
32001	Bacterial Viability and Gram Stain Kit
30027	Viability/Cytotoxicity Assay Kit for Bacterial Live and Dead Cells

Please visit our website at [www.biotium.com](http://www.biotium.com) for information on our life science research products.

GelRed™ and its uses are covered by granted US patents. EvaGreen® is a registered trademark of Biotium, Inc. EvaGreen® dye and applications are covered under granted and pending US and international patents. Cheetah™ Taq and its uses are covered by a granted US patent. PMA™ and PMAxx™ are trademarks of Biotium, Inc. PMAxx™ and its uses are covered by granted and/or pending US patents. Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.