

Product Information

PMA™ Real-Time PCR Bacterial Viability Kit - *E. coli* (*uidA*)

Catalog Numbers:

31050 (kit containing PMA)

31050-X (kit containing PMAxx)

Unit Size: 1 kit (200 PCR reactions)

Kit Contents

| Component | 31050 | 31050-X |
|--|------------|------------|
| 40019: PMA™ dye, 20 mM in H ₂ O | 1 X 100 uL | |
| 40069: PMAxx™ dye, 20 mM in H ₂ O | | 1 X 100 uL |
| 31038: PMA Enhancer for Gram Negative Bacteria, 5X Solution | 1 X 16 mL | 1 X 16 mL |
| 31003A: Fast EvaGreen® qPCR Master Mix | 2 X 1 mL | 2 X 1 mL |
| 99939: 10X ROX Reference Dye | 1 X 1 mL | 1 X 1 mL |
| 31050A: <i>uidA</i> primer mix, 5 uM each primer For: 5'-CGGTGATATCGTCCACCCAG-3' Rev: 5'-TGGATCGCGAAAACGTGTGA-3' | 1 X 400 uL | 1 X 400 uL |

Storage and Handling

Store kit at -20 °C. After first thaw, PMA Enhancer should be stored at 4°C. Store PMA, PMAxx and Fast EvaGreen Master Mix protected from light. Protect PMA and PMAxx from light during use. Components are stable for at least 6 months when stored as recommended. Kit components are stable for several freeze/thaw cycles.

Spectral Properties

PMA and PMAxx: $\lambda_{abs} = 464$ nm (before photolysis);

$\lambda_{abs}/\lambda_{em} = \sim 510/\sim 610$ nm (following photolysis and reaction with DNA/RNA)

EvaGreen: $\lambda_{abs} = 471$ nm (without DNA)

$\lambda_{abs}/\lambda_{em} = 500/530$ nm (with DNA)

Product Description

Viability-PCR kits are designed for selective detection of viable bacteria by real-time PCR. Each kit contains a viability dye (PMA™ or PMAxx™), Fast EvaGreen® qPCR Master Mix, and PCR primers for detection of a specific strain of bacteria.

This kit contains primers for amplification within the *E. coli uidA* gene, with reagents sufficient to treat 80 bacterial cultures and perform 200 PCR reactions. The number of samples that can be treated with PMA or PMAxx using the kit may vary depending on sample type.

PMA is a photoreactive DNA binding dye developed by Biotium. It is cell membrane-impermeable and so selectively binds to DNA from dead cells with compromised membrane integrity, while leaving DNA from viable cells intact. Upon photolysis, the dye forms a stable covalent bond, resulting in permanent DNA modification. PMA inhibits PCR amplification of modified DNA templates by a combination of removal of modified DNA during purification and inhibition of template amplification by DNA polymerases (1). Thus the dye is useful in the selective detection of viable pathogenic cells by real-time qPCR (Figure 1).

PMAxx was developed by Biotium as an improved version of our popular PMA dye. In experiments with laboratory bacterial strains, PMAxx increases the difference between live and dead a further 3-7 Ct compared to PMA. Therefore viability PCR with PMAxx is more effective at discriminating between live and dead bacteria. Because PMAxx works the same way as PMA, it can directly replace PMA in your current PMA-PCR protocol.

PMA Enhancer for Gram Negative Bacteria is designed for use with PMA or PMAxx. When Enhancer is added to gram negative bacteria before treatment with PMA, dead cell DNA levels are further decreased, and thus live-dead cell discrimination is improved. The amount of improvement varies depending on such factors as the bacterial strain and the way that the bacteria were killed. PMA Enhancer gives the most improvement when bacteria are dead but their membranes are not completely disrupted, as occurs after mild heat treatment.

Fast EvaGreen Master Mix is a ready-to-use hot-start mix for qPCR and DNA melt curve analysis of PCR amplicons. It is formulated for qPCR using a fast cycling protocol, but can also be used for qPCR using regular cycling protocols. Fast EvaGreen Master Mix contains Cheetah™ Taq, Biotium's fast-activating chemically-modified hot-start Taq polymerase, which is particularly suitable for fast PCR cycling protocols.

Escherichia coli is a commonly used laboratory bacteria, and some strains can cause digestive illness. PCR to amplify the gene *uidA* has been published and shown to be highly specific for *E. coli* (2). The primers provided in the kit have been validated at Biotium for real-time qPCR using EvaGreen Master Mix (Figures 2-4). An example of PMAxx-PCR in *E. coli* is shown in Figure 5.




Figure 1. The cell membrane-impermeable PMA dye selectively and covalently modifies DNA from dead bacteria with compromised membranes. Subsequent PCR amplification of PMA-modified DNA templates is inhibited, allowing selective quantitation of viable bacteria.

Quick guide protocol

(Detailed protocol on following page)

1. Aliquot 400 uL cell culture or sample into tubes. If desired, prepare live and dead cell controls.
2. Add 100 uL 5X PMA Enhancer to a final concentration of 1X. See detailed protocol for more information.
3. Working in dim light, add 25-50 uM PMA or PMAxx viability dye to tubes. Include no-dye controls.
4. Incubate for 10 min, rocking, protected from light.
5. Expose samples to light to crosslink dye to DNA. We recommend 15 min in the PMA-Lite.
6. Isolate DNA using a commercial kit or other protocol.
7. Set up qPCR reactions, using 2 uL of each isolated DNA sample as templates. Do not normalize the DNA concentrations.
8. Compare the amount of total and live-cell-derived DNA in your sample by calculating the dCt ($dCt = Ct$ with viability dye - Ct without viability dye). See detailed protocol for more information.

Detailed protocol for treating gram-negative bacteria with PMA or PMAxx plus Enhancer for qPCR

The following is a protocol for treating cultured laboratory strains of gram-negative bacteria with PMA or PMAxx. Treatment of complex biological or environmental samples such as feces or soil may require optimization of sample dilution for dye and light treatment. PMA Enhancer generally improves the activity of PMA and PMAxx on gram-negative bacteria, but has a detrimental effect on gram-positive bacteria. However, you may want to test whether it is beneficial in your assay of interest. If both gram-negative and gram-positive bacteria are to be treated in one sample, Enhancer should not be used.

- Inoculate an appropriate media broth with bacteria (volume is dependent on size of experiment).
- Shake cultures at 200 RPM at 37°C overnight.
- Continuing culturing bacteria until the OD₆₀₀ of the culture is approximately 1.
- For dead cell control samples, heat inactivate bacteria at 58°C for 3 hours or 95°C for 5 min. To confirm killing of bacteria, plate 10 uL of heat inactivated bacteria on the appropriate media plate, and 10 uL of a 1:100 dilution of control bacteria on another plate. Place the plate at 37°C and check for colony growth after 24-48 hours.
- Pipette 400 uL aliquots of bacterial culture into clear microcentrifuge tubes.
- [Optional]: Add 100 uL of 5X Enhancer to each tube, for a 1X final Enhancer concentration.
- Working quickly and in low light, thaw the 20 mM PMA or PMAxx stock and prepare a working stock by diluting to 5 mM in water. If using Enhancer, add 2.5 uL of working stock to each tube for a final concentration of 25 uM. If Enhancer is not used, we recommend adding 4 uL of PMA working stock to 400 uL of sample for a final concentration of 50 uM. 25 uM PMAxx should be sufficient for most assays.
- Incubate tubes in the dark for 10 minutes at room temperature. Flick tubes occasionally to mix, or incubate on a rocker covered with aluminum foil.
- Expose samples to light to cross-link PMA or PMAxx to DNA.
 - For best results, we recommend that the photo-crosslinking be carried out on Biotium's PMA-Lite LED Photolysis Device (see next page for more information). 15 min exposure should be sufficient for complete PMA or PMAxx activation.
 - Commercial halogen lamps (>600 W) for home use have been employed for photoactivating PMA in some publications, though results have not been consistent due to inevitable variation in the set-up configurations. If you decide to use a halogen lamp, we recommend that you lay tubes on a block of ice set 20 cm from the light source. The ice block should be in a clear tray with a piece of aluminum foil under the clear tray to reflect the light upward. Set the lamp so that the light source is pointing directly downward onto the samples. Expose samples to light for 5-15 min.
- Pellet cells by centrifuging at 5,000 x g for 10 minutes.
- Extract genomic DNA using a standard protocol or commercially available kit. Use an appropriate protocol or kit for DNA extraction from complex biological or environmental samples (e.g., feces or soil).

- Perform qPCR using the uidA primers included with this kit to detect *E. coli*. See reaction setup and fast cycling parameters below.
- Data analysis: Compare the amount of total and live-cell-derived DNA in your sample by calculating the dCt (dCt = Ct with viability dye - Ct without viability dye). The dCt of a control sample of killed cells can be calculated to determine the maximum inhibition that can be achieved by PMA or PMAxx in your sample, and the dCt of control live cells can be calculated to control for false negatives that may arise from dye getting into live cells.

Note 1: Part of the proposed mechanism of action of PMA is the removal of PMA-bound DNA from samples via precipitation; therefore the amount of input DNA in each sample should not be normalized between samples. Instead, PCR should be performed using equal volumes of gDNA eluate from each sample. For a positive control, 1 ng of live cell gDNA per reaction should be sufficient for achieving good signal. For gDNA extracted from bacterial cultures using a commercial extraction kit, 2 uL of eluted DNA can be used as a starting point for optimization.

PCR Reaction Setup

Add reaction components to each PCR tube or well according to the table below:

| Reaction component | Amount per 20 uL reaction | Final concentration |
|-----------------------------|---------------------------|---------------------|
| 2X Fast EvaGreen Master Mix | 10 uL | 1X |
| uidA primer mix, 5 uM | 2 uL | 0.5 uM each |
| Template | x uL See Notes 1&2 | See Note 1 |
| ROX | Optional | See Table 1 |
| dH ₂ O | Add to 20 uL | |

Note 2: Template volume should not exceed 10% of final reaction volume.

Fast-cycling parameters for uidA real-time PCR on *E. coli* gDNA

| | |
|-------------------------------------|----------------|
| Hold | |
| 95 °C for 2-10 minutes (see Note 3) | |
| Cycling | |
| 95 °C for 5 seconds | Cycle 40 times |
| 64 °C for 30 seconds (acquire data) | |
| Melt | |
| 57 °C to 99 °C | |

Note 3 - Activation of Cheetah™ Taq DNA Polymerase requires only 2 minutes at 95 °C, but genomic DNA can take longer to fully denature. If you observe high background fluorescence during initial amplification cycles, try increasing the hold time.

Table 1. Recommended ROX Concentration for PCR Instruments

| PCR Instrument | Recommended Rox Concentration | Amount of 10X ROX per 20 uL reaction |
|---|-------------------------------|---|
| BioRad: iCycler, MyiQ, MiQ 2, iQ 5, CFX-96, CFX-384, MJ Opticon, Option2, Chromo4, MiniOpticon Qiagen: Rotor-Gene Q, Rotor-Gene3000, Rotor-Gene 6000 Eppendorf: Mastercycler realplex Illumina: Eco RealTime PCR System Cepheid: SmartCycler Roche: LightCycler 480, LightCycler 2.0 | No ROX | None |
| ABI: 7500, 7500 Fast, Viia 7 Stratagene: MX4000P, MX3000P, MX3005P | Low ROX 0.05-0.1X final | Dilute 10X ROX 1:10 with dH ₂ O to obtain 1X ROX; add 1 to 2 uL of 1X ROX per 20 uL reaction |
| ABI: 5700, 7000, 7300, 7700, 7900, 7900HT, 7900HT Fast, StepOne, StepOne plus | High ROX 1X final | 2 uL of 10X ROX per 20 uL reaction |




Figure 2. Fast EvaGreen® qPCR Master Mix was used to amplify a fragment of *uidA* from 10 ng, 1 ng, 100 pg, 10 pg, or 1 pg of *E. coli* gDNA (ATCC). The real-time PCR was performed on a RotorGeneQ (Qiagen).




Figure 3. Melt curve analysis of the *uidA* real-time PCR product generated in Figure 2, from 1 ng of *E. coli* gDNA input.




Figure 4. Reaction product from real-time PCR amplification of *uidA* (500 bp fragment) from 1 ng *E. coli* gDNA input. Biotium's 100 bp DNA ladder was run in the first lane. The 1% agarose 1X TBE gel was post-stained with 3X GelRed in water and imaged on a UVP GelDoc-iT using UV illumination and an ethidium bromide filter (3 second exposure).

Light sources for photoactivation

Biotium offers the PMA-Lite™ LED Photolysis Device for light-induced cross-linking of PMA to dsDNA. The PMA-Lite™ LED Photolysis Device is a thermally-stable blue LED light source that provides even illumination to all samples. It contains a cooling unit to prevent sample overheating as well as several timer settings to allow for precisely timed light treatment.



PMA™ Real-Time PCR Bacterial Viability Kit - *E. coli* (*uidA*)
PSF006




Figure 5. Live and heat-killed *E. coli* were treated with 25 μ M PMAxx, followed by photoactivation. Real-time PCR was performed with primers that amplify a region of the *uidA* gene of *E. coli* and Fast EvaGreen Master Mix, on a RotorGeneQ (Qiagen). PMAxx reduced the dead cell signal more strongly than PMA, but had no effect on the live cell signal.

References

- Nocker A., et al. Comparison of propidium monoazide with ethidium monoazide for differentiation of live vs. dead bacteria by selective removal of DNA from dead cells. *J. Microbiol. Meth.* 67(2), 310-320 (2006).
- Maheux AF, et al. Analytical comparison of nine PCR primer sets designed to detect the presence of *Escherichia coli*/*Shigella* in water samples. *Water Res.* (43) 3019-3028 (2009).

Related Products

| Catalog number | Product |
|----------------|---|
| E90002 | PMA-Lite™ LED Photolysis Device |
| 40013 | PMA™ dye, 1 mg |
| 40019 | PMA™ dye, 20 mM in dH ₂ O, 100 μ L |
| 40069 | PMAxx™ dye, 20 mM in dH ₂ O, 100 μ L |
| 31038 | PMA Enhancer for Gram Negative Bacteria, 5X Solution |
| 31033 | PMA Real-Time PCR Bacterial Viability Kit - <i>Salmonella enterica</i> (<i>invA</i>) |
| 31034 | PMA Real-Time PCR Bacterial Viability Kit - <i>Mycobacterium tuberculosis</i> (<i>groEL2</i>) |
| 31035 | PMA Real-Time PCR Bacterial Viability Kit - <i>Staphylococcus aureus</i> (<i>nuc</i>) |
| 31036 | PMA Real-Time PCR Bacterial Viability Kit - <i>Staphylococcus aureus</i> (<i>mecA</i>) |
| 31037 | PMA Real-Time PCR Bacterial Viability Kit - <i>E. coli</i> O157:H7 (Z3276) |
| 31051 | PMA Real-Time PCR Bacterial Viability Kit - <i>Listeria monocytogenes</i> (<i>hly</i>) |
| 31053 | PMA Real-Time PCR Bacterial Viability Kit - <i>Legionella pneumophila</i> (<i>mip</i>) |
| 31003 | Fast EvaGreen® qPCR Master Mix (200 rxn), 2 x 1 mL |
| 31022 | Ready-to-Use 1 kb DNA Ladder, 150 applications (1.5 mL) |
| 31032 | Ready-to-Use 100 bp DNA Ladder, 150 applications (1.5 mL) |
| 41003 | GelRed™ Nucleic Acid Gel Stain, 10,000X in water, 0.5 mL |
| 32000-1 | Live Bacterial Gram Stain Kit |
| 32001 | Bacterial Viability and Gram Stain Kit |
| 30027 | Viability/Cytotoxicity Assay Kit for Bacterial Live and Dead Cells |

Please visit our website at www.biotium.com for information on our life science research products.

GelRed™ and its uses are covered by granted US patents. EvaGreen® is a registered trademark of Biotium, Inc. EvaGreen® dye and applications are covered under granted and pending US and international patents. Cheeta™ Taq and its uses are covered by a granted US patent. PMA™ and PMAxx™ are trademarks of Biotium, Inc. PMAxx™ and its uses are covered by granted and/or pending US patents. Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.