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IVD

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Defibrinated horse blood (BH-F)

IVD

Culture of fastidious bacteria. Determination of hemolysis.

SUMMARY AND EXPLANATION

Defibrinated horse blood provides a rich nutritive additive suited to the growth of most fastidious bacterial species regardless of their metabolism.

It reveals hemolysis which is a basic criterion for orienting bacterial identification (1).

PRINCIPLE

Horse blood is used to supplement agar culture media (such as Columbia agar, Trypcase soy agar, etc.) in order to prepare fresh blood agars (3).

CONTENT OF THE KIT

REF 55 832	1 x 100 ml bottle
REF 55 833	1 x 20 ml bottle

COMPOSITION

Defibrinated horse blood without preservative or anticoagulant

REAGENTS REQUIRED BUT NOT PROVIDED

- Columbia agar (Ref. 41 244 or 51 026).
- Trypcase soy agar (Ref. 41 466, 41 467 or 51 044).

WARNINGS AND PRECAUTIONS

- **For *in vitro* diagnostic use only.**
- **For professional use only.**
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- All specimens, microbial cultures and inoculated products should be considered infectious and handled appropriately. Aseptic technique and usual precautions for handling the bacterial group studied should be observed throughout this procedure. Refer to "CLSI/NCCLS M29-A, *Protection of Laboratory Workers from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids, and Tissue; Approved Guideline - Current Revision*". For additional information on handling precautions, refer to "Biosafety in Microbiological and Biomedical Laboratories – CDC/NIH – Latest edition", or the current regulations in the country of use.
- Do not use reagents past the expiry date.
- Do not use bottles which show signs of contamination or hemolysis.
- Before use, make sure that the tamper-proof seal of the bottle stopper is intact.

- The performance data presented were obtained using the procedure indicated in this package insert. Any change or modification in the procedure may affect the results.

STORAGE CONDITIONS

- **Store the bottles at 2-8°C until the expiry date.**
- After opening, the bottles can be stored for 1 week at 2-8°C

INSTRUCTIONS FOR USE

1. **Allow the bottle of blood to come to room temperature.**
2. Shake gently to mix.
3. Add 5% blood to the previously liquefied agar bottle maintained at 45-50 °C. The user is responsible for choosing the blood concentration according to the intended use of the medium and in accordance with current standards.
4. Mix carefully avoiding the formation of bubbles.
5. Dispense into Petri dishes (18-20 ml per plate).
6. After reconstitution and cooling of the medium, keep the plates at 2-8°C.

READING AND INTERPRETATION

Refer to the package insert of the agar used.

QUALITY CONTROL

Protocol:

The control is performed using blood added to Trypcase soy agar.

The nutrient capacity of the medium can be tested using the following strains (incubation in CO₂-enriched atmosphere):

- *Streptococcus pyogenes* ATCC® 19615.
- *Streptococcus pneumoniae* ATCC® 6305.

Range of expected results:

Strain	Results at 33-37°C	
<i>Streptococcus pyogenes</i> ATCC® 19615	Growth after 24 hours	β hemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6305		α hemolysis

Note:

It is the responsibility of the user to perform Quality Control taking into consideration the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature and environment, etc.).

PERFORMANCE

Performance was evaluated at 37°C using Columbia and Trypcase soy agars.

Forty strains were tested: 39 bacterial strains (*Streptococci*, *Streptococcus pneumoniae*, *Listeria*, other Gram (+) and Gram (-) bacteria and anaerobic bacteria) and 1 yeast (*Candida*).

Nutrient capacity:

All the tested strains grew within 24-48 hours on both blood agars.

Hemolysis:

The following types of hemolyses are observed after 24 hours of incubation:

- β hemolysis for the 5 strains of *Streptococcus pyogenes*, for 2 out of the 4 strains of streptococci D / enterococci and for the 4 strains of *Listeria* (*L. monocytogenes*, *L. seeligeri* and *L. ivanovii*).
Note: β hemolysis is not characteristic of these species: certain staphylococci and streptococci (*S. agalactiae* and group C) have also produced β hemolysis.
- α hemolysis for the 5 strains of *Streptococcus pneumoniae*.

WASTE DISPOSAL

Dispose of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

LITERATURE REFERENCES

1. DELMAS P., FRENEY J. – Les streptocoques. – *Lyon Pharm.*, 1989, vol. 40, n°5, p. 353-369.
2. FACKLAM R.R., PADULA J.F., MORTHAM E.C. and al. – Presumptive identification of group A, B, and D streptococci on agar plate media. – *J. Clin. Microbiol.*, 1979, vol. 9, n°6, p. 665-672.
3. FINEGOLD S.M., MARTIN W.J., SCOTT E.G. – *Bailey and Scott's diagnostic microbiology* – MOSBY, 1978.
4. FLANDROIS J.P., CHOMARAT M. – *Bactériologie médicale pratique* - MEDSI / Mac GRAW-HILL, 1989 – ISBN 2-86439-161-9.
5. MURRAY P.R., BARON E.J., PFALLER M.A. and al. – *Manual of Clinical Microbiology* - 6th Ed.- ASM Press, 1995 – ISBN 1-55581-086-1.
6. RODLOFF A.C., APPELBAUM P.C., ZABRANSKY R.J. – *Cumitech 5A. Practical anaerobic bacteriology*. - American Society for Microbiology, 1991 – ISBN 1-55581-C05A.
7. WILKINS T.D., CHALGREN S.L., LIMENEZ-ULATE F. and al. – Inhibition of *Bacteroides fragilis* on blood agar plates and reversal of inhibition by added hemin. – *J. Clin. Microbiol.*, 1976, vol. 3, n°3, p. 359-363.

INDEX OF SYMBOLS

Symbol	Meaning
	GB : Catalogue number US : Catalog number
	In Vitro Diagnostic Medical Device
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use

WARRANTY

bioMérieux disclaims all warranties, express or implied, including any implied warranties of MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE. bioMérieux shall not be liable for any incidental or consequential damages. IN NO EVENT SHALL BIOMERIEUX'S LIABILITY TO CUSTOMER UNDER ANY CLAIM EXCEED A REFUND OF THE AMOUNT PAID TO BIOMERIEUX FOR THE PRODUCT OR SERVICE WHICH IS THE SUBJECT OF THE CLAIM.

ATCC is a used, pending and/or registered trademark belonging to American Type Culture Collection.

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bioMérieux and the blue logo are used, pending and/or registered trademarks belonging to bioMérieux SA or one of its subsidiaries.

69280 Marcy-l'Etoile / France
Tel. 33 (0)4 78 87 20 00
Fax 33 (0)4 78 87 20 90
<http://www.biomerieux.com>

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