

Consumables - Storage

HPLC Magic Box™

**HPLC Magic Box :
Stop Dead Volume**
Sure-Fit capillary,
Uptisur™ connectors,
PEEK tubing, Upti-Cut,
wrenches. (23 parts)



BZ3050

HPLC Magic Box : Stainless Steel Capillaries (1/16")

Perfectly perpendicular, burr-free
end cuts, precision lengths
5, 10, 15cm with 0.13, 0.17,
0.25mm (5 units of each).
(15 parts)



BZ3300

HPLC Magic Box PEEK Capillaries (1/16")

Perfectly perpendicular, burr-free
end cuts, precision lengths
5, 10, 20cm with 0.13, 0.17,
0.25mm (5 units of each).
(15 parts)




BZ3310



Consumables - Storage

GC Magic Box™

GC Magic Box : "Gas line"
Copper tubing 1/8", Tubing Cut,
Reamer, Nuts & Back and Front
Ferrules, Wrenches. (38 parts)



BZ1580



BZ1870

GC Magic Box : Brass fittings 1/4", 1/8"
Nuts, Front and Back Ferrules, Unions, Tees,
Caps. (182 parts).