


Consumables - Syringes

SGE Syringes

Needle Selection



Cone Tip

minimizes septum coring for GC autosampler use



Bevel Tip

standard general purpose needle tip style for dispensing liquids
ideal for GC analysis with manual injection




Dual gauge

narrow gauge at front end for on column injection
wide gauge to increase strength for autosampler use



LC Tip

square tip needle prevents damage to the valve's rotor seal and stator face of the injector
suitable for the pipetting of liquids



Dome Tip

rounded and polished tip to facilitate septum penetration



Side Hole/ Dome Tip

liquid sample is filled and dispensed through the side hole
recommended for large volume GC injection



Teflon coated Tip

ideal for dispensing aqueous samples without wetting the outside of the needle



Sheathed/Bevel Tip

same as the bevel tip, except a 0.7 mm OD, re-enforcing sheath is silver soldered over the needle to within 10 mm of the tip. Sheath provides exceptional strength to the needle

Consumables - Syringes

SGE Syringes

Syringe Selection



Fixed Needle Syringe

Easy to use
Economical



Removable Needle Syringe

Easily replace the needle in case of bending
Able to change to a different length or outer diameter needle with the same syringe



MicroVolume Syringe with Plunger-in-needle

Plunger extends to needle tip
Sample is contained only in the needle, i.e no glass contact
Ideal for dispensing very small liquid volumes



MicroVolume Syringe with Standard Plunger

Stainless steel plunger individually fitted to its own barrel
Plunger not replaceable
Industry standard syringe



Gas Tight Syringe

Teflon tip plunger
Suitable for both liquid and gas samples



Plunger Protection

Extended protection tube from top of the flange helps prevent plunger bending during injection



Guided Plunger

Extended barrel guides plunger during injection
Suitable for technician and student use



SuperfleX™ Flexible Plunger

Made from Titanium /Nickel
Indestructible plunger
Ideal for student use

Consumables - Syringes

SGE Syringes

GC autosampler syringes

Agilent Technologies

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty	P/N	Qty
--------	--------	-------	----	-----------	-----	-----	-----	-----	-----	-----

Fixed needle

5 µl	42 mm	26	0.47 mm	Cone	001800	1				
5 µl (M)	42 mm	23	0.63 mm	Cone	001810	1	001814	6		
10 µl	42 mm	26	0.47 mm	Cone	002800	1	002804	6		
10 µl (M)	42 mm	23	0.63 mm	Cone	002810	1	002814	6	002813	25
10 µl Gas Tight (M)	42 mm	23	0.63 mm	Cone	002812	1				

Fixed needle - dual gauge

5 µl	42 mm	23-26	0.63 /0.47 mm	Cone	001821	1	001822	6		
10 µl	42 mm	23-26	0.63 /0.47 mm	Cone	002821	1	002822	6	002824	25
10 µl Gas Tight	42 mm	23-26	0.63 /0.47 mm	Cone	002826	1	002827	6		
25 µl Gas Tight	42 mm	23-26	0.63 /0.47 mm	Cone	003668	1				
50 µl Gas Tight	42 mm	23-26	0.63 /0.47 mm	Cone	004668	1				
100 µl Gas Tight	42 mm	23-26	0.63 /0.47 mm	Cone	005668	1				

Fixed needle - SuperfleX™ - Flexible plungers

10 µl	42 mm	26	0.47 mm	Cone	002830	1	002831	6		
10 µl	42 mm	23	0.63 mm	Cone	002835	1	002838	6		

Removable needle

0.5 µl	42 mm	26	0.47 mm	Cone	000400	1	033708	2		
0.5 µl (M)	42 mm	23	0.63 mm	Cone	000410	1	033715	2		
1 µl	42 mm	23	0.63 mm	Cone	000610	1	037715	2		
5 µl	42 mm	26	0.47 mm	Cone	001805	1	036710	2		
5 µl (M)	42 mm	23	0.63 mm	Cone	001815	1	036720	2		
10 µl	42 mm	26	0.47 mm	Cone	002805	1	037715	2		
10 µl (M)	42 mm	23	0.63 mm	Cone	002815	1	037717	2		
10 µl Gas Tight	42 mm	26	0.47 mm	Cone	002817	1	037715	2		
10 µl Gas Tight	42 mm	23	0.63 mm	Cone	002818	1	038717	2		
25 µl (M)	42 mm	23	0.63 mm	Cone	003665	1	038717	2		
50 µl	42 mm	23	0.63 mm	Cone	004665	1	037717	1		
100 µl	42 mm	23	0.63 mm	Cone	005665	1	037717	1		
250 µl	42 mm	23	0.63 mm	Cone	006665	1	038717	1		

Replacement plunger in needle kit


Removable needle - Dual gauge

0.5 µl	42 mm	23-26	0.63 /0.47 mm	Cone	000415	1	033730	1		
5 µl	42 mm	23-26	0.63 /0.47 mm	Cone	001825	1	033730	1		
10 µl	42 mm	23-26	0.63 /0.47 mm	Cone	002825	1	033730	1		
10 µl Gas Tight	42 mm	23-26	0.63 /0.47 mm	Cone	002829	1	033730	1		

Replacement plunger in needle kit

Removable needle - SuperfleX™ - Flexible plungers

10 µl	42 mm	26	0.47 mm	Cone	002844	1	037715	1		
10 µl	42 mm	23	0.63 mm	Cone	002845	1	037717	1		



(M) : suitable w/Merlin Microseal™ injector

Consumables - Syringes


SGE Syringes

CTC/Leap

Volume	Length	Gauge	OD	Tip style	P/N	Qty	P/N	Qty
Fixed needle								
5 µl	50 mm	23	0.63 mm	Cone	001981	1		
5 µl	50 mm	26	0.47 mm	Cone	001982	1		
10 µl	50 mm	23	0.63 mm	Cone	002981	1		
10 µl	50 mm	26	0.47 mm	Cone	002980	1	002986	6
10 µl Gas Tight	50 mm	23	0.63 mm	Cone	002987	1		
25 µl	50 mm	26	0.47 mm	Cone	003980	1		
25 µl Gas Tight	50 mm	23	0.63 mm	Cone	003987	1		
100 µl Gas Tight	50 mm	23	0.63 mm	Cone	005335	1		
Removable needle					Replacement needle			
0.5 µl	50 mm	26	0.47 mm	Cone	000490	1	037785	2
0.5 µl	50 mm	23	0.63 mm	Cone	000492	1	037787	2
10 µl	50 mm	26	0.47 mm	Cone	002982	1	037785	1
10 µl	50 mm	23	0.63 mm	Cone	002984	1	033772	1
10 µl Gas Tight	50 mm	26	0.47 mm	Cone	002985	1	038749	2
25 µl Gas Tight	50 mm	26	0.47 mm	Cone	003988	1	033770	1
100 µl Gas Tight	50 mm	26	0.47 mm	Cone	005333	1	038250	2

PerkinElmer

Volume	Length	Gauge	OD	Typ Style	P/N	Qty	P/N	Qty
AutoSystem TM & Clarus 500								
Fixed needle								
5 µl	70 mm	26	0.47 mm	Cone	001953	1		
5 µl	70 mm	23	0.63 mm	Cone	001954	1		
5 µl Gas Tight	70 mm	26	0.47 mm	Cone	001955	1		
5 µl Gas Tight	70 mm	23	0.63 mm	Cone	001957	1		
5 µl	70 mm	23	0.63 mm	Cone	004670	1		
Removable needle					Replacement needle			
0.5 µl	70 mm	26	0.47 mm	Cone	000475	1	033750	1
5 µl	70 mm	23	0.63 mm	Cone	000478	1	033765	1



Consumables - Syringes


SGE Syringes

Shimadzu

Volume	Length	Gauge	OD	Tip style	P/N	Qty	P/N	Qty
AOC 14, 17, 20 & 20 i								
Fixed needle								
5 µl	42	26	0.47	Cone	001987	1		
5 µl	42	23	0.63	Cone	001988	1		
10 µl	42	26	0.47	Cone	002887	1		
Removable needle					Replacement needle			
0.5 µl	50	26	0.47	Cone	000435	1	033732	1
0.5 µl	42	26	0.47	Cone	000440	1	033738	1
0.5 µl	42	23	0.63	Cone	000445	1	033745	1
10 µl	50	26	0.47	Cone	002885	1	033732	1
10 µl	42	26	0.47	Cone	002897	1	033732	1
10 µl	42	23	0.63	Cone	002898	1	033745	1
10 µl Gas Tight	42	23	0.63	Cone	002902	1	033745	1

Thermo Finnigan

Volume	Length	Gauge	OD	Typ style	P/N	Qty	P/N	Qty
AS200, AS800 & AS2000								
Fixed needle								
5 µl	50 mm	23	0.63 mm	Cone	001981	1		
10 µl	50 mm	26	0.47 mm	Cone	002980	1	002986	6
10 µl	70 mm	26	0.47 mm	Cone	002978	1		
10 µl	80 mm	26	0.47 mm	Cone	002992	1		
Removable needle					Replacement needle			
10 µl	80 mm	26	0.47 mm	Cone	002993	1	031535	1



Consumables - Syringes

SGE Syringes

Varian

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty	P/N	Qty
Removable needle										
10 µl	50 mm	25	0.50 mm	Dome	002926	1				
8035, 8100 & 8200 Series										
Fixed needle										
10 µl	53 mm	25	0.50 mm	Side Hole	002923	1				
Removable needle										
							Needle		Replacement plunger	
1 µl	51 mm	26	0.47 mm	Cone	000655	1	034720	1		
10 µl	53 mm	25	0.50 mm	Side Hole	002924	1	037777	1	031218	1
10 µl	53 mm	25	0.50 mm	Side Hole	005921	1	038745	1	031824	1
8400										
Fixed needle										
10 µl	50 mm	26	0.47 mm	Bevel	002950	1				



002924

Consumables - Syringes


SGE Syringes

MicroVolume syringes (manual)

0.5 µl - 5 µl

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty	
Plunger in needle 6.5 mm OD barrel					Replacement plunger in needle kit				
0.5 µl	50 mm	23	0.63 mm	Cone	000300	1	033010	1	
0.5 µl	50 mm	23	0.63 mm	Bevel	000301	1	033011	1	
0.5 µl	50 mm	26	0.47 mm	Cone	000303	1	033012	1	
Plunger in needle 8 mm OD barrel					Replacement plunger in needle kit				
0.5 µl	70 mm	23	0.63 mm	Cone	000310	1	033057	1	
0.5 µl	70 mm	23	0.63 mm	Bevel	000311	1	033060	1	
1 µl	50 mm	23	0.63 mm	Cone	000500	1	034055	1	
1 µl	50 mm	23	0.63 mm	Bevel	000501	1	034056	1	
1 µl	70 mm	23	0.63 mm	Cone	000505	1	034057	1	
1 µl	70 mm	23	0.63 mm	Bevel	000506	1	034060	1	
1 µl	115 mm	23	0.63 mm	Cone	000510	1	034059	1	
1 µl	70 mm	26	0.47 mm	Cone	000570	1	034610	1	
5 µl	50 mm	23	0.63 mm	Cone	000800	1	035055	1	
5 µl	50 mm	23	0.63 mm	Bevel	000801	1	035056	1	
5 µl	70 mm	23	0.63 mm	Cone	000802	1	035057	1	
5 µl	70 mm	23	0.63 mm	Bevel	000803	1	035058	1	
5 µl	115 mm	23	0.63 mm	Cone	000804	1	035059	1	
5 µl	70 mm	23	0.63 mm	Nozzle	000805	1			
Plunger in needle 6.5mm OD barrel with repeating adaptor					Replacement plunger in needle kit				
0.5 µl	50 mm	23	0.63 mm	Cone	000350	1	033010	1	
0.5 µl	50 mm	26	0.47 mm	Cone	000353	1	033012	1	
Plunger in needle 6.5 mm OD barrel					Replacement plunger in needle kit				
0.5 µl	70 mm	23	0.63 mm	Cone	000355	1	033057	1	
1 µl	70 mm	23	0.63 mm	Cone	000550	1	034055	1	
1 µl	70 mm	23	0.63 mm	Cone	000553	1	034057	1	
5 µl	50 mm	23	0.63 mm	Cone	000850	1	035055	1	
5 µl	70 mm	23	0.63 mm	Cone	000852	1	035057	1	

Features



Plunger in needle 5 µl syringe with nozzle tip (P/N 000805)
 Nozzle tip with 0.12 mm opening and cone shape for better septum penetration.
 Full sample displacement
 Plunger extended to the tip of the needle - Zero dead volume
 No sample carry over
 Improved sample volume reproducibility

Consumables - Syringes


SGE Syringes

5 µl - 10 µl

Volume	Length	Gauge	OD	Tip style	P/N	Qty	P/N	Qty	P/N	Qty
Fixed needle Plunger protection syringes										
5 µl (M)	42 mm	23	0.63 mm	Cone	001839	1				
5 µl	50 mm	26	0.47 mm	Bevel	001000	1				
5 µl	70 mm	26	0.47 mm	Bevel	001002	1				
10 µl (M)	42 mm	23	0.63 mm	Cone	002839	1				
10 µl	50 mm	26	0.47 mm	Bevel	002000	1	002030	6	002033	25
10 µl	70 mm	26	0.47 mm	Bevel	002003	1				
10 µl	50 mm	26	0.47 mm	Cone	002005	1				
10 µl	115 mm	26	0.47 mm	Bevel	002007	1				
Fixed needle - Superflex™ plunger										
5 µl	50 mm	26	0.47 mm	Bevel	001100	1				
5 µl	50 mm	26	0.47 mm	Cone	001105	1				
10 µl	50 mm	26	0.47 mm	Bevel	002100	1	002130	6		
10 µl	70 mm	26	0.47 mm	Bevel	002101	1				
10 µl	50 mm	26	0.47 mm	Cone	002105	1	002133	6		
10 µl	50 mm	23	0.63 mm	Cone	002108	1	002135	6		
Fixed needle Guided plunger										
5 µl	50 mm	26	0.47 mm	Bevel	001400	1				
10 µl	50 mm	26	0.47 mm	Bevel	002400	1				
Fixed needle 6.5 mm OD barrel with repeating adaptor										
10 µl	50 mm	26	0.47 mm	Bevel	002040	1				
Removable needle Plunger Protection										
5 µl	50 mm	26	0.47 mm	Bevel	001050	1	036110	5		
10 µl	50 mm	26	0.47 mm	Bevel	002050	1	037110	5		
10 µl	70 mm	26	0.47 mm	Bevel	002055	1	037130	5		

Suitable for use with Merlin Microseal injector

Plunger protection.
prevents plunger bending



Injection volume :
never inject less than 10 %
of syringe capacity.
E.g. using a 10 µl
syringe, the sample should
be 1 µl or more to ensure
accuracy

Consumables - Syringes

SGE Syringes

5 µl - 10 µl

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty
Removable needle SuperfleX™ plunger								
5 µl	50 mm	26	0.47 mm	Bevel	001150	1	036110	5
10 µl	50 mm	26	0.47 mm	Bevel	002150	1	037110	5
Removable needle Guided plunger								
5 µl	50 mm	26	0.47 mm	Bevel	001450	1	036110	5
5 µl	70 mm	26	0.47 mm	Bevel	001455	1	036130	5
10 µl	50 mm	26	0.47 mm	Bevel	002450	1	037110	5
10 µl	70 mm	26	0.47 mm	Bevel	002453	1	037130	5
Removable needle Shout (Half scale) guided plunger								
5 µl	50 mm	26	0.47 mm	Bevel	001495	1	037110	5
Removable needle - 6.5 mm OD with repeating adaptor								
5 µl	50 mm	26	0.47 mm	Bevel	001090	1	036110	5
10 µl	50 mm	26	0.47 mm	Bevel	002090	1	037110	5
Removable needle Guided plunger with repeating adaptor								
5 µl	50 mm	26	0.47 mm	Bevel	001490	1	036110	5
10 µl	50 mm	26	0.47 mm	Bevel	002490	1	037110	5
10 µl	70 mm	26	0.47 mm	Bevel	002493	1	037130	5

25 µl - 500 µl

Volume	Length	Gauge	OD	Tip style	P/N	Qty	P/N	Qty
Fixed needle standard plunger								
25 µl	50 mm	25	0.50 mm	Bevel	003000	1		
50 µl	50 mm	25	0.50 mm	Bevel	004000	1		
100 µl	50 mm	25	0.50 mm	Bevel	005000	1		
250 µl	50 mm	25	0.50 mm	Bevel	006000	1		
500 µl	50 mm	25	0.50 mm	Bevel	007000	1		
Removable needle								
25 µl	50 mm	25	0.50 mm	Bevel	003050	1	038110	5
50 µl	50 mm	25	0.50 mm	Bevel	004050	1	038110	5
100 µl	50 mm	25	0.50 mm	Bevel	005050	1	038110	5
250 µl	50 mm	25	0.50 mm	Bevel	006050	1	038110	5
500 µl	50 mm	25	0.50 mm	Bevel	007050	1	038110	5
Removable needle 8 mm OD with repeating adaptor								
25 µl	50 mm	25	0.50 mm	Bevel	003090	1	038110	5
50 µl	50 mm	25	0.50 mm	Bevel	004090	1	038110	5
100 µl	50 mm	25	0.50 mm	Bevel	005090	1	038110	5
500 µl	50 mm	25	0.50 mm	Bevel	007090	1	038110	5



Repeating Adaptor.
For reproducible dispensed volumes.
Removes the guess work of lining up
the plunger.

Consumables - Syringes

SGE Syringes

Gas Tight syringes (manual)

10 µl - 500 µl

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty	P/N	Qty
Fixed needle										
10 µl	50 mm	26	0.47 mm	Bevel	002200	1	002202	6		
10 µl	70 mm	26	0.47 mm	Bevel	002208	1				
25 µl	50 mm	25	0.50 mm	Bevel	003200	1				
50 µl	50 mm	25	0.50 mm	Bevel	004200	1				
100 µl	50 mm	25	0.50 mm	Bevel	005200	1				
250 µl	50 mm	25	0.50 mm	Bevel	006200	1				
500 µl	50 mm	25	0.50 mm	Bevel	007200	1				



Fixed needle luer tip						Replacement plunger				
50 µl	-	-	-	-	004229	1			031820	1
100 µl	-	-	-	-	005229	1			031825	1
250 µl	-	-	-	-	006229	1			031830	1
500 µl	-	-	-	-	007229	1			031835	1



Fixed needle luer lock						Replacement plunger				
50 µl	-	-	-	-	004230	1			031820	1
100 µl	-	-	-	-	005230	1			031825	1
250 µl	-	-	-	-	006230	1			031830	1
500 µl	-	-	-	-	007230	1			031835	1

Fixed needle luer lock with plunger stop										
50 µl	-	-	-	-	004232	1				
100 µl	-	-	-	-	005232	1				
250 µl	-	-	-	-	006232	1				
500 µl	-	-	-	-	007232	1				



Removable needle						Replacement needle		Replacement plunger		
10 µl	50 mm	26	0.47 mm	Bevel	002250	1	037110	1	031811	1
25 µl	70 mm	25	0.50 mm	Bevel	003250	1	038110	1	031815	1
50 µl	50 mm	25	0.50 mm	Bevel	004250	1	038110	1	031820	1
100 µl	50 mm	25	0.50 mm	Bevel	005250	1	038110	1	031825	1
250 µl	50 mm	25	0.50 mm	Bevel	006250	1	038110	1	031830	1
500 µl	50 mm	25	0.50 mm	Bevel	007250	1	038110	1	031835	1



Removable needle guided plunger						Replacement needle		Replacement plunger		
10 µl	50 mm	26	0.47 mm	Bevel	002455	1	037110	1	031805	1

Removable needle guided plunger										
10 µl	75 mm	-	0.23 mm	Flat	002500	1				

Consumables - Syringes

SGE Syringes

LC autosampler syringes

Agilent Technologies

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty
1090 & 1100							Replacement plunger	
25 µl	-	-	-	-	003670	1	0318181	1
50 µl	-	-	-	-	006670	1	031829	1



CTC

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty
Fixed needle							Replacement plunger	
10 µl	2"mm	22	0.028"mm	LC	002991	1	031120	2
25 µl	2"mm	22	0.028"mm	LC	003984	1		1
50 µl	2"mm	22	0.028"mm	LC	004810	1	031821	1
100 µl	2"mm	22	0.028"mm	LC	005331	1	0318261	1



Removable needle							Replacement plunger	
25 µl	2"mm	22	0.028"mm	LC	003985	1	031816	1
100 µl	2"mm	22	0.028"mm	LC	005330	1	031826	1
250 µl	2"mm	22	0.028"mm	LC	006330	1	031831	1

Hitachi

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty
							Replacement plunger	
500 µl	-	-	-	-	007660	1	031837	1



PerkinElmer

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty
							Replacement plunger	
50 µl	-	-	-	-	004995	1	0318221	1
100 µl	-	-	-	-	005990	1	0318271	1
250 µl	-	-	-	-	006995	1	031833	1
500 µl	-	-	-	-	007995	1	0318381	1
1 ml	-	-	-	-	008185	1	0318441	1



Consumables - Syringes

SGE Syringes

Thermo Finnigan

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty
250 μ l	-	-	-	-	006660	1	Replacement plunger 031833	1
500 μ l	-	-	-	-	007680	1	031838	1
1 ml	-	-	-	-	008180	1	031844	1
2.5 ml	-	-	-	-	008660	1	031853	1
5 ml	-	-	-	-	008780	1	031856	1

Waters WISP®

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty
25 μ l	-	-	-	-	003990	1	Replacement plunger 031819	1
250 μ l	-	-	-	-	006690	1	031834	1



Consumables - Syringes

SGE Syringes

LC manual syringes

Beckman / Altex, Rheodyne, SSI, Valco Valves

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty	
Fixed needle						Replacement plunger			
5 µl	2" mm	22	0.028" mm	LC	001301	1			
10 µl	2" mm	22	0.028" mm	LC	002301	1			
10 µl	2" mm	22	0.028" mm	LC	002335	1	031810	2	
25 µl	2" mm	22	0.028" mm	LC	003300	1			
50 µl	2" mm	22	0.028" mm	LC	004300	1			
100 µl	2" mm	22	0.028" mm	LC	005300	1			
250 µl	2" mm	22	0.028" mm	LC	006300	1			
500 µl	2" mm	22	0.028" mm	LC	007300	1			

Fixed needle - SuperfleX™ plunger

10 µl	2" mm	22	0.028" mm	LC	002300*			
-------	-------	----	-----------	----	---------	--	--	--

Removable needle

						Replacement needle		Replacement plunger	
5 µl	2" mm	22	0.028"mm	LC	001310	1	035250	5	
10 µl	2" mm	22	0.028" mm	LC	002310	1	037250	5	
10 µl	2" mm	22	0.028" mm	LC	002313	1	037250	5	031811 2
25 µl	2" mm	22	0.028" mm	LC	003310	1	038250	5	
25 µl	2" mm	22	0.028" mm	LC	003312	1	038250	5	031815 1
50 µl	2" mm	22	0.028" mm	LC	004310	1	038250	5	
50 µl	2" mm	22	0.028" mm	LC	004312	1	038250	5	031820 1
100 µl	2" mm	22	0.028" mm	LC	005310	1	038250	5	
100 µl	2" mm	22	0.028" mm	LC	005312	1	038250	5	031825 1
250 µl	2" mm	22	0.028" mm	LC	006310	1	038250	5	
250 µl	2" mm	22	0.028" mm	LC	006312	1	038250	5	031830 1
500 µl	2" mm	22	0.028" mm	LC	007310	1	038250	5	
500 µl	2" mm	22	0.028" mm	LC	007312	1	038250	5	031835 1
1 ml	2" mm	22	0.028" mm	LC	008105	1	039250	5	031842 1
2.5 ml	2" mm	22	0.028" mm	LC	008505	1	039250	5	031852 1

Removable needle - SuperfleX™ plunger

10 µl	2" mm	22	0.028" mm	LC	002350			
-------	-------	----	-----------	----	--------	--	--	--

Removable needle with repeating adaptor

						Replacement needle			
10 µl	2" mm	22	0.028" mm	LC	002345	1	037250	5	
25 µl	2" mm	22	0.028" mm	LC	003320	1	038250	5	
50 µl	2" mm	22	0.028" mm	LC	004320	1	038250	5	
100 µl	2" mm	22	0.028" mm	LC	005320	1	038250	5	
250 µl	2" mm	22	0.028" mm	LC	006320	1	038250	5	
500 µl	2" mm	22	0.028" mm	LC	007320	1	038250	5	



Consumables - Syringes

SGE Syringes - LC manual syringes

Waters U6K Valves

Volume	Length	Gauge	OD	Tip Style	P/N	Qty	P/N	Qty
Removable needle with repeating adaptor							Replacement needle	
25 µl	50 mm	25	0.50 mm	LC	003375		038260	5
50 µl	50 mm	25	0.50 mm	LC	004375		038260	5
100 µl	50 mm	25	0.50 mm	LC	005375		038260	5
250 µl	50 mm	25	0.50 mm	LC	006375		038260	5
500 µl	50 mm	25	0.50 mm	LC	007375		038260	5

