



## Certificate of Analysis / Data Sheet

**Product:** Alpha-Amylase

**Catalog No:** BTNM-0010

**Size:** 25g

**Lot No:** N/A

**Release Date:**

**Retest Date:** 3 years from release date

**Storage:** Room Temp.

### Specifications

**CAS Number:** N/A

**Appearance:** 9000-90-2

| Test            | Specification   | Result                   |
|-----------------|---|--------------------------|
| Source          | <i>Bacillus subtilis</i>  | <i>Bacillus subtilis</i> |
| Identity Test   | PASS  | PASS                     |
| Activity        | ≥ 50,000 BAU/g  | PASS                     |
| Unit Definition | One unit will liberate 1.0mg of Maltose from Starch in three minutes. | -----                    |

This document has been produced electronically and is valid without signature

**Date of Print:** February 8, 2022

**G-Biosciences/Geno Technology, Inc.**

Rev 02/08/22 DB

