

CERTIFICATE OF QUALITY

Product: BETA-GALACTOSIDASE, E. COLI
Grade: CONJUGATION GRADE

Code: 0641

Appearance: White flaky powder

Mol. Formula: Not Applicable

Mol. Weight: Not Applicable

Description:

<u>Assay:</u>	<u>Min</u>	<u>Max</u>	<u>Units</u>
Activity (Protein) @37C	600		u/mg
Aggregation		3.0	%

Storage: Store frozen. Keep dry. Warm to room temperature before opening.