
BCIP/INT solution

IHC / Blotting Alkaline Phosphatase substrate

Product Number UP096051, 100ml

The Uptima BCIP/INT solution substrate is an one component chromogenic substrate for immobilized Alkaline Phosphatase (AP). The BCIP/INT is an excellent choice for *in situ* hybridizations, when tissues will be counterstained with hematoxylin. Its brick-red end product is clearly visible against a blue hematoxylin background, avoiding the misleading results that can arise from observing NBT formazan under these conditions.

This substrate is ready to use, straight from the bottle, with no dilutions or mixing.

Sensitivity Sensitivity is excellent — in immunoblotting assays, BCIP/INT solution will detect sub-nanogram levels of mouse IgG.

Storage Stored cold (2-8°C). (L)
The BCIP/INT solution substrate is stable 1 year from date purchase, and 6 months after opening

Directions for use

Protocol for Western Blot

1. Wash the blot to remove unbound antibodies.
Note: *Insufficient washes may lead to background.*
Remove the excess of buffer
2. Add ready-to-use BCIP/INT solution to cover well the blot.
Note: *Insufficient volume may cause background on some area of the blot.*
3. Incubate 10-15 min at room temperature or until desired stain development.
4. Rinse the blot with water. Allow the blot to dry to the air.

Other information

Related / associated products and documents

See [BioSciences Innovations catalogue](#) and [e-search tool](#).
[page A352+](#)

saturation agents for immunoassays
HRP labeled secondary antibodies

Streptavidin-AP [#UP39588](#)

[SeaBlock](#) , BSA solution [#UP900100](#)
Anti Human IgG – AP [#UP893090](#)
Anti Mouse IgG – AP [#UP764810](#)
Anti Rabbit IgG – AP [#UP659711](#)

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. Uptima is not liable for any damage resulting from handling or contact with this product. rev. : 103E-B06VE

Contact your local distributor



Uptima, powered by

213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

Uptima@interchim.com