



## Dimethyl Formamide (DMF)

### Product Description

<b>Name :</b>	<b>Dimethyl Formamide (DMF), High purity Grade</b>
<b>Catalog Number</b>	12137K, 1L 12137L, 2,5L
<b>CAS:</b>	68-12-2
<b>Molecular Weight :</b>	115 g/mol
<b>Storage:</b>	Room temperature / z

### Introduction

Dimethyl Formamide (DMF) is used in Peptide synthesis having a very low amine content suitable for HPLC , GC analytical purposes.

### Specifications

Exp date: 2 years from DOM

<u>Test</u>	<u>Specification</u>
Purity (%):	>99,95
Water (%):	<0,030
Dimethylamine (ppm)	<10
Formic Acid (ppm)	<3
Methanol (%):	<0,0050

### Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.  
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

[Order on-line](#) or [Contact](#) your local distributor

**Disclaimer :** Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses.  
Uptima is not liable for any damage resulting from handling or contact with this product.

Rev.S08vb