

CA-14345I

Nitro-Blue Tetrazolium Chloride

Product Information

Chemical name : Nitro-Blue Tetrazolium Chloride (NBT)


Cat. Number : UP14345I, 500g

Mol. Formula : C₄₀H₃₀N₁₀O₆C₁₂

Mol. Weight : 817.65 g/mol

CAS : 298-83-9

Storage : Store +2-8°C. Protect from light.
-20°C for long term ^(M)



Test	Specification
------	---------------

Appearance :	Yellow crystalline powder
---------------------	---------------------------

Em (λ max, MetOH) :	≥ 54000
-----------------------------	---------

λ max, MetOH :	255-260 nm
-----------------------	------------

Loss on drying :	≤ 10 %
-------------------------	--------

Solubility: (1%, water):	soluble
---------------------------------	---------

For any information, please contact / Pour toute information complémentaire, s'adresser à :
Uptima -Interchim, Division BioSciences uptima@interchim.com; Hotline : +33 (0)4 70 03 73 06

Contact your local distributor

uptima@interchim.com

Uptima, powered by

 213 Avenue J.F. Kennedy - BP 1140
 93103 Montluçon Cedex - France
 Tél. 04 70 03 88 55 - Fax 04 70 03 82 60