



211 bis av. J.F. Kennedy - B.P. N° 1140
03103 Montluçon Cedex - France
☎(33) (0)4 70 03 88 55 (+)
Fax (33) (0)4 70 03 82 60
interbiotech@interchim.com



ISO 9001:2000 CERTIFIED

CERTIFICATE OF QUALITY

Product: LDH PIG HEART
Grade: REAGENT GRADE
Code: 0253
Appearance: Clear to slightly cloudy liquid
Mol. Formula: Not Applicable
Mol. Weight: Not Applicable

Description:

<u>Assay:</u>	<u>Min</u>	<u>Max</u>	<u>Units</u>	<u>Disposition</u>
Concentration	5000		U/ml	PASS
GOT		0.04	%	PASS
GPT		0.04	%	PASS
MDH		0.01	%	PASS
Specific Activity	300		U/mg	PASS

Storage: Store cold. Keep dry. Warm to room temperature before opening.

0253RG

This signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature:

Title: QC MANAGER Date: 1/30/2009