
Specification Sheet / Fiche de spécification

UP36859 Bovine Albumin standard

Product Information / Informations Produit

Part Number / référence : UP36859A, 10 x 1ml UP36859D, 30ml

Product / Produit : Bovine Albumin Serum

Form: Solution at 2mg/ml in distilled water, 0.05% NaN₃
0.22 µm filtrated

Applications: Standard for quantification of protein by BC Assay, MicroBC Assay, COO Assay

Storage / Conservation: +4°C protect from light(L) – best use before 6 months
(store at -20°C for longer period – avoid freeze-thaw cycles)
(may be stored at room temperature for use within 2 weeks)

Quality control / Contrôle qualité

Lot :

Test	Specification / Spécification
Appearance / Aspect	Limpid Solution, free of particle
Total proteins / Protéines totales (colorimetric BC Assay)	2 +/- 0.08 mg/ml
Abs. (OD280nm)	1.30-1.45

We certify that the bovine raw material (plasma) used for the manufacture of this lot, is of french origin.
Nous garantissons que la matière première d'origine bovine (plasma) utilisée pour la fabrication de ce lot est d'origine française.
A veterinary and sanitary certificate is available on inquire.

For any information, please ask :
Pour toute information complémentaire, s'adresser à :
Uptima / Interchim; Hotline : +33(0)4 70 03 73 06
Uptima@interchim.com ; www.interchim.com



Visa Quality Uptima

Date :