

FT-80146A

interchim®

INTERFINE CHEMICALS ANALYTICAL SCIENCES BIOCHROMATOGRAPHY LIFE SCIENCES

+33 4 70 03 73 06 interbiotech@interchim.com

Big Endothelin-1 22-38 Human

Product Description

Catalog #: 80146A, 0.5mg
Name: **Big Endothelin-1 22-38 Human (Big ET-1)**
 MW : 1806.1
 $C_{80}H_{125}N_{23}O_{25}$
 Val-Asn-Thr-Pro-Glu-His-Val-Val-Pro-Tyr-Gly-Leu-Gly-Ser-Pro-Arg-Ser
Stability : Solution is stable for 1-2 weeks when stored at -20°C. Prepare solution just before use.
Storage: -20°C
 For Research Use Only

Directions for use

References

- **Ahlborg G. et al.**, Central and regional hemodynamic effects during infusion of Big endothelin-1 in healthy humans, *J Appl Physiol*, Jun 1996; 80: 1921 - 1927. [Abstract](#)
- **Itoh Y. et al.**, Cloning and sequence analysis of cDNA encoding the precursor of a human endothelium-derived vasoconstrictor peptide, endothelin: identity of human and porcine endothelin, *FEBS Lett.* 1988 Apr 25;231(2):440-4. [Abstract](#)
- **Klaar U. et al.**, Prognostic value of serial B-type natriuretic peptide measurement in asymptomatic organic mitral regurgitation, *Eur J Heart Fail*, Feb 2011; 13: 163 - 169. [Article](#)

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.
 Please inquire for higher quantities (availability, shipment conditions).

Please contact InterBioTech – Interchim for any other information
 Hotline : +33(0)4 70 03 73 06 – Interbiotech@interchim.com

Disclaimer : Materials are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. Interchim® is not liable for any damage resulting from handling or contact with this product.