

## TG-SDS buffer 40X Powder

*Used for Western blotting and gel electrophoresis*

### Product Description

<b>Name:</b>	<b>40x Tris Glycine SDS Powder (1 sachet makes 40L)</b>
	One sachet contains enough material to formulate 1liter of 40X Tris Glycine SDS (TGSDS) solution or 40 liters of 1x TGSDS (working strength) solution.
<b>Cat.number:</b>	<b>865984, 40L (1L à 40X)</b>
<b>Purity:</b>	Tris >99.9% Glycine >99.9% SDS >99.9%
<b>Storage:</b>	Room temperature (powder) +2 to 8°C (solution)

### Applications

Tris-Glycine SDS Running Buffer (1X) is formulated for separation of proteins, in their denatured state, for use with Tris-Glycine gels. The sachet is a convenient and quick method for preparing the correct formulation of buffer.

### Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.

For any information, please ask : Biosciences / Interchim; Hotline : +33(0)4 70 03 88 55

[Order on-line](#) or [Contact](#) your local distributor.

**Disclaimer :** Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product.