

TG-SDS Buffer 10X

Used for SDS-PAGE (protein polyacrylamide gel electrophoresis)

Product Information

Name:	TG-SDS 10X (0.25M Tris, 1.92M Glycine, 1% SDS)
Chemical description:	When reconstituted with water, TG-SDS 1X final concentrations are 0.025M Tris, 0.192M Glycine, 0.1% SDS, pH 8.30-8.70.
Cat.number:	UP914959, 1L to make 10L UP91495E, 5L to make 50L
Form:	Liquid - Dissolve with distilled water
Storage:	Room temperature (Z)
Applications:	TG-SDS is used in various laboratory techniques, including electrophoresis...

Other information

For any information, please ask : Biosciences / Interchim; Hotline : +33(0)4 70 03 88 55
[Order on-line](#) or [Contact](#) your local distributor.