



## HRP Donkey Anti-Mouse IgG (H+L) (cross adsorbed)

### Product information

<b>Name</b>	HRP Donkey Anti-Mouse IgG (H+L) min X Bov,Ck,Gt,GP,Hms,Hrs,Hu,Rb,Rat,Shp Sr Prot)
<b>Cat. number</b>	UP98432E
<b>Concentration</b>	1mg/ml
<b>Quantity</b>	50 µg
<b>Physical state</b>	Freeze dried powder
<b>Buffer</b>	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA , Protease IgG free
<b>Preservative</b>	0.1% (v/v) Kathon CG
<b>Reconstitution</b>	Rehydrate with 55 µl distilled water. Centrifuge to remove any particulates.
<b>Working solution</b>	1:20-1:2,000 for most applications
<b>Storage</b>	Store freeze-dried product at 2-8°C. For long term storage after reconstitution, dilute with 50% glycerol and store at -20°C.
<b>Expiry date</b>	One year from date of reconstitution
<b>Specificity</b>	<p><b>Based on IEP, this antibody reacts with:</b></p> <p>Heavy (g) chains on mouse IgG Light chains on all mouse Immunoglobulins</p> <p><b>Based on IEP no reactivity is observed to</b></p> <p>Non-immunoglobulin mouse serum immunoglobulins IgG from bovine, chicken, goat, guinea Pig, hamster, horse, human, rabbit, rat and sheep.</p>

**For research use only**