



Small, but convenient!

# Mini Cell Strainers

For more information : [http://www.funakoshi.co.jp/exports\\_contents/80727](http://www.funakoshi.co.jp/exports_contents/80727)



Small and convenient cell strainers, which fit for 1.5 / 2.0 ml tubes, 5 ml sampling tubes and 15 ml conical tubes.

## Mini Cell Strainer (for 5 ml sampling tubes and 15 ml conical tubes)



5 ml tube    15 ml tube    5ml tube



25  $\mu$ m    40  $\mu$ m    70  $\mu$ m



## Features

- ◆ Facilitate speedy cell separation from primary tissues or cell clumps
- ◆ Slits on the rim avoid overflow and help speedy filtration
- ◆ Fit for 5 ml sampling and 15 ml conical tubes
- ◆ Color-coded according to mesh pore sizes (25  $\mu$ m : White, 40  $\mu$ m : Blue, 70  $\mu$ m : Clear)
- ◆ Materials : Frame ; PP, Mesh ; Nylon
- ◆ EOG sterilized and individually wrapped
- ◆ Made in Japan

# Mini Cell Strainer II (for 1.5 / 2.0 ml or 5 ml (φ12 × 75 mm) tubes)



1.5 ml tube



2 ml tube



5 ml tube  
( φ12×75 mm )



## Features

- ◆ For prefiltration of samples before flow cytometry or mass spectrometry
- ◆ Slits on the rim avoid overflow and help speedy filtration
- ◆ Fit for 1.5 / 2.0 ml tubes (φ12 mm) and 5 ml (φ12 x 75 mm) sampling tubes
- ◆ EOG sterilized and individually wrapped
- ◆ Made in Japan

## Product Information

### ■ Mini Cell Strainer

[ Manufacturer : HIT ]

Product Name	Size	Code
Mini Cell Strainer, 25 μm, White	100 pieces	HT-AMS-12502
Mini Cell Strainer, 40 μm, Blue	100 pieces	HT-AMS-14002
Mini Cell Strainer, 70 μm, Clear	100 pieces	HT-AMS-17002

### ■ Mini Cell Strainer II

[ Manufacturer : HIT ]

Product Name	Size	Code
Mini Cell Strainer II, 40 μm, Blue	100 pieces	HT-AMS-04002
Mini Cell Strainer II, 40 μm, Blue with 5 ml (φ12 x 75 mm) tube set (100 pieces each)	1 set	HT-AMS-04008

#### NOTE

\* All products shown here are for research use only, not for diagnostic use.  
\* Products might be changed for improvement without notice.

\* All brand and product names are registered as trademarks of their respective companies.  
\* Please contact your local distributors for orders, quote request and any inquiry.

Your Local Distributor



funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,  
Tokyo 113-0033, JAPAN  
Phone: +81-3-5684-6296  
Fax: +81-3-5684-6297  
Email: export@funakoshi.co.jp