



Tris/HCl, 1.0M Solution, pH 7.5

Product Description

Name :	Tris/Hcl 1M solution pH 7.5
Catalog Number :	Molecular Biology Grade AYO67A, 100ml AYO67B, 500ml AYO67C, 1L
Structure :	CAS: [1185-53-1]

Storage: Room temperature (z)

Directions for use

Description

1M solution of Tris in high purity dH₂O, adjusted to pH 7.5. Solution is 0.2 µm filtered and dispensed into durable, square Nalgene bottles.

Contaminant Testing:

DNase (endo): Passes test
 DNase (exo): Passes test
 RNase: Passes test
 Protease: Passes test

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.
 Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

[Order on-line](#) or [Contact](#) your local distributor

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses.
 Uptima is not liable for any damage resulting from handling or contact with this product.