

DBCO-PEG4-vc-PAB-MMAE

Cat. No: ADCsR-053CL

PRODUCT INFORMATION

Product Overview

This set is made by MMAE conjugated to DBCO-PEG4-vc-PAB linker. Monomethyl auristatin E (MMAE), a potent tubulin inhibitor, is a toxin payload in antibody drug conjugate.


Product Type

Linker-Toxin SET

Classification

Click chemistry Reactive LTK

Structure



Molecular Weight

1658.1

Purity

95%

For research use only. Not intended for pharmaceutical development, diagnostic, therapeutic or any in vivo use. No products from Creative Biolabs may be resold, modified for resale or used to manufacture commercial products without prior written approval from Creative Biolabs.