


ProChambers™ MultiWell Chambers

- Uses standard slides
- Easy loading by multi-channel pipettes or automated Robots
- Utilizes, leak free, compression fit gasket
- Embossed with letter and number panels for easy identification for each subarray
- Durable construction
- Will not break slides


16 Well ProChamber Microarray System



- 16 (7mm x 7mm) printing surfaces
- rectangular openings

Assembly includes lower carrier tray, upper structure plate, slide, printing Template and gasket (patent pending)

24 Well ProChamber Microarray System



- 24 (6 mm x 7mm) printing surfaces
- rectangular openings

Assembly includes lower carrier tray, upper structure plate, slide, printing Template and gasket (patent pending)