

High Capacity Streptavidin Magnetic Beads

Product No.:	1497
Product name:	High Capacity Streptavidin Magnetic Beads, 1 μ m
Loading capacity:	≥ 12 nmol free biotin per milligram of beads
Formulation:	10 mg/mL suspension in water
Preservative:	0.05% sodium azide
Support:	Highly uniformed supermagnetic particles
Coating:	Highly crosslinked, proprietary (non-styrene) polymer matrix
Mean bead diameter:	0.97 μ m
Surface area per mg:	5.7×10^9 μ m ² /mg
Density:	1.33 g/cm ³
Storage:	Upon receipt store at 2-8°C. Do not freeze
Shipping:	Product shipped at ambient temperature