




H10E

Interchim BioSciences innovation proposes a complete range of other products for electrophoresis analysis

Electrophoresis chart



See below [info & ordering](#) of these **Protein Electrophoresis great tools** and the [Protein Electrophoresis catalog](#). See also the [Genomics catalog](#) (electrophoresis of nucleic acids).

Electrophoresis tools overview

GebaGel electrophoresis system

An economic, easy-to-use horizontal electrophoresis system with resolutive & cost-effective pre-cast gels

Product	Comb color	Cat.#	Unit
10% Protein GeBaGel /Proteins	Blue	BI9690	8u
12% Protein GeBaGel /Proteins	Green	BI9700	8u
4-12% Protein GeBaGel /Proteins	Violet	BI9710	8u
8-16% Protein GeBaGel /Proteins	Red	BI9720	8u
17% Protein GeBaGel /Peptides	Black	BI9730	8u
GEBARUNNER Electrophoresis system		BI9740	1u

Also available with 1, 2, and 15 wells

[More information](#)



NEXT gel technology

The new generation of ultra-resolutive gels, ultra-rapid, ready in ONE min !


Description	Fine resolution range	cat.# / 100ml	cat.#/500ml (1)
NEXT Gel™ 5%	50-200 KDa	GS4270	GS4271
NEXT Gel™ 7.5%	20-100 KDa	GS4280	GS4281
NEXT Gel™ 10%	10-70 KDa	BG6290	BG6291
NEXT Gel™ 12.5%	10-50 KDa	GS4290	GS4291
NEXT Gel™ 15%	5-40 KDa	GS4300	GS4301
NEXT Gel™ Running Buffer, 20X		GS4310	GS4311
LP-NEXT Gel™	1-3 000KDa	BI6150 / 1kit (3)	

T08860(1) Each package of 500ml contains reagents for 50 mini and 20 regular gels.

(2) This ready to dilute buffer will achieve high resolution on a wide range of protein fragment sizes.

(3) This Kit includes : NEXT Gel™ sample buffer 4X ,NEXT Gel™ running buffer, 20X Agarose High Resolution Protein (HRP)

[More information](#)



[Prices and technical sheets on-line](#)

Staining for proteins in gel

The best of stains, whatever you work with coomassie, silver or by fluorescence:

Description	cat.#
CooBlue FX Instant Protein Stain Quick staining, safe/no MetOH or acid, 20ng sensitivity	UPG4562A, 500ml UPG4562B, 4.5L
ProSave™ Protein Gel stain Ultrafast staining (7min!), reversible and compatible/any appl.	BP7121, kit
SilverBullit Very sensitive and convenient, compatible with MS	T08860, 1 kit (qsp 2.5L)
LavaPurple Protein Gel AND Blot stain fluorescent staining, the most sensitive (50pg); identifies more proteins by LC/MS, very economic (9€/gel)!	67433A, kit 20 minigels 67433B, kit 100 minigels 67433C, kit 20 minigels

[More information](#)

[Prices and technical sheets on-line](#)

Preparative Electrophoresis: OFF Gel

Ideal to fractionate proteins according their pI from complex samples

Resolution 0.1-06 pH- load µg to mg

Increase considerably the number of proteins identified by LC/MS

[More information](#)

3100 OFFGEL system

CF3900-G3100AA , 1u



Other products/documents

[Protein Electrophoresis catalog](#)

[Products HighLights Overview](#)

Information inquire

Reply by Fax : +33 (0) 4 70 03 82 60 or email at interbiotech@interchim.com

I wish to receive the complete documentation about: _____

Name: _____ 2nd name: _____ Position: _____

Company/Institute: _____ Service, Lab: _____

Adress: _____

Zip code: _____ Town: _____

Tel _____ Fax _____ Email: _____