

Interchim Innovations

Interbiotech - BioScience Innovations

K03E

■ Laboratory Labware Organizers

Drawer/Ganizers® and Refrig/Arrangers®

Drawer/Ganizers® Spacers are designed to maximize laboratory working space while increasing safety. The durable, light-weight polystyrene Spacers fit in standard size lab drawers with separate compartment designations that organize and protect supplies commonly found in the lab, as well as reduce breaks and spills that pose potential hazards in the laboratory setting.

- Protect supplies
- Conserve space
- Increase handling efficiency
- Transport
- Promote safety

Go to [Drawer/Ganizer® Organizer Trays](#)
[Refrig/Arranger® Storage Trays](#)



Drawer/Ganizer® Portable Organizer Trays

Spectrum introduces new Drawer/Ganizer Portable Trays molded from 1/8 inch thick polystyrene plastic with easy gripping folded edges for transporting assorted lab items and mobilizing experimental activities. The deeper compartments are better suited than most other lab trays for holding larger quantities as well as serving as secondary containment for spills. Like the original Drawer/Ganizer Drawer Spacers, these Portable Organizer Trays are designed to fit standard lab drawers, refrigerators and bench tops to organize and protect supplies, conserve valuable space and promote lab safety. Spectrum offers 8 functional configurations (1 tray/pkg).

- Organize & Protect Lab Supplies
- Designed for Portability & Storage
- Stronger & More Durable
- Promotes Lab Safety
- Provides Secondary Containment



9 Portable Organizer Tray Configurations:

Pipette / Column Tray

Two sizes for sorting 4 or 6 types of pipettes, chromatography columns or long instruments in a protected parallel position.



Stopper / Fitting Tray

Great for separating large quantities of different types of small items. 2 sizes available with 12 or 16 equal compartments.



Parts / Chrom Tray

Various size compartments for easy to find storage of lab & chromatography supplies like utensils, valves, connectors, etc.



Micropipettor Tray

2 sizes designed for storing up to 7 micropipettors, pipettes and/or glass syringes as well as 2 types of pipette tips.



Clamp / Instrument Tray

Two tray sizes available, each with 3 compartment sizes for storing clamps, hardware, instruments and other odd-shaped items.



Utility / Equipment Tray

Ideal for the separation of bulky and/or fragile items that can otherwise be damaged. 2 sizes available, each with 4 spaces.



Grid / Vial Tray

Secures and protects up to 99 vials or fragile items in upright position. 1 tray size with 9 rows of 11 spaces, 1 inch deep.



Beaker / Fleaker Tray

Round depressions with spout cut-outs for secure storage of 50-600 ml beakers. Great protection for drying glassware.



Ordering Information: Drawer/Ganizer Portable Organizer Trays (1/pkg)

Part No.	Tray Configuration	Size	Dimensions (cm)	No. of Comp.	Comp. Dim. (cm)
155510	Pipette / Column	Large	49 x 44 x 5	6	42 x 7
155585		Small	34 x 44 x 5	4	"
155495	Stopper / Fittings	Large	49 x 44 x 7	16	11 x 10
155570		Small	34 x 44 x 7	12	10 x 10
155450	Parts /	Large	49 x 44 x 7	15	various
155645	Chromatography	Small	34 x 44 x 7	7	"
155430	Micropipettor	Large	49 x 44 x 7	7, 2	27 x 6, 23 x 23
155675		Small	34 x 44 x 7	5, 2	27 x 6, 16 x 14
155525	Clamp / Instrument	Large	49 x 44 x 7	7	various
155600		Small	34 x 44 x 7	5	"
155465	Utility / Equipment	Large	49 x 44 x 7	4	23 x 20
155630		Small	34 x 44 x 7		20 x 15
155665	Grid / Vial	Small	34 x 44 x 6	70	3.1 x 3.1
155480	Beaker / Fleaker	Large	49 x 44 x 5	27	5.6, 7.5, 9.4 dia.

[Prices and technical sheets on line](#)

Refrig/Arranger[®] Storage Trays

Refrig/Arranger Trays are the ideal way to conserve valuable space in laboratory refrigerators and freezers as well as organize laboratory bench tops and cabinets. Made of high-impact polystyrene, Refrig/Arranger Trays are also well suited for freezer storage (-90°C).

These lightweight, sturdy Storage Trays are perfect for safe and convenient transportation and refrigerated storage of glass vials or small bottles.

Refrig/Arranger Trays are available in channel and grid designs with a raised front edge designed for labeling and preventing spills.

- Efficient Use of Valuable Space
- Protects Glass & Plastic Containers
- Transport from Bench-top to Fridge
- Secondary Containment for Spills
- Freezer Storage as low as -90 C
- 3 Channel Tray Sizes
- 2 Grid Tray Sizes



Ordering Information: Refrig/Arranger[®] Storage Trays

Part No.	Tray Configuration	Size	Dimensions (cm)	No. of Comp.	Comp. Dim. (cm)	Qty
183552	Channel Tray	Large	34.3 x 40	7	4.8 (wide)	1
183550		Medium		10	3.2 (wide)	
183547		Small		13	1.9 (wide)	
183554	Grid Tray	Large	34.3 x 40	60	3.2 x 3.2	1
183556		Small		104	1.9 x 1.9	

[Prices and technical sheets on line](#)

Related products/documents

[Products HighLights Overview](#)

[BioSciences innovation catalogue](#)

Information inquire

Reply by Fax : +33 (0) 4 70 03 82 60 or email at interbiotech@interchim.com

I wish to receive the complete documentation about: _____

Name: _____ 2nd name: _____ Position: _____

Company/Institute: _____ Service, Lab: _____

Adress: _____

Zip code: _____ Town: _____

Tel _____ Fax _____ Email: _____