

**ImmunoBioScience CORP (IBS)****UNIVERSAL ANTIBODY DILUTION BUFFER™**

|                        |            |         |              |
|------------------------|------------|---------|--------------|
| <b>Catalog number:</b> | AR-6562-02 | 100 ml  | Ready to use |
|                        | AR-6562-03 | 250 ml  | Ready to use |
|                        | AR-6562-05 | 1000 ml | Ready to use |

**Description:** Universal antibody dilution buffer is ready to use for dilution of antibodies in all immunoassays (Immunofluorescence, IHC, ELISA, WB). This buffer does not contain any mammalian proteins, phosphate, sodium azide or mercury preservative and can be used for dilution of all antibodies, including peroxidase, and antibodies to phosphoproteins. Not suitable for dilution of antibodies to S100 proteins.

**Reagent:** Buffer ready to use.

**Storage:** 2-8°C.

**Application:** Dilute antibodies, primary or secondary as recommended in the protocol by the manufacturer of these products. The diluted antibodies can be stored at 2-8°C for several weeks depending upon the stability of antibodies.

-----

**Limitation and warranty:** Our warranty is limited to the actual price paid for the product. We are not liable for any property damage, personnel injury, time, effort or economic loss due to our product.

**MSDS:** Avoid skin and eye contact with all laboratory products. Use appropriate lab. gear, lab coat, gloves and safety glasses. Do not ingest any lab. products. This product is not approved for administration in human or animals.

-----

Universal antibody dilution buffer is a trade mark of ImmunoBioScience corp.  
For research use only, not for use in diagnostic or therapeutic procedures or in vivo use