

FT-DY8950



## 1% w/v Sodium Azide Solution

### Product Description

<b>Name :</b>	<b>1% w/v Sodium Azide Solution</b>
<b>Catalog Number :</b>	DY8950, 100ml
<b>Structure :</b>	N <sub>3</sub> Na CAS [26628-22-8] MW : 65,01 g/mole

**Storage:** 2-8°C

Sodium Azide solution is widely used as a preservative and deemed to be less hazardous when used as a solution. Suitable for Gel Filtration Chromatography , Affinity Chromatography, Anti-body preservation , Molecular biology, Cell Biology and General Biological use.

1% w/v Sodium Azide Solution is generally used and diluted to 0.01% w/v strength with deionised RO purified water, it is also 0.2micron filtered.

### Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>. Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

[Order on-line](#) or [Contact](#) your local distributor

**Disclaimer :** Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product.