


FT-EQH960



dNTP Set

Product Description

Name :	dNTP Set, 100 mM
Catalog Number :	FP-EQH960, 4 x 250 µl
Component :	dATP, 100 mM, FP-EQH96a, 250 µl dCTP, 100 mM, FP-EQH96b, 250 µl dGTP, 100 mM, FP-EQH96c, 250 µl dTTP, 100 mM, FP-EQH96ad 250 µl

Storage: -20°C. Product is stable for at least 1 year from date of receipt when stored as recommended.

Highly purified dNTPs (>99% by HPLC) in purified water. Each deoxynucleotide solution (dATP, dCTP, dGTP, dTTP) is provided in a separate vial at 100 mM. The dNTPs are suitable for use in PCR, cDNA synthesis, sequencing, and other molecular biology applications

Technical and scientific information

Related products

- Cheetah HoStart Taq DNA Polymerase, DU0570
- Uptima UptiTherm DNA Polymerase, UPS53921
- EvaGreen Dye 20X Water, BI179T
- Forget-Me-Not qPCR Master Mix, 1K9000
- Fast EvaGreen qPCR MasterMix, DV722E
- HotStart Polymerase Modification Kit, IWU520
- FluoProbes 490 - dUTP, FP-JO2800
- FluoProbes 594 - dUTP, FP-1H0840
- FluoProbes 647 - dUTP, FP-1H0860
- Biotin-11-dUTP, 1mM, FP-AM539A
- GelRed Nucleic Acid Gel Stain, 10000X in water, BY1740
- GelGreen Nucleic Acid Gel Stain, 10000X in water, CJ2730
- DNA Gel Extraction Kit, 1B9354
- 100bp DNA Ladder, S54811
- 1kb DNA Ladder, S54801
- TBE Buffer 10X, 86510A

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.