

HIGH LOAD STREPTAVIDIN-POLYHRP CONJUGATE (85R-202)

High Load Streptavidin Poly-HRP Conjugate for use in immunoassays

OVERVIEW

Synonyms	Streptavidin-PolyHRP, Streptavidin-Poly-HRP, Streptavidin Poly HRP, Streptavidin-PolyHRP Conjugate (High Load)
Applications	ELISA

SPECIFICATIONS

Tag/Conjugate	Poly-HRP
Form & Buffer	Supplied in universal IHC Diluent/Blocker/Stabilizer buffer (Cat No: 85R-1021) + 50% glycerol

USAGE & ASSAY INFORMATION

Usage Recommendations	Optimal conditions to be determined by the end user
------------------------------	---

STORAGE & SAFETY

Storage	Store at 4 deg C for short term storage or aliquot and store at -20 deg C for longer term. Do not freeze below -25 deg C
----------------	--

GENERAL INFORMATION

Biological Significance	Streptavidin Poly-HRP Conjugate is streptavidin biotin-binding protein that is conjugated with polymers of horseradish peroxidase, enabling signal amplification and detection of biotinylated antibodies.
--------------------------------	--